

Drawing Index

These sheets are a document set and should not be separated. Electrical information and references are contained on all sheets.

SITE READINESS

C1

EQUIPMENT LAYOUT

A1

(Equipment locations, heat loads, component weights, environmental specs)

STRUCTURAL LAYOUT

S1

(Structural support/mounting locations for floor/wall/ceiling, wall support elevations)

STRUCTURAL DETAILS

S2 THRU S3

(Floor and Ceiling loading information)

ELECTRICAL LAYOUT

E1

(Contractor supplied wiring, interconnect methods, junction point locations and descriptions)

ELECTRICAL SPECIFICATIONS

E2

(Maximum wiring run lengths, interconnect diagram, system power specifications)

ELECTRICAL DETAILS

E3 THRU E4

EQUIPMENT DETAILS

D1 THRU D3

These drawings indicate the placement and interconnection of the listed equipment components. These drawings are not construction or site preparation drawings. Customer remains ultimately responsible for preparing the site to accommodate the operation of such equipment in compliance with GE Healthcare's written specifications and all applicable federal, state, and/or local requirements.

\* REQUIRED REFERENCE \*

Innova IGS BiPlane  
Pre Installation Manual  
5435414-1-1EN

A mandatory component of this drawing set is the GE Healthcare Pre Installation manual. Failure to reference the Pre Installation manual will result in incomplete documentation required for site design and preparation.

Pre Installation documents for GE Healthcare products can be accessed on the web at:

www.gehealthcare.com/siteplanning

GE Healthcare



Interventional  
Site Planning

CUSTOMER ACCEPTANCE



imagination at work

Customer Site Readiness  
Requirements

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare Installation Project Manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare Installation Project Manager can supply a reference list of rigging contractors.
- New construction requires the following; 1. Secure area for equipment, 2. Power for drills and other test equipment, 3. Capability for image analysis, 4. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- Contact a radiation physicist or consultant to specify radiation containment requirements.

GE Equipment Delivery  
Requirements

The items on the GE Healthcare Site Readiness Checklist are REQUIRED to facilitate equipment delivery to the IS site. Equipment will not be delivered if these requirements are not satisfied.

GE Healthcare Site Readiness Checklist Rev 19

Before using this document ensure you have the latest Rev from MyWorkshop on DOC0422752

GEHC Global Order #: \_\_\_\_\_ Customer: \_\_\_\_\_  
GEHC PMI: \_\_\_\_\_ FE / Installer: \_\_\_\_\_

The customer is responsible for proper site preparation regardless of any GEHC measurements/inspections/assessments.

	Storage is lean ready?	PMI is lean ready?	FE is lean ready?	Comments if "N", enter comments or action plan
<b>GEHC Minimum Requirements</b>				
1 <b>MR Magnet Delivery Requirements:</b> Ensure cryogen venting system is available for magnet connection as defined by GEHC Pre-Installation Manual (PIM) requirements, exhaust fan system is installed and operational, 480V power, and chilled water supply is available 24x7 that meets system cooling requirements. External connectivity is available for magnet monitoring and phone service is available during delivery. Surface mount vibratmat installed where required. Magnet room final flooring is in place.				
2 <b>MR RF Screen Room Requirements:</b> RF Screen Room is tested with copy of Test Report, emailed to <a href="mailto:5435414@ge.com">5435414@ge.com</a> , that it is compliant with GEHC specifications. Dock Bolt and magnet anchors ( if applicable) installed using 2 part anchor. For HDx systems, blower box mount bolts installed by RF vendor using 2 part anchors				
3 <b>State Regulatory Requirements:</b> Facility registration number provided for states of <u>IL, KY, HI, RI, SC, TX</u> . X-ray shielding plan and state acknowledgment letter provided to installer for <u>AR, DC, NC, SC, CO</u> is <u>WA</u> . <b>Site Drawing Requirements:</b> Final version of equipment network and antenna, installation drawings (including red lined versions) verified to match actual room and has been provided to installer.				
4 <b>Surface Penetration Requirements:</b> Customer/Contractor scheduled to provide required drilling or cutting into floors, ceilings, and walls; OR surface penetration permit available and posted in the room when GEHC will perform the work.				
5 <b>Pre-Delivery Route Requirements:</b> The equipment delivery route from the truck to the final destination within the facility has been reviewed with all key stakeholders to safely meet the minimum requirements for equipment access, and all communications/notifications have occurred. Arrangements have been made for special handling (elevator, rigging, floor protection, fork lift, rollback truck, etc).				
6 <b>Finished Room Requirements:</b> Rooms that will contain equipment, including storage areas not in scan suite, are dust free. Provisions taken to maintain a dust free room. Precautions must be taken to prevent dust from entering rooms containing equipment when construction is incomplete in adjacent areas. All walls primed (final coat not needed on Day 1). Shielding, doors, and windows are to be installed. No contractor work being done during or after the installation that will cause dust in the installation areas or potential equipment damage. Room security to prevent unauthorized access and theft has been discussed with customer. The customer is aware of these security issues, implications and responsibility. For Storage: Room must meet PIM requirements for storage.				
7 <b>Electrical Requirements:</b> Lockable (LOTO) Main Disconnect Panel (MDP) is installed per GE guidelines and system power is available. Conduits, electrical cable ducting/dividers/cable trays, and access flooring is installed in proper location and height. Surface floor duct and load-side wires can be installed at time of system installation. Validate outlet location and requirements meet specifications for device/equipment.				
8 <b>HVAC Requirements:</b> The HVAC/Chilled Water systems designed to maintain the environment per spec/PIM is at running state and appears to provide the desired environmental conditions including location of vents, temperature and humidity for system operation.				
9 <b>Flooring Requirements:</b> Floor is clean and prepared for final floor covering. Floor levelness/flatness is measured and within tolerance, and there are no visible defects per GEHC specifications. Confirm customer anchoring plan aligns with designed floor thickness. Final flooring installed where required for network racks.				
10 <b>Ceiling Requirements:</b> Unistrut (or equivalent) location, levelness and spacing is measured (or vendor confirmed) and consistent with the requirement of the installation drawings. Ensure unistrut and rails are not used as mounting surfaces. Ceiling grid is installed. Permanent lighting is installed and operational. HVAC diffusers are installed and connected to ductwork. Ceiling tiles installed per PMI discretion.				
11 <b>Staging Requirements:</b> Space has been identified to support the active installation process only. This area meets PIM/project book requirements. Storage space has been identified, if needed. This secured space would be used to store equipment indefinitely. If offsite, transportation plan has been developed at customer expense. This space must meet PIM requirements.				
12 <b>Network Connectivity:</b> Hardware for network connectivity/network drop is in place prior to delivery with specified network firewall configuration where required. Site Surveys for wireless mobile XR units have been completed.				
13 <b>Medical Gases Requirements:</b> Systems (hard piped or portable) in place to allow testing and calibration of equipment (anesthesia, including ventilation).				

GE Healthcare

Healthcare Project Implementation – Design Center  
Milwaukee, Wisconsin  
Copyright 2009 General Electric Company – Proprietary to GE

SHEET TITLE: SITE READINESS

MODALITY TYPE: INNOVA IGS 630 BIPLANE

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS, IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST GEHC PIM. THE USER OF THIS PLAN SHALL BE RESPONSIBLE FOR ANY ACTUAL CONSTRUCTION ERRORS. GE HEALTHCARE SHALL NOT BE HELD RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE/ROOM: IR BP 1Z107

JL MCCLELLAN  
MEMORIAL VA HOSPITAL  
LITTLE ROCK, ARKANSAS

PROJECT	REVISION
142509	01
DATE:	21.Jul.14
DRAWN BY:	LLM
CHECKED BY:	LLM
GON NO:	4222033
GON DT:	08.Aug.14

REVISION HISTORY:

SHEET

C1

This drawing is based on Sketch No.: FloorPlan-X-FP1

PIM R2

RQ – 145731



GE EQUIPMENT LISTING							EQUIPMENT CROSS REFERENCE CHART	
EQUIPMENT ON ORDER FROM GE HEALTHCARE, INSTALLED BY GE HEALTHCARE, PER GON 4222033 DATED 08.Aug.14							P = PREAPPROVAL	
NOTE: LOCAL CONDITIONS MAY DICTATE THAT ITEMS IDENTIFIED IN THIS CATEGORY BE INSTALLED BY OTHERS.							C = CALCULATIONS/ PENDING APPROVAL	
							S = SPECIFICATIONS ONLY	
ITEM NO.	QUANTITY ORDERED	REFER TO SHEET "D"						
		ITEM DESCRIPTION (* = EXISTING/REINSTALL)	WEIGHT	HEAT OUTPUT (PER HOUR)	DETAIL NO.	STRC PLAN	ELEC PLAN	
(1)	2	LONGITUDINAL STATIONARY RAIL FOR XT SUSPENSION	68 lbs			B20 078		C
(2)	1	AW WORKSTATION	81 lbs	1201 btu	M1013AW C7619D B5030	---		C
(3)	1	MEDRAD MARK V INJECTOR ON PEDESTAL	90 lbs	320 btu				C
(4)	1	INNOVA IQ TABLE	1750 lbs	614 btu	B8162	B50	LW5	C
(5)	1	COUNTERBALANCED EYE AND THYROID SHIELD WITH R96 LAMP (TRACK NOT ON ORDER)	143 lbs		B5031E	49N 31F	LMP	S
(6)	1	LATERAL DETECTOR CHILLER	33 lbs	709 btu	B5150A	-	DC	-
(7)	1	AP DETECTOR CHILLER	33 lbs	709 btu	B5150A	-	DC	-
(8)	1	LATERAL COOLIX 4100 WATER CHILLER	264 lbs	11737 btu	B-IGS03 B-IGS04	-	CHLR	C
(9)	1	LATERAL COOLIX 4100 AUTOTRANSFORMER	66 lbs	153 btu	B-IGS05 B-IGS04	-	AT	-
(10)	1	AP COOLIX 4100 WATER CHILLER	264 lbs	11737 btu	B-IGS03 B-IGS04	-	CHLR	C
(11)	1	AP COOLIX 4100 AUTOTRANSFORMER	66 lbs	153 btu	B-IGS05	-	AT	-
(12)	1	LARGE DISPLAY MONITOR ON SINGLE MONITOR SUSPENSION 7 ft. x 9 in. INBOARD BRIDGE (MOUNT TWO GE MONITORS ON BACKSIDE OF LP MONITOR)	784 lbs	1706 btu	B2004 B2015	-	LDM WBM1	
(13)	1	LARGE DISPLAY MONITOR CABINET	253 lbs	3412 btu	B2014	-	LDC	C
(14)	1	3 KVA UPS CABINET (LARGE DISPLAY SUBSYSTEM OPTION)	99 lbs	546 btu	B2016	-	UPS3	C
(15)	1	XR BUZZER (LOCATED ABOVE CEILING)	2 lbs		B5150H	-	XR8	-
(16)	1	OPERATORS CONSOLE	22 lbs	546 btu	C7502 B5050D C7619D	-	WBC1	C
(17)	1	INNOVA POSITIONER (REFERENCE TABLE BASE-PLATE DETAIL FOR FLOOR MOUNTING INFORMATION)	1653 lbs	2416 btu	B5150D B5150E B5150F B5150G B5050E B5050F B5050G B5050H B5050N	---	LC1	C
(18)	1	UPS CABINET	1170 lbs	4061 btu	E4502SC	-	UPS	-
(19)	1	3 KVA UPS CABINET	81 lbs	546 btu		-	UPS1	-
(20)	1	LC/LP CABINET (C2)	621 lbs	4570 btu		-	C2	-
(21)	1	AP FRONTAL CABINET (C1)	888 lbs	4413 btu	B0558B	.	C1	-
(22)	1	LATERAL CABINET (C3)	703 lbs	2945 btu		-	C3	-
(23)	1	LATERAL POSITIONER BRIDGE MOUNT ASSEMBLY MOUNTED FROM CEILING SUPPORTS	1421 lbs	4126 btu	B5150B B5150C B5050K B5050L B5050M B5050N	-	LP4	C
(24)	2	LONGITUDINAL STATIONARY RAIL FOR LATERAL GENTRY INNOVA POSITIONER	68 lbs			B20 083		C

6	VITALING SPEAKER							
51	1 CABLE DRAPE RAIL FOR LP POSITIONER						B20 043	-
52	1 VITALING CONSOLE						-	-
53	1 VITALING MICROPHONE <ONE ON MONITOR BRIDGE IN EXAM ROOM>					B0573	-	-
54	1 VITALING MICROPHONE <ONE ON COUNTERTOP IN CONTROL ROOM>						-	-
55	1 INNOVA MAIN DISCONNECT. REFERENCE JUNCTION POINT 'A' ON SHEET E1 FOR DETAILED DESCRIPTION.	899 lbs	2215 btu	E4502B	-	PDB		

8'-10"

26'-5"

13'-0"

4'-2 5/8"

2'-10"

10'-0"

1'-5"

13'-10"

17'-10"

22'-6"

20'-4"

13'-3"

10'-11"

CL. TABLE PIVOT

ISO CENTER

LOCATED ELSEWHERE

70

18

IR BP 1Z107

EQUIP

CORRIDOR

CONTROL

SCRUB

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

40

41

42

43

44

45

46

47

48

49

50

51

52

53

54

55

56

57

58

59

60

61

62

63

64

65

66

67

68

69

70

71

72

73

74

75

76

77

78

79

80

81

82

83

84

85

86

87

88

89


90

UNIT CL


THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

## TYPICAL WALL SUPPORT ELEVATIONS

S120




S100



SCALE: 1/4" = 1'-0"

## STRUCTURAL LAYOUT

REQUIRED CEILING HEIGHT = 9'-4" = /-0.2"



NO CEILING MOUNTED ITEMS SUCH AS LIGHTS, SPRINKLER HEADS, EXHAUST FANS, ECT. CAN BE PLACED BETWEEN LP POSITIONER UNISTRUT (HATCHED AREA)

## STRUCTURAL SUPPORT METHODS

## CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
1	>>COMPONENTS BELOW CEILING<< CABLE DRAPE RAIL, UNISTRUT CAT. NO. CPGESS OR EQUIVALENT, * TO ORDER, CALL UNISTRUT WISCONSIN AT 262-796-8710.
2	UNISTRUT OR EQUIVALENT SUPPORT IN CEILING FOR FASTENING CEILING SUPPORTED EQUIPMENT. SUPPORTS TO RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL. RUN WALL TO WALL, BE PARALLEL, SQUARE, AND IN THE SAME HORIZONTAL PLANE, FLUSH WITH THE FINISHED CEILING. RAILS ARE MOUNTED TO THESE SUPPORTS EVERY 8'-2" AND REQUIRE 430 LBS. (597 LBS. IN SEISMIC REGIONS) PER BOLT LOAD. METHODS OF SUPPORT THAT PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION.
3	AREA RADIATION SHIELD TRACK MOUNTED TO GRIDDED CEILING UNISTRUT
4	HATCHED AREA INDICATES MONITOR BRIDGE BEARING BLOCK PATH.
5	STATIONARY RAILS ATTACHED TO UNISTRUT GRID IN CEILING.
6	HATCHED AREA INDICATES LP POSITIONER BEARING BLOCK PATH.
7	MOUNT XR BUZZER BRACKET ON WALL, ABOVE CEILING
8	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S120, FOR MAIN DISCONNECT CONTROL.
9	AREA OCCUPIED BY GE SUPPLIED TABLE BASE
10	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S100, FOR ATLAS CABINET.
11	AREA OCCUPIED BY GE SUPPLIED POSITIONER BASEPLATE

## STRUCTURAL NOTES

- ALL STEEL WORK AND PARTS NECESSARY TO SUPPORT CEILING MOUNTED TUBE HANGER OR OTHER EQUIPMENT ARE TO BE SUPPLIED BY THE CUSTOMER OR HIS CONTRACTORS. THE UNISTRUT OR EQUIVALENT STRUCTURE SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, RUN WALL TO WALL, BE PARALLEL, SQUARE AND IN THE SAME HORIZONTAL PLANE FLUSH WITH FINISHED CEILING. THE SYSTEM IS TO BE CROSS BRACED VERTICALLY, HORIZONTALLY AND DIAGONALLY TO ALLOW NO MOVEMENT AND A MAXIMUM OF 1,58mm (1/16") DEFLECTION. CLOSURE STRIPS SHALL BE PROVIDED FOR AREAS OF UNISTRUT EXPOSED AND WITHOUT MOUNTING UNITS.
- METHODS OF SUPPORT FOR THE STEELWORK THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE CONCRETE OR MASONRY ANCHORS IN DIRECT TENSION.
- ALL UNITS THAT ARE WALL MOUNTED OR WALL SUPPORTED ARE TO BE PROVIDED WITH SUPPORTS WHERE NECESSARY. WALL SUPPORTS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS. SEE PLAN AND DETAIL SHEETS FOR SUGGESTED LOCATIONS AND MOUNTING HOLE LOCATIONS.
- ALL CEILING MOUNTED FIXTURES, AIR VENTS, SPRINKLERS, ETC. TO BE FLUSH MOUNTED, OR SHALL NOT EXTEND MORE THAN 6,35mm (1/4") BELOW THE FINISHED CEILING.
- CONTROL WALLS WITH TUBE HANGER PASSAGE ABOVE SHALL BE CONSTRUCTED TO 2130mm (7'-0") HIGH.
- FLOOR SLABS ON WHICH EQUIPMENT IS TO BE INSTALLED MUST BE LEVEL TO 3,17mm (1/8") in 3050mm (10'-0")
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.
- CUSTOMERS CONTRACTOR MUST PROVIDE ALL PENETRATIONS IN POST TENSION FLOORS.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL ANY NON-STANDARD ANCHORING. DOCUMENTS FOR STANDARD ANCHORING METHODS ARE INCLUDED WITH GE EQUIPMENT DRAWINGS FOR GEOGRAPHIC AREAS THAT REQUIRE SUCH DOCUMENTATION.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL HARDWARE FOR "THROUGH THE FLOOR" ANCHORING AND/OR ANY BRACING UNDER ACCESS FLOORS. THIS CONTRACTOR MUST ALSO PROVIDE FLOOR DRILLING THAT CANNOT BE COMPLETED BECAUSE OF AN OBSTRUCTION ENCOUNTERED WHILE DRILLING BY THE GE INSTALLER SUCH AS REBAR ETC.
- IT IS THE CUSTOMER'S RESPONSIBILITY TO PERFORM ANY FLOOR OR WALL PENETRATIONS THAT MAY BE REQUIRED. THE CUSTOMER IS ALSO RESPONSIBLE FOR ENSURING THAT NO SUBSURFACE UTILITIES (E.G., ELECTRICAL OR ANY OTHER FORM OF WIRING, CONDUITS, PIPING, DUCT WORK OR STRUCTURAL SUPPORTS (I.E. POST TENSION CABLES OR REBAR)) WILL INTERFERE OR COME IN CONTACT WITH SUBSURFACE PENETRATION OPERATIONS (E.G. DRILLING AND INSTALLATION OF ANCHORS/SCREWS) PERFORMED DURING THE INSTALLATION PROCESS. TO ENSURE WORKER SAFETY, GE INSTALLERS WILL PERFORM SURFACE PENETRATION OPERATIONS ONLY AFTER THE CUSTOMER'S VALIDATION AND COMPLETION OF THE "GE SURFACE PENETRATION PERMIT"

GE Project Manager: JOHN COOPER  
Telephone: 913-221-9439

THE GE HPV TECHNICAL SUPPORT GROUP IS AN ADDITIONAL RESOURCE THAT CAN PROVIDE ANSWERS FOR GENERAL GE PRODUCT SITING QUESTIONS AND CAN BE REACHED AT (877)-305-9677 OR MAILTO:HPVTECHS@ge.com

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

GE Healthcare

Healthcare Project Implementation - Design Center  
Milwaukee, Wisconsin

SHEET TITLE: STRUCTURAL LAYOUT

MODALITY TYPE: INNOVA ICS 630 BIPLANE

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE ACTUAL CONSTRUCTION. HOWEVER, THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: ROOM: IR BP 1Z107

JL MCCLELLAN  
MEMORIAL VA HOSPITAL  
LITTLE ROCK, ARKANSAS

PROJECT	REVISION
142509	01
DATE:	21.Jul.14
DRAWN BY:	LLM
CHECKED BY:	LLM
GON NO:	4222033
GON DT:	08.Aug.14





REVISION HISTORY:

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SHEET


S1

472-353






**FLOOR MOUNTING : INNOVA 2100-3100-4100 (UNITY)/OMEGA V LONG TABLE (WITH IQ TILT TABLE BASEPLATE) INSTALLATION (TEMPLATE NO. 2360133)**

**B5049N**  
REV. DATE: 06/04/09



**WARNING!! THE RELATIONSHIP BETWEEN THE TABLE BASE AND THE POSITIONER BASEPLATE IS CRITICAL.**

PRIOR TO DRILLING MOUNTING HOLES CONTACT LOCAL GE HEALTHCARE INSTALLATION PROJECT MANAGER OR LEAD FIELD ENGINEER TO VERIFY THAT THE PROPER FULL SIZE FLOOR MOUNTING TEMPLATE IS USED.



**GE Healthcare**

Healthcare Project Implementation - Design Center  
Milwaukee, Wisconsin

**STRUCTURAL DETAILS**

SHEET TITLE: INNOVA ICS 630 BIPLANE

MODALITY TYPE: INNOVA ICS 630 BIPLANE

PROJECT TITLE: ROOM: IR BP 1Z107

JL MCCLELLAN  
MEMORIAL VA HOSPITAL  
LITTLE ROCK, ARKANSAS


PROJECT: 142509  
REVISION: 01

DATE: 21.Jul.14  
DRAWN BY: LLM  
CHECKED BY: LLM  
GON NO: 4222033  
GON DT: 08.Aug.14

REVISION HISTORY:

SHEET  
S2

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED






SCALE: 1/4" = 1'-0"

# ELECTRICAL PLAN

REQUIRED CEILING HEIGHT = 9'-4" = +/- 0.2"

## JUNCTION POINT DESCRIPTIONS



**CONDUIT RUNS: INNOVA IGS/ PLUS BIPLANE**

**CNDS. REQ'D. FOR BASE SYSTEM (LATERAL PLANE) (CONDUITS ARE LOCATED ABOVE CEILING)**

REV DATE: 10/30/08

(1)	LP4	TO	C3	FOUR 4" CNDS. USABLE CABLE LENGTH UP TO 42 FT.
(2)	WC2	TO	WC1	ONE EMPTY 3" CND. (FOR WATER LINES) USABLE CABLE LENGTH UP TO 68 FT.

NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS

**CNDS. REQ'D. FOR BASE SYSTEM (AP PLANE) (CONDUITS ARE LOCATED BELOW FLOOR)**

REV DATE: 10/30/08

(3)	LC1	TO	C1/C2	FOUR 4" CNDS. USABLE CABLE LENGTH UP TO 60 FT.
(4)	LC1	TO	LU5	ONE 4" & ONE 2" CND. USABLE CABLE LENGTH UP TO 13 FT.
(6)	WBC1	TO	C1/C2	ONE 3 1/2" & TWO 2 1/2" CNDS. USABLE CABLE LENGTH UP TO 60 FT.

NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS

**CONDUITS REQUIRED FROM POINT "XRLC" (CONDUITS ARE LOCATED ABOVE CEILING)**

REV DATE: 10/30/08

(7)	XRL1	TO	PDB	ONE 1/2" CND.
(9)	XRLC	TO	PDB	ONE 1/2" CND.
(10)	XRLC	TO	120-V 1Ø POWER	CND. AS REQ'D

**CONDUITS REQUIRED FOR REMOTE "20KVA UPS" (CONDUITS ABOVE CEILING OR BELOW FLOOR)**

REV DATE: 10/01/08

(29)	UPS	TO	UIB	ONE 2 1/2" CND. AND ONE 1" CND. USABLE CABLE LENGTH 70 FT.
(30)	UPS	TO	C1	ONE 2 1/2" CND. USABLE CABLE LENGTH 70 FT.
(31)	UIB	TO	C1	CABLES RUN IN DUCT USABLE CABLE LENGTH 15 FT.

NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS

**CONDUITS REQUIRED FROM POINT "PDB" (CONDUITS ABOVE CEILING OR BELOW FLOOR)**

REV DATE: 10/30/08

(18)	PDB	TO	UPS1	EXTERNALLY CONNECTED
(19)	PDB	TO	UPS	TWO CNDS. AS REQ'D. USABLE CABLE LENGTH 70 FT.
(20)	PDB	TO	RDS1	ONE 1/2" CND.
(21)	PDB	TO	RDS2	ONE 1/2" CND.
(22)	PDB	TO	C1	TWO 2 1/2" CNDS. FOR TWO CUSTOMER SUPPLIED POWER/ GROUND RUNS (JEDI/ CHLR) (AND GE SUPPLIED WIRES) CABLE LENGTH 19 FT.
(23)	PDB	TO	C1	ONE 1 1/2" CND. FOR TWO SUPPLIED SIGNAL CABLES CABLE LENGTH 19 FT.
(24)	PDB	TO	C1	ONE 1 1/2" CND. FOR 230-VGE SUPPLIED CABLES CABLE LENGTH 19 FT.
(25)	PDB	TO	C2	ONE CND. AS REQ'D. FOR ONE CUSTOMER SUPPLIED POWER/ GROUND RUN CABLE LENGTH 19 FT.
(26)	PDB	TO	C2	ONE 1 1/2" CND. FOR SIGNAL CABLES (RML1, XRL1, XRLC)
(27)	PDB	TO	C3	TWO 2 1/2" CNDS. FOR TWO CUSTOMER SUPPLIED POWER/ GROUND RUNS (AND GE SUPPLIED WIRES) CABLE LENGTH 19 FT.
(28)	PDB	TO	LU5	(TABLE POWER) RUN IN DUCT/ CONDUIT SYSTEM (IF CANNOT RUN IN CND./ DUCT SYSTEM, THEN RUN ONE ADDITIONAL 2" CND.)
(29)	PDB	TO	PDB1	CONDUIT AS REQUIRED
(30)	PDB1	TO	480-V 3Ø POWER	CONDUIT AS REQUIRED
(31)	PDB	TO	IE	(INJECTOR POWER) CONSULT MFG. (RUN IN DUCT/ CONDUIT SYSTEM)

NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS



**JUNCTION POINT NOTES**

- ALL JUNCTION BOXES, CONDUIT, DUCT, DUCT DIVIDERS, SWITCHES, CIRCUIT BREAKERS, ETC., ARE TO BE SUPPLIED AND INSTALLED BY CUSTOMERS ELECTRICAL CONTRACTOR.
- CONDUIT AND DUCT RUNS SHALL HAVE SWEEP RADIUS BENDS
- CONDUITS AND DUCT ABOVE CEILING OR BELOW FINISHED FLOOR MUST BE INSTALLED AS NEAR TO CEILING OR FLOOR AS POSSIBLE TO REDUCE RUN LENGTH.
- CEILING MOUNTED JUNCTION BOXES ILLUSTRATED ON THIS PLAN MUST BE INSTALLED FLUSH WITH FINISHED CEILING.
- ALL DUCTWORK MUST MEET THE FOLLOWING REQUIREMENTS:
  - DUCTWORK SHALL BE METAL WITH DIVIDERS AND HAVE REMOVABLE, ACCESSIBLE COVERS.
  - DUCTWORK SHALL BE CERTIFIED/RATED FOR ELECTRICAL POWER PURPOSES.
  - DUCTWORK SHALL BE ELECTRICALLY AND MECHANICALLY BONDED TOGETHER IN AN APPROVED MANNER.
  - PVC AS A SUBSTITUTE MUST BE USED IN ACCORDANCE WITH ALL LOCAL AND NATIONAL CODES.
- ALL OPENINGS IN ACCESS FLOORING ARE TO BE CUT OUT AND FINISHED OFF WITH GROMMET MATERIAL BY THE CUSTOMERS CONTRACTOR.
- GENERAL CONTRACTOR TO INSERT PULL CORDS FOR ALL CABLE RUN CONDUITS BETWEEN THE EQUIPMENT ROOM AND THE OPERATORS CONTROL ROOM.
- 10 FOOT PIGTAILS AT ALL JUNCTION POINTS.
- ALL WIRING MUST BE THIN OR TFFN STRANDED COPPER THERMOPLASTIC 600 VOLT OR EQUIVALENT INSULATION. ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.
- GROUNDING IS CRITICAL TO EQUIPMENT FUNCTION AND PATIENT SAFETY. SITE MUST CONFORM TO WIRING SPECIFICATIONS SHOWN ON THIS PLAN.

**ELECTRICAL OUTLET LEGEND**

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS. HEIGHT ABOVE FLOOR DETERMINED BY LOCAL CODES UNLESS OTHERWISE SPECIFIED.

⊕	DUPLEX HOSPITAL GRADE, DEDICATED WALL OUTLET 120-V, SINGLE PHASE POWER
⊕	DEDICATED TELEPHONE LINES (SEE ELECTRICAL DETAIL ELEC-1 OR ELEC-67)
⊕	NETWORK OUTLET (SEE ELECTRICAL DETAILS ELEC-83 AND ELEC-84 OR ELEC-87)
⊕	5-15R NEMA RECEPTACLE, DEDICATED OUTLET 120-V, SINGLE PHASE POWER
⊕	DUPLEX HOSPITAL GRADE, DEDICATED OUTLET 120-V EMERGENCY, SINGLE PHASE POWER, 15A
⊕	NEMA 15-30R RECEPTACLE, DEDICATED OUTLET 120-V, SINGLE PHASE POWER
⊕	6-GANG HOSPITAL GRADE, DEDICATED WALL OUTLET 115-V, SINGLE PHASE POWER



**CONDUITS REQUIRED FROM POINT "XRLC" (CONDUITS ARE LOCATED ABOVE CEILING)**

REV DATE: 10/30/08

(7)	XRL1	TO	PDB	ONE 1/2" CND.
(9)	XRLC	TO	PDB	ONE 1/2" CND.
(10)	XRLC	TO	120-V 1Ø POWER	CND. AS REQ'D

**CONDUITS REQUIRED FOR REMOTE "20KVA UPS" (CONDUITS ABOVE CEILING OR BELOW FLOOR)**

REV DATE: 10/01/08

(29)	UPS	TO	UIB	ONE 2 1/2" CND. AND ONE 1" CND. USABLE CABLE LENGTH 70 FT.
(30)	UPS	TO	C1	ONE 2 1/2" CND. USABLE CABLE LENGTH 70 FT.
(31)	UIB	TO	C1	CABLES RUN IN DUCT USABLE CABLE LENGTH 15 FT.



NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS

**CONDUITS REQUIRED FROM POINT "PDB" (CONDUITS ABOVE CEILING OR BELOW FLOOR)**

REV DATE: 10/30/08

(18)	PDB	TO	UPS1	EXTERNALLY CONNECTED
(19)	PDB	TO	UPS	TWO CNDS. AS REQ'D. USABLE CABLE LENGTH 70 FT.
(20)	PDB	TO	RDS1	ONE 1/2" CND.
(21)	PDB	TO	RDS2	ONE 1/2" CND.
(22)	PDB	TO	C1	TWO 2 1/2" CNDS. FOR TWO CUSTOMER SUPPLIED POWER/ GROUND RUNS (JEDI/ CHLR) (AND GE SUPPLIED WIRES) CABLE LENGTH 19 FT.
(23)	PDB	TO	C1	ONE 1 1/2" CND. FOR TWO SUPPLIED SIGNAL CABLES CABLE LENGTH 19 FT.
(24)	PDB	TO	C1	ONE 1 1/2" CND. FOR 230-VGE SUPPLIED CABLES CABLE LENGTH 19 FT.
(25)	PDB	TO	C2	ONE CND. AS REQ'D. FOR ONE CUSTOMER SUPPLIED POWER/ GROUND RUN CABLE LENGTH 19 FT.
(26)	PDB	TO	C2	ONE 1 1/2" CND. FOR SIGNAL CABLES (RML1, XRL1, XRLC)
(27)	PDB	TO	C3	TWO 2 1/2" CNDS. FOR TWO CUSTOMER SUPPLIED POWER/ GROUND RUNS (AND GE SUPPLIED WIRES) CABLE LENGTH 19 FT.
(28)	PDB	TO	LU5	(TABLE POWER) RUN IN DUCT/ CONDUIT SYSTEM (IF CANNOT RUN IN CND./ DUCT SYSTEM, THEN RUN ONE ADDITIONAL 2" CND.)
(29)	PDB	TO	PDB1	CONDUIT AS REQUIRED
(30)	PDB1	TO	480-V 3Ø POWER	CONDUIT AS REQUIRED
(31)	PDB	TO	IE	(INJECTOR POWER) CONSULT MFG. (RUN IN DUCT/ CONDUIT SYSTEM)

NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS



**CONDUITS REQUIRED FROM POINT "XRLC" (CONDUITS ARE LOCATED ABOVE CEILING)**

REV DATE: 10/30/08

(7)	XRL1	TO	PDB	ONE 1/2" CND.
(9)	XRLC	TO	PDB	ONE 1/2" CND.
(10)	XRLC	TO	120-V 1Ø POWER	CND. AS REQ'D

**CONDUITS REQUIRED FOR REMOTE "20KVA UPS" (CONDUITS ABOVE CEILING OR BELOW FLOOR)**

REV DATE: 10/01/08

(29)	UPS	TO	UIB	ONE 2 1/2" CND. AND ONE 1" CND. USABLE CABLE LENGTH 70 FT.
(30)	UPS	TO	C1	ONE 2 1/2" CND. USABLE CABLE LENGTH 70 FT.
(31)	UIB	TO	C1	CABLES RUN IN DUCT USABLE CABLE LENGTH 15 FT.



NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS

**CONDUITS REQUIRED FROM POINT "PDB" (CONDUITS ABOVE CEILING OR BELOW FLOOR)**

REV DATE: 10/30/08

(18)	PDB	TO	UPS1	EXTERNALLY CONNECTED
(19)	PDB	TO	UPS	TWO CNDS. AS REQ'D. USABLE CABLE LENGTH 70 FT.
(20)	PDB	TO	RDS1	ONE 1/2" CND.
(21)	PDB	TO	RDS2	ONE 1/2" CND.
(22)	PDB	TO	C1	TWO 2 1/2" CNDS. FOR TWO CUSTOMER SUPPLIED POWER/ GROUND RUNS (JEDI/ CHLR) (AND GE SUPPLIED WIRES) CABLE LENGTH 19 FT.
(23)	PDB	TO	C1	ONE 1 1/2" CND. FOR TWO SUPPLIED SIGNAL CABLES CABLE LENGTH 19 FT.
(24)	PDB	TO	C1	ONE 1 1/2" CND. FOR 230-VGE SUPPLIED CABLES CABLE LENGTH 19 FT.
(25)	PDB	TO	C2	ONE CND. AS REQ'D. FOR ONE CUSTOMER SUPPLIED POWER/ GROUND RUN CABLE LENGTH 19 FT.
(26)	PDB	TO	C2	ONE 1 1/2" CND. FOR SIGNAL CABLES (RML1, XRL1, XRLC)
(27)	PDB	TO	C3	TWO 2 1/2" CNDS. FOR TWO CUSTOMER SUPPLIED POWER/ GROUND RUNS (AND GE SUPPLIED WIRES) CABLE LENGTH 19 FT.
(28)	PDB	TO	LU5	(TABLE POWER) RUN IN DUCT/ CONDUIT SYSTEM (IF CANNOT RUN IN CND./ DUCT SYSTEM, THEN RUN ONE ADDITIONAL 2" CND.)
(29)	PDB	TO	PDB1	CONDUIT AS REQUIRED
(30)	PDB1	TO	480-V 3Ø POWER	CONDUIT AS REQUIRED
(31)	PDB	TO	IE	(INJECTOR POWER) CONSULT MFG. (RUN IN DUCT/ CONDUIT SYSTEM)

NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS



**CONDUITS REQUIRED FROM POINT "XRLC" (CONDUITS ARE LOCATED ABOVE CEILING)**

REV DATE: 10/30/08

(7)	XRL1	TO	PDB	ONE 1/2" CND.
(9)	XRLC	TO	PDB	ONE 1/2" CND.
(10)	XRLC	TO	120-V 1Ø POWER	CND. AS REQ'D

**CONDUITS REQUIRED FOR REMOTE "20KVA UPS" (CONDUITS ABOVE CEILING OR BELOW FLOOR)**

REV DATE: 10/01/08

(29)	UPS	TO	UIB	ONE 2 1/2" CND. AND ONE 1" CND. USABLE CABLE LENGTH 70 FT.
(30)	UPS	TO	C1	ONE 2 1/2" CND. USABLE CABLE LENGTH 70 FT.
(31)	UIB	TO	C1	CABLES RUN IN DUCT USABLE CABLE LENGTH 15 FT.


NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS

**CONDUITS REQUIRED FROM POINT "PDB" (CONDUITS ABOVE CEILING OR BELOW FLOOR)**

REV DATE: 10/30/08

(18)	PDB	TO	UPS1	EXTERNALLY CONNECTED
(19)	PDB	TO	UPS	TWO CNDS. AS REQ'D. USABLE CABLE LENGTH 70 FT.
(20)	PDB	TO	RDS1	ONE 1/2" CND.
(21)	PDB	TO	RDS2	ONE 1/2" CND.
(22)	PDB	TO	C1	TWO 2 1/2" CNDS. FOR TWO CUSTOMER SUPPLIED POWER/ GROUND RUNS (JEDI/ CHLR) (AND GE SUPPLIED WIRES) CABLE LENGTH 19 FT.
(23)	PDB	TO	C1	ONE 1 1/2" CND. FOR TWO SUPPLIED SIGNAL CABLES CABLE LENGTH 19 FT.
(24)	PDB	TO	C1	ONE 1 1/2" CND. FOR 230-VGE SUPPLIED CABLES CABLE LENGTH 19 FT.
(25)	PDB	TO	C2	ONE CND. AS REQ'D. FOR ONE CUSTOMER SUPPLIED POWER/ GROUND RUN CABLE LENGTH 19 FT.
(26)	PDB	TO	C2	ONE 1 1/2" CND. FOR SIGNAL CABLES (RML1, XRL1, XRLC)
(27)	PDB	TO	C3	TWO 2 1/2" CNDS. FOR TWO CUSTOMER SUPPLIED POWER/ GROUND RUNS (AND GE SUPPLIED WIRES) CABLE LENGTH 19 FT.
(28)	PDB	TO	LU5	(TABLE POWER) RUN IN DUCT/ CONDUIT SYSTEM (IF CANNOT RUN IN CND./ DUCT SYSTEM, THEN RUN ONE ADDITIONAL 2" CND.)
(29)	PDB	TO	PDB1	CONDUIT AS REQUIRED
(30)	PDB1	TO	480-V 3Ø POWER	CONDUIT AS REQUIRED
(31)	PDB	TO	IE	(INJECTOR POWER) CONSULT MFG. (RUN IN DUCT/ CONDUIT SYSTEM)

NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS



POINT	DESCRIPTION	QTY.	HARDWARE	DETAIL NO., SHT. E3
AT	COOLIX 4100 AUTOTRANSFORMER	1	EXTERNALLY CONNECTED TO 'CHLR' (WATER CHILLER)	ELEC-5
AT	COOLIX 4100 AUTOTRANSFORMER	1	EXTERNALLY CONNECTED TO 'CHLR' (WATER CHILLER)	ELEC-6
C1	AP FRONTAL CABINET	1	32 IN. OF GROMMET MATERIAL FOR AN 8" X 8 IN. OPENING IN DUCT COVER	ELEC-5
C2	LC/LP CABINET	1	32 IN. OF GROMMET MATERIAL FOR AN 8" X 8 IN. OPENING IN DUCT COVER	ELEC-6
C3	LATERAL CABINET	1	32 IN. OF GROMMET MATERIAL FOR AN 8" X 8 IN. OPENING IN DUCT COVER	ELEC-5
CHLR	COOLIX 4100 WATER CHILLER	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-5
CHLR	COOLIX 4100 WATER CHILLER	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-6
DC	LATERAL DETECTOR CHILLER	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-5
DC	AP DETECTOR CHILLER	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-6
LC1	INNOVA LC	1	24 X 24 X 12 IN. BOX SUITABLE LENGTH OF 8 IN. DIA. THREADED CONDUIT OR PIPE DIA. LOCKNUTS 1 IN. DIA. LOCKNUT SUPPLIED FITTING 12 X 12 X 6 IN. BOX 6 IN. DIA. BUSHING	ELEC-100
LDC	LARGE DISPLAY CABINET	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-6
LDM	LARGE DISPLAY MONITOR	1	COVERPLATE 12 IN. DIA. CHASE NIPPLE 18 X 18 X 6 IN. FLUSH CEILING BOX 3/4 IN. DIA. CHASE NIPPLE	ELEC-8
LMP	SURGICAL LAMP	1	4 X 4 X 6 IN. BOX 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
LP4	LATERAL POSITIONER	1	COVERPLATE 12 IN. DIA. CHASE NIPPLE 18 X 18 X 6 IN. FLUSH CEILING BOX 3/4 IN. DIA. CHASE NIPPLE	ELEC-8
LU5	OMEGA TABLE	1	COVERPLATE 12 IN. DIA. CHASE NIPPLE 18 X 18 X 6 IN. FLUSH CEILING BOX 3/4 IN. DIA. CHASE NIPPLE	ELEC-48
PDB	MAIN DISCONNECT	1	150-AMP PANE INCLUDED IN ORDER	ELEC-143
PDB1	LOCAL SERVICE DISCONNECT	1	150-AMP LOCAL SERVICE DISCONNECT (CUSTOMER SUPPLIED)	ELEC-16
RDS1	EMERGENCY OFF	1	PROVIDE A SINGLE GANG, 2 1/8 IN. DEEP, FLUSH MTD. WALL BOX.	ELEC-16
RDS2	EMERGENCY OFF	1	PROVIDE A SINGLE GANG, 2 1/8 IN. DEEP, FLUSH MTD. WALL BOX.	ELEC-8
UPS	UPS CABINET	1	COVERPLATE 12 IN. DIA. CHASE NIPPLE 18 X 18 X 6 IN. FLUSH CEILING BOX 3/4 IN. DIA. CHASE NIPPLE	ELEC-8
UPS1	3 KVA UPS	1	EXTERNALLY CONNECTED TO PDB	ELEC-5
UPS3	3 KVA UPS (LD SUBSYSTEM)	1	EXTERNALLY CONNECTED TO LARGE DISPLAY CABINET - 'LDC'	ELEC-6
WBC1	OPERATORS CONSOLE	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-8
WBM1	TV MONITOR	2	1/2 IN. DIA. CHASE NIPPLE 18 X 18 X 6 IN. FLUSH CEILING BOX 3/4 IN. DIA. CHASE NIPPLE	ELEC-8
WC1	WATER CHILLER HOSE OUTLET	1	6 X 6 X 6 IN. FLUSH CEILING BOX 3/4 IN. DIA. CHASE NIPPLE	ELEC-9
WC	WATER CHILLER HOSE OUTLET	1	3 IN. CONDUIT STUBBED 2 IN. ABOVE FLOOR	ELEC-8
WC2	WATER CHILLER HOSE OUTLET	1	COVERPLATE 12 IN. DIA. CHASE NIPPLE 18 X 18 X 6 IN. FLUSH CEILING BOX 3/4 IN. DIA. CHASE NIPPLE	ELEC-8
XRB	XR BUZZER (LOCATED ABOVE CEILING)	1	SINGLE GANG BOX	ELEC-8
XRL1	WARNING LIGHT	1	COVERPLATE 12 IN. DIA. CHASE NIPPLE 18 X 18 X 6 IN. FLUSH CEILING BOX 3/4 IN. DIA. CHASE NIPPLE	ELEC-146
XRLC	WARNING LIGHT (ROOM LIGHT CONTROL OR EQUIVALENT MAX 24V CONTROLLER)	1	RW120-10-RL WARNING LIGHT CONTROL LIGHT CONTROL OR EQUIVALENT MAX 24V CONTROLLER.	


**CONTRACTOR SUPPLIED AND INSTALLED WIRING**

ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS.

WIRE RUN, FROM - TO	QUANTITY, WIRE SIZE/COLOR
<30> 3 PHASE > PDB1	3-BLACK, 1-WHITE, 1-GREEN (REFER TO FEEDER TABLE)
<29> PDB1 > PDB	3-BLACK, 1-WHITE, 1-GREEN (REFER TO FEEDER TABLE)
<22> PDB > C1 <JEDI>	3-ND. 1 BLACK, 1-ND. 1 GREEN
<27> PDB > C3 <JEDI>	3-ND. 1 BLACK, 1-ND. 1 GREEN
<26> PDB > C2	3-ND. 8 BLACK, 1-ND. 8 GREEN
<19> PDB > UPS	6-ND. 6 BLACK, 2-ND. 6 WHITE, 2-ND. 4 GREEN
<22> PDB > AT <AP>	3-ND. 10 BLACK, 1-ND. 10 GREEN
<22> PDB > AT <LAT>	3-ND. 10 BLACK, 1-ND. 10 GREEN
<20> PDB > RDS1	2-ND. 14 BLACK, 2-ND. 14 WHITE, 1-ND. 14 GREEN

INTERCONNECT DIAGRAM

TYPICAL VIEWS



EQUIPMENT DESCRIPTIONS

ITEM	DESCRIPTION	WEIGHT (lb)	HEAT DISSIPATION (btu)	DRAWING DESIGNATOR
1	- XR BUZZER	2		XRB
2	- ATLAS CABINET C2	630	4570	C2
3	- ATLAS CABINET C1	890	4413	C1
4	- ATLAS CABINET C3	705	2945	C3
5	- DETECTOR CONDITIONER (FRONTAL)	33	709	DC
6	- WATER CHILLER (FRONTAL)	447	18723	CHLR
7	- DETECTOR CONDITIONER (LATERAL)	33	709	DC
8	- WATER CHILLER (LATERAL)	447	16320	CHLR
9	- 20kva UPS CABINET	1170	4061	UPS
10	- 3kva UPS CABINET	81	546	UPS1
11	- TV CEILING SUSPENSION (8 MONITOR)	630	1638	WBM1
12	- INNOVA LC POSITIONER	1653	2416	LC1
13	- INNOVA LP POSITIONER	1679	4126	LP4
14	- OMEGA V LONG TABLE	1750	614	LU5
15	- VCM OPERATOR CONSOLE	22	546	WBC1
16	- ROOM LIGHTS			RML1
17	- XRAY WARNING LAMP			XRLL1
18	- XRAY WARNING LAMP CONTROLLER			XRLLC
19	- RDS1 PUSHBUTTON			RDS1
20	- RDS2 PUSHBUTTON			RDS2
21	- PDB MAIN DISCONNECT	899	2215	PDB
22	- LOTO DISCONNECT BREAKER			PDB1

OPTIONS

ITEM	DESCRIPTION	WEIGHT (lb)	HEAT DISSIPATION (btu)	DRAWING DESIGNATOR
23	- BOLUS CHASE HANDSWITCH	2		WBBC
24	- ADVANTAGE WINDOWS WORKSTATION	81	1201	AW
25	- IVUS VOLCANO CONSOLE	68	1631	IVUS
26	- IVUS VOLCANO COLOR PRINTER	X		CP
27	- INJECTOR HEAD	15		IH
28	- INJECTOR ELECTRONICS	37	320	IE
29	- REMOTE CONTROL FOR INJECTOR	4		IEC
30	- LAMP (RADIATION SHIELD TRACK)	143		LMP
31	- LD CABINET	254	3412	LDC
32	- LD MONITOR	784	1706	LDM
33	- MICRO PACE	X		MP
34	- MACH 3 TRANSFORMER	70	X	MST
35	- SKYTRON LAMP	50	341	SL
36	- PHYSIO. PRINTER	X	309	
37	- 150 KVA UPS	2160	31802	UPS
38	- UPS BATTERY CABINET	3529	X	
39	- MAIN BYPASS PANEL	350	X	MBP
40	- 3kva UPS CABINET	81	546	UPS3
41	- AUTO TRANSFORMER	99	239	AT

UPDATED: 12/21/10

INNOVA BIPLANE SYSTEMS

REV. DATE: 10/22/07

VOLTAGE PRIMARY SOURCE IS REQUIRED FOR ALL INSTALLATIONS.  
RANGE OF LINE VOLTAGES  
NOMINAL LINE VOLTAGE OF 360 TO 480, 3 PHASE, 50 OR 60 Hz

REQUIRED POWER SUPPLY: WYE DISTRIBUTION

MAXIMUM DAILY VOLTAGE VARIATION MUST FALL WITHIN ONE OF THE RANGES IN TABLE A.

TABLE A  
ALLOWABLE  
INPUT  
VOLTAGES/  
CURRENT  
DEMAND

NOMINAL VOLTAGE	NORMAL RANGE ±10 PERCENT	CURRENT (AMPS)	
		MAX MOMENTARY	CONTINUOUS
360	324-396	289	32
380	342-418	274	31
400	360-440	260	29
420	378-462	248	28
440	396-484	236	26
460	414-506	226	25
480	432-528	217	24

ALL CALCULATIONS BASED UPON NOMINAL VOLTAGE

NOTE LOW LINE CONDITIONS MAY INHIBIT SOME HIGH KVP TECHNIQUES.  
THE GENERATOR AUTOMATICALLY ESTABLISHES THESE INHIBITS  
BASED ON ACTUAL LINE CONDITIONS AND SYSTEM REGULATION.

PHASE--  
BALANCE.

PHASE-TO-PHASE VOLTAGES MUST BE WITHIN +2 PERCENT  
OF THE LOWEST PHASE-TO-PHASE VOLTAGE. MAXIMUM ALLOWABLE  
TRANSIENT VOLTAGE EXCURSIONS ARE 2.5 PERCENT OF RATED  
LINE VOLTAGE AT A MAXIMUM DURATION OF 5 CYCLES AND  
FREQUENCY OF 10 TIMES PER HOUR.

POWER  
DEMAND

CONTINUOUS POWER DEMAND = 20KVA. (MAX DEMAND = 171 KVA)

TABLE B  
MAXIMUM  
MOMENTARY  
POWER  
DEMAND.

DEMAND	INNOVA JEDI
kVa * POWER FACTOR AT	180 0.9
mA	1250
kVp	80

\* DEMAND INCLUDES POWER FOR ENTIRE ADVANTX SYSTEM.  
LINE VOLTAGE REGULATION AT MAXIMUM POWER DEMAND  
MUST BE LESS THAN OR EQUAL TO 6 PERCENT.

DISTRI-  
BUTION  
TRANS-  
FORMER FOR A SINGLE UNIT INSTALLATION, THE MINIMUM TRANSFORMER SIZE  
IS 225 KVA.

ELECTRICAL NOTES

- NOTE 1: ALL WIRES SPECIFIED SHALL BE COPPER STRANDED, FLEXIBLE, THERMO-PLASTIC, COLOR CODED, CUT 10 FOOT LONG AT OUTLET BOXES, DUCT TERMINATION POINTS OR STUBBED CONDUIT ENDS.  
ALL CONDUCTORS, POWER, SIGNAL AND GROUND, MUST BE RUN IN A CONDUIT OR DUCT SYSTEM. ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS. WIRE RUNS MUST BE CONTINUOUS COPPER STRANDED AND FREE FROM SPLICES. **ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.**
- NOTE 2: WIRE SIZES GIVEN ARE FOR USE OF EQUIPMENT. LARGER SIZES MAY BE REQUIRED BY LOCAL CODES.
- NOTE 3: IT IS RECOMMENDED THAT ALL WIRES BE COLOR CODED, AS REQUIRED IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 4: CONDUIT SIZES SHALL BE VERIFIED BY THE ARCHITECT, ELECTRICAL ENGINEER OR CONTRACTOR, IN ACCORDANCE WITH LOCAL OR NATIONAL CODES.
- NOTE 5: CONVENIENCE OUTLETS ARE NOT ILLUSTRATED. THEIR NUMBER AND LOCATION ARE TO BE SPECIFIED BY OTHERS. LOCATE AT LEAST ONE CONVENIENCE OUTLET CLOSE TO THE SYSTEM CONTROL, THE POWER DISTRIBUTION UNIT AND ONE ON EACH WALL OF THE PROCEDURE ROOM. USE HOSPITAL APPROVED OUTLET OR EQUIVALENT.
- NOTE 6: GENERAL ROOM ILLUMINATION IS NOT ILLUSTRATED. CAUTION SHOULD BE TAKEN TO AVOID EXCESSIVE HEAT FROM OVERHEAD SPOTLIGHTS. DAMAGE CAN OCCUR TO CEILING MOUNTING COMPONENTS AND WIRING IF HIGH WATTAGE BULBS ARE USED. RECOMMEND LOW WATTAGE BULBS NO HIGHER THAN 75 WATTS AND USE DIMMER CONTROLS (EXCEPT MR). DO NOT MOUNT LIGHTS DIRECTLY ABOVE AREAS WHERE CEILING MOUNTED ACCESSORIES WILL BE PARKED.
- NOTE 7: **ROUTING OF CABLE DUCTWORK, CONDUITS, ETC., MUST RUN DIRECT AS POSSIBLE OTHERWISE MAY RESULT IN THE NEED FOR GREATER THAN STANDARD CABLE LENGTHS (REFER TO THE INTERCONNECTION DIAGRAM FOR MAXIMUM USABLE LENGTHS POINT TO POINT).**
- NOTE 8: CONDUIT TURNS TO HAVE LARGE, SWEEPING BENDS WITH MINIMUM RADIUS IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 9: A SPECIAL GROUNDING SYSTEM IS REQUIRED IN ALL PROCEDURE ROOMS BY SOME NATIONAL AND LOCAL CODES. IT IS RECOMMENDED IN AREAS WHERE PATIENTS MIGHT BE EXAMINED OR TREATED UNDER PRESENT, FUTURE, OR EMERGENCY CONDITIONS. CONSULT THE GOVERNING ELECTRICAL CODE AND CONFER WITH APPROPRIATE CUSTOMER ADMINISTRATIVE PERSONNEL TO DETERMINE THE AREAS REQUIRING THIS TYPE OF GROUNDING SYSTEM.
- NOTE 10: THE MAXIMUM POINT TO POINT DISTANCES ILLUSTRATED ON THIS DRAWING MUST NOT BE EXCEEDED.
- NOTE 11: PHYSICAL CONNECTION OF PRIMARY POWER TO GE EQUIPMENT IS TO BE MADE BY CUSTOMERS ELECTRICAL CONTRACTOR WITH THE SUPERVISION OF A GE REPRESENTATIVE. THE GE REPRESENTATIVE WOULD BE REQUIRED TO IDENTIFY THE PHYSICAL CONNECTION LOCATION, AND INSURE PROPER HANDLING OF GE EQUIPMENT.
- NOTE 12: GEHC CONDUCTS POWER AUDITS TO VERIFY QUALITY OF POWER BEING DELIVERED TO THE SYSTEM. THE CUSTOMER'S ELECTRICAL CONTRACTOR IS REQUIRED TO BE AVAILABLE TO SUPPORT THIS ACTIVITY.

DIAGRAM KEY

--- CUSTOMER/CONTRACTOR SUPPLIED WIRING. ROUTE IN  
ADEQUATE CONDUIT OR RACEWAY.  
— GE FURNISHED CABLE RUNS. ROUTE IN EMPTY  
CONDUIT OR RACEWAY.

REV. DATE: 05.Mar.14

GE Healthcare

SHEET TITLE: ELECTRICAL SPECIFICATIONS  
MODALITY TYPE: INNOVA ICS 630 BIPLANE

PROJECT TITLE: ROOM: IR BP 1Z107  
JL MCCLELLAN  
MEMORIAL VA HOSPITAL  
LITTLE ROCK, ARKANSAS

PROJECT	REVISION
142509	01
DATE:	21.Jul.14
DRAWN BY:	LLM
CHECKED BY:	LLM
GON NO:	4222033
GON DT:	08.Aug.14

REVISION HISTORY:

SHEET

E2



ELECTRICAL DETAIL  
TABLE INTERCONNECT DETAIL, UNDER FLOOR

ELEC-134  
REV. DATE: 18.Jul.14

OMEGA TABLE  
4" [102mm] PIPE THROUGH  
FLOOR TO CABLE ACCESS

ONLY ONE PIPE THROUGH  
FLOOR REQUIRED

TILTING TABLE  
4" [102mm] PIPE THROUGH  
FLOOR TO CABLE ACCESS

FLUSH MOUNTED  
FLOOR PLATE

POSITIONER SIDE

(1) 2" [51mm] AND  
(1) 4" [102mm]  
CONDUIT FROM POSITIONER  
\*\*\* OR \*\*\*  
DUCTWORK AS SHOWN ON E1

6" x 6" x 16" BOX  
[152mm x 152mm x 406mm]

NOTE: PIPE, JUNCTION BOX AND DUCT or CONDUIT ARE TO BE SUPPLIED AND  
INSTALLED BY CUSTOMER or CUSTOMER'S CONTRACTOR.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL  
VERTICAL WALL DUCT (TYPICAL)

ELEC-6  
REV. DATE: 03/19/04

REFER TO CHART FOR MINIMUM DIVIDER REQUIREMENT  
LOCAL CODES MAY REQUIRE ADDITIONAL DIVIDERS

ELECTRICAL DUCT

DUCT WIDTH

EQUAL

EQUAL

REMOVABLE  
DUCT COVER

GROMMETED  
OPENING

RUBBER GROMMET

COVER PLATE  
TO BE REMOVABLE

ELECTRICAL CONTRACTOR TO FURNISH  
AND INSTALL SCREWS AS SHOWN

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" [610mm]	2
18" [457mm]	2
10" [254mm]	2
6" [152mm]	1
4" [102mm]	1

REMOVABLE SECTION  
OF WALL DUCT

REMOVABLE  
DUCT COVER

DETAIL NOT TO SCALE

ELECTRICAL DETAIL  
J.B. / WALL DUCT DETAIL (TYPICAL)

ELEC-2  
REV. DATE: 09/30/94

REMOVABLE COVER

JUNCTION BOX ABOVE CEILING

PARTITION

CNDS. ABOVE CEILING

REMOVABLE COVER

FINISHED CEILING

DETAIL NOT TO SCALE

ELECTRICAL DETAIL  
BOX WITH COVERPLATE AND NETWORK JACK

ELEC-83  
REV. DATE: 10/06/98

BOX

NETWORK JACK

COVERPLATE

DETAIL NOT TO SCALE

ELECTRICAL DETAIL  
NETWORK CONNECTION (TYPICAL)

ELEC-84  
REV. DATE: 03/06/04

LOCAL AREA NETWORK

FINISHED CEILING

1/2" CONDUIT FROM J.B.  
TO ABOVE FINISHED CEILING.

TO BE DETERMINED

FINISHED FLOOR

COVERPLATE WITH NETWORK  
RECEPTACLE

SINGLE GANG J.B.

FOR NUCLEAR SYSTEMS A DIRECT NETWORK  
CONNECTION IS TO BE MADE BETWEEN THE  
SYSTEM AND THE REVIEW WORKSTATION.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL  
INSITE CONNECTION (TYPICAL)

ELEC-1  
REV. DATE: 04/24/02

ONE OF THE FOLLOWING TWO SELECTIONS MUST BE INSTALLED AT THE LOCATION SHOWN ON THE  
ELECTRICAL PLAN (SHEET E1) FOR GE INSITE CONNECTION BASED UPON SYSTEM CONFIGURATION.

A) ONE INTERNET ACCESSIBLE VIRTUAL PRIVATE NETWORK (VPN) CONNECTION WITH A STATIC IP  
ADDRESS, AND ONE TELEPHONE LINE - DEDICATED-DIRECT-DIALING, VOICE GRADE.

OR

B) TWO TELEPHONE LINES - ONE DEDICATED DIRECT-DISTANCE-DIALING, VOICE GRADE AND  
ONE A DEDICATED DATA LINE.

FINISHED CEILING

1" CONDUIT FROM J.B.  
TO ABOVE FINISHED CEILING.

TO BE DETERMINED

FINISHED FLOOR

SINGLE GANG J.B.

COVERPLATE WITH TWO  
TELEPHONE RECEPTACLES  
OR  
ONE TELEPHONE RECEPTACLE AND  
ONE NETWORK RECEPTACLE

ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER OR THEIR CONTRACTOR.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL  
HORIZONTAL WALL DUCT (TYPICAL)

ELEC-5  
REV. DATE: 03/19/04

TYPICAL WALL DUCT

REMOVABLE  
DUCT COVER

FINISHED FLOOR

GROMMETED  
OPENING

REMOVABLE SECTION OF  
WALL DUCT COVER

REFER TO CHART FOR MINIMUM DIVIDER REQUIREMENT  
LOCAL CODES MAY REQUIRE ADDITIONAL DIVIDERS

ELECTRICAL DUCT

RUBBER GROMMET

COVER PLATE  
TO BE REMOVABLE

DUCT WIDTH

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" [610mm]	2
18" [457mm]	2
10" [254mm]	2
6" [152mm]	1
4" [102mm]	1

DETAIL NOT TO SCALE

ELECTRICAL DETAIL  
HORIZONTAL WALL DUCT (TYPICAL)

ELEC-5A  
REV. DATE: 06/16/08

TYPICAL WALL DUCT

REMOVABLE  
DUCT COVER

FINISHED FLOOR

REMOVABLE SECTION OF  
WALL DUCT COVER

REFER TO CHART FOR MINIMUM DIVIDER REQUIREMENT  
LOCAL CODES MAY REQUIRE ADDITIONAL DIVIDERS

ELECTRICAL DUCT

2" x 4" OPENING  
CUT INTO TOP OF DUCT  
FOR 12" OF GROMMETED  
MATERIAL

DUCT WIDTH

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" [610mm]	2
18" [457mm]	2
10" [254mm]	2
6" [152mm]	1
4" [102mm]	1

DETAIL NOT TO SCALE

ELECTRICAL DETAIL  
BOX WITH COVERPLATE (TYPICAL)

ELEC-8  
REV. DATE: 09/30/94

OUTLET BOX

HARDWARE

DETAIL NOT TO SCALE

ELECTRICAL DETAIL  
CONDUITS THRU-FLOOR (TYPICAL)

ELEC-9  
REV. DATE: 08/08/94

FINISHED FLOOR

HARDWARE

1.5" (38 mm) TYP.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL  
X-RAY WARNING LIGHT & ROOM LIGHT CONTROL PANEL

ELEC-146  
REV. DATE: 05.SEP.12

24 VAC POWERED  
FROM GE IMAGE SYSTEM  
ON SIGNAL IN PDB

MAXIMUM  
24 VAC  
48 WATTS

RW120-10-RL  
X-RAY ROOM WARNING  
LIGHT / ROOM LIGHTING  
CONTROL PANEL

INPUT TO PANEL &  
WARNING LIGHT POWER

120-VAC 20A MAX.

X-RAY WARNING LIGHT-2A MAX.

ROOM LIGHT POWER  
120 VAC, 277V, 480V, 600V

ROOM LIGHTS

FROM GE ROOM  
LIGHT SIGNAL IN PDB  
(DRY CONTACTS)  
OR  
REMOTE FOOT SWITCH  
WITH 100 FT. CABLE

MAXIMUM  
24 VAC  
48 WATTS

SUPPLIED BY  
THIS PANEL

X-RAY WARNING LIGHT  
OR ROOM LIGHT ARE NOT  
PART OF THIS CAT. NO.

THE RW120-10-RL IS RECOMMENDED IF 120 VAC "X-RAY ON" WARNING LIGHT AND  
ROOM LIGHT CONTROL ARE UTILIZED (FOR USE WITH BIPLANE SYSTEM)

THE E4500AK IS RECOMMENDED IF 120 VAC "X-RAY ON" WARNING LIGHT ONLY

24 VAC POWERED  
FROM GE IMAGE SYSTEM  
ON SIGNAL IN PDB

MAXIMUM  
24 VAC  
48 WATTS

E4500AK  
X-RAY ROOM WARNING  
LIGHT CONTROL PANEL

INPUT TO PANEL FOR  
"X-RAY IN USE" LIGHTS

120-VAC 20A MAX.

X-RAY WARNING LIGHT-2A MAX.

X-RAY IN USE WARNING LIGHT  
IS NOT PART OF THIS  
CAT. NO.

E4500AK NOTE:  
IF 24VAC X-RAY  
IN USE WARNING  
LIGHTS ARE UTILIZED,  
THIS RELAY PANEL  
IS NOT REQUIRED.

CONTROL PANEL CAN BE LOCATED ABOVE THE CEILING NEAR THE WARNING LIGHT  
UNLESS SPECIFIED ON SHEET A1 AS BEING INCLUDED ON EQUIPMENT ORDER,  
ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER'S CONTRACTOR

ELECTRICAL DETAIL  
EMERGENCY OFF BUTTON

ELEC-16  
REV. DATE: 05/14/09

PLAN VIEW

FRONT VIEW

OFF

BACK VIEW

SINGLE GANG BOX  
SUPPLIED BY CONTRACTOR

2 1/2" [64mm]

DETAIL NOT TO SCALE

GE Healthcare

Healthcare Project Implementation - Design Center  
Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL DETAILS

MODALITY TYPE: INNOVA ICS 630 BIPLANE

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT  
AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS  
IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO DETAILS  
AND SPECIFICATIONS OF THE MANUFACTURER. THE USER SHALL BE RESPONSIBLE FOR  
OBTAINING THE LATEST EDITIONS OF THE MANUFACTURER'S LITERATURE AND FOR  
ACTUAL CONSTRUCTION PURPOSES. THE COMPANY CANNOT ACCEPT  
RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: ROOM: IR BP 1Z107

JL MCCLELLAN  
MEMORIAL VA HOSPITAL  
LITTLE ROCK, ARKANSAS

PROJECT  
142509

REVISION  
01

DATE: 21.Jul.14

DRAWN BY: LLM

CHECKED BY: LLM

GON NO: 4222033

GON DT: 08.Aug.14

REVISION HISTORY:

SHEET  
E3




THIS drawing is based on Sketch No.: FloorPlan-X-FP1

PIM R2

RQ - 145731

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED





This drawing is based on Sketch No.: FloorPlan-X-FP1

PIM R2

RQ – 145731

PROJECT TITLE: ROOM: IR BP 1Z107  
JL MCCLELLAN  
MEMORIAL VA HOSPITAL  
LITTLE ROCK, ARKANSAS


PROJECT	REVISION
142509	01
DATE:	21.Jul.14
DRAWN BY:	LLM
CHECKED BY:	LLM
GON NO:	4222033
GON DT:	08.Aug.14

REVISION HISTORY:

SHEET  
E4

SHEET TITLE: ELECTRICAL DETAILS  
MODALITY TYPE: INNOVA ICS 630 BIPLANE













THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST EDITIONS OF THE NATIONAL ELECTRICAL CODE, THE NATIONAL FIRE PROTECTION ASSOCIATION (NFPA) CODES, AND THE COMPANY'S DESIGN STANDARDS. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.



GE Healthcare

Healthcare Project Implementation – Design Center

Milwaukee, Wisconsin



GE Healthcare

Healthcare Project Implementation - Design Center

Minneapolis, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS

MODALITY TYPE: INNOVA ICS 630 BIPLANE

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN. EVERY EFFORT HAS BEEN MADE TO CONFORM TO DETAILS OF THE MANUFACTURER'S LATEST CATALOGS. THE USER SHALL BE RESPONSIBLE FOR VERIFYING THE ACTUAL CONSTRUCTION DIMENSIONS. GE HEALTHCARE SHALL NOT BE HELD RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: ROOM: IR BP 1Z107

JL MCCLELLAN

MEMORIAL VA HOSPITAL

LITTLE ROCK, ARKANSAS

PROJECT	REVISION
142509	01
DATE:	21.Jul.14
DRAWN BY:	LLM
CHECKED BY:	LLM
GON NO:	4222033
GON DT:	08.Aug.14

REVISION HISTORY:

SHEET

D1




## REV. DATE: 15.MAR.12



DETAIL NOT TO SCALE

## REV. DATE: 17.Dec.12



DETAIL NOT TO SCALE

## REV. DATE: 12/07/09



## REV. 00: 10/30/08



NOTE: XR-BUZZER BRACKET IS MOUNTED ON WALL, ABOVE CEILING. PLACE SPEAKER ABOVE GRILLED CEILING TILE FOR SOUND PENETRATION.

## REV. DATE: 08/06/07




## REV. DATE: 10/25/10



## REV. DATE: 08/26/08



## REV. DATE: 26.JUL.12



DETAIL NOT TO SCALE

## REV. DATE: 05/10/05



## REV. DATE: 01/04/07



## REV. DATE: 06.SEP.12



## REV. DATE: 06.SEP.12



THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

This drawing is based on Sketch No.: FloorPlan-X-FP1

PIM R2

RQ - 145731

SHEET  
D2

SHEET TITLE: EQUIPMENT DETAILS  
MODALITY TYPE: INNOVA IGS 630 BIPLANE

PROJECT TITLE: ROOM: IR BP 1Z107  
JL MCCLELLAN  
MEMORIAL VA HOSPITAL  
LITTLE ROCK, ARKANSAS

PROJECT	REVISION
142509	01
DATE: 21.Jul.14	
DRAWN BY: LLM	
CHECKED BY: LLM	
GON NO: 4222033	
GON DT: 08.Aug.14	

REVISION HISTORY:

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

SHEET  
D2

**GE Healthcare**



Healthcare Project Implementation – Design Center  
Milwaukee, Wisconsin

