

**Utility P:osition 1**  
(Marked on floor P1)

**S-Lab I/O Power**  
1-230 VAC 60 Hz 1-phase, 15 Amp dedicated / Isolated circuit.  
Receptacle Nema #6-15R.  
**Buck / Boost Transformer must be used on circuit.**

**Centrifuge Power**  
1-230 VAC 60 Hz 1-phase, 15 Amp dedicated / Isolated circuit.  
Receptacle Nema #6-15R.  
**Buck / Boost Transformer must be used on circuit.**

**Accessory Power**  
1-120 VAC 60 Hz 1-phase, 15 Amp dedicated circuit  
Receptacle Duplex Nema #5-15R

**Also in this Position**  
RJ-45 Ethernet Connection LAB LIS  
Air Supply Line from Compressor 1/2" copper tubing run with 3/8" female threaded fitting, shut off valve, moisture trap, & pressure gauge.

**Utility Position ?**  
**Dimension Vista Power**  
1- 208-220VAC 60hz 1-Phase, 50Amp  
Dedicated Circuit 6 AWG Wire  
Required.  
Receptacle Hubbel CS 8269 or equivalent

**Accessory Power**  
1- 120 VAC 60 Hz 1-phase, 15 Amp  
dedicated circuit Receptacle Duplex  
Nema #5-15R

**Vista Water Supply**  
**Feed Water Connection**  
Tap Water 1/2" NPT Male Fitting  
See Vista Specs for details

**Utility Position ?**  
**Dimension Vista Power**  
1- 208-220VAC 60hz 1-Phase, 50Amp  
Dedicated Circuit 6 AWG Wire  
Required.  
Receptacle Hubbel CS 8269 or equivalent

**Accessory Power**  
1- 120 VAC 60 Hz 1-phase, 15 Amp  
dedicated circuit Receptacle Duplex  
Nema #5-15R


**Vista Water Supply**  
**Feed Water Connection**  
Tap Water 1/2" NPT Male Fitting  
See Vista Specs for details

**Utility Position ??**  
(Location to be determined)  
**Jun-Air Air Compressor Power**  
1-230 VAC 60 Hz, 1-phase, 20 Amp  
dedicated / Isolated circuit.  
Receptacle Nema #6-20R  
**Buck / Boost Transformer must be used on circuit.**

**Also in this Position**  
Air Supply Line 1/2" copper tubing run to power position 1 with 3/8" female threaded fitting.

**Floor Drain**  
(Location to be determined)

**Considerations:**  
Instrument Waste  
Possible Trip Hazard



**CONCEPTUAL DRAWING ONLY,  
NOT TO SCALE**

EXACT POWER LOCATIONS WILL BE MARKED ON THE LABORATORY FLOOR.

For any Questions or concerns Contact:

Dennis A. Clement  
Customer Project Specialist  
Laboratory Automation  
Siemens Healthcare Diagnostics  
Cellular: 707-592-3398  
Voice: 800-948-3233, opt 2, ext.

**SIEMENS**

**Vista Water Supply**  
(Location to be determined)

**Considerations:**  
Instrument Location