# 

# Healthcare www.healthcare.philips.com

# **Final Site Preparation Support Document**

The equipment components shown in this drawing package are based on the current proposed purchase and are subject to change if modifications are made to the configuration.

		Revision History Note for Architects and/or Contractors: If revisions are listed, these drawings must be thoroughly reviewed so that all changes can be incorporated into your project		
Rev.	Date	Revision Descriptions	Ву	Section A
А	5/7/2014	AL/A1: Added 25kVA UPS in equipment closet per PM request.	LP	Section A General N
В	9/10/2014	Created Final Site Preparation Support Document. Updated equipment per OA# 6600224430.010000.	LP	Equipmer
С	11/14/2014	AL/A1/AD6/SL/S1/EL/E1/E3/E4/ED4: Updated 25kVA UPS information and added knife switch, "SWC".	LP	Equipme Transpor Equipme
				Section S
				Support Support Support Support Support Support Electrical Electrical Raceway Electrical
				Remote Check Li

lotes nt Lege Plan Detai t Deta

### Support Plan

lotes egen 'lan -'lan -Details

### Electrical Plan

Notes Leger Plan · & Cor Detai

Service ----

rights

₹

2012.

ics N.V.

### Table of Contents

### **Equipment Plan**

	AN
end	AL
n	A1
iils	AD1
ails	AD2 - AD6

	SN
d	SL
Floor & Wall	
Ceiling	S2
s SD1 - S	D3

s	EN
nd	EL
	E1
nduit Information	E2 - E4
ls	ED1 - ED4

e Network	N1-N2
	CHK

Project	Allura FD20 Ceiling		VA Oklahoma Citv	Oklahoma City, OK	-Room 1 1st Floor	THE DRAWINGS AND RELATED INSTRUCTIONS PROVIDED BY PHILIPS ARE ACCEPTABLE FOR USE BY THE HOSPITAL'S ARCHITECT OR ENGINEER TO USE FOR THE DEVELOPMENT OF
Philips Contacts	Project Manager: John Wright	Contact Number: (214) 704-8619	Email: john.wright@philips.com		Drawn By: Laura Phillips	PHILIPS ARE ACCEPTABLE FOR USE BY THE HOSPITAL
Project Details	Drawing Number	N-WES140162 C	Date Drawn: 11/14/2014	Quote: 1-ZMC9DA Rev. 2	Order: 6600224430.010000	AND RELATED INSTRUCTIONS PROVIDED BY
		C	;1			THE DRAWINGS /
			8	3.20	J.14	

### **General Specifications**

### 1. Responsibility

The customer shall be solely responsible, at its expense for preparation of site, including any required structural alterations. The site preparation shall be in accordance with plans and specifications provided by Philips. Compliance with all safety electrical and building codes relevant to the equipment and its installation is the sole responsibility of customer. The customer shall advise Philips of conditions at or near the site which could adversely affect the carrying out of the installation work and shall ensure that such conditions are corrected and that the site is fully prepared and available to Philips before the installation work is due to begin. The customer shall provide all necessary plumbing, carpentry work, or conduit wiring required to attach and install products ready for use.

### 2. Permits

Customer shall obtain all permits and licenses required by federal, state/provincial or local authorities in connection with the construction, installation and operation of the products and related rules, regulations, shall bear any expense in obtaining same or in complying with any ordinances and statutes.

### 3. Radiation Protection

The customer or his contractor, at his own expense, shall obtain the service of a licensed radiation physicist to specify radiation protection. (X-Ray Tube output 125 KVp max.)

### 4. Asbestos and Other Toxic Substances

Philips assumes no hazardous waste (i.e., PCB's in existing transformers) exists at the site. If any hazardous material is found, it shall be the sole responsibility of the customer to properly remove and dispose of this material at its expense. Any delays caused in the project for this special handling shall result in Philips time period for completion being extended by like period of time. Philips assumes that no asbestos material is involved in this project in any ceilings, walls or floors. If any asbestos material is found anywhere on the site, it shall be the customer's sole responsibility to properly remove and/or make safe this condition, at the customer's sole expense.

### 5. Labor

In the event local labor conditions make it impossible or undesirable to use Philips' regular employees for such installation and connection, such work shall be performed by laborers supplied by the customer, or by an independent contractor chosen by the customer at the customer's expense, and in such case, Philips agrees to furnish adequate engineering supervision for proper completion of the installation.

### 6. Schedule

2012.

The general contractor should provide Philips with a schedule of work to assist in the coordination of delivery of Philips supplied products which are to be installed by the contractor and delivery of the primary equipment.

### 7. Extended Installation or Turnkey Work by Philips

Any room preparation requirements for Philips equipment indicated on these drawings is the responsibility of the customer. If an extended installation or turnkey contract exists between Philips and the customer for room preparation work required by the equipment represented on these drawings, some of the responsibilities of the customer as depicted in these drawings may be assumed by Philips. In the event of a conflict between the work described in the turnkey contract workscope and these drawings, the turnkey contract workscope shall govern.

### 8. Infection Control and Interim Life Safety Measures

Compliance with all Infection Control and Interim Life Safety Measures shall be the sole responsibility of the customer. The customer shall provide all means and methods necessary for compliance with Infection Control (IC) and Interim Life Safety Measures (ILSM) in connection with the construction and installation/operation of the products shown herein and shall bear any expenses related to same.

(14.0)

### **Minimum Site Preparation Requirements**

A smooth efficient installation is vital to Philips and their customers. Understanding what the minimum site preparation requirements are will help achieve this goal. The following list clearly defines the requirements which must be fulfilled before the installation can begin.

1. Walls to be painted or covered, baseboards installed, floors to be tiled and/or covered, ceiling shall have grid tiles and lighting fixtures installed and operational.

2. Doors and windows, especially radiation protection barriers, installed and finished with locksets operational.

3. All electrical convenience, conduit, raceway, knockouts, cable openings, chase nipples, and junction boxes installed and operational.

- 4. Incoming mains power operational and connected to room x-ray breaker.
- 5. 120V convenience outlets operational.

6. All support structure correctly installed. All channels, pipes, beams and/or other supporting devices should be level, parallel, and free of lateral or longitudinal movements.

- 7. All contractor supplied cables pulled and terminated
- 8. A dust-free environment in and around the procedure room.

9. All HVAC (heating, ventilating and air conditioning) installed and operational as per specifications.

10. Architectural features such as computer floor, wood floor, casework, bulkheads, installed and finished. When technical cabinets are installed in a closet with doors, it is suggested that the customer install a temperature alarm in the event of an air conditional failure.

11. All plumbing installed and finished.

12. Philips does not install or connect developing tanks, automatic processors or associated equipment, built in illuminators, cassette pass boxes, loading benches and cabinets, lead protective screens, panels or lead glass window and frame. This is to be done by the customer/contractor.

13. Clear door openings for moving equipment into the building must be 42" (1067mm) W x 82" (2083mm) H min, 48" (1219mm) W x 82" (2083mm) H rec., Or larger contingent on an 8'-0" (2438mm) corridor width.

14. Countertop is 30" (765mm) for seated height and 36" (915mm) for standing height.

### Note

Once Philips has moved equipment into the suite and started the installation, the contractor shall schedule his work around the Philips installation team on site. It is suggested that a telephone be provided in the room to receive telephone calls. This would alleviate facility staff from answering calls for Philips personnel.

### Remote Service Diagnostics

Medical imaging equipment to be installed by Philips Medical is equipped with a service diagnostic feature which allows for remote and on site service diagnostics. To establish this feature, a RJ45 type ethernet 10/100/1000 Mbit network connector must be installed as shown on plan. Access to customer's network via their remote access server is needed for Remote Service Network (RSN) connectivity. All cost with this feature are the responsibility of the customer.

(12.0)

Temperature
Temperature gradient
Humidity (non-condensing Humidity shall be stable within

		_
Equi	pment	Room

\*Average heat emission during clinic Data applicable for basic system: Large monitor + 4 x small monitor in 1 workstation + 2 x small monitor in

Add 1706 BTU/hr for additional large Add 273 BTU/hr for additional small Add 1024 BTU/hr for additional work

Equipment's designed airflow is from handling in the rack cabinet equipme

	Electi Ma
Power Output:	100kW
Supply Configuration:	3 phase, ider bonding conc
Nominal Line Voltage:	480 VAC, 60
Branch Power Requireme	ent: 225 k\
Circuit Breaker:	3 phase, Typ

### Remote

The control of customer lighting must demonstrated on Sheet ED3. Lighti

HVAC Re	equirement for Gen	eral Equipment Locations				V	
	Opera	ation					
Temp	perature	50°F (10°C) to 86°F (30°C)					
Temperat	ture gradient	Max. 1°F / Minute (0.5°C / Minute	)				
	on-condensing) e stable within 10%	20% to 80%					
Exan	n Room	*6483 BTU/hr					
Equipm	nent Room	*8189 BTU/hr					
cable for basi hitor + 4 x sm ion + 2 x sma BTU/hr for ado BTU/hr for ado	n during clinical use ic system: all monitor in Monitor Cei all monitor in Control Roo dditional large monitor ditional small monitor dditional workstation	5 1					
	airflow is from bottom to to binet equipment area acc	op and front to back. Please design the ordingly.	air				
			(14.0)				
	Electrical Red Mains 40E			Project Allura FD20 Ceiling	0	a City	loor
out:	100kW			20 (		Ö	st
figuration:	3 phase, identical 3 wir bonding conductor, delt	e power and isolated unity ground with a (preferred) or wye		ct ra FD	- - 5	VA Oklahoma Oklahoma Citv	-Room 1 1s
ne Voltage:	480 VAC, 60 Hz			Project Allura		<b>A</b> N N	Å Å
ver Requirem	nent: 225 kVA						
aker:	3 phase, Type D 125 A	with long-time delay and shunt trip	(14.2)	n Wright	(214) 704-8619	hilips.com	sd
	Remote Control o	of Room Lighting		r: Joh	r. (21	ght@p	a Philli
	r lighting must incorporate	e an electrical isolation system such as s the responsibility of the customer.	(12.0)	Philips Contacts Project Manager: John Wright	Contact Number:	Email: john.wright@philips.com	Drawn By: Laura Phillips
				Project Details Drawing Number	N-WES140162 C	Date Drawn: 11/14/2014 Quote: 1-ZMC9DA Rev. 2	Order: 6600224430.010000
					A	N	

THE DRAWINGS AND RELATED

Γ			
F	A	A Furr	nished and installed by Phi nished by customer/contract
_	C	3 Furr C Inst	nished by customer/contract alled by customer/contract nished by Philips and instal
	E	E Exis F Futu	sting
	G	G Opti	onal item furnished by Phi
	1	Г	— Equip
		$\downarrow$	
A	، (s	SP	Clea Stand
A	. (M;	ISA	Angio Diagnost 7 w
A		VE)	Certeray iX Generat
A			Peripheral 40E Cab
A			Mains 40E Cabinet
A			Image 40E Cabinet
A			Viewing/Control
A			Documentation Box
			(Final location to be
	<u> </u>		and/or local Philips Exam Room Auxilia
			58" LCD Monitor Su
A			
A		~	Video Connection B
A			Video Connection B
A	. (М.	IAV	Mavig Ceiling Track
A	. (т	TR)	3MC Surgical Light M LED 3MC Transfe
			UPS Cabinet - 25 k
			Universal Power Co
			Remote Status Pan
			Knife Switch
			Dose Aware - Base
			One LCD Monitor C
		<u>v</u> 2	
_			

copyright

of the

ž

ion in whole or in part is

Rep eserved

onics N.V. 2012. All rights

Elect Philips inklijke |

0

### Equipment Legend

ilips ctor and installed by customer/contractor alled by contractor

ilips oment Designation Detail Sheet \_\_\_\_ Weight Heat Load Description (lbs) (btu/hr) 2557 1706 AD2 ith Pivot and Tilt 1693 205 AD2 320 2971 AD3 tor Cabinet 441 2049 AD3 inet 826 5464 AD3 441 1877 AD3 126 567 AD3 176 0 AD4 - Mounted on Wheels coordinated with customer Service) 1.7 AD4 ary Box 7 603 1020 AD4 uspension Box 11 34 AD4 11 34 AD4 Зох 350 AD5 w/ Radiation Shield and M LED 167 AD5 ormer 17 998 AD6 kVA 11564 1020 ontroller - 25 kVA AD6 nel (for UPS) 12 50 AD6 22 - AD6 85 AD5 Station 3.2 239 AD5 68 Carriage

Ъ Ш 뿓 <sup>1</sup>OR Project Allura FD20 Ceiling VA Oklahoma City Oklahoma City, OK -Room 1 1st Floor ARCHITECT OR ENGINEER TO Philips Contacts Project Manager: John Wright Contact Number: (214) 704-8619 Email: john.wright@philips.com ΗH В≺ Drawn By: Laura Phillips ILIPS ARE ACCEPTABLE FOR USE PHILIPS ARE 
 Project Details

 Drawing Number


 Drawing Number

 N-WES140162

 Date Drawn: 11/14/2014

 Quote: 1-ZMC9DA Rev. 2

 Order: 6600224430.010000
 THE DRAWINGS AND RELATED INSTRUCTIONS CONSTRUCTION DOCUMENTS. AL






Source	Location	Displayed
Physio	Control	FlexVision
IH	Control	FlexVision
MB	Exam	TV2

	Project Details	Philips Contacts	Project
	Drawing Number	Project Manager: John Wright	Allura FD20 Ceiling
	N-WES140162 C	Contact Number: (214) 704-8619	
	Date Drawn: 11/14/2014	Email: john.wright@philips.com	VA Oklahoma Citv
-	Quote: 1-ZMC9DA Rev. 2		Oklahoma City, OK
-	Order: 6600224430.010000	Drawn By: Laura Phillips	-Room 1 1st Floor


# THE DRAWINGS AND RELATED CONSTRUCTION DOCUMENTS.

### Detail - Clea Ceiling (C-ARM) Transport Details




		Transport Possibilities		
	Crate	Pallet	Kick Wheels Wide	Kick Wheels Small
Height	77.95" (1980mm)	76.22" (1936mm)	69.02" (1753mm)	77.76" (1975mm)
Weight	2050 lbs (930 kg)	1940 lbs (880 kg)	2061 lbs (935 kg)	1764 lbs (800 kg)

### Detail - Clea Ceiling (L-ARM) Transport Details




			Transport Possibilities			
	Crate	Pallet	Klick Wheels Wide	Klick Wheels Wide Elevator	Klick Wheels Small	Klick Wheels Small Elevator
Height	57.09" (1450mm)	54.80" (1392mm)	49.25" (1251mm)	79.53" (2020mm)	49.25" (1251mm)	79.53" (2020mm)
Weight	2094 lbs (950 kg)	1973 lbs (895 kg)	1896 lbs (860 kg)	1896 lbs (860 kg)	1896 lbs (860 kg)	1896 lbs (860 kg)


	Project Allura FD20 Ceiling	VA Oklahoma City Oklahoma City, OK	Order: 6600224430.010000 Drawn By: Laura Philips - ROOM 1 1St Floor - Proom 1 St Floor - Proom
	Philips Contacts Project Manager: John Wright Contact Number: (214) 704-8619	Email: john.wright@philips.com	Drawn By: Laura Phillips
	Project Details Drawing Number N-WES140162 C	Date Drawn: 11/14/2014 Quote: 1-ZMC9DA Rev. 2	Order: 6600224430.010000 THE DRAWINGS AND RELATED INSTRUCTIONS PROVIDED BY
(12.0)	A	D1	




0

	Philips Contacts	Project
Drawing Number N-WES140162 C	Project Manager: John Wright Contact Number: (214) 704-8619	Allura FD20 Ceiling
Date Drawn: 11/14/2014 Quote: 1-ZMC9DA Rev. 2	Email: john.wright@philips.com	VA Oklahoma City
Order: 6600224430.010000	Drawn By: Laura Phillips	-Room 1 1st Floor










The CRC Cover must be attached to the the wall box.

Acoustic noise level: <= 65 dB(A) @ 1 meter in front of the rack and 1 meter high (1 meter = 39.37")

	Peripheral 4	(14.0) IOE Cabinet
MP	Weight	Heat Dissipation
	441 lbs	2049 btu/hr







ME

2012.

cs N.V.

Acoustic noise level: <= 55 dB(A) @ 1 meter in front of the rack and 1 meter high (1 meter = 39.37")


	(14.0)
Certeray iX Gene	rator 40E Cabinet
Weight	Heat Dissipation
320 lbs	2971 btu/hr


2 CHITECT OR ENG Å USE PTABLE FOR THE DRAWINGS AND RELATED INS CONSTRUCTION DOCUMENTS.

Drawn By: Laura Phillips

AD3

Project Details







	58" LCD Mon	itor Suspension
ти	Weight	Heat Dissipation
	603 lbs	1020 btu/hr

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED.


Ш 뿥 <sup>1</sup>OR ARCHITECT OR ENGINEER TO ВҮ ТНЕ PHILIPS ARE ACCEPTABLE FOR USE Å THE DRAWINGS AND RELATED INSTRUCTIONS CONSTRUCTION DOCUMENTS.

Ъ






$11\frac{3}{4}$ $298\text{mm}$ $Front$ $298\text{mm}$ $Front$ $5ide$	
(12.0)         Dose Aware - Base Station         DBS       Weight       Heat Dissipation         3.2 lbs       85 btu/hr	Project Allura FD20 Ceiling VA Oklahoma City Oklahoma City, OK -Room 1 1st Floor
	Project Allura FD20 Ceiling VA Oklahoma City, OK -Room 1 1st Floor
$\begin{array}{c} 2^{2}-7,\frac{3}{3},\\ \hline 1,2\\ \hline 0,2\\ \hline 0,3\\ \hline 0,2\\ \hline 0,2$	Project DetailsPrilips ContactsDrawing NumberPhilips ContactsDrawing NumberProject Manager: John WrightDrawing NumberProject Manager: John WrightDate Drawn: 11/14/2014Project Manager: John WrightQuote: 1-ZMC9DA Rev. 2Contact Number: (214) 704-8619Order: 6600224430.010000Drawn By: Laura Phillips
$\begin{array}{c c} \hline 1'-9\frac{1}{2}" & \hline 2'-5\frac{3}{8}" \\ \hline [546mm] & \hline [746mm] \\ \hline Front & \underline{Side} \\ \hline \hline \hline 1'-9\frac{1}{2}" & \hline 2'-5\frac{3}{8}" \\ \hline \hline [746mm] & \hline [746mm] \\ \hline \hline \hline 12.0 \\ \hline $	Project Details Drawing Number N-WES140162 C Date Drawn: 11/14/2014 Quote: 1-ZMC9DA Rev. 2 Order: 6600224430.010000
	AD5



$11\frac{3}{4}^{-}$	
(12.0) Dose Aware - Base Station DBS Weight Heat Dissipation 3.2 lbs 85 btu/hr	Project Details     Project Details       Project Details     Project Details       Drawing Number     Project Manager: John Wright       N-WES140162 C     Project Manager: John Wright       Date Drawn: 11/14/2014     Project Manager: John Wright       Quote: 1-ZMC9DA Rev. 2     Contact Number: (214) 704-8619       Order: 6600224430.010000     Drawn By: Laura Philips
$\begin{array}{c} 2^{1} \cdot 7^{2} \cdot$	Philips Contacts Project Manager: John Wright Contact Number: (214) 704-8619 Email: john.wright@philips.com Drawn By: Laura Philips
1'-9 <sup>1</sup> / <sub>2</sub> "     2'-5 <sup>3</sup> / <sub>8</sub> "       [546mm]     [746mm]       Front     Side       (12.0)       One LCD Monitor Carriage       TV2     Weight       Heat Dissipation       68 lbs     239 btu/hr	Project Details Drawing Number N-WES140162 C Date Drawn: 11/14/2014 Quote: 1-ZMC9DA Rev. 2 Order: 6600224430.010000
	AD5 8.20.14

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED.



		(14.0)
	25 kVA UP	S with UPC
	Weight	Heat Dissipation
UPS	998 lbs	11564 btu/hr
UPC	1020 lbs	11504 Dlu/III
RSP	12 lbs	50 btu/hr
SWC	22 lbs	- btu/hr

Project Details	Philips Contacts	Project	
Drawing Number N-WES140162 C	Project Manager: John Wright Contact Number: (214) 704-8619	Allura FD20 Ceiling	
Date Drawn: 11/14/2014	Email: john.wright@philips.com		
Order: 6600224430.010000	Drawn By: Laura Phillips		

### **Equipment Support Information**

### 1. General

The customer shall be solely responsible, at its expense, for preparation of the site, including any required structural alterations. The site preparation shall be in accordance with this plan and specifications, the architectural/construction drawings and in compliance with all safety and building codes. The customer shall be solely responsible for obtaining all construction permits from jurisdictional authority.

### 2. Equipment Anchorage

Philips provides, with this plan and specifications, information relative to equipment size, weight, shape, anchoring hole locations and forces which may be exerted on anchoring fasteners. The customer shall be solely responsible, through the engineer of record for the building, to provide on the architectural/construction drawings, information regarding the approved method of equipment anchoring to floors, wall and/or ceiling of the building. Any anchorage test required by local authority shall be the customer's responsibility. Stud type anchor bolts should not be specified as they hinder equipment removal for service. Consult with Philips service prior to specifying anchor methods. Philips equipment must be electrically isolated from anchorage.

### 3. Floor Loading and Surface

Philips provides, with this plan and specifications, information relative to size, weight and shape of floor mounted equipment. The customer shall be solely responsible, through the engineer of record for the building, to provide on the architectural/construction drawings confirmation of the structural adequacy of the floor upon which the equipment will be placed. Any load test required by local authority, shall be the customer's responsibility.

The floor surface upon which Philips equipment is to be placed/anchored shall be flat and level to within plus or minus 1/16" (2mm) over a length of 39" (1m).

### 4. Ceiling Support Apparatus

a. Philips provides, with this plan and specifications, information relative to size, weight and shape of ceiling supported equipment. The customer shall be solely responsible, through the engineer of record for the building, to provide on the architectural/construction drawings, information regarding the approved method of structural support apparatus, fasteners and anchorage to which Philips will attach equipment. Any anchorage and/or load test required by local authority shall be the customer's responsibility. Philips equipment must be electrically isolated from anchorage.

b. Contractor to clearly mark Philips equipment longitudinal centerline on bottom of each structural support.

c. The structural support apparatus surface to which Philips equipment is to be attached, shall have horizontal equipment attachment surfaces parallel, square and level to within plus or minus  $\frac{1}{16}$ " (2mm) per entire span.

d. Any drilling and/or tapping of holes required to attach Philips equipment to the structural support apparatus shall be the responsibility of the customer.

e. Fasteners/anchors (i.e., bolts, spring nuts, lock and flat washers) and strip closures shall be provided by the customer.

### 5. Lighting

Lighting fixtures shall be placed in such a position that they are not obscured by equipment or its movement, nor shall they interfere with Philips ceiling rails and equipment movement or otherwise adversely affect the equipment. Such lighting fixture locations shall be the sole responsibility of the customer.

### 6. Ceiling Obstructions

There shall be no obstructions that project below the finished ceiling in the area covered by ceiling suspended equipment travel.

### 7. Seismic Anchorage (For Seismic Zones Only)

All seismic anchorage hardware, including brackets, backing plates, bolts, etc., shall be supplied and installed by the customer/contractor unless otherwise specified within the support legend on this sheet. Installation of electronic cabinets to meet seismic anchorage requirements must be accomplished using flush mounted expansion type anchor/bolt systems to facilitate the removal of a cabinet for maintenance. Do not use threaded rod/adhesive anchor systems. Consult with Philips regarding any anchor system issues. Philips equipment must be electrically isolated from anchorage.

### 8. Floor Obstructions/ Floor Coverings

There shall be no obstructions on the floor (sliding door tracks, etc.) in front of the Philips technical cabinets. Floor must be clear to allow cabinets to be pulled away from the wall for service. Contractor to verify with Philips the preferred floor covering installation method.

### 9. Safety Factors

In a worst case situation the dynamic bolt force of a floor or ceiling must be multiplied by factor 4. (static bolt force of the ceiling must be multiplied by factor 8). All safety factors are included in the bearing force values in sheet SD1.

### 10. Stiffness Requirements of Ceiling

Stiffness: 10,000,000 Newton/meter - 57.1 klb/in Stiffness: 20,000,000 Newtonmeter/Rad - 177,014 (klb in)/Rad The maximum deflection on the Philips rails must not exceed 0.04" (1mm) caused by the static load (weight) of the ceiling stand

### 11. Vibration

The maximal allowed external frequency that will not destroy the image quality of our equipment is:

a. 0 Hz till 20 Hz (frequency area of our equipment) - Displacement amplitude is smaller than 0.005mm


b. Greater than 20 Hz - Displacement amplitude is smaller than 0.01mm

	Project Details	Philips Contacts	Project
	Drawing Number	Project Manager: John Wright	Allura FD20 Ceiling
S	N-WES140162 C	Contact Number: (214) 704-8619	)
N	Date Drawn: 11/14/2014	Email: john.wright@philips.com	VA Oklahoma Citv
]	Quote: 1-ZMC9DA Rev. 2	-	Oklahoma Citv. OK
	Order: 6600224430.010000	Drawn By: Laura Phillips	-Room 1 1st Floor
THE DRAWINGS AND RELATED CONSTRUCTION DOCUMENTS.	AND RELATED INSTRUCTIONS PROVIDED BY P DOCUMENTS.	THE DRAWINGS AND RELATED INSTRUCTIONS PROVIDED BY PHILIPS ARE ACCEPTABLE FOR USE BY THE HOSPITAL'S ARCHITECT OR ENGINEER TO USE FOR THE CONSTRUCTION DOCUMENTS.	S ARCHITECT OR ENGINEER TO USE FOR THE

(14.0)

Dec Od for Elson & Wall Oursport Loueut	Floor & Wall Support Legend	
<ul> <li>See S1 for Floor &amp; Wall Support Layout</li> <li>Notes: <ol> <li>Anchors for items that are installed/anchored by customer/contractor shall be provided by customer/contractor.</li> <li>Anchors for items that are installed/anchored by Philips shall be provided by Philips. If customer's engineering documents specify anchors other than those listed in this document, the anchors shall be provided by customer/contractor and installed by Philips.</li> <li>In all instances, the wall and/or floor support are the sole responsibility of the customer/contractor. The customer's architect/engineer of record shall specify wall and/or floor support sufficient for the bolt forces shown on the details.</li> </ol> </li> </ul>	A Furnished and installed/anchored by Philips (exceptions may exist, see Note 2) B Furnished and installed by customer/contractor C Installed/anchored by customer/contractor D Furnished by Philips and installed/anchored by contractor E Existing F Future G Optional Item Number Detail Sheet Description B F1 A F1 Anchors in wall for Control Room Connection Box (CY) SD3	HLP
	D       F2       AD7 Adaptation Plate       SD1         C       F3       Anchors in wall for Dose Aware Base Station       SD2         C       F4       Anchors in wall for SWC       -	
See S2 for Ceiling Support Layout	Ceiling Support Legend         A Furnished and installed by Philips         B Furnished by customer/contractor and installed by customer/contractor         C Installed by customer/contractor         D Furnished by Philips and installed by contractor         E Existing         F Future         G Optional	Project Allura FD20 Ceiling VA Oklahoma City, OK Oklahoma City, OK -Room 1 1st Floor
	Item Number       Detail Sheet         V       V         A       C1       2 - Philips Clea Rails         A       C2       2 - Philips Monitor Equipment Rails         B       C3       Unistrut (P1001 or equal) - Bottom of Unistrut 0" to $\frac{1}{4}$ " (6mm) Below Finished Ceiling         A       C4       Mavig Ceiling Track	Philips Contacts Project Manager: John Wright Contact Number: (214) 704-8619 Email: john.wright@philips.com Drawn By: Laura Phillips
		Project Details Drawing Number N-WES140162 C Date Drawn: 11/14/2014 Quote: 1-ZMC9DA Rev. 2 Order: 6600224430.010000

THE DRAWINGS AND RELATED CONSTRUCTION DOCUMENTS.




# Floor & Wall Support Layout

Required Unistrut Height: 9' - 6  $\frac{3}{16}$ ", + $\frac{3}{8}$ " / -0 (2900mm, +10mm / -0) Unistrut height measured from finished floor to bottom of Unistrut.



0

	DHLIPS
	Project Details         Philips Contacts         Project           Drawing Number         Project Details         Project           Drawing Number         Project Manager: John Wright         Project           N-WES140162 C         Contact Number: (214) 704-8619         Project           Date Drawn: 11/14/2014         Email: john.wright@philips.com         VA Oklahoma City, OK           Order: 6600224430.010000         Drawn By: Laura Philips         VA Oklahoma City, OK
	Philips Contacts Project Manager: John Wright Contact Number: (214) 704-8619 Email: john.wright@philips.com Drawn By: Laura Phillips
	Project Details Drawing Number N-WES140162 C Date Drawn: 11/14/2014 Quote: 1-ZMC9DA Rev. 2 Order: 6600224430.010000 ND RELATED INSTRUCTIONS PROVIDED BY1
Refer to Floor/Wall Support Legend - Sheet SL	S1
	 8.20.14



# **Ceiling Support Layout**


Required Unistrut Height: 9' - 6  $\frac{3}{16}$ ", + $\frac{3}{8}$ " / -0 (2900mm, +10mm / -0) Unistrut height measured from finished floor to bottom of Unistrut.




0

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED.

	Project	Allura FD20 Ceiling	•	VA Oklahoma City	Oklahoma City. OK	-Room 1 1st Floor	THE DRAWINGS AND RELATED INSTRUCTIONS PROVIDED BY PHILIPS ARE ACCEPTABLE FOR USE BY THE HOSPITAL'S ARCHITECT OR ENGINEER TO USE FOR THE DEVELOPMENT OF
	Philips Contacts	Project Manager: John Wright	Contact Number: (214) 704-8619	Email: john.wright@philips.com		Drawn By: Laura Phillips	PHILIPS ARE ACCEPTABLE FOR USE BY THE HOSPITAL
	Proiect Details	Drawing Number	N-WES140162 C	Date Drawn: 11/14/2014	Quote: 1-ZMC9DA Rev. 2	Order: 6600224430.010000	ND RELATED INSTRUCTIONS PROVIDED BY
Refer to Ceiling Support Legend - Sheet SL				52			HE DRAWINGS AN
 				8	3.20	0.14	」 ⊨ à





2 OR ENG β EPTABLE FOR USE THE DRAWINGS AND RELATED INS CONSTRUCTION DOCUMENTS.



ABLE FOR THE DRAWINGS AND RELATED CONSTRUCTION DOCUMENTS.

Equipment	Option	Anchor Style (provided by Philips)	Anchor Size (provided by Philips)	Qty.	Support Size & Material (provided & installed by customer/contractor)
Mavig Ceiling Track	A	Bolts, flat washer, lock washer, spring nuts	A307 Grade or ASME Grade 5 Bolts: $\frac{3}{8}$ " (10mm) x 2" (50mm) L Spring Nuts: $\frac{3}{8}$ " (10mm)	8	Unistrut
	A	Round Phillips Head Self Drilling Screws	#10-16 x 1 ½" (38mm) L	3	Drywall with minimum 20 gauge Steel backing
Control Room Connection Box (CY)	В	SPAX Multipurpose flat head screw	#10 x 1 ½" (38mm) L	3	Drywall with minimum 20 gauge Steel backing
	С	Toggler Snaptoggle and (round head screws)	#BA and (#10-24 x 2 <sup>1</sup> / <sub>2</sub> " (63.5mm) L)	3	Minimum $rac{5}{8}$ (16mm) Drywall





THE DRAWINGS AND RELATED INS CONSTRUCTION DOCUMENTS.


Emergency Power	General Electrical Information
<ul> <li>Philips does not require equipment to be on emergency power. If the customer deems it necessary for the equipment to be supplied with emergency power, the following specifications must be applied:</li> <li>The Mains 40E cabinet feeding an Allura Xper system will have an absolute peak surge current of &lt;380A.</li> <li>The transfer switch must be double actuator type with a minimum time delay of 400 milliseconds in both directions (utility to emergency - emergency to utility). This time is required to allow filters to dissipate their stored energy before a different mains voltage is applied. Russelectric type RMTD, Asco Series 7000 delayed transition transfer switch or equivalent is recommended.</li> <li>To reduce the emergency power generator load demand, Philips equipment can be put into a lower power mode of operation by the connection of a potential free closure from the transfer switch. This potential free, normally open contact, has to be rated for 24VDC/100mA. For Philips cardio/vascular Allura equipment, the two wires from this contact have to be routed to the equipment area and connected to the System Coordinator cabinet (MA).</li> </ul>	<ol> <li>General         The customer shall be solely responsible, at its expense, for preparation of the site, including any required electrical alter with this plan and specifications, the architectural/construction drawings and in compliance with all safety and electrical obtaining all electrical permits from jurisdictional authority.     </li> <li>Materials and Labor         The customer shall be solely responsible, at its expense, to provide and install all electrical ducts, boxes, conduit, cable herein.     </li> <li>Electrical Ducts and Boxes         Electrical ducts and boxes shall be accessible and have removable covers. Floor ducts and boxes shall have watertight separate channels by metal dividers, separately specified herein, to separate wiring and/or cables into groups as follow protective earth wiring (PE). Group B: Output power wiring with associated protective earth wiring (PE). Group C: signal high-voltage cables, the use of 90 deg. ells is not acceptable. On ceiling duct and wall duct use 45 deg. bends at all cor tunnels supplied and installed by contractor to maintain separation of cables.     </li> <li>Conduit         Conduit point - to - point runs shall be as direct as possible. Empty conduit runs used for cables may require pull boxes or cord shall be installed in each conduit run. All conduits which enter duct prior to their termination point must maintair over tunnels, or conduit supplied and installed by contractor from entrance into duct to exit from duct. Do not use flex or conduit supplied and installed by contractor form entrance into duct to exit from duct. Do not use flex or conduit supplied and installed by contractor from entrance into duct to exit from duct. Do not use flex or conduit supplied and installed by contractor from entrance into duct to exit from duct.     </li> </ol>
Electrical Requirement Notes for Systems with Mains 40E Cabinet	5. Conductors
Electrical power distribution at the facility shall comply with:	All conductors, separately specified, shall be 75°C stranded copper, rung out and marked.
Utilization voltages per ANSI C84.1 - 1982 range A.	<ol> <li>Disconnecting Means</li> <li>A disconnecting means shall be provided as separately specified.</li> </ol>
Voltage to be supplied is 3 phase, delta.	7. Warning Lights and Door Switches
Phase conductors to be sized for instantaneous voltage drop per NEC 517 - 73 and Philips recommendations.	"X-ray on" warning lights and x-ray termination door switches should be provided at all entrances to x-ray rooms as requ
Metal conduit shall not be used as the equipment ground conductor.	<ol> <li>Dimmer Switches</li> <li>X-ray room lights should be provided with dimmer switches.</li> </ol>
The Philips system uses an isolated ground scheme grounding only the Allura system per clause 250.96B of the NEC. The raceway from the X-ray breaker (CB) to the Mains 40E Cabinet shall be supplemented by an internal insulated equipment grounding conductor installed in accordance with clause 250.146(D) of the NEC.	
ANSI / NFPA 70 - National Electrical Code	Electrical Notes
Article 250 - Grounding Article 517 - Healthcare Facilities ANSI / NFPA 99 - Healthcare Facilities NEMA standard XR9 - Power Supply Guideline for X-ray Machines	1. The contractor will supply & install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.
Article 250 - Grounding Article 517 - Healthcare Facilities ANSI / NFPA 99 - Healthcare Facilities NEMA standard XR9 - Power Supply Guideline for X-ray Machines	1. The contractor will supply & install all breakers, shunt trip and incoming power to the breakers. The exact location
Article 250 - Grounding Article 517 - Healthcare Facilities ANSI / NFPA 99 - Healthcare Facilities NEMA standard XR9 - Power Supply Guideline for X-ray Machines <b>Power Quality Guidelines</b>	<ol> <li>The contractor will supply &amp; install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.</li> <li>The contractor shall supply &amp; install all pull boxes, raceways, conduit runs, stainless steel covers, etc. Conduit/rac entire length. A Greenlee pull string/measuring tape (part no. 435, or equivalent) shall be provided with conduit runs.</li> <li>All pre - terminated, cut to length cables, will be supplied and installed by Philips. All cables to the breakers, will be</li> </ol>
Article 250 - Grounding Article 517 - Healthcare Facilities ANSI / NFPA 99 - Healthcare Facilities NEMA standard XR9 - Power Supply Guideline for X-ray Machines	<ol> <li>The contractor will supply &amp; install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.</li> <li>The contractor shall supply &amp; install all pull boxes, raceways, conduit runs, stainless steel covers, etc. Conduit/rac entire length. A Greenlee pull string/measuring tape (part no. 435, or equivalent) shall be provided with conduit runs.</li> </ol>
Article 250 - Grounding Article 517 - Healthcare Facilities ANSI / NFPA 99 - Healthcare Facilities NEMA standard XR9 - Power Supply Guideline for X-ray Machines <b>Power Quality Guidelines</b> 1. Power supplied to medical imaging equipment must be separate from power feeds to air conditioning, elevators, outdoor lighting, and other frequently switched or motorized loads. Such loads can cause waveform	<ol> <li>The contractor will supply &amp; install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.</li> <li>The contractor shall supply &amp; install all pull boxes, raceways, conduit runs, stainless steel covers, etc. Conduit/rac entire length. A Greenlee pull string/measuring tape (part no. 435, or equivalent) shall be provided with conduit runs.</li> <li>All pre - terminated, cut to length cables, will be supplied and installed by Philips. All cables to the breakers, will be local arrangements.</li> <li>Provide and install 50mm diameter chase nipples between adjacent wall boxes.</li> <li>Electrical raceway shall be installed with removable covers. The raceway should be accessible for the entire lengtl an adequate number of access hatches should be supplied to enable installation of cabling. Approved conduits may be</li> </ol>
Article 250 - Grounding Article 517 - Healthcare Facilities ANSI / NFPA 99 - Healthcare Facilities NEMA standard XR9 - Power Supply Guideline for X-ray Machines <b>Power Quality Guidelines</b> 1. Power supplied to medical imaging equipment must be separate from power feeds to air conditioning, elevators, outdoor lighting, and other frequently switched or motorized loads. Such loads can cause waveform distortion and voltage fluctuations that can hinder high quality imaging. 2. Equipment that utilizes the facility power system to transmit control signals (especially clock systems) may	<ol> <li>The contractor will supply &amp; install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.</li> <li>The contractor shall supply &amp; install all pull boxes, raceways, conduit runs, stainless steel covers, etc. Conduit/rac entire length. A Greenlee pull string/measuring tape (part no. 435, or equivalent) shall be provided with conduit runs.</li> <li>All pre - terminated, cut to length cables, will be supplied and installed by Philips. All cables to the breakers, will be local arrangements.</li> <li>Provide and install 50mm diameter chase nipples between adjacent wall boxes.</li> <li>Electrical raceway shall be installed with removable covers. The raceway should be accessible for the entire length an adequate number of access hatches should be supplied to enable installation of cabling. Approved conduits may be manner that will not allow cables to fall out of the raceway when the covers are removed. In most cases, this will require removable from the top. Raceway system as illustrated on this drawing are based upon length of furnished cables. Any</li> </ol>
Article 250 - Grounding Article 517 - Healthcare Facilities ANSI / NFPA 99 - Healthcare Facilities NEMA standard XR9 - Power Supply Guideline for X-ray Machines <b>Power Quality Guidelines</b> 1. Power supplied to medical imaging equipment must be separate from power feeds to air conditioning, elevators, outdoor lighting, and other frequently switched or motorized loads. Such loads can cause waveform distortion and voltage fluctuations that can hinder high quality imaging. 2. Equipment that utilizes the facility power system to transmit control signals (especially clock systems) may interfere with medical imaging equipment, thus requiring special filtering.	<ol> <li>The contractor will supply &amp; install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.</li> <li>The contractor shall supply &amp; install all pull boxes, raceways, conduit runs, stainless steel covers, etc. Conduit/rac entire length. A Greenlee pull string/measuring tape (part no. 435, or equivalent) shall be provided with conduit runs.</li> <li>All pre - terminated, cut to length cables, will be supplied and installed by Philips. All cables to the breakers, will be local arrangements.</li> <li>Provide and install 50mm diameter chase nipples between adjacent wall boxes.</li> <li>Electrical raceway shall be installed with removable covers. The raceway should be accessible for the entire length an adequate number of access hatches should be supplied to enable installation of cabling. Approved conduits may be manner that will not allow cables to fall out of the raceway when the covers are removed. In most cases, this will require removable from the top. Raceway system as illustrated on this drawing are based upon length of furnished cables. Any maximum allowable length of furnished cables. Conduit or raceway above - ceiling must be kept as near to finished ceiling the supplied ceiling must be kept as near to finished ceiling the supplied ceiling must be kept as near to finished ceiling maximum allowable length of furnished cables.</li> </ol>
<ul> <li>Article 250 - Grounding Article 517 - Healthcare Facilities</li> <li>ANSI / NFPA 99 - Healthcare Facilities</li> <li>NEMA standard XR9 - Power Supply Guideline for X-ray Machines</li> <li>Power Quality Guidelines</li> <li>1. Power supplied to medical imaging equipment must be separate from power feeds to air conditioning, elevators, outdoor lighting, and other frequently switched or motorized loads. Such loads can cause waveform distortion and voltage fluctuations that can hinder high quality imaging.</li> <li>2. Equipment that utilizes the facility power system to transmit control signals (especially clock systems) may interfere with medical imaging equipment, thus requiring special filtering.</li> <li>3. The following devices provide a high impedance, nonlinear voltage source, which may affect image quality: Static UPS systems, Series filters, Power conditioners, and Voltage regulators.</li> <li>Do not install such devices at the mains supply to medical imaging equipment without consulting Philips</li> </ul>	<ol> <li>The contractor will supply &amp; install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.</li> <li>The contractor shall supply &amp; install all pull boxes, raceways, conduit runs, stainless steel covers, etc. Conduit/rac entire length. A Greenlee pull string/measuring tape (part no. 435, or equivalent) shall be provided with conduit runs.</li> <li>All pre - terminated, cut to length cables, will be supplied and installed by Philips. All cables to the breakers, will be local arrangements.</li> <li>Provide and install 50mm diameter chase nipples between adjacent wall boxes.</li> <li>Electrical raceway shall be installed with removable covers. The raceway should be accessible for the entire length an adequate number of access hatches should be supplied to enable installation of cabling. Approved conduits may be manner that will not allow cables to fall out of the raceway when the covers are removed. In most cases, this will require removable from the top. Raceway system as illustrated on this drawing are based upon length of furnished cables. Any maximum allowable length of furnished cables. Conduit or raceway above - ceiling must be kept as near to finished ceile.</li> <li>Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or National contractor.</li> </ol>
<ul> <li>Article 250 - Grounding Article 517 - Healthcare Facilities</li> <li>ANSI / NFPA 99 - Healthcare Facilities</li> <li>NEMA standard XR9 - Power Supply Guideline for X-ray Machines</li> <li>Power Quality Guidelines</li> <li>1. Power supplied to medical imaging equipment must be separate from power feeds to air conditioning, elevators, outdoor lighting, and other frequently switched or motorized loads. Such loads can cause waveform distortion and voltage fluctuations that can hinder high quality imaging.</li> <li>2. Equipment that utilizes the facility power system to transmit control signals (especially clock systems) may interfere with medical imaging equipment, thus requiring special filtering.</li> <li>3. The following devices provide a high impedance, nonlinear voltage source, which may affect image quality: Static UPS systems, Series filters, Power conditioners, and Voltage regulators.</li> <li>Do not install such devices at the mains supply to medical imaging equipment without consulting Philips installation or service personnel.</li> </ul>	<ol> <li>The contractor will supply &amp; install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.</li> <li>The contractor shall supply &amp; install all pull boxes, raceways, conduit runs, stainless steel covers, etc. Conduit/rac entire length. A Greenlee pull string/measuring tape (part no. 435, or equivalent) shall be provided with conduit runs.</li> <li>All pre - terminated, cut to length cables, will be supplied and installed by Philips. All cables to the breakers, will be local arrangements.</li> <li>Provide and install 50mm diameter chase nipples between adjacent wall boxes.</li> <li>Electrical raceway shall be installed with removable covers. The raceway should be accessible for the entire length an adequate number of access hatches should be supplied to enable installation of cabling. Approved conduits may be manner that will not allow cables to fall out of the raceway when the covers are removed. In most cases, this will require removable from the top. Raceway system as illustrated on this drawing are based upon length of furnished cables. Any maximum allowable length of furnished cables. Conduit or raceway above - ceiling must be kept as near to finished ceiles.</li> <li>Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or Nationa 7. Convenience outlets are not illustrated. Their number and location are to be specified by the customer/architect.</li> </ol>
<ul> <li>Article 250 - Grounding Article 517 - Healthcare Facilities</li> <li>ANSI / NFPA 99 - Healthcare Facilities</li> <li>NEMA standard XR9 - Power Supply Guideline for X-ray Machines</li> <li>Power Quality Guidelines</li> <li>1. Power supplied to medical imaging equipment must be separate from power feeds to air conditioning, elevators, outdoor lighting, and other frequently switched or motorized loads. Such loads can cause waveform distortion and voltage fluctuations that can hinder high quality imaging.</li> <li>2. Equipment that utilizes the facility power system to transmit control signals (especially clock systems) may interfere with medical imaging equipment, thus requiring special filtering.</li> <li>3. The following devices provide a high impedance, nonlinear voltage source, which may affect image quality: Static UPS systems, Series filters, Power conditioners, and Voltage regulators.</li> <li>Do not install such devices at the mains supply to medical imaging equipment without consulting Philips installation or service personnel.</li> <li>4. Line impedance is the combined resistance and inductance of the electrical system and includes the impedance of the power source, the facility distribution system, and all phase conductors between the source and</li> </ul>	<ol> <li>The contractor will supply &amp; install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.</li> <li>The contractor shall supply &amp; install all pull boxes, raceways, conduit runs, stainless steel covers, etc. Conduit/rac entire length. A Greenlee pull string/measuring tape (part no. 435, or equivalent) shall be provided with conduit runs.</li> <li>All pre - terminated, cut to length cables, will be supplied and installed by Philips. All cables to the breakers, will be local arrangements.</li> <li>Provide and install 50mm diameter chase nipples between adjacent wall boxes.</li> <li>Electrical raceway shall be installed with removable covers. The raceway should be accessible for the entire length an adequate number of access hatches should be supplied to enable installation of cabling. Approved conduits may be manner that will not allow cables to fall out of the raceway when the covers are removed. In most cases, this will require removable from the top. Raceway system as illustrated on this drawing are based upon length of furnished cables. Any maximum allowable length of furnished cables. Conduit or raceway above - ceiling must be kept as near to finished ceile.</li> <li>Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or National contractor.</li> </ol>
<ul> <li>Article 250 - Grounding Article 517 - Healthcare Facilities</li> <li>ANSI / NFPA 99 - Healthcare Facilities</li> <li>NEMA standard XR9 - Power Supply Guideline for X-ray Machines</li> <li>Power Quality Guidelines</li> <li>1. Power supplied to medical imaging equipment must be separate from power feeds to air conditioning, elevators, outdoor lighting, and other frequently switched or motorized loads. Such loads can cause waveform distortion and voltage fluctuations that can hinder high quality imaging.</li> <li>2. Equipment that utilizes the facility power system to transmit control signals (especially clock systems) may interfere with medical imaging equipment, thus requiring special filtering.</li> <li>3. The following devices provide a high impedance, nonlinear voltage source, which may affect image quality: Static UPS systems, Series filters, Power conditioners, and Voltage regulators.</li> <li>Do not install such devices at the mains supply to medical imaging equipment without consulting Philips installation or service personnel.</li> <li>4. Line impedance is the combined resistance and inductance of the electrical system and includes the impedance of the power source, the facility distribution system, and all phase conductors between the source and the imaging equipment. Philips publishes recommended conductor sizes based on equipment power requirements, acceptable voltage drops, and assumptions about the facility source impedance. The minimum conductor size is based on the total line impedance and NEC requirements. Unless impedance calculations are</li> </ul>	<ol> <li>The contractor will supply &amp; install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.</li> <li>The contractor shall supply &amp; install all pull boxes, raceways, conduit runs, stainless steel covers, etc. Conduit/race entire length. A Greenlee pull string/measuring tape (part no. 435, or equivalent) shall be provided with conduit runs.</li> <li>All pre - terminated, cut to length cables, will be supplied and installed by Philips. All cables to the breakers, will be local arrangements.</li> <li>Provide and install 50mm diameter chase nipples between adjacent wall boxes.</li> <li>Electrical raceway shall be installed with removable covers. The raceway should be accessible for the entire length an adequate number of access hatches should be supplied to enable installation of cabling. Approved conduits may be manner that will not allow cables to fall out of the raceway when the covers are removed. In most cases, this will requir removable from the top. Raceway system as illustrated on this drawing are based upon length of furnished cables. Any maximum allowable length of furnished cables. Conduit or raceway above - ceiling must be kept as near to finished ceiles.</li> <li>Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or Nationa 7. Convenience outlets are not illustrated. Their number and location are to be specified by the customer/architect.</li> <li>Electrical contractor shall install ground bond wires at conduit openings within wall boxes as required by national and set of the contractor.</li> </ol>
<ul> <li>Article 250 - Grounding Article 517 - Healthcare Facilities</li> <li>ANSI / NFPA 99 - Healthcare Facilities</li> <li>NEMA standard XR9 - Power Supply Guideline for X-ray Machines</li> <li>Power Quality Guidelines</li> <li>1. Power supplied to medical imaging equipment must be separate from power feeds to air conditioning, elevators, outdoor lighting, and other frequently switched or motorized loads. Such loads can cause waveform distortion and voltage fluctuations that can hinder high quality imaging.</li> <li>2. Equipment that utilizes the facility power system to transmit control signals (especially clock systems) may interfere with medical imaging equipment, thus requiring special filtering.</li> <li>3. The following devices provide a high impedance, nonlinear voltage source, which may affect image quality: Static UPS systems, Series filters, Power conditioners, and Voltage regulators.</li> <li>Do not install such devices at the mains supply to medical imaging equipment without consulting Philips installation or service personnel.</li> <li>4. Line impedance is the combined resistance and inductance of the electrical system and includes the impedance of the power source, the facility distribution system, and all phase conductors between the source and the imaging equipment. Philips publishes recommended conductor sizes based on equipment power requirements, acceptable voltage drops, and assumptions about the facility source impedance. The minimum</li> </ul>	<ol> <li>The contractor will supply &amp; install all breakers, shunt trip and incoming power to the breakers. The exact location the architect or contractor.</li> <li>The contractor shall supply &amp; install all pull boxes, raceways, conduit runs, stainless steel covers, etc. Conduit/racentire length. A Greenlee pull string/measuring tape (part no. 435, or equivalent) shall be provided with conduit runs.</li> <li>All pre - terminated, cut to length cables, will be supplied and installed by Philips. All cables to the breakers, will be local arrangements.</li> <li>Provide and install 50mm diameter chase nipples between adjacent wall boxes.</li> <li>Electrical raceway shall be installed with removable covers. The raceway should be accessible for the entire length an adequate number of access hatches should be supplied to enable installation of cabling. Approved conduits may be manner that will not allow cables to fall out of the raceway when the covers are removed. In most cases, this will requir removable from the top. Raceway system as illustrated on this drawing are based upon length of furnished cables. Conduit or raceway above - ceiling must be kept as near to finished ceiles. Conduit sizes shall be verified by the architect, electrical engineer or contractor, in accordance with local or Nationa 7. Convenience outlets are not illustrated. Their number and location are to be specified by the customer/architect.</li> <li>Electrical contractor shall install ground bond wires at conduit openings within wall boxes as required by national ar shall be installed in such a way to prevent the inadvertent contact with the installed Philips equipment to maintain Philip</li> <li>Install an insulated stranded ground wire per feeder/conductor size from the Main Disconnect (CB) to the ERB and</li> </ol>

0

I alterations. The site preparation shall be in accordance rical codes, the customer shall be solely responsible for ables, wires, fittings, bushing, etc., As separately specified tight covers. Ducts shall be divided into as many as four llows: Group A: incoming power wiring with associated gnal and/or data wiring and/or cables. Group D: X-Ray I corners. All intersecting points in duct to have cross over						
required by code. (12.0)		FD20 Ceiling	oma City	City, OK	lst Floor	HOSPITAL'S ARCHITECT OR ENGINEER TO USE FOR THE DEVELOPMENT OF
tion of the breakers and shunt trips will be determined by	Project	Allura FD	VA Oklahoma	Oklahoma	-Room 1 1	AL'S ARCHITECT OR E
t/raceways must be free from burrs and sharp edges over its ns.		t 619	р Б Е			ш
vill be supplied and installed by the contractor, subject to	acts	Project Manager: John Wright Contact Number: 7214) 704-8619	john.wright@philips.col		ura Phillips	ABLE FOR USE BY TH
ength. In case of non - accessible floors, walls and ceilings, y be substituted. All raceways will be designed in a equire above - ceiling raceway to be installed with the covers Any changes in routing of raceway system could exceed d ceiling as possible.	<b>Philips Contacts</b>	Project Mana	Email: john.w		Drawn By: Laura Phillips	HILIPS ARE ACCEPT
tional Electrical Codes, whichever govern.					0	ED BY P
ct. aal and local electrical codes. Ground bond wires and lugs hilips isolated ground scheme and maintain patient safety. and from the ERB to the Mains 40E Cabinet (per NEC ring, etc.	Project Details	Drawing Number N-WES140162 C	Date Drawn: 11/14/2014	Quote: 1-ZMC9DA Rev. 2	Order: 6600224430.010000	DRAWINGS AND RELATED INSTRUCTIONS PROVIDED BY PHILIPS ARE ACCEPTABLE FOR USE BY TH
(14.1)		E	EN			E DRAWINGS A

		Electrical Legend					Electrical Legend		
	B Furni C Insta	re				B Furni C Instal	JIE		
		Item Number     Detail Sheet				<u> </u>	Item Number     Detail Sheet		
	$\bigvee$	Description	↓ [			$\downarrow$	Description	$\neg \downarrow  $	
	СВ	480V, 3 phase, Type D 125 A circuit breaker with long-time delay and shunt trip (e.g. Square D HDL36125 or equivalent). Run power from breaker to "MA", leaving an 8' (2440mm) tail at "MA". See Sheet "ED1" for power quality requirements. Location per local code or owner requirements. (Not shown on plan) Shunt Trip (emergency off) - Large mushroom-head button on remote control station with contacts to operate	ED1	E	з _/	Nì	RJ45 type Ethernet 10/100/1000 Mbit network connector with access to customer's network. Locate within 10' (3050mm) of network card. Network fiber optic and Ethernet cabling, connectors, wall boxes, patch panels, etc. are the responsibility of the purchaser. Philips assumes no responsibility for procurement, installation, or maintenance of these components.		
	ST	feature of "CB" (if required by local code or owner, and mandatory for VA and D.O.D installations). (Not shown on plan)		E	з /	N2	RJ45 type Ethernet 10/100/1000 Mbit network connector. Access to customer's network via their remote access server is needed for Remote Service Network (RSN) connectivity.	N1	
,	GE	Local building steel (i.e. structural steel, cold water pipe > 2" (50mm), ground rod). (Not shown on plan)		E	з (	$\oplus_{\mathbf{s}}$	120V/20A dedicated duplex outlet for service in the equipment room. (Not shown on plan)		
5	ERB	Equi-Potential Reference Bar mounted in a 12" (305mm) W x 12" (305mm) H x 4" (105mm) D pull box with hinged cover, surface mounted to the bottom of "WR2" when possible.	ED2	E	з (	φ	120V/20A dedicated duplex outlet for each of the Video Connection Boxes. Verify electrical requirements for customer provided equipment.		
		19 $\frac{1}{4}$ " (490mm) W x 67" (1705mm) H x 4" (105mm) D flanged-edge terminal wall box, surface mounted 75" (1905mm) A.F.F. to top of box. General contractor to cut top and/or bottom of box as required.	ED3	E	(VB		Iocal Philips Service. LCD on cart must remain within 25 (7.62m) of VB9.		bu
' }-	CY /B1	Grommet opening on "WR3". Approximate location shown is recommended and may be changed - verify relocation with local Philips Service.		F	VB.				FD20 Ceiling
	MSA	10" (255mm) W x 10" (255mm) L x 6" (155mm) D floor box, flush mounted with underside of AD7 Adaptation Plate.				H BS	4" (105mm) W x 4" (105mm) H x 4" (105mm) D pull box with removable screw-type cover plate, flush mounted. Exact height to be determined. Location shown is recommended and may be changed - verify relocation with local Philips Service.		Project Allura Fl
5	SP	18" (460mm) W x 18" (460mm) L x 6" (155mm) D ceiling box, flush mounted with removable screw-type cover plate. Provide one 3" (80mm) diameter knockout.		C			4" (105 mm)W x 4" (105 mm)L x 2 1/2" (65 mm)D box for the transformer used for the M LED 3MC Surgical Light. Location to be determined locally. Recommended location is above the ceiling near the surgical light. Verify location with local Philips Service. Hardwire 115V/20A hospital power to "TR".	AD5	
5	TV	18" (460mm) W x 18" (460mm) L x 6" (155mm) D ceiling box, flush mounted with removable screw-type cover plate. Provide a $2\frac{1}{2}$ " (65mm) round cutout (Two $2\frac{1}{2}$ " (65mm) round cutouts are required for systems with two monitor carriages - verify with local Philips Service).		E	з (	мз	4" (105 mm)W x 4" (105 mm)L x 2 1/2" (65 mm)D ceiling box, flush mounted with removable screw-type cover plate. Recommended location is near or above the cable spooler.		John Wright
6	VR1 WR2	10" (255mm) W x 4" (105mm) D wall raceway, surface mounted with removable screw-type cover plate. "WR1" is at finished floor. "WR2" is at 75" (1905mm) A.F.F. to bottom of raceway.	ED3	E		⊕ 	120V/20A dedicated quadplex outlet for "TV2".		s Contacts t Manager: John Wrig
5	WR3	10" (255mm) W x 4" (105mm) D wall raceway, surface mounted with removable screw-type cover plate. "WR3" is at finished floor. "WR3" may need to be cut at the location of the "CY" connection box.	ED3				UPS - 25 kVA. Universal Power Controller - 25 kVA.	ED4 ED4	Philips Contacts Project Manager:
	РНУ	Stub up point for physiological monitoring cables. Run conduit to customer's physiological console location. Contact manufacturer for power requirements, etc.					Remote Status Panel (wall mounted in the control area) - 4" (105mm) W x 4" (105mm) H x 4" (105mm) D pull box with removable screw-type cover plate, flush mounted. Exact height to be determined. Location shown is		
5	ATY	Auxiliary Box - 6" (155mm) W x 6" (155mm) H x 4" (105mm) D wall box, flush mounted 70" (1780mm) A.F.F. to the bottom of the box with removable screw-type cover plate. Height and location shown are recommended and may be changed - verify height and relocation with local Philips Service.		E		RSP)	recommended and may be changed - verify relocation with local Philips Service.	ED4	( _ !
	WL	Warning Light - Provide a surface or flush mounted light fixture above door to indicate when X-ray is on, if required by local code or physicist of record. (Not shown on plan)	ED3			Т 	Knife switch. On/off 3 phase switch rate 480V 100 A with auxiliary contacts rated 1 phase 120V 4 A. Surface		it Details Ing Number
		Door Switch - 120V/5A switch limited to open when door is open. Mount in upper corner on strike side of main entry door(s) (Cooper no. 1665 or equivalent), if required by local code or physicist of record. See Sheet "ED3" diagram for connection details. (Not shown on plan)	ED3		o (s'	w9	mounted 50" (1270mm) above finished floor to bottom of box.	ED4	Project [ Drawing
							See E1 - E4 sheets for conduit and raceway requirements.		E

THE DRAWINGS AND RELATED CONSTRUCTION DOCUMENTS.



# **Electrical Layout**


Required Unistrut Height: 9' - 6  $\frac{3}{16}$ ", + $\frac{3}{8}$ " / -0 (2900mm, +10mm / -0) Unistrut height measured from finished floor to bottom of Unistrut.



0

ф

Image: Dispect Details       Perfect Details       Perfect Details         Image: Dispect Details       Project Details       Project Details		<b>DHLIPS</b>	
Image: Markage in the propertion of the propert of the propertion of the propertice of th		Project Allura FD20 Ceiling VA Oklahoma City Oklahoma City, OK -Room 1 1st Floor	S ARCHITECT OR ENGINEER TO USE FOR THE DEVELOPMENT OF
HE DRAWINGS AND RELATED INSTRUCTIONS PROVIDED BY		John Wright (214) 704-8619 t@philips.com	PHILIPS ARE ACCEPTABLE FOR USE BY THE HOSPITAL
and Raceway/Conduit - Sheet E2-E4	Refer to Electrical Logand Shoot El	Project Details Drawing Number N-WES140162 C Date Drawn: 11/14/2014 Quote: 1-ZMC9DA Rev. 2 Order: 6600224430.010000	33 AND RELATED INSTRUCTIONS PROVIDED BY F ON DOCUMENTS.
	and Raceway/Conduit - Sheet E2-E4	E1	THE DRAWING





Project Details	Philips Contacts	Project	
Drawing Number N-WES140162 C	Project Manager: John Wright Contact Number: (214) 704-8619	Allura FD20 Ceiling	
Date Drawn: 11/14/2014 Ounter 1-7MC9DA Rev 2	Email: john.wright@philips.com	VA Oklahoma City	
Order: 6600224430.010000	Drawn By: Laura Phillips	-Room 1 1st Floor	

Note: The use of 90 degree ells is not acceptable. Use 45 degree bends at all raceway corners. For conduit runs, use the minimum bending radius specific to the conduit diameter. The use of crossover tunnels at all applicable locations is required. The above mentioned recommendations will help to ensure the integrity of the cables and fiber optic runs.

- \* Countertop Height Guide: 30" (765mm) for standard seated height. 36" (915mm) for standard standing height.
- \* Ensure that the wall junction boxes are mounted perpendicular to the

\* Verify exact ceiling height of Equipment and Control Room Area. \* Architect to coordinate with end users/technicians to determine final placement of control desk components prior to installation in order to avoid rework. Architect to coordinate with Philips Project Manager to reflect final placement on Philips drawings.

Construction         Construction         Norman         Special Requirements         Special Requirements         Special Requirements         Construction	Conduit supp Conduit supp Conduits and Conduit exist Conduit exist Conduit exist	blied/installe blied/installe d cables su ting - cable ting - cable ting - cable	must take n must have a ad by contractor pplied and insta s supplied and i s supplied and i rify with local Ph	- Philips cat - Philips cat lled by contra nstalled by P nilips and ins nstalled by c	t route poir bles installed b bles installed b actor Philips talled by contri- ontractor	y Philips y contractor	* P Power (AC) D Power (DC) G Ground S Signal H High Tension C Cooling Hose A Air Supply Hose		A C B C C C D C E C F C	conduit sup conduit sup conduits an conduit exis conduit exis conduit exis	blied/installed blied/installed d cables sup ting - cables ting - cables ting - cables	d by contractor d by contractor plied and instal supplied and ir	- Philips cab - Philips cab lled by contra nstalled by P nilips and inst nstalled by co	t route point Ig. bles installed by bles installed by actor hilips talled by contra- ontractor	Philips contractor	s P Power (AC) D Power (DC) G Ground * S Signal H High Tension C Cooling Hose A Air Supply Hose	
Novem         Col         1         P         Per NEC. Per NEC.         See conductanguand size share.         C         3         Gas         1         S         2         2         -         Finduar spectra (Paine Monitoring). Velop Windual Pripes Bookshare)         Velop Windual Science (Paine Monitoring). Velop Windual Science (Paine Moninoring). Velop Windual Science (				Туре	Conduit	Conduit	Special Requirements		1	_			Туре	Conduit	Conduit	Special Requirements	
GB:         GB: <td></td> <td></td> <td>1</td> <td>Р</td> <td>Per N.E.C</td> <td>. Per N.E.C.</td> <td>See conductor/ground size chart.</td> <td>С</td> <td></td> <td>MSA</td> <td><b>РНУ</b></td> <td>1</td> <td>S</td> <td>2 <u>1</u>"</td> <td>-</td> <td></td> <td></td>			1	Р	Per N.E.C	. Per N.E.C.	See conductor/ground size chart.	С		MSA	<b>РНУ</b>	1	S	2 <u>1</u> "	-		
Gas         Gas         T         P         2         50'           Gas         Gas         Trad         -         -         -         -         For Injects A anding Box Paint Monitoring, Video           Gas         Gas         Trad         -         -         -         -         For Injects A anding Box Paint Monitoring, Video           Gas         Trad         -         -         -         -         -         -         -         -         -         For Injects A anding Box Paint Monitoring, Video           Keek         Gas         1         P         2         55'         See Sheet *ED2' for details.         A         35'         (Wb)         1         S         1'         B22'           Crit         Using State         1         S         1'         S2'         4''         Tube Cooling Horse.         A         4''         4'''         S'''         1''''         S'''         S'''         S'''         S''''         S''''         S''''         S''''         S''''         S''''         S''''         S''''         S'''''         S''''         S'''''         S'''''         S'''''         S'''''         S'''''         S''''''         S'''''''         S'''''         S'''''''''			1	Р	2 <sup>1</sup> / <sub>2</sub> "	Per N.E.C.		G	32	PHY		1	S	2"	33'		
Ref       A       1       P       2       6         Ref       A       35       Ref       Ref <td></td> <td>  ∕st ∕</td> <td>1</td> <td>Р</td> <td></td> <td>50'</td> <td></td> <td>G</td> <td>33</td> <td></td> <td>Third</td> <td>-</td> <td>-</td> <td>-</td> <td>-</td> <td></td> <td></td>		∕st ∕	1	Р		50'		G	33		Third	-	-	-	-		
Ref       Norm       I       P <td><math>\rightarrow</math></td> <td><math>  \succ  </math></td> <td>1</td> <td>Р</td> <td><u>3</u>"</td> <td>6'</td> <td></td> <td>G</td> <td></td> <td>Third</td> <td></td> <td>-</td> <td>G</td> <td>-</td> <td>-</td> <td>For Injector, Auxiliary Box, Patient Monitoring, Video</td> <td></td>	$\rightarrow$	$  \succ  $	1	Р	<u>3</u> "	6'		G		Third		-	G	-	-	For Injector, Auxiliary Box, Patient Monitoring, Video	
With With With With 1       P       2*       65       A       37       Visit With 1       S       1*       S       1*       S2       1*       S2       1*       S2*       1*       S2* <th< td=""><td><math>\mid \succ \mid</math></td><td>Room</td><td>1</td><td>Р</td><td><u>3</u>"</td><td>-</td><td>See Sheet "ED2" for details.</td><td>-</td><td></td><td>1.7</td><td>  -&gt;=-&lt; - </td><td>- · · · <u></u> · · ·   1</td><td></td><td> ·····   1"</td><td>82'</td><td></td><td></td></th<>	$\mid \succ \mid$	Room	1	Р	<u>3</u> "	-	See Sheet "ED2" for details.	-		1.7	->=-< -	- · · · <u></u> · · ·   1		·····   1"	82'		
ATT       BD       1       S       2       55'       41'	$  \cdot \rangle \rightarrow \langle \cdot  $	7-1-1-	1	· · · P	<u>3</u> "	· · · · _ 55'		A	36	VB2	MB	1	S	1"	82'		
NMA       1       S       2 <sup>1</sup> / <sub>2</sub> 44       Tube Cooling Hoses.       A       33       VMA       VMA       1       S       1'       B2'         GP       VMC       1       S       1'       B2'       C       1 <sup>1</sup> / <sub>2</sub> 44       Tube Cooling Hoses.       A       43       66'       NMB       1       S       1'       B2'         GP       VMC       1       P/G       1 <sup>1</sup> / <sub>2</sub> 39'       High Tension Cables.       A       44       VMD       1       S       1'       B2'         GP       VMC       1       P/G       2 <sup>1</sup> / <sub>2</sub> 39'       High Tension Cables.       A       44       VMD       1       S       1'       B2'       91'         GP       VMC       1       S       2 <sup>1</sup> / <sub>2</sub> 39'       MMD       1       S       1'       S2'       91'         GP       VMC       1       S       2 <sup>1</sup> / <sub>2</sub> 39'       MMD       1       S       2 <sup>1</sup> / <sub>2</sub> 91'         GP       VMD       1       S       2 <sup>1</sup> / <sub>2</sub> 39'       91'       91'       5       91'         GP       VMD       1       S       <	$\mid \succ \mid$	$  \succ  $	1	S	<u>3</u> "	55'		A	37	$\mid$	$  \succ  $	1	S	1"	82'		
VT       TN       1       S       1/2       65       65       1/2       67       1/2 <th< td=""><td>ATY</td><td><math>  \succ  </math></td><td>1</td><td>S</td><td>2 <sup>1</sup>/<sub>2</sub>"</td><td>41'</td><td></td><td>A</td><td>38</td><td><math>\mid</math></td><td><math>  \succ  </math></td><td>1</td><td>S</td><td>1"</td><td>82'</td><td></td><td></td></th<>	ATY	$  \succ  $	1	S	2 <sup>1</sup> / <sub>2</sub> "	41'		A	38	$\mid$	$  \succ  $	1	S	1"	82'		
SP         ME         1         P/G         1 <sup>1</sup> / <sub>2</sub> 39'           SP         ME         1         S         1 <sup>1</sup> / <sub>4</sub> 39'         High Tension Cables.         A         41         WB         1         S         1 <sup>1</sup> / <sub>4</sub> 82'           SP         ME         1         H         2 <sup>2</sup> / <sub>2</sub> 37'         High Tension Cables.         A         44         WB         1         S         1 <sup>1</sup> / <sub>4</sub> 82'           SP         ME         1         B         2 <sup>1</sup> / <sub>4</sub> 39'         39'         A         44         WB         1         S         1 <sup>2</sup> / <sub>4</sub> 82'           SP         ME         1         S         2 <sup>1</sup> / <sub>4</sub> 39'         39'         A         44         WB         1         S         2 <sup>1</sup> / <sub>4</sub> 91'           SP         MA         1         S         2 <sup>1</sup> / <sub>4</sub> 91'         A         44         WB         OCY         1         S         2 <sup>1</sup> / <sub>4</sub> 91'           MSD         MA         1         S         2 <sup>1</sup> / <sub>4</sub> 91'         A         48         WB'         1         S         2 <sup>1</sup> / <sub>4</sub> 91'           MSD	$  \rightarrow  $	$  \succ  $	1	S		65'		A		$\mid$	$  \succ  $	1	S	1"	82'		bu
SP       AE       1       PIG       1 <sup>1</sup> / <sub>2</sub> 39         SP       ME       1       S       1 <sup>1</sup> / <sub>2</sub> 39       A       41       (BP)       (MB)       1       S       1 <sup>1</sup> / <sub>2</sub> 82         SP       ME       1       H       2 <sup>1</sup> / <sub>2</sub> 37       High Tension Cables.       A       43       (BP)       (MB)       1       S       1 <sup>1</sup> / <sub>2</sub> 82         SP       MIP       1       PIG       2 <sup>1</sup> / <sub>2</sub> 37       High Tension Cables.       A       44       480       (BP)       1       S       2 <sup>1</sup> / <sub>2</sub> 91'         SP       MIP       1       S       2 <sup>1</sup> / <sub>2</sub> 37       91'       A       44       46       (BB)       (CY)       1       S       2 <sup>1</sup> / <sub>2</sub> 91'         SP       MA       1       S       2 <sup>1</sup> / <sub>2</sub> 37       91'       A       45       6/B       CY)       1       S       2 <sup>1</sup> / <sub>2</sub> 91'         SP       MA       1       S       2 <sup>1</sup> / <sub>2</sub> 39'       91'       A       46       CP)       1       S       2 <sup>1</sup> / <sub>2</sub> 91'         SP       MA       1       S <td>·</td> <td></td> <td>2</td> <td> · ·</td> <td>1 <u>1</u></td> <td>+ · · · · -</td> <td>Tube Cooling Hoses.</td> <td>A</td> <td></td> <td><math>\mid</math></td> <td><math>  \succ  </math></td> <td>1</td> <td></td> <td>1"</td> <td></td> <td></td> <td>Ceilli</td>	·		2	· ·	1 <u>1</u>	+ · · · · -	Tube Cooling Hoses.	A		$\mid$	$  \succ  $	1		1"			Ceilli
A       VIC       I       S       I	$\mid \succ \mid$	$  \succ  $	1		-					$\mid$	$  \succ  $	1		1"			20 0
SP       ME       1       H       2½       37       High Tension Cables.       A       43       VB2       V       1       S       ½       91'	$\mid \succ \mid$	$\rightarrow$	1		_					$\mid$	$  \succ  $	1	-	1"			Ĝ
SP       MP       1       P/G       2*       39'         SP       MP       1       S       2*       39'       A       44       VB3       C/V       1       S       3*       91'         SP       MP       1       S       2*       39'       A       44       VB3       C/V       1       S       3*       91'         SP       MA       1       S       2*       39'       A       44       VB3       C/V       1       S       3*       91'         SP       MA       1       S       2*       39'       91'       A       46       VB3       C/V       1       S       3*       91'         MSA       MP       1       S       2*       39'       91'       A       48       VB7       C/V       1       S       3*       91'       91'         MSA       MP       1       S       2*       39'       91'       A       48       VB7       C/V       1       S       2*       91'       91'         MSA       MP       1       S       2*       55'       Monitor on Cart.       For N LED 3MC Surgical Light.<	$\mid \succ \mid$	$  \succ  $	1				High Tension Cables	-		$  \rangle \rightarrow \langle \rangle$	->=-< -	- · · · <u>-</u> · · · 1	· ·	· · · · · · · · · · · · · · · · · · ·	+ · · · · +		iject I <b>ura</b>
SP       MP       1       S       2'       39'       A       45       VBA       CC       1       S       3'       91''       <	$\mid \succ \mid$	$  \succ$	1							$\mid$	$  \succ  $	1		2 <u>1</u> "			Pro All
SP       (MA       1       S       2"       37	$\mid \succ \mid$	$  \succ  $	1							$\mid$	$  \succ  $	1	-	2 <u>1</u> "			
MSA       MP       1       S       2 ½       39'       A       49       VB8       CV       1       S       ½       91'       Image: Constraint of the constraint	$\mid \succ \mid$	$  \rightarrow  $	1	-	-					$\mid$	$  \succ  $	1	-	2 <u>1</u> "	•		t
MSA       MP       1       S       2 $\frac{1}{2^4}$ 39'       A       49       VB8       Cr/       1       S $\frac{1}{2^4}$ 91'       Image: Construction of the construlicon of	·   · )—( · ·	->	- · · · · · · · · · · · · · · · · · · ·	· ·	· · · · _	+ · · · · -	- · · – · – · · – · · – · · – · · – · · – · · – · · – · · – · · – · · – · · – · · – · · – · · – ·			$\rightarrow$	$  \succ  $	1					Wrig
MSA       MP       1       S       2 ½       39'       A       49       VB8       CV       1       S       ½       91'       Image: Constraint of the constraint	$\rightarrow$	$  \succ$	1							$\mid$	$  \succ  $	1					lohn
IV       IV <th< td=""><td><math>  \rightarrow  </math></td><td><math>  \succ</math></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td><math>\mid</math></td><td><math>  \succ  </math></td><td>4</td><td></td><td>_</td><td></td><td></td><td><b>Contacts</b> Manager: John Wright</td></th<>	$  \rightarrow  $	$  \succ$								$\mid$	$  \succ  $	4		_			<b>Contacts</b> Manager: John Wright
IV       WIP       I       S       2       S2       S2       IV       I	·   · )—( · ·	->		· ·	· · · · _	+ · · · · -	-··-·-	-		$  \rangle \rightarrow \langle v  $		- · · ·	· ·	··· — · · —	91		t <b>onta</b> anag
IV       IV <th< td=""><td><math>\mid</math> <math>\rightarrow</math> <math>\mid</math></td><td><math>\mid \succ</math></td><td></td><td></td><td></td><td></td><td></td><td>-</td><td>· · · —</td><td><math>  \cdot \rangle \rightarrow \langle \cdot  </math></td><td>  <del>- &gt;=-&lt;</del>; -+</td><td>- · · ·</td><td> · ·</td><td> · · · ·</td><td>+ · · - · · +</td><td>For M LED 3MC Surgical Light</td><td>ips C ect M</td></th<>	$\mid$ $\rightarrow$ $\mid$	$\mid \succ$						-	· · · —	$  \cdot \rangle \rightarrow \langle \cdot  $	<del>- &gt;=-&lt;</del> ; -+	- · · ·	· ·	· · · ·	+ · · - · · +	For M LED 3MC Surgical Light	ips C ect M
IV       IV <th< td=""><td><math>\mid \succ \mid</math></td><td><math>\mid \succ</math></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td><math>\mid</math></td><td><math>  \succ  </math></td><td>1</td><td></td><td></td><td>-  </td><td></td><td>Proj</td></th<>	$\mid \succ \mid$	$\mid \succ$								$\mid$	$  \succ  $	1			-		Proj
TV       WM       1       S       3"       65'       For Intercom.       A       54       VB9       MB       1       S       1 <sup>1</sup> / <sub>2</sub> "       55'       Monitor on Cart.         CY       MP       1       S       2"       50'       For Intercom.       A       54       VB9       MB       1       S       1 <sup>1</sup> / <sub>2</sub> "       55'       Monitor on Cart.         CY       MA       1       P/G       1 <sup>1</sup> / <sub>2</sub> "       55'       ST       1       P       3"       150'         CY       MA       1       S       2 <sup>1</sup> / <sub>2</sub> "       55'       B       56       MA       UPC       1       P       3"       150'         MA       WM       1       S       2 <sup>1</sup> / <sub>2</sub> "       55'       S5'       B       56       MA       UPC       1       P       1"       45'       9"       9"       9"       9"       1"       45'       9"       9"       9"       9"       9"       9"       1"       45'       9"       9"       9"       9"       9"       1"       45'       9"       9"       9"       9"       9"       9"       9"       9"       9"       9"       9"	$\mid$ $\rightarrow$ $\mid$	$  \succ$					For Flow/Ginn M	-	· · ·	$  \rangle \rightarrow \langle \rangle$	->=		· ·	· · _ · · _			
CY       MP       1       S       2"       50'         CY       MA       1       P/G       1 ½"       55'         CY       MA       1       S       2½"       55'         MA       1       S       2½"       55'       B       56       MA       UPC       1       P       3"       150'         MA       1       S       2½"       55'       B       56       MA       UPC       1       P       1"       45'         MA       WM       1       S       1"       82'       For optional equipment (IE. Physio Monitor/ Slave       B       59       MA       UPC       1       P       1"       45'	$\mid \succ \mid$	$\mid \succ$								$\rightarrow$	$  \succ  $	1	-				
CY       MA       1       P/G       1 $\frac{1}{2}$ "       55'         CY       MA       1       S       2 $\frac{1}{2}$ "       55'       B       56       MA       UPC       1       P       1"       45'         MA       WM       1       S       1"       82'       For optional equipment (IE. Physio Monitor/Slave       B       59       MA       UPC       1       P       1"       45'	·	->	1	· ·	···	+ · · · · -		-	· · ·	$  \cdot \rangle \rightarrow \langle \cdot  $	∣ →→→∢ →	1 - · · · <u>-</u> · · ·	· ·	· · · ·	+ · · · · +		
$ \langle IV \rangle \langle WRS   2   S   1^{+}$	$  \rightarrow  $	$  \succ$	1							$\mid$		1	-				er s
$ \langle IV \rangle \langle WR \rangle  2   S   1 = 1$	$\rightarrow$	$\mid \succ$	1							$\mid$	$  \rightarrow  $	1					<b>etail:</b> Jumb
$ \langle IV \rangle \langle WRS   2   S   1^{+}$	$  \rightarrow  $	$\mid \succ$	1							$\mid$	$  \rightarrow  $	1	-				sct D 'ing ∧
$ \langle IV \rangle \langle WRS \rangle  /  S   1 =   $	·   · )—{ · ·	->	1   · · · <u>−</u> · ·	· ·	··· ··-	82'	For ontional equipment /IF. Physic Monitor/ Slave			$\mid$	$  \succ  $	1			45'		<b>Proj</b> € Draw
$MSA$ $WR3$ 2 S $1\frac{1}{2}$ - For future options (Patient Monitoring). Verify with local B 60 $MA$ $SWO$ 1 P 1" 45'	$  \rightarrow  $	$\rightarrow$	2	S		-	Monitor)	В		$\mid$	$  \succ  $	1	-	1"	45'		

c <u>e</u>. ₹ 2012. N.< ۵. THE DRAWINGS AND RELATED CONSTRUCTION DOCUMENTS.

[ 1. 2.	All c	onduit	t runs m	nust take n	nost direc		nduit Requ eneral Not		1. 2.	All con	luit runs	must take n	nost direc		nduit Required eneral Notes to point.		
A B C D E F	Conduit Conduit Conduit Conduit Conduit Conduit	supplied supplied s and ca existing existing existing	d/installed d/installed ables supp g - cables s g - cables s g - cables s	by contractor by contractor lied and insta supplied and i	- Philips cal - Philips cal lled by contr nstalled by F nilips and ins nstalled by c	oles installed by oles installed by actor 'hilips talled by contractor	Philips contractor	<ul> <li>P Power (AC)</li> <li>D Power (DC)</li> <li>G Ground</li> <li>* Signal</li> <li>H High Tension</li> <li>C Cooling Hose</li> <li>A Air Supply Hose</li> </ul>	A B C D E F	Conduit sup Conduit sup Conduits an Conduit exis Conduit exis Conduit exis	blied/installe blied/installe d cables su ting - cables ting - cables ting - cables	ed by contractor ed by contractor oplied and insta s supplied and i	r - Philips cat r - Philips cat alled by contra installed by P hilips and ins installed by c	oles installed by oles installed by actor hilips talled by contra- pontractor	Philips contractor	* P Power (AC) D Power (DC) G Ground * Signal H High Tension C Cooling Hose A Air Supply Hose	
Run	Cor			Conduit Quantity	Cable Type (*)	Minimum Conduit Size	Maximum Conduit Length	Special Requirements	Rur		t To	Conduit Quantity	Cable Type (*)	Minimum Conduit Size	Maximum Conduit Length	Special Requirements	
No. 61	sv	$\neg +$		1	P	2"	15'	Flex Conduit	/ No					Size			-
62		$\prec$ $\mid$ )	swo	1	P	2"	15	Flex Conduit									
63		$\prec$ $\mid$ )	swo	1	P	2"	15'	Flex Conduit									
64		v v Sol (	UPS	1	Р	1"	15'										
65		vc) (	UPS	1	Р	1"	15'										
66		v v v	UPS	1	Р	1"	15'										
67		» م	UPS	1	Р	1"	15'										
68	R	SP) (	UPS	1	Р	1 <u>1</u> "	250'	Remote Status Panel.									
																	Proiect
																	Philips Contacts
																	Project Details

 Quote:
 1-ZMC9DA Rev. 2
 Oklahoma City, OK

 Order:
 6600224430.010000
 Drawn By: Laura Phillips
 Oklahoma City, OK

 THE DRAWINGS AND RELATED INSTRUCTIONS PROVIDED BY PHILIPS ARE ACCEPTABLE FOR USE BY THE HOSPITAL'S ARCHITECT OR ENGINEER TO USE FOR THE DEVELOPMENT OF CONSTRUCTION DOCUMENTS.
 OKLAHOMA CITY, OK


Pov	ver Quality Requirements (Mains 40E Cabinet)
Power Output	100kW
Supply Configuration	3 phase, identical 3 wire power and isolated unity ground with bonding conductor, delta (preferred) or wye
Nominal Line Voltage	480 VAC, 60 Hz
Line Voltage Variation	Voltage Variations are never to exceed $\pm 10\%$ when measured using 10 minute mean RMS values with a measurement window of 1 week. At least 95% of all measured 10 minute mean RMS values shall be within $\pm 5\%$ of the configured nominal voltage.
Line Voltage Balance	2% maximum of nominal voltage between phases
Frequency Variation	± 1.0 Hz
Voltage Surges	To 110% of steady-state voltage 100 msecs. Maximum duration, 6 per hour max.
Voltage Sags	To 90% of steady-state voltage 100 msecs. Maximum duration, 6 per hour max.
Line Impulses	1000 VPK above phase-neutral RMS absolute maximum. No more than 1 impulse per hour to exceed 500 VPK.
Neutral-Ground Voltage	2.0 volts maximum RMS value
Neutral-Ground Impulses	No more than 1 per hour that exceeds 25 volts and 1 Mjoule
High Frequency Noise	3.0 volts steady-state maximum. Over 3.0 volts permitted for 100 msec. maximum, 1 per hour max.
Grounded Conductor Impedance	0.1 Ohms @ 60 hz. maximum

Branch Power	225 kVA
Max. Standby Current	8 A @ 3mA, 100 kVP continuous
Circuit Breaker (CB)	3 phase, Type D 125 A with long-time delay and shunt
Recommended conductor/ground sizes for 1% impe	responsible for calculating conductor/ground sizes. dance of supply conductors to circuit breaker (CB). pper conductors:
Nominal Line Voltage (in VAC) (60 Hz)	480
1/0 AWG	76.92ft
2/0 AWG	96.74ft
3/0 AWG	121.95ft
4/0 AWG	155.34ft
250 KCM	181.82ft
300 KCM	217.98ft
400 KCM	294.12ft
Max. Instantaneous Power (1000mA @ 100 kVP)	249 KVA
Max. Inst. Current @ CB (RMS value over half-cycle)	300 A
Max. Phase-phase impedance @ CRC	0.465 Ω
Max. Load Voltage Drop @ CB (RMS value over half-cycle)	139.5 V
Output Voltage Mains 40E Cabinet	480 VAC ± 10%
Max. Inst. Current @ Mains 40E Cabinet output (RM value over half-cycle)	300 A
Max Phase-phase impedance @ Mains 40E Cabinet CRC input terminal	0.545 Ω
Max. Load Voltage Drop @ Mains 40E Cabinet output	163.5 V

Ω¥

### Branch Circuit and Wire Gauge Requirements (Mains 40E Cabinet)

Project	Allura FD20 Ceiling			
Philips Contacts	Project Manager: John Wright	Contact Number: (214) 704-8619	Email: John.wright@philips.com	
Project Details	Drawing Number	N-WES140162 C	Outote: 1-ZMC9DA Rev 2	Drder: 6600224430.010000



copyrigl

the


or in

whole ion in

g éd.


All rights

nics N.V. 2012.



rough Divisions o scale) ed by metal barriers into four sections: iated PE. ociated PE. stands. s.		HLIPS
High-Voltage Signal, Data and Video (if not in conduit) are placed in the appropriate trough and les from one division cross cables from must be continuous from the beginning. teel dividers grounded to building ground. estraints in all troughs.		
	(12.0)	Project DetailsPhilips ContactsProjectDrawing NumberDrawing NumberProject Manager: John WrightProjectDrawing NumberProject Manager: John WrightProjectN-WES140162 CContact Number: (214) 704-8619Allura FD20 CeilingDate Drawn: 11/14/2014Email: john.wright@philips.comVA Oklahoma City, OKOrder: 6600224430.010000Drawn By: Laura PhillipsPoklahoma City, OK
n be configured by n off the desired its everytime es or the exam be set to turn by the Philips rify with customer om lights are		Philips Contacts Project Manager: John Wright Contact Number: (214) 704-8619 Email: john.wright@philips.com Drawn By: Laura Phillips
Use" indicator e energized ng exposures.		Project Details Drawing Number N-WES140162 C Date Drawn: 11/14/2014 Quote: 1-ZMC9DA Rev. 2 Order: 6600224430.010000
	(14.0)	ED3 8.20.14

**JCTION DOCUMENTS** 



THE DRAWINGS AND RELATED INS CONSTRUCTION DOCUMENTS.

## Philips Healthcare Remote Services Network (RSN)


Secure broadband connection required for Philips remote technical support, diagnostics, and applications assistance

### Broadband Site-to-Site Connectivity (Preferred)

This connectivity method is designed for customers who prefer a connection from the RSN Data Center to the Health Care Facility (HCF) utilizing their existing VPN equipment.

### **Connectivity Details:**

- A Site-to-Site connection from the RSN data center's Cisco router will be established to the HCF's VPN concentrator.
- The VPN Tunnel will be an IPSEC, 3DES encrypted Tunnel using IKE as standard, but alternative standards are also available, such as AES, MD5, SHA, Security Association lifetime and Encryption Mode.
- Every system that we will be servicing remotely will have a static NAT IP that we configure on the RSN Data center side.



### Action Required by Hospital

- Review and approve connection details.
- Complete appropriate Site Checklist.

- Configure and allow Site-to-Site access prior to setting up connectivity depending on the access criteria that the HCF decides to implement (ex: Source IP filtering, destination IP

filtering, NAT assignment, etc.).

- Route traffic from within the hospital network with destination addresses 192.68.48.0/22 to the designed IP provided by Philips.

### Broadband Router Installed at Health Care Facility


This connectivity method is designed for customers who have a dedicated high speed connection for Philips equipment.

### **Connectivity Details:**

- An RSN Cisco 1711 or 1712 router will be preconfigured and installed at the HCF by Philips in conjunction with the HCF IT representative.
- The VPN Tunnel will be an IPSEC, 3DES encrypted Tunnel using IKE and will be established from the RSN-DC and terminated at the RSN Router on-site.
- One to One NAT is used to limit access to Philips equipment only.
- Router Config and IP auditing is enabled for Customer IT to view via website 24/7. - Dedicated DSL connections are also supported.

### **Option 1: Parallel to HCF Firewall Connectivity Method**

This connectivity method is designed for customers who prefer a Philips RSN Router installed on site utilizing all the security features provided and managed by Philips.




### Action Required by Hospital:

- Assign a fixed public IP Address from the ISP to be configured on the Philips router. This is the DOTTED link on the picture connected to the firewall.
- Assign a Back end IP for the Philips router on the Hospital Network.
- Complete appropriate Site Checklist.

- Route traffic from within the hospital network with destination addresses 192.68.48.0/22 to internal Philips router Ethernet interface. This is the DASHED line connected to the firewall.


This connectivity method is designed for customers who prefer a Philips RSN Router installed equipment.



### Action Required by Hospital:

- Assign a fixed public IP Address from the ISP to be configured on the Philips router. This is the DOTTED link on the picture connected to the firewall. - Assign a Back end IP for the Philips router on the Hospital Network. - Complete appropriate Site Checklist. - Route traffic from within the hospital network with destination addresses 192.68.48.0/22 to internal Philips router Ethernet interface. This is the DASHED line connected to the firewall. - Configure and allow on the firewall on the DASHED line interface access between the IP address allocated by the hospital to the Philips internal Ethernet router interface and the target modality IP address.

Option 3: Router Installed Inside the HCF's DZM and existing, or new DMZ, allowing access to Philips equipment.



### Action Required by Hospital:

- Assign a fixed public IP Address from the ISP to be configured on the Philips router. This is the DOTTED link on the picture connected to the firewall. - Assign a Back end IP for the Philips router on the Hospital Network. - Complete appropriate Site Checklist. - Route traffic from within the hospital network with destination addresses 192.68.48.0/22 to internal Philips router Ethernet interface. This is the DASHED line connected to the firewall. - Configure and allow on the firewall on the DASHED line interface IPSec protocol communication by opening protocol 500, 50, 51, 47 and port 23 + TACACS. Traffic should be between external IP Address located on the Philips router and the RSN Data center IP address 192.68.48/24 and IP address AOSN TACAS. - Configure and allow on the firewall on the DASHED line interface access between the IP address allocated by the hospital to the Philips internal Ethernet router interface and the target modality IP address.

### **Option 2: Back End Connected to the HCF Firewall Connectivity Method**

on site by setting up an IP-Based policy allowing access thru existing HCF Firewall to Philips

This connectivity method is designed for customers who prefer the RSN Router installed inside



(12.0)

# System Network Information MPORTANT NOTE: It is the customer's responsibility to coordinate with the local Philips Engineer to provide ALL required network information and install ALL required network cabling & drops according to Philips specifications PRIOR to the scheduled installation start date. Failure to do so may delay system installation and jeopardize the customer hand over date.

Allura	IP Sec [ ]yes [ ]no	XperIM	IP Sec [	]yes [ ]	no			Time Synchronizatio	on
Physical Location:			Location	n 1	Location 2	2	Location3	Physical Location:	
Hostname:		Physical Location:						Server Name:	
MAC Address:		Hostname:						RIS	Physical Location
IP Address		MAC Address:							Basic Local RIS
letmask:		IP Address						Hostname:	
Gateway:		Netmask:						IP Address:	
AE Title:		Gateway:						AE Title:	
Port Number (5101):		AE Title:						Max PDU Size:	16384 or
XtraVision	IP Sec [ ]yes [ ]no	Port Number (3010):						Port Number:	
Physical Location:		Remote Software In	stallation (R	RPS)				Secure Node:	
Hostname:		Enable Distribution:		[	]yes [	] no		Encryption:	
MAC Address:		Enable Installation:		[	]yes [	] no		Certificate Name:	
? Address		Dicom Printer	1	-	-			PPSM IHE	
letmask:			Location 1	Locati	on 2 Lo	ocation3	Location 4	Compatible:	
Bateway:		Physical Location:						Time Synchronizatio	on
E Title XtraVision:								Allura Xper:	20/21(ftp), 80(htt
ort Number (3110):		Hostname:							9903(fsf.net)
E Title for X-Ray Mod:		IP Address						Allura CV20:	20/21(ftp), 80(htt
for X-Ray Modality:		AE Title:						XtraVision:	20/21(ftp), 80(htt 5900(vnc), 9905
P Navigator	IP Sec [ ]yes [ ]no	Port Number :						EP Navigator (R3):	20/21(ftp), 443(h
nysical Location:		PACS	Physical Lo		01 /				20/21(ftp), 80(htt
ostname:			Store/ Import 1	Store/ Import 2	Store/ Export	Query Retriev		EP Cockpit (R1.2):	9903(fsf.net)
AC Address:		Hostname:						CX50:	
Address		IP Address							
tmask:		AE Title:						Xper IM:	
eway:		Port Number :						View Forum	
Title:		PACS	Physical Lo	cation:					
t Number:			Store/ Import 1	Store/ Import 2	Store/ Export	Query, Retriev		Hospital Network	M2M Server
ew Forum	IP Sec [ ]yes [ ]no	Hostname:							(PRS)
sical Location:		IP Address						Scheme (https):	
stname:		AE Title:						IP Address (192.68.49.50):	
C Address:		Port Number :						Portnumber (443):	
ddress		Audit Trail						Use Proxy Server:	[]yes []no
mask:		Physical Location:						IP Address	
ateway:		Hostname:							
Title:		IP Address						Port Number:	
		AE Title:						User Name:	
Port Number:		Port Number :						Password:	

Locatio	n:	
cal RIS	WLM	MPPS
	[ ]yes [ ]no	[]yes []no
	[]yes []no	[]yes []no
		[]yes []no

, 80(http), 443(https), 5900(vnc),

, 80(http), 4440(fsf)

, 80(http), 443(https), 5660(ist/ice), ), 9905(lots)

, 443(https), 5660(ist/ice), 9055(lots)

١.	80(http).	443(https),	5900(vnc).
,	oo(incp),	110(inapo),	0000(1110),

rver )	Proxy	ePO Server (PRS)
] no		

	Project Details	Philips Contacts	Project	
	Drawing Number	Project Manager: John Wright	Allura FD20 Ceiling	
N	N-WES140162 C	Contact Number: (214) 704-8619		
2	Date Drawn: 11/14/2014	Email: john.wright@philips.com	VA Oklahoma Citv	
)	Quote: 1-ZMC9DA Rev. 2		Oklahoma Citv. OK	
	Order: 6600224430.010000	Drawn By: Laura Phillips	-Room 1 1st Floor	

(13.0)

### Instructions

he ъ

part is I or in whole .⊑

Å ēd.

onics N.V. 2012. All rights

© Koninklijke Philips Electr

This form is to be used by Project Manager, Contractor and Service Engineer.

Information is used to develop and determine site ready date.

Items listed are go/no go items for delivery unless noted as delay only items.

Items identified with \*\*\* as delayed items must be completed after hours or on weekend. These items cannot be accomplished while installation is in progress. Also, these items must be completed within two days of installation start or they may stop installation.

Site Readiness Checklist
Modality:
Order:
Site Name:
Location:
Contact Name:
Contact Phone Number
Customer site preparation verified in general against the Philips final planning drawings.
☐ Walls finished including painting.
Doors installed.
☐ Floor leveled according to Philips drawings and specifications.
☐ Floors are tiled/covered finished. Flooring is covered with protective covering (scratch protection).
Ceiling lights installed.
Cable conduit and ductwork installed and clean. Position checked. Duct covers in place but not finally closed. Cable opening are clear, without sharp edges. Pull strings in conduit. Installation per Philips specifications.
HVAC environmental equipment installed and working according to Philips specifications.
Ceiling installation completed.
Electrical preparation according to Philips specifications.
All network cabling, drops installed according to Philips specifications (including hardcopy cameras).
All pre-cabling identified on Philips drawings has been installed.
Pre-move survey completed - Delivery route identified.
Lead glass installed ***.
X-Ray warning lights installed ***.
Dedicated phone line for modem use***.
Room has been cleaned ***.
Cabinets and casework installed (with insulation and building steel) according Philips specifications***.
RSN survey completed and submitted
Philips RSN Champion contacted.
Approved for Delivery
Project Manager
Service Engineer
HIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

### Items Specific for the Cardio/Vascular Modality

Date

Date

Unistrut installed and level according to Philips specifications.

Floor plates installed and level according to Philips specifications.

All cover plates have holes punched and nipples required and bushings installed.

Emergency power requirements installed according to Philips specifications.

Building steel ground installed to ECB section of ERB.

Non-Philips provided room electrical equipment grounds installed to PCB middle section of ERB.

Conduit lengths measured according to Philips specifications. Note: Specifications is from source box to destination box (not just conduit run length).

Routing of ductwork and conduits must be installed according to Philips specifications.

	Project Details	Philips Contacts	Project	
С	Drawing Number	Project Manager: John Wright	Allura FD20 Ceiling	
; <b> </b>	N-WES140162 C	Contact Number: (214) 704-8619		
<b>  </b>	Date Drawn: 11/14/2014	Email: john.wright@philips.com	VA Oklahoma Citv	
<b>K</b>	Quote: 1-ZMC9DA Rev. 2		Oklahoma City, OK	
) 14	Order: 6600224430.010000	Drawn By: Laura Phillips	-Room 1 1st Floor	
THE DRAWING: CONSTRUCTIO	THE DRAWINGS AND RELATED INSTRUCTIONS PROVIDED BY CONSTRUCTION DOCLIMENTS	PHILIPS ARE ACCEPTABLE FOR USE BY THE HOSPITAL	THE DRAWINGS AND RELATED INSTRUCTIONS PROVIDED BY PHILIPS ARE ACCEPTABLE FOR USE BY THE HOSPITAL'S ARCHITECT OR ENGINEER TO USE FOR THE DEVELOPMENT OF CONSTRUCTION DOCUMENTS.	