


THIS LAYOUT IS NOT REFLECTIVE OF THE LATEST SALES QUOTE AVAILABLE AT THE TIME OF DRAWING. IT HAS BEEN MODIFIED TO REFLECT ANTICIPATED CHANGES TO YOUR CONFIGURATION. IF REQUESTED CHANGES ARE NOT IMPLEMENTED A REVISION WILL BE REQUIRED. PLEASE ENSURE THAT ANY CONFIGURATION CHANGES ARE COMPLETED PRIOR TO REQUESTING FINAL DRAWINGS.


XF (TRANSFORMER) MAY BE LOCATED IN THE CEILING BY THE CUSTOMER/CONTRACTOR.

SYMBIA INTEVO SYSTEM REQUIRES A MINIMUM FINISHED CEILING HEIGHT OF 8'-0".

COLLIMATOR CART PLACEMENT TO BE DETERMINE BY THE CUSTOMER WHEN LOCATED OUTSIDE THE SCANNER ROOM.



HEIGHT OF WINDOW TO BE COORDINATED WITH COUNTERTOP/DESK HEIGHT.



B
C

SCALE: 3/16" = 1'-0"

PROJECT MANAGER: TIM SADRAKULA
PHONE: (845) 270-3659
EMAIL: TIMOTHY.SADRAKULA@SIEMENS.COM

PROVIDENCE VAMC

830 CHALKSTONE AVENUE, PROVIDENCE, RI 02908

SIEMENS

C 02/15/18 REDRAWN WALL BACKGROUND PER CAD 'VAMC PROV REV BASE 020518.DWG.'

SHEET OF
1 2

DATE:
05/25/17


DRAWN BY:
M. COPLEY

REF:
30193891

SCANNER ROOM 284 - SYMBIA INTEVO 16
REFERENCE DOCUMENT - NOT FOR CONSTRUCTION

1702183

EQUIPMENT LEGEND

NO	DESCRIPTION	SMS SYM	WEIGHT (LBS)	BTU/HR TO AIR	DIMENSIONS (INCHES)			REMARKS
					W	D	H	
①	SYMBIA INTEVO 16 GANTRY W/ COLLIMATORS	ⓑ	7,413	24,574	93	84 1/2	90 1/2	6,826 BTU ON STANDBY. WORST CASE WEIGHT 7,963 LBS. WITH (2) HIGH ENERGY COLLIMATORS AT 275 LBS. EACH
②	FRONT PHS	ⓔ	2,512	—	31 1/8	97 1/2	23 3/16	MAXIMUM HEIGHT 41"
③	UNDER THE FLOOR PHS CABLE STANDARD	ⓔ	—	—	—	—	—	UNDER FLOOR
④	PHS STANDARD PIVOT 	ⓔ	—	—	—	—	—	STANDARD PIVOT — 33 DEGREES
⑤	REAR PHS WITH SNAC	ⓔ	415.2	—	—	—	—	ON FLOOR
⑥	PATIENT BOOM SWING ARM	ⓔ	—	—	—	—	—	
⑦	LINE CONNECTION BOX	Ⓛⓑ	245	1,365	29 1/2	11 3/4	32	ON FLOOR
⑧	UPS FOR SPECT	ⓈⓅ	124	1,024	10	28 3/8	17 7/8	ON FLOOR
⑨	TRANSFORMER REQUIRED	ⓍⓉ	—	—	—	—	—	CUSTOMER SUPPLIED PRIOR TO INSTALLATION. SEE POWER SCHEDULE
⑩	COLLIMATOR CART (EMPTY)	ⓔ	400	—	47 3/8	32 5/8	47 1/2	WORST CASE 1330 LBS. WITH 1 SET HE AND 1 SET ME
⑪	IMAGE CONSTRUCTION SYSTEM FOR SYNGO MI (ACQUISITION) WORKPLACE	ⓈⓈ	50	2,389	8	22	18	OFF FLOOR.
⑫	IMAGE RECONSTRUCTION SYSTEM FOR SYNGO MI (ACQUISITION) WORKPLACE	ⓈⓈ	47	*	8	22	18	OFF FLOOR. *BTU'S INCL. WITH ICS
⑬	CONTROL AND KEYBOARD	ⓔ	—	—	—	—	—	ON CUSTOMER'S COUNTER
⑭	18" MONITOR	ⓔ	31	—	18 3/8	2 5/8	14 13/16	ON CUSTOMER'S COUNTER
⑮	SYNGO MI (ACQUISITION) WORKPLACE UPS FOR IMS STANDARD COMPONENT	ⓈⓈ	71	—	5	19	17 1/4	OFF FLOOR.
⑯	DEDICATED RECONSTRUCTION SYSTEM WORKPLACE — STANDARD INTEVO 16	ⓈⓈ	55	*	8	20 3/4	18	OFF FLOOR. *BTU'S INCL. WITH ICS
⑰	UPS FOR E.SOFT/C.CAM	ⓔ	45	—	6 3/4	17 1/2	8 1/2	ON FLOOR

ARCHITECTURAL NOTES

- 1) ALL PRELIMINARY EQUIPMENT LAYOUTS SUBMITTED BY SIEMENS HEALTHCARE ARE BASED ON THE RECOMMENDED SPACE NECESSARY FOR THE OPERATION AND SERVICEABILITY OF THE EQUIPMENT BEING PROPOSED. SIEMENS WILL NOT SUBMIT AN EQUIPMENT LAYOUT THAT IS NOT IN THE BEST INTEREST OF BOTH THE CUSTOMER AND SIEMENS. ALL EQUIPMENT LAYOUTS ARE BASED EITHER ON AN ACTUAL SITE SURVEY OR ARCHITECTURAL DRAWINGS SUPPLIED TO SIEMENS. SIEMENS WILL NOT BE RESPONSIBLE FOR ANY ALTERATIONS THAT ENCR OACH WITHIN DESIGNATED SAFETY AND SERVICE CLEARANCE ZONES AS INDICATED ON DRAWINGS (I.E., PIPE CHASES, VENTILATION DUCTS, CASEWORK, AND SOFFITS, ETC.) MADE BY THE CUSTOMER OR REQUIRED BY A CUSTOMER'S ARCHITECTURAL FIRM ONCE PRELIMINARY DRAWINGS HAVE BEEN SUBMITTED AND APPROVED. DO NOT ALTER ANY SPECIFICATIONS AND/OR DIMENSIONS WITHOUT CONTACTING AND RECEIVING WRITTEN CONFIRMATION FROM SIEMENS PROJECT MANAGER.
- 2) SIEMENS HEALTHCARE IS NOT AN ARCHITECTURAL OR ENGINEERING FIRM. DRAWINGS SUPPLIED BY SIEMENS ARE NOT CONSTRUCTION DRAWINGS. THEREFORE, THESE DRAWINGS ARE TO BE USED ONLY FOR INFORMATION TO COMPLEMENT ACTUAL CONSTRUCTION DRAWINGS AVAILABLE FROM A CUSTOMER APPOINTED ARCHITECTURAL REPRESENTATIVE OR A CUSTOMER'S ENGINEERING DESIGN GROUP. THE CUSTOMER'S ARCHITECT AND GENERAL CONTRACTOR SHALL BE ULTIMATELY RESPONSIBLE FOR COMPLIANCE WITH ALL APPLICABLE CODES AND PROFESSIONAL DESIGN REQUIREMENTS INCLUDING OSHA/NEC SAFETY CLEARANCE REQUIREMENTS IN ADDITION TO SIEMENS—REQUIRED SAFETY/SERVICE CLEARANCES SHOWN.
- 3) THE CUSTOMER IS RESPONSIBLE FOR ALL ROOM AND AREA PREPARATION COSTS, PROFESSIONAL FEES, PERMITS, REPORTS, AND INSPECTION FEES.
- 4) EQUIPMENT WARRANTIES, EXPRESSED OR IMPLIED ON THE PART OF SIEMENS SHALL BE CONTINGENT UPON STRICT COMPLIANCE WITH THE ARCHITECTURAL, STRUCTURAL, ELECTRICAL, MECHANICAL AND RECOMMENDATIONS AND REQUIREMENTS CONTAINED IN THESE DRAWINGS, UNLESS SPECIFIED OTHERWISE.
- 5) ALL DIMENSIONS SHOWN ARE FROM FINISHED SURFACES UNLESS SPECIFIED OTHERWISE.
- 6) THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST. ACTUAL PROTECTION REQUIREMENTS SHALL BE SPECIFIED BY A REGISTERED RADIATION PHYSICIST AT CUSTOMER'S ENGAGEMENT AND EXPENSE. RESPONSIBILITY FOR ALL INFORMATION AS TO THE ROOM LOCATION, USE, AND NUMBER OF ANTICIPATED EXAMINATIONS TO BE PERFORMED PER TIME PERIOD SHALL BE PROVIDED TO THE PHYSICIST BY THE CUSTOMER. THE CUSTOMER SHALL FURTHER TAKE ALL RESPONSIBILITY IN THE COMMUNICATION AND COORDINATION OF ACTIVITIES OF THE RADIATION PHYSICIST AND THE ARCHITECTURAL REPRESENTATIVE.
- 7) SIEMENS HEALTHCARE SHALL BE RESPONSIBLE FOR SIEMENS EQUIPMENT INSTALLATION, CALIBRATION, CONNECTION AND INSTALLATION OF SIEMENS PROVIDED CABLES. THE CUSTOMER/ELECTRICAL CONTRACTOR IS RESPONSIBLE FOR TERMINATIONS OF CUSTOMER/ELECTRICAL CONTRACTOR—SUPPLIED CABLES TO SIEMENS EQUIPMENT. IN THE EVENT THAT SPECIFIC TRADE RULES OR LICENSE REQUIREMENTS PROHIBIT THIS, THE CUSTOMER SHALL INITIATE THE SERVICES OF APPROVED OTHER CONTRACTORS AND PAY FOR SELECTED, APPROVED PARTIES TO PERFORM THIS WORK WITH SUPERVISION PROVIDED BY SIEMENS. CALIBRATION WHEN ACCOMPLISHED OUTSIDE OF NORMAL INSTALLATION SEQUENCES DUE TO CONTRACTOR OR TRADE RULE ACTIONS OR REQUIREMENTS SHALL BE SUPPORTED BY, CHARGED TO, AND ACCEPTED BY THE CUSTOMER AS AN ADDITIONAL INSTALLATION EXPENSE.
- 8) THE CUSTOMER SHALL COORDINATE WITH SIEMENS PROJECT MANAGER THE LOCATIONS AND TRAVEL OF ALL ANCILLARY EQUIPMENT TO BE CEILING OR WALL MOUNTED (I.E.: O.R. LIGHTS, MEDICAL GAS COLUMNS, PHYSIOLOGICAL MONITORING INJECTORS, CRT PLATFORMS, SPRINKLER HEADS, SMOKE DETECTORS, ELECTRICAL OUTLETS, HVAC GRILLES, SPEAKERS, AND GENERAL ROOM LIGHTING, ETC.).
- 9) THE GENERAL CONTRACTOR/CUSTOMER SHALL BE RESPONSIBLE FOR ALL FINAL PAINT, TOUCH-UP AND ANY COSMETIC OR TRIM WORK WHICH NEEDS TO BE OR IS REQUIRED TO BE COMPLETED AFTER THE INSTALLATION OF THE SIEMENS EQUIPMENT AND ANY ASSOCIATED SUPPORT APPARATUS.

FOR ADDITIONAL PLANNING RELATED QUESTIONS, PLEASE REFER TO TYPICAL DRAWING SET OR SIEMENS PROJECT MANAGEMENT.

- RESTRICTIONS:
- 1) SERVICE CLEARANCE HAS NOT BEEN MAINTAINED IN SOME AREAS AS INDICATED ON THE PLAN.
 - 2) THIS DIMENSION IS NOT ADEQUATE SPACE FOR A THROUGH WALKING PATH AS INDICATED ON THE PLAN.

WARNINGS:

THERE IS NO VIEW DOWN THE BORE OF THE GANTRY FROM THE CONTROL ROOM AREA.

Authorization:

I have reviewed the below referenced plan(s) and understand and agree with the proposed room configuration and equipment arrangement. By signing below, I agree that I have the authority to sign on behalf of the customer. I have been fully informed by Siemens about the limitations (if applicable) that exist regarding the proposed equipment installation. Siemens has explained that the listed limitations may interfere with certain functions of the equipment. We agree to assume responsibility to inform employees, agents, contractors, users or any other persons handling this equipment of these limitations and will assure that the equipment is handled in compliance with these parameters.

Signature: _____ Date: _____
Name Printed: _____ Title: _____

PROJECT MANAGER: TIM SADRAKULA
PHONE: (845) 270-3659
EMAIL: TIMOTHY.SADRAKULA@SIEMENS.COM

PROVIDENCE VAMC

830 CHALKSTONE AVENUE, PROVIDENCE, RI 02908

C 02/15/18 REDRAWN WALL BACKGROUND PER CAD 'VAMC PROV' REV BASE 020518.DWG.' SHEET OF 2 2 DATE: 05/25/17 DRAWN BY: M. COPLEY REF: 30193891

SCANNER ROOM 284 - SYMBIA INTEVO 16
REFERENCE DOCUMENT - NOT FOR CONSTRUCTION

SIEMENS

1702183