

1. EXISTING AIR TERMINAL UNIT TO REMAIN.
2. EXISTING DUCTWORK TO REMAIN.
3. EXISTING AIR DEVICE TO REMAIN.
4. REMOVE EXISTING DUCTWORK AND ASSOCIATED DAMPERS, HANGERS, ETC.
5. REMOVE EXISTING DUCTWORK BACK TO THIS POINT FOR RECONNECTION IN NEW WORK.
6. REMOVE EXISTING THERMOSTAT AND ASSOCIATED WIRING, ETC.
7. REMOVE EXISTING AIR DEVICE AND ASSOCIATED DUCTWORK, HANGERS, ETC.
8. REMOVE EXISTING AIR DEVICE FOR RELOCATION IN NEW WORK.


[illegible]

Heapy Engineering
MEP Design Technology Planning Commissioning Energy

Nationally Recognized Leader in Sustainability

1400 W Dorothy Lane, Dayton, OH 45409-1310
Ph 937-224-0861 Fax 937-224-5777 www.heapy.com

HEAPY PROJECT No.: 2015-04016 FIRM LICENSE No.: 01528



ARCHITECT/ENGINEERS:

pfb architects
Cincinnati • Chicago

PFB Architects, Inc.
9481 Kerwood Road
Cincinnati, Ohio 45242
(513) 861 3200

Approved: Project Director

Date	10/30/2015
------	------------

Checked	DLE
---------	-----

Project No.	
VA Project No.	539-CSI-201 3585.00
Building Number 01	
Drawing Number MD-102	
Dwg. of	



Department of
Veterans Affairs

1/23/2015 1:37:14 PM

Scale: 1/4" = 1'-0"

Drawing Title

SECOND FLOOR PLAN HVAC NEW
WORK

Approved: Project Director

Project Title	MRI SITE PREP
---------------	---------------

Location	Cincinnati, Ohio
----------	------------------

Date

Checked

Drawn

Project No.	
VA Project No.	539-CSI-201 3585.00

Drawing Number
M-100


M

Office of
Construction
and Facilities
Management

 Department of
Veterans Affairs

CONSULTANTS:

Heapy Engineering
MEP Design Technology Planning Commissioning Energy
Nationally Recognized Leader in Sustainability
1400 W Dorothy Lane, Dayton, OH 45409-1310
PH 937-224-0861 Fax 937-224-5777 www.heapy.com
HEAPY PROJECT No.: 2015-04016 FIRM LICENSE No.: 01528



ARCHITECT/ENGINEERS:

pflc architects
Cincinnati • Chicago

PFB Architects, Inc.
9461 Kenwood Road
Cincinnati, Ohio 4524
(513) 861 3200

FULLY SPRINKLERED

GENERAL NOTES

A. REFER TO SHEET M-100 FOR INDEX, LEGEND AND ADDITIONAL GENERAL NOTES.

NOTES

- EXISTING EXHAUST FAN TO REMAIN.
- EXISTING DUCTWORK TO REMAIN.
- EXISTING CONDENSING UNIT TO REMAIN.
- EXISTING AIR HANDLING UNIT TO REMAIN.
- EXISTING PIPE PENETRATION ROOF CURB TO REMAIN.
- EXISTING PIPING TO REMAIN.
- PROVIDE ROOF EQUIPMENT SUPPORT CURB. REFER TO DETAIL ON SHEET M-501.
- RECEIVE AND INSTALL NEW KKT KRAUS MRI CHILLER FURNISHED BY SIEMENS.
- 2" ST. GLYCOL CHILLED WATER UP FROM SECOND FLOOR THRU EXISTING PIPING PENETRATION ROOF CURB. REFER TO SHEET M-102 FOR CONTINUATION.
- PROVIDE ALUMINUM JACKET ON ALL EXTERIOR PIPING AS SPECIFIED.
- PROVIDE QUENCH VENT PER DETAIL ON SHEET M-106.
- SHEAR WALL BEING PROVIDED UNDER VA PROJECT #539-328. COORDINATE WITH COR.
- COLUMN BEING PROVIDED UNDER VA PROJECT #539-328. COORDINATE WITH COR.
- AIR HANDLING UNIT BEING PROVIDED UNDER VA PROJECT #539-328. COORDINATE WITH COR.
- EMERGENCY EXHAUST DUCT UP FROM SECOND FLOOR. PROVIDE DUCT PENETRATION ROOF CURB. REFER TO SHEET M-102 FOR CONTINUATION.
- TURN DUCT DOWN AND CONNECT TO PLENUM CASING PROVIDED UNDER VA PROJECT #539-328. IF PLENUM CASING NOT PROVIDED IN TIME, TURN DUCT DOWN AND TERMINATE WITH BRIDGECRACK.
- PROVIDE GRAVITY BACKDRAFT DAMPER FULL SIZE OF EXHAUST DUCT. INSULATE DUCT BETWEEN DAMPER AND LOUVERED PLENUM.
- AIR HANDLING UNIT BEING PROVIDED UNDER VA PROJECT #539-15-102. COORDINATE WITH COR.
- PROVIDE ALUMINUM DUCTWORK.


NOTE: THIRD FLOOR IS CURRENTLY A ROOF. VA PROJECT #539-328 (IN CONSTRUCTION) IS ADDING FIFTH FLOOR AND TURNING THIRD FLOOR INTO A TWO STORY MECHANICAL ROOM. STRUCTURE SHOULD BE INSTALLED PRIOR TO CONSTRUCTION STARTING FOR THIS PROJECT.

THIRD FLOOR PLAN - NEW WORK

Scale: 1/4" = 1'-0"

CONSULTANTS:

Heapy Engineering
MEP Design Technology Planning Commissioning Energy
Nationally Recognized Leader in Sustainability
1400 W Dorothy Lane, Dayton, OH 45409-1310
Ph 937-224-0861 Fax 937-224-5777 www.heapy.com
HEAPY PROJECT No.: 2015-04016 FIRM LICENSE No.: 01528



ARCHITECT/ENGINEERS:

pfb architects
Cincinnati • Chicago
PFB Architects, Inc.
9481 Kenwood Road
Cincinnati, Ohio 45242
(513) 861 3200

Drawing Title

THIRD FLOOR PLAN HVAC NEW WORK

Approved: Project Director

Project Title

MRI SITE PREP

Location

Cincinnati, Ohio

Date

10/30/2015

Checked

DLE

Drawn

WJS

Project No.

539-CSI-201

Building Number

01

Drawing Number

M-103


Dwg. of

Office of
Construction
and Facilities
Management



[illegible]


NOTE: FIFTH FLOOR CONSTRUCTION IS BEING PROVIDED UNDER VA PROJECT #539-328 SIMULTANEOUSLY WITH THIS PROJECT. STRUCTURE SHOULD BE INSTALLED PRIOR TO CONSTRUCTION STARTING FOR THIS PROJECT.



 **ROOF PLAN - NEW WORK**
Scale: 1/4" = 1'-0"

[illegible]

CONSULTANTS:



ARCHITECT/ENGINEERS:

pfb architects
Cincinnati • Chicago

PFB Architects, Inc.
9461 Kerwood Road
Cincinnati, Ohio 45242
(513) 861 3200

	Drawing Title
--	---------------

ROOF PLAN HVAC NEW WORK

Approved: Project Director

Project Title

MRI SITE PREP

Location	
----------	--

Cinci

Checked

Checked

Drawn

Project No.	
VA Project No.	539-CSI-201 3585.00

01


Drawing Number

M-106

Dwg. of


Office of
Construction
and Facilities
Management






- ① 18 GAUGE GALVANIZED STEEL COUNTER-FLASHING.
- ② WELDED 14 GAUGE EQUIPMENT SUPPORT CURB, MEETING ASTM A-446, 525, 526 AND 527 REQUIREMENTS, WITH WELDED CORNERS WITH SEAMS JOINED BY CONTINUOUS WELDS. CURB SHALL BE INTERNALLY REINFORCED WITH BULKHEADS AND SPREADERS, 24" ON CENTER TO MEET LOAD RATINGS OF EQUIPMENT. CURB TO EXTEND 6" BEYOND EQUIPMENT. REFER TO FLOOR PLANS FOR HEIGHT.
- ③ SECURE CURB TO ROOF WITH EXPANSION BOLTS (CONCRETE ROOF) OR RUST RESISTANT BOLTS (METAL DECK AND BAR JOIST ROOF), 12" O.C.
- GENERAL NOTE:
1. THIS DETAIL IS NOT INTENDED FOR ROOFTOP AHU SUPPORT. REFER TO STRUCTURAL DRAWINGS FOR REQUIREMENTS OF ROOFTOP AHU SUPPORT.

EQUIPMENT/DUCT SUPPORT ROOF CURB




QUENCH VENT THRU ROOF

SCALE: NONE




- NOTES:
1. ALL VANE ELBOWS SHALL BE CONSTRUCTED AND INSTALLED AS DETAILED BY SMACNA.
2. WHEN W1 DOES NOT EQUAL W2, VANE SHALL BE SINGLE THICKNESS VANE TYPE REGARDLESS OF W DIMENSION.
3. ALL SINGLE THICKNESS VANES SHALL HAVE A 2" RADIUS, 1 1/2" MAXIMUM SPACE BETWEEN VANES AND A 3/4" TRAILING EDGE.
4. WHEN W EQUALS W2 AND W1 IS GREATER THAN 20", VANES SHALL BE DOUBLE VANE TYPE.

DUCTWORK SQUARE VANE ELBOWS




- NOTES:
1. BRANCH DUCT TAKE-OFF WITH MANUAL DAMPER.
2. BRANCH DUCT SIZES, UNLESS NOTED ON PLANS ARE TO BE SIZED AS FOLLOWS:
- 100 CFM AND LESS - 6" DIA.
 - 101 CFM TO 250 CFM - 8" DIA.
 - 251 CFM TO 400 CFM - 10" DIA.
 - 401 CFM TO 700 CFM - 12" DIA.


RETURN OR EXHAUST GRILLE/REGISTER CONNECTION



FLEXIBLE AIR DUCT CONNECTOR



EXHAUST OR RETURN BRANCH DUCTWORK




TYPICAL DUCTWORK TRANSITION WITH EQUIPMENT MOUNTED IN DUCT, PLAN OR SIDE VIEW.

TYPICAL DUCTWORK TRANSITION, PLAN OR SIDE VIEW.


NOTE: UNLESS OTHERWISE INDICATED ON PLANS, MAXIMUM ANGLES SHOWN SHALL APPLY.

DUCTWORK TRANSITIONS




- NOTES:
1. THE INTERIOR SURFACE OF ALL RADIUS ELBOWS SHALL BE MADE ROUND.
2. ALL STANDARD RADIUS ELBOWS CAN BE SUBSTITUTED WITH SHORT RADIUS ELBOWS. ALL SHORT RADIUS ELBOWS SHALL HAVE VANES. VANES SHALL BE CONSTRUCTED, SUPPORTED AND FASTENED AS RECOMMENDED BY SMACNA.

DUCTWORK RADIUS ELBOWS




0.5 INCH WG TO 2 INCHES WG (LOW PRESSURE) DUCT RISER SUPPORT



2 INCHES WG TO 4 INCHES WG (HIGH AND MEDIUM PRESSURE) DUCT RISER SUPPORT

NOTE: ALL DUCT WORK RISERS WHICH ARE RUN EXPOSED, SUCH AS THRU ATTIC FLOORS AND FAN ROOM FLOORS SHALL BE PROVIDED WITH A 3" HIGH CONCRETE CURB AROUND OPENING FOR DUCT.


DUCT RISER SUPPORTS



HANGER STRAPS OR RODS			
MAX. DUCT DIA. - IN.	QUANTITY/SIZE IN.	MAX. LOAD LBS.	MAX. SPACING IN.
26	ONE 1 x 22 GA. STRAP	260	144
36	ONE 1 x 18 GA. STRAP	420	144
50	ONE 1 x 16 GA. STRAP	700	144
60	TWO 3/8 DIA. RODS	1320	144
84	TWO 1/2 DIA. RODS	2500	144

NOTE: TABULATED DATA FROM SMACNA ALLOWS FOR DUCT REINFORCING AND INSULATION, BUT NO EXTERNAL LOAD.

ROUND DUCT HANGERS




SECURING HANGER RODS IN CONCRETE

NEW MRI CONTROLS

1. EMERGENCY EXHAUST FAN
- 1.1 PROVIDE OXYGEN SENSOR ON WALL OF NEW MRI SCANNING ROOM. COORDINATE LOCATION WITH SIEMENS AND COR.
- 1.2 EXHAUST FAN EF-A2001B SHALL BE STARTED AND STOPPED BY THE DCP OR REMOTELY AT THE EOC. H-O-A SWITCH SHALL BE KEPT IN THE "AUTO" POSITION. "HAND" AND "OFF" SHALL BE USED ONLY FOR MAINTENANCE.
- 1.3 EXHAUST FAN EF-A2001B SHALL OPERATE IF OXYGEN LEVEL OF MRI SCANNING ROOM DROPS BELOW 20.0% (ADJUSTABLE). FAN SHALL ALSO OPERATE IF ALARM IS PRESENT FROM MRI CONTROL PANEL.
- 1.4 AUTOMATIC DAMPER IN RETURN DUCT SHALL CLOSE IF EF-A2001B OPERATES.
2. MRI CONTROL PANEL INTERFACE
- 2.1 EXTEND ALARM SIGNAL FROM MRI CONTROL PANEL TO DDC SYSTEM.

EXISTING MRI CONTROLS

1. EMERGENCY EXHAUST FAN
- 1.1 PROVIDE OXYGEN SENSOR ON WALL OF EXISTING MRI SCANNING ROOM. COORDINATE LOCATION WITH SIEMENS AND COR.
- 1.2 EXHAUST FAN EF-5-2 OPERATE IF OXYGEN LEVEL OF MRI SCANNING ROOM DROPS BELOW 20.0% (ADJUSTABLE). FAN SHALL ALSO OPERATE IF ALARM IS PRESENT FROM MRI CONTROL PANEL.
2. MRI CONTROL PANEL INTERFACE
- 2.1 EXTEND ALARM SIGNAL FROM MRI CONTROL PANEL TO DDC SYSTEM.



- 1.1 IF THE SPACE TEMPERATURE IS BELOW SETPOINT, THE BOX DAMPER SHALL BE AT THE LISTED DEAD BAND MINIMUM CFM. IF SPACE TEMPERATURE RISES ABOVE THE COOLING SETPOINT, THE CONTROL SHALL MODULATE THE BOX DAMPER BETWEEN THE LISTED DEAD BAND MINIMUM CFM AND THE LISTED COOLING MAXIMUM CFM TO SATISFY THE SPACE COOLING SETPOINT. BOX CONTROLS SHALL REVERSE ACTION DURING "WARM-UP" CYCLES.

AIR TERMINAL UNIT CONTROLS

11/23/2015 1:36:56 PM

Revisions

Date

CONSULTANTS:

Heapy Engineering

MEP Design Technology Planning Commissioning Energy

Nationally Recognized Leader in Sustainability

1400 W Dorothy Lane, Dayton, OH 45409-1310
Ph 937-224-0861 Fax 937-224-5777 www.heapy.com

HEAPY PROJECT No.: 2015-04016 FIRM LICENSE No.: 01528

ARCHITECT/ENGINEERS:

pfb architects

Cincinnati • Chicago

PFB Architects, Inc.
9461 Kenwood Road
Cincinnati, Ohio 45242
(513) 861-3200

Drawing Title

DETAILS

Approved: Project Director

Project Title

MRI SITE PREP

Location

Cincinnati, Ohio

Date

10/30/2015

Checked

DLE

Drawn

WJS

Project No.

VA Project No. 539-CSI-201 3585.00

Building Number

01

Drawing Number


M-501

Dwg. of

Office of Construction and Facilities Management

Department of Veterans Affairs

MECHANICAL PLAN



ENVIRONMENTAL REQUIREMENTS

- AIR CONDITIONING IS TO PROVIDE A TEMPERATURE OF 70°F ±5°F IN THE CONTROL & EQUIPMENT ROOMS 65°F-71°F IN EXAM ROOM. RELATIVE HUMIDITY OF 40-60% (NON-CONDENSING) IS REQUIRED EXAMINATION ROOM AND 40-60% (NON-CONDENSING) IN ALL OTHER AREAS WHERE SIEMENS EQUIPMENT IS INSTALLED. THESE CONDITIONS ARE TO BE MET AT ALL TIMES; 24 HOURS A DAY, 7 DAYS A WEEK.
- A DEDICATED AIR CONDITIONING AND HUMIDIFICATION SYSTEM IS RECOMMENDED FOR THE EXAM ROOM. A MINIMUM AIR EXCHANGE RATE OF 6 TIMES PER HOUR FOR THE EXAM ROOM IS REQUIRED. IT IS RECOMMENDED TO INSTALL A FRESH AIR SYSTEM WITH 30%-50% FRESH AIR INTAKE.
- THE HEAT INTO THE EXAM ROOM IS LESS THAN 10,236 BTU/HR. THE HEAT INTO THE EQUIPMENT ROOM IS LESS THAN 3,412 BTU/HR. THIS HEAT DISSIPATION IS FROM THE SIEMENS EQUIPMENT ONLY. AUXILIARY SUPPORT EQUIPMENT (e.g. UPS) AND LIGHTING MUST BE CONSIDERED FOR TOTAL HEAT LOADS.
- IT IS IMPORTANT FOR FRESH AIR INTAKE SYSTEMS TO EXHAUST AIR DIRECTLY OUT OF THE BUILDING. THE EXHAUST AIR MUST NOT BE DEFLECTED INTO ANOTHER ROOM. THE MAGNET ROOM EXHAUST AIR SHOULD BE INSTALLED AT LEAST 6'-6" ABOVE FINISHED FLOOR.
- THE AIR INTAKE OF THE AIR CONDITIONING SYSTEM MUST NOT BE LOCATED IN THE VICINITY OF THE QUENCH VENT EXHAUST.
- IF THE INPUT DRAWS UPON AIR FROM OUTSIDE THE BUILDING, IT IS RECOMMENDED TO INSTALL AN ON-SITE FILTER TO REMOVE DUST PARTICLES GREATER THAN 10 MICRONS.
- DO NOT LOCATE ANY HVAC DIFFUSERS ABOVE THE MAGNET. THERE SHALL NOT BE AIR BLOWING DIRECTLY ON THE MAGNET.

CHILLED WATER SUPPLY


A CHILLED WATER SUPPLY IS REQUIRED TO THE MRI SYSTEM 24 HOURS A DAY, YEAR ROUND FOR THE COLD HEAD AND GRADIENT SYSTEMS. THIS CAN BE PROVIDED BY A DEDICATED KRAUS ECO CHILLER AND INTERFACE PANEL. CHILLED WATER CAN ALTERNATIVELY BE SUPPLIED BY OTHER MEANS IN COMBINATION WITH A SEPARATOR CABINET PROVIDED BY SIEMENS.

THE PIPE SIZE BETWEEN THE KRAUS CHILLER AND INTERFACE PANEL, OR BETWEEN THE WATER SUPPLY AND SEP MUST BE 2 INCH UP TO 82 FEET, 2-1/2 INCH UP TO 148 FEET, CONSULT KKT KRAUS FOR LONGER PIPE. PERMISSIBLE MATERIALS THAT CAN BE USED FOR THE PIPING ARE: STAINLESS STEEL (VQA, V4A), NON-FERROUS METAL (COPPER, BRASS), SYNTHETIC MATERIAL, PLASTICS, BRAZING SOLDER, HARD SOLDER, OR FITTING SOLDER TYPE 3 AND 4. THERE ARE MATERIALS THAT MAY CAUSE DAMAGE TO THE COOLING SYSTEM AND CANNOT BE USED, THESE MATERIALS ARE ALUMINUM, IRON, CARBON STEEL, ZINC, ZINC PLATED STEEL, OR STANDARD STEEL PIPES.

THESE REQUIREMENTS ARE REQUIRED FOR NEW INSTALLATIONS, IF EXISTING WATER PIPES COMPLY WITH SIEMENS WATER SPECIFICATIONS, THEY DO NOT NEED TO BE REPLACED.


NORMAL TAP WATER MUST BE AVAILABLE FOR FILLING THE SECONDARY WATER CIRCUIT. THERE SHALL BE A HOSE BIB LOCATED WITHIN 65' OF THE SEP, EPC OR THE KRAUS CHILLER.

THE SUPPLY AND RETURN CHILLED WATER PIPES MUST BE LABELED. THE LOCATION OF THE LABELS MUST BE AT ALL CONNECTION AND REFILLING POINTS AND MUST CONTAIN FLOW DIRECTION AND CONTENTS.



CEILING HEIGHTS

MAGNET EXAMINATION ROOM: 7'-11" MINIMUM
EQUIPMENT ROOM: 7'-3" MINIMUM WITH RESTRICTION
ALL ANCILLARY AREAS: 6'-11" MINIMUM



NOTES:

ALL PIPING AND PLUMBING FIXTURES SHALL BE FURNISHED, INSTALLED, PRESSURE TESTED AND CHARGED BY THE MECHANICAL CONTRACTOR PRIOR TO THE DELIVERY AND INSTALLATION OF THE SIEMENS SUPPLIED AND INSTALLED EQUIPMENT UNLESS SPECIFIED OTHERWISE.

AT THE HIGHEST POINT OF THE WATER SUPPLY PIPE FROM THE KRAUS CHILLER AN AUTOMATIC DEAERATION DEVICE (AIR VENT) WITH BALL VALVE MUST BE INSTALLED BY THE MECHANICAL CONTRACTOR.

SYSTEM MUST BE PROVEN TO BE FREE FROM LEAKAGE.

THE MECHANICAL ENGINEER OF RECORD SHALL BE ULTIMATELY RESPONSIBLE FOR THE SITE SPECIFIC DESIGN AND SPECIFICATION OF THE MECHANICAL AND PIPING SYSTEMS AS SHOWN AND SHALL BE IN ACCORDANCE WITH ALL APPLICABLE LOCAL, STATE AND NATIONAL CODES. ALL WORK SHALL BE PERFORMED BY THE MECHANICAL CONTRACTOR AND SHALL BE SUBJECT TO COMPLIANCE WITH ALL APPLICABLE LOCAL, STATE AND NATIONAL CODES.

THE SUPPLY AND RETURN PIPES FROM THE CHILLED WATER SUPPLY TO THE IFP MUST BE LABELED TO SHOW FLOW DIRECTION AND CONTENT (WATER/GLYCOL).

A TAP WATER SUPPLY MUST BE AVAILABLE WITHIN 45' OF THE IFP AND CHILLER CONNECTION FOR FILLING THE CIRCUIT.

1 PIPING SCHEMATIC FOR CHILLED WATER-KRAUS ECO CHILLER

CHILLED WATER REQUIREMENTS

FLOW RATE:	31.7 GPM
WATER TEMPERATURE:	66.2°F - 71.6°F
HEAT DISSIPATION TO WATER	204,910 BTU/HR
WATER PRESSURE	MAXIMUM 87 PSI
PRESSURE LOSS	NA
CHILLED WATER ACIDITY RANGE	6 pH TO 8 pH
CHILLED WATER HARDNESS	<250 ppm CALCIUM CARBONATE
CHLORINE GAS CONCENTRATION	<200 ppm
FILTRATION	700 µm

FOR INSTALLATION OF A KKT CHILLER, IT IS THE RESPONSIBILITY OF THE MECHANICAL CONTRACTOR TO PROVIDE A MIXTURE OF WATER WITH 35%-38% ETHYLENE GLYCOL. PRIOR TO CHILLER START UP DO NOT USE PROPYLENE GLYCOL OR AUTOMOTIVE ANTI-FREEZE.

THE AMOUNT OF THE MIXTURE MUST FILL THE CHILLER, MR SYSTEM AND PIPING (SUPPLY AND RETURN), SEE EXAMPLES BELOW.

(1) GALLON OF UNDILUTED GLYCOL, OR (2) GALLONS OF WATER/GLYCOL MIXTURE MUST REMAIN ON SITE FOR USE AFTER START UP.

PIPE DIAMETER	TOTAL LENGTH	MIXTURE VOLUME	GLYCOL NEEDED
2"	100'	31.3 GALLONS	11.9 GALLONS
2"	200'	47.6 GALLONS	18.1 GALLONS
2.5"	100'	40.5 GALLONS	15.4 GALLONS
2.5"	200'	66.0 GALLONS	25.1 GALLONS

MIXTURE VOLUME = 3.14 x (PIPE RADIUS)² x PIPE LENGTH + 15 GALLONS. GLYCOL AMOUNT = 35-38% OF MIXTURE VOLUME.

MECHANICAL NOTES

- THE AIR H.V.A.C. SYSTEM MUST OPERATE FOR A MINIMUM OF 48 CONSECUTIVE HOURS PRIOR TO THE DELIVERY OF THE EQUIPMENT.
- THE FILTERS MUST BE CHANGED IMMEDIATELY PRIOR TO THE DELIVERY OF THE EQUIPMENT.
- SIEMENS REQUIRES THE USE OF A DEDICATED H.V.A.C. SYSTEM FOR THE EQUIPMENT ROOM TO BE LOCATED, SIZED AND SPECIFIED BY THE MECHANICAL ENGINEER OF RECORD AND TO BE SUPPLIED AND INSTALLED BY THE MECHANICAL CONTRACTOR.
- SIEMENS RECOMMENDS THAT THE CUSTOMER PROVIDE AND INSTALL AN OXYGEN MONITORING SYSTEM WITH VISUAL AND AUDIBLE ALARMS TO INDICATE WHEN THE OXYGEN CONTAINED IN AMBIENT AIR FALLS BELOW PRE-PROGRAMMED SAFETY LEVELS WITH THE SENSOR TO BE LOCATED IN THE SCAN ROOM IN THE AREA DESIGNATED FOR CRYOGEN FILLING.
- THE SIEMENS ACTIVE SHIELDED MAGNET RECIRCULATES LIQUID HELIUM, ELIMINATING THE NEED FOR A DEDICATED CRYOGEN STORAGE AREA. THE RECIRCULATING SYSTEM SIGNIFICANTLY REDUCES THE HELIUM "BOIL OFF". THE MAGNET WILL REQUIRE OCCASIONAL FILLING. A DELIVERY ROUTE FOR CRYOGEN DEWARMS MUST BE ESTABLISHED. A MINIMUM 36" CLEARANCE IS REQUIRED.

FIRE CONTROL NOTES

- SIEMENS HAS NO SPECIFIC REQUIREMENT FOR FIRE PROTECTION. FIRE PROTECTION REQUIREMENTS SHALL BE IN ACCORDANCE WITH LOCAL CODES AND CUSTOMER'S INSURANCE REQUIREMENTS. ALL FIRE PROTECTION SYSTEMS SHALL BE DEFINED BY THE ARCHITECT OF RECORD WITH DESIGN, SPECIFICATION AND DETAILING OF THE FIRE PROTECTION SYSTEM BY THE MECHANICAL ENGINEER OF RECORD IN ACCORDANCE WITH SIEMENS GUIDELINES AS STATED HEREIN. THE ELECTRONIC EQUIPMENT OF THE MR SYSTEMS WILL BE DAMAGED BY WATER. REDUCTION OR ELIMINATION OF WATER USED FOR FIRE SUPPRESSION WILL REDUCE POTENTIAL WATER DAMAGE. PRE-ACTION WATER CATCHERS OR OTHER METHODS CAN REDUCE OR ELIMINATE WATER. REFER TO YOUR FIRE PROTECTION PROFESSIONAL.
- THE USE OF SMOKE DETECTORS INSIDE OF THE MR EXAMINATION ROOM IS NOT RECOMMENDED. SMOKE DETECTORS, BY DESIGN, CAN GENERATE NOISE THAT MAY INTERFERE WITH THE MR EXAMINATION AND CAUSE IMAGE ARTIFACTS. IF THE USE OF A SMOKE DETECTOR IN THE EXAMINATION ROOM IS MANDATED BY LOCAL REQUIREMENTS, SPECIAL NOISE TESTS MUST BE PERFORMED BY SIEMENS SERVICE AFTER THE MRI IS OPERATIONAL. MRI EQUIPMENT PERFORMANCE PROBLEMS DUE TO SMOKE DETECTORS ARE THE RESPONSIBILITY OF THE CUSTOMER AND ARE NOT COVERED UNDER WARRANTY OR SERVICE AGREEMENT.
- ALL MATERIAL USED INSIDE THE MAGNET ROOM SHALL BE NON-MAGNETIC.
- ALL PENETRATIONS IN THE RF CABIN/SHIELD SHALL BE THROUGH A WAVE GUIDE TO BE EQUIPPED WITH A SIEMENS APPROVED DIELECTRIC COUPLER ON BOTH ENDS OF THE WAVE GUIDE. ALL WAVE GUIDES SHALL BE DESIGNED, DETAILED AND SPECIFIED BY THE RF CABIN/SHIELD CONTRACTOR WITH ALL LOCATIONS TO BE DETERMINED BY THE ARCHITECT AND MECHANICAL ENGINEER OF RECORD TO BE ESTABLISHED IN A PRE-PLANNING MEETING PRIOR TO THE DESIGN, SPECIFICATION, AND FABRICATION OF THE RF CABIN/SHIELD.
- EACH ELECTRICAL PENETRATION OF THE RF CABIN/SHIELD FOR ELECTRICAL SERVICING OF THE FIRE PROTECTION SYSTEM SHALL BE THROUGH AN RF FILTER TO BE SUPPLIED BY THE RF SHIELD CONTRACTOR WITH FILTER LOCATIONS TO BE DETERMINED BY THE ARCHITECT AND THE ELECTRICAL ENGINEER OF RECORD TO BE ESTABLISHED IN A PRE-PLANNING MEETING PRIOR TO THE DESIGN, SPECIFICATION AND FABRICATION OF THE RF CABIN/SHIELD.
- IT IS PERMISSIBLE TO RUN "BLACK PIPE" UP TO THE DIELECTRIC COUPLER ON THE OUTSIDE OF THE RF SHIELD.
- THERE MUST BE NO GROUND CONNECTIONS MADE DURING THE INSTALLATION OF EITHER THE PIPING OR ELECTRICAL FOR THE FIRE PROTECTION SYSTEM.
- THE USE OF HALON IS NOT ACCEPTABLE.
- THE LOCATION OF FIRE CONTROL SYSTEM COMPONENTS SHALL BE COORDINATED THROUGH THE ARCHITECT OF RECORD WITH ALL LOCATIONS TO BE COORDINATED WITH SIEMENS EQUIPMENT LOCATIONS AS SHOWN ON THE 1/4" SCALE EQUIPMENT LOCATION PLAN.
- THE FIRE CONTROL CONTRACTOR SHALL VERIFY EQUIPMENT MOUNTING PROCEDURES ON ANY WALLS CONTAINING RF SHIELDING WITH THE SIEMENS PROJECT MANAGER PRIOR TO THE COMMENCEMENT OF WORK.

COMPRESSOR LINE INSULATION

COMPRESSOR LINES RUNNING FROM THE COMPRESSOR (OR SEP CABINET) TO THE MAGNET ARE INSULATED BY SIEMENS. ADDITIONAL INSULATION (ARMAFLEX OR EQUIVALENT) FOR NOISE REDUCTION (CHIRPING) MAY BE REQUIRED. ADDITIONAL INSULATION NOT PROVIDED BY SIEMENS.

ATTENTION:

THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.

THIS SET OF PLANS REPRESENTS A COMPLETE SET OF DETAILS AND SHOULD NOT BE SEPARATED.

IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.

ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.

THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

PROJECT MANAGER: TRENT CHILDRESS
TEL: (317) 341-1094
EXT: 1000
FAX: 1000
EMAIL: trent.childress@siemens.com

SIEMENS
CINCINNATI VAMC
3200 VINE STREET, CINCINNATI, OH 45220
MRI SUITE 2024 - MAGNETOM SKYRA W/MOBILE TABLE


PROJECT #: **1400888**
SHEET #: **M-101**

SHEET 9 OF 10
DRAWN BY: R. SUTHERS
DATE: 11/03/15

SCALE: AS NOTED
REF: #30182993

CONSULTANTS:

Heapy Engineering
MEP Design Technology Planning Commissioning Energy
Nationally Recognized Leader in Sustainability
1400 W Dorothy Lane, Dayton, OH 45409-1310
Ph 937-224-0861 Fax 937-224-5777 www.heapy.com
HEAPY PROJECT No.: 2015-04016 FIRM LICENSE No.: 01528



architects
Cincinnati • Chicago
PFB Architects, Inc.
9401 Kenwood Road
Cincinnati, Ohio 45242
(513) 861 3200

Drawing Title
SIEMENS SUPPLEMENT MRI

Approved: Project Director

Project Title
MRI SITE PREP

Location
Cincinnati, Ohio

Date
10/30/2015

Checked
DLE

Drawn
WJS

Project No.
VA Project No. 539-CSI-201
3585.00








Building Number
01

Drawing Number
M-502

Dwg. of

Office of
Construction
and Facilities
Management





CRYOGEN NOTES

1) "CRYOGENS" IS A TERM USED TO IDENTIFY THE REFRIGERANT USED TO MAKE THE MAGNET "SUPER-CONDUCTING". IN THIS APPLICATION, LIQUID AND GASEOUS HELIUM. SPECIAL CARE MUST BE TAKEN DURING THE TRANSFER OF THE MAGNET WITH CRYOGENS AND NORMAL EXHAUST OF CRYOGENS FROM THE SYSTEM. ASIDE FROM THE OBVIOUS DANGER OF FREEZING, HELIUM GAS WILL ALSO DISPLACE THE OXYGEN IN THE ROOM. THE INSTALLATION OF AN APPROVED TOXGARD MONITORING SYSTEM IS RECOMMENDED.

2) THERE SHALL BE A TRANSPORT ROUTE FOR DELIVERY OF CRYOGENS TO THE EXAM ROOM. SPECIAL VESSELS CALLED DEWARs ARE USED TO TRANSPORT HELIUM. A 250 LITER DEWAR WEIGHS 335 POUNDS AND HAS A 32" DIAMETER. A 500 LITER IS 540 POUNDS, AND IS 42" IN DIAMETER.

3) HELIUM GAS CYLINDERS MAY BE USED DURING THE INITIAL FILLING OF HELIUM INTO THE MAGNET. THE FACILITY IN WHICH THESE MAY BE USED NEEDS TO HAVE THE ABILITY TO TEMPORARILY STORE AND SECURE THESE CYLINDERS THAT WILL PREVENT THEM FROM INADVERTENTLY FALLING OVER.

4) OUTSIDE VENTING OF THE HELIUM IS TO BE PROVIDED BY MEANS OF A VENT PIPE OF NON-MAGNETIC MATERIAL CALLED A QUENCH VENT.

REV 0

HELIUM CONTENT	
LITERS AT 100%	1,190
TYPICAL BOIL OFF RATE	0.0 L/HR
TYPICAL REFILL INTERVAL	10 YEARS
FOR TYPICAL CLINICAL USE, DEPENDING ON SEQUENCES AND OPERATING TIME.	
WITHOUT THE COLD HEAD RUNNING THE LIQUID HELIUM WILL BOIL OFF FROM 98% TO 0% IN APPROX. 29 DAYS. THE LOSS DURING SHIPPING IS APPROX. 3.3% PER DAY.	

- QUENCH VENT NOTES**
- 1) IN THE EVENT OF A QUENCH, THE THERMAL ENERGY DISSIPATED CAUSES AN EXTREMELY RAPID BOIL OFF OF THE LIQUID HELIUM. THE SYSTEM MUST BE CAPABLE OF VENTING THE LARGE VOLUME OF GAS GENERATED AT THE APPROXIMATE EXPANSION RATIO OF 1:700 FROM LIQUID AT 4.2°K TO ROOM TEMPERATURE GAS. THE EXHAUST SYSTEM IS CRITICAL FOR THE SAFE OPERATION OF THE MAGNET. THE DATA IN THIS DOCUMENT MUST BE FOLLOWED. SINCE HELIUM VENTED IN A QUENCH IS AN ASPHYXIANT & AN EXTREMELY COLD GAS, THE QUENCH TUBE MUST ALWAYS END AT A POINT WHERE ACCESS BY PEOPLE IS NOT POSSIBLE. QUENCH TUBE PLANNING MUST ONLY BE DONE BY QUALIFIED PERSONNEL. IT IS THE OWNER'S RESPONSIBILITY TO ENSURE THAT THE QUENCH TUBE IS MAINTAINED IN AN OPERABLE STATE.
- 2) IF THE QUENCH VENT IS NOT CONFIGURED CORRECTLY THERE IS A RISK OF DANGER THAT MAY LEAD TO DEATH OR SERIOUS INJURY AND CAN RESULT IN STRUCTURAL DAMAGE. THE EXHAUST MUST NOT BE VENTED IN AN ENCLOSED SPACE. THE OPERATOR OF THE SYSTEM MUST PREPARE AN EMERGENCY PLAN IN THE EVENT OF A QUENCH.
- 3) THE QUENCH TUBE CONSISTS OF STRAIGHT, HYDRAULICALLY SMOOTH SECTIONS, BENDS UP TO 90° AND A DIFFUSER, IF REQUIRED. THE END OF THE TUBE MUST BE PROTECTED FROM RAIN, SNOW, AND FOREIGN OBJECTS. ROUND SECTIONS ONLY, NO SQUARE SECTIONS.
- 4) THE SIEMENS MAGNET HAS A QUENCH VALVE ASSEMBLY FOR CONNECTION TO THE TUBE LOCATED AT THE TOP LEFT SIDE OF THE MAGNET (SEE MAGNET ELEVATION). THE MECHANICAL CONTRACTOR WILL SUPPLY AND INSTALL A QUENCH VENT TUBE WITH CAP. TO BE NON-MAGNETIC STAINLESS STEEL (≥22 GAUGE RECOMMENDED). GRADES AISI304, 309, 316, OR 321 ONLY. THERMAL CONDITIONS MAY CAUSE THE TUBE TO CONTRACT UP TO 3mm/METER SO A STAINLESS STEEL BELLOWS OR FLEXIBLE SECTION MUST BE INSTALLED A MINIMUM OF EVERY 32'-9" NOT TO EXCEED 2% OF THE OVERALL LENGTH. THE QUENCH TUBE MAY ALSO BE MADE OF ALUMINUM. EXTRUDED TUBE ALUMINUM GRADES 6063 AND 6062 ONLY MUST BE USED. ROLLED AND WELDED TUBE FROM SHEET ALUMINUM GRADE 5083 ONLY MUST BE USED. THE WALL SECTIONS OF ALUMINUM TUBE MUST BE A MINIMUM 14 GAUGE. THERMAL CONTRACTION OF 4.5 mm/METER MUST BE CONSIDERED FOR ALUMINUM QUENCH TUBES. THE MOVEMENT OF THE BELLOWS MUST BE RESTRICTED TO PREVENT EXCESSIVE EXPANSION DUE TO PRESSURE. THE WEIGHT OF THE TUBE MUST BE SUPPORTED BY THE BUILDING AND BE FLEXIBLE ENOUGH TO ALLOW MOVEMENT FROM THERMAL CONTRACTION. THE WALL EXIT SHOULD ALSO BE FLEXIBLE.
- 5) THE MAXIMUM INTERNAL PRESSURE IS CALCULATED AT 1.45 PSI. THE MAXIMUM PRESSURE SHOULD BE ENGINEERED FOR 6.5 PSI.
- 6) USE THE QUENCH VENT CALCULATOR PROVIDED BY SIEMENS TO DESIGN A QUENCH VENT THAT MEETS DESIGN REQUIREMENTS FOR DIAMETER, LENGTH, NUMBER OF ELBOWS AND PRESSURE DROP. ALL BENDS MUST BE SMOOTH WALLED AND HAVE A CENTERLINE TO INTERNAL PIPE DIAMETER RATIO OF 1.5 TO 5.0. EXPANSIONS TO PIPE DIAMETER CAN BE DONE WITH A DIFFUSER. ONLY ROUND TUBE SECTIONS MAY BE USED. RECTANGULAR SECTIONS ARE NOT ALLOWED.
- 7) THERE MUST BE A 12-19 INCH FLEXIBLE SECTION OF PIPE FOR CONNECTION TO THE QUENCH VALVE AT THE MAGNET WITH AN INSIDE DIAMETER GREATER THAN 4" (1.5T) OR 6" (3.0T) AND ABLE TO WITHSTAND 6.5 PSI.
- 8) SECTIONS OF THE PIPE CAN ONLY BE JOINED BY WELDING OR BOLTED FLANGES WITH FIBER GASKETS. ROTARY FLANGES ARE PERMITTED, VEE CLAMPED FLANGES MAY NOT BE USED.
- 9) THE PROTECTION AT THE END OF THE TUBE SHALL BE 3/8" WIRE MESH COVERING AN AREA AT LEAST 2.5 TIMES THE CROSS SECTION AREA OF THE QUENCH PIPE.
- 10) WHERE THE QUENCH TUBE EXITS THROUGH A FLAT ROOF, THE OUTLET MUST BE ABOVE A LEVEL WHERE WATER COULD ENTER THE DUCT. THE ROOF DRAINS BECOME BLOCKED. IN THE CASE OF A HORIZONTAL EXIT THROUGH A WALL, THE OUTLET SHALL BE ANGLED DOWNWARD NOT LESS THAN 1 PIPE DIAMETER TO PREVENT RAIN INGRESS. THE EXIT SHALL BE LOCATED ABOVE THE LEVEL OF DRIFTING SNOW.
- 11) WHERE THE QUENCH TUBE EXITS VERTICALLY, A RAIN COVER MUST ALSO BE FITTED WITH THE DIAMETER TO BE TWO TIMES THE DIAMETER OF THE QUENCH TUBE. THE CLEARANCE BETWEEN THE RAIN GUARD AND THE MESH SHALL 2 TIMES THE DIAMETER OF THE TUBE. A DEFLECTOR PLATE SHALL BE WELDED TO THE TUBE WHERE IT EXITS THE ROOF TO PREVENT HELIUM FROM RE-ENTERING THE BUILDING. THE DEFLECTOR SHALL BE AT LEAST THE DIAMETER OF THE RAIN GUARD AND LOCATED TWO PIPE DIAMETERS ABOVE THE ROOF AND TWO PIPE DIAMETERS BELOW THE RAIN GUARD.
- DURING A QUENCH THE HELIUM GAS EXITING THE QUENCH PIPE MAY BE AT TEMPERATURES OF LESS THAN -400°F. DUE TO THIS TEMPERATURE ROOFING MATERIALS OR ITEMS AROUND THE VENT EXIT MAY BE ADVERSELY AFFECTED. CONSIDERATION OF MATERIALS AND ITEMS PLACED NEAR THE VENT EXIT SHOULD BE TAKEN INTO ACCOUNT SO DAMAGE DOES NOT OCCUR.
- 12) TO AVOID INJURY FROM COLD BURNS AND ASPHYXIATION ACCESS TO THE QUENCH VENT MUST BE RESTRICTED BY 9'-11" ON EACH SIDE AND BELOW, AND 19'-9" ABOVE WITH WARNING SIGNS. THE EXIT MUST NOT BE LOCATED WHERE HELIUM GAS COULD BE DRAWN INTO AN AIR INLET OR OPEN WINDOW. A WARNING MUST BE PLACED NEAR THE QUENCH VENT OUTLET.
- 13) THE QUENCH TUBE MUST HAVE MINIMUM 1" INSULATION FOR THE FULL LENGTH. WITHIN THE RF ROOM THERE SHOULD BE A 1" LAYER OF MINERAL FIBER INSULATION WITH A VAPOR BARRIER AND 1" CLASS 0 OR CLASS AP ARMAFLEX. OUTDOOR PIPES MUST BE WEATHERPROOF. THE TUBE MUST HAVE A WARNING POSTED ALONG ITS ENTIRE LENGTH FOR EXTREMELY COLD HELIUM GAS - AUTHORIZED PERSONNEL ONLY. THE INSULATION MUST NOT TOUCH THE MAGNET COVERS. TO AVOID RF DISTURBANCES THE INSULATION MUST NOT MAKE ELECTRICAL CONTACT WITH THE WAVEGUIDE.
- 14) GALVANIC SEPARATION MUST BE PROVIDED BETWEEN THE MAGNET, THE QUENCH VENT, THE RF ROOM, AND THE BUILDING. TWO SEPARATIONS ARE REQUIRED USING STAINLESS STEEL BOLTS, INSULATING BUSHES AND LOCKING NUTS. NO OTHER DESIGNS ARE PERMITTED FOR SAFETY.
- 15) THE DESIGN AND CONSTRUCTION OF THE QUENCH PIPE MUST BE DOCUMENTED WITH DRAWINGS AND CALCULATIONS THAT ARE KEPT WITH INSTALLATION DOCUMENTS. IT MUST COMPLY WITH THE REQUIREMENTS IN THIS DOCUMENT BEFORE BEING CONNECTED TO THE MAGNET.
- REV 2

ATTENTION:

THIS DRAWING IS DESIGNED TO CONFORM TO FEATURES AND EQUIPMENT REQUIREMENTS PRESENTED AT THE TIME OF THEIR PREPARATION. SINCE BOTH THESE FACTORS ARE SUBJECT TO DESIGN MODIFICATION, THEY ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.

IT IS RECOMMENDED THAT THE SIEMENS DRAWINGS BE INCORPORATED WITH THE CONSTRUCTION DOCUMENTS FOR REFERENCE.

ALL DIMENSIONS SHOWN ON THIS DRAWING ARE FROM FINISHED SURFACES.

THIS DRAWING DOES NOT PROVIDE RADIATION SHIELDING REQUIREMENTS FOR X-RAY AND ASSOCIATED EQUIPMENT. THE CUSTOMER IS RESPONSIBLE FOR CONSULTING WITH A REGISTERED RADIATION PHYSICIST TO SPECIFY RADIATION PROTECTION.

SYM	DATE	DESCRIPTION
Δ	11/03/15	CHANGED NO CENTER LOCATION PER CUSTOMER
Δ	10/22/15	R101(RS) DATED 09/29/15 APPROVED BY CUSTOMER FOR FINALS

PROJECT MANAGER: TRENT CHILDRESS
TEL: (317) 341-1094 EXT: _____
FAX: _____
EMAIL: trent.childress@siemens.com

SIEMENS CINCINNATI VAMC

3200 VINE STREET, CINCINNATI, OH 45220
MRI SUITE 2024 - MAGNETOM SKYRA W/MOBILE TABLE

PROJECT #: **1400888**

THE USE OR REPRODUCTION OF THIS TITLE BLOCK WITHOUT SIEMENS' AUTHORIZATION WILL RESULT IN PROSECUTION UNDER FULL EXTENT OF THE LAW.

ALL RIGHTS ARE RESERVED.

SCALE: AS NOTED REF: #30182993

DATE: 11/03/15

SHEET: **M-501**

DRAWN BY: R. SUTHERS

DATE: 11/03/15

CONSULTANTS:


Heapy Engineering

MEP Design Technology Planning Commissioning Energy

Nationally Recognized Leader in Sustainability

1400 W Dorothy Lane, Dayton, OH 45409-1310
Ph 937-224-0861 Fax 937-224-5777 www.heapy.com

HEAPY PROJECT No.: 2015-04016 FIRM LICENSE No.: 01528



ARCHITECT/ENGINEERS:

pfb architects

Cincinnati • Chicago

PFB Architects, Inc.
8401 Kenwood Road
Cincinnati, Ohio 45242
(513) 861 3200

Drawing Title

SIEMENS SUPPLEMENT MRI

Approved: Project Director

Project Title

MRI SITE PREP

Location

Cincinnati, Ohio

Date

10/30/2015

Checked

DLE

Drawn

WJS

Project No.

539-CSI-201

3585.00

Building Number

01

Drawing Number

M-503

Dwg. of

Office of Construction and Facilities Management

Department of Veterans Affairs

11/23/2015 1:37:05 PM