
Item No. Qty

Description

1 1

Optima CT660 System

The Optima CT660 is GE's latest generation intelligent CT system. It is a scalable 32 to 128* slice platform including advanced innovations from our Discovery Series (TM). This means that Optima CT660 can grow as your clinical needs expand. With the Optima CT660 you get fast, high-quality acquisitions at optimized dose for patients young and old, large and small, across a wide spectrum of procedures: angiography, brain, chest, abdomen, orthopedic, and more.

* Overlapped Reconstruction is required to obtain additional reconstructed slices in one axial rotation.

Key Features:

- Exclusive V-Res (TM) Detector technology providing 20mm of 0.625mm or 40mm of 1.25mm acquisitions
- Volara* XT Digital DAS (Data Acquisition System): The Volara* XT digital DAS for faster sampling and improved image performance and reduced artifacts
- Fast coverage speed of 110mm/sec
- Diode technology providing true 32 channel acquisition and is ready for future expansion
- Full 360 degree rotation in 0.4 (pediatric), 0.5, 0.6, 0.7, 0.8, 0.9, 1.0 and 2.0 (axial) seconds, ensuring short breath holds, comfortable exams and flexibility to customize protocols for unique patient needs with minimal coverage impact
- Routine thin slice scanning, as thin as 0.625mm or 1.25mm optimizing the use of thinner images for sagittal, coronal, oblique, and volume image presentation and review
- Highly efficient compact geometry design delivering optimum performance of the x-ray tube and generator
- Image decomposition to:
 - Retrospective thin images from data sets where thicker images were initially reconstructed
 - Facilitates more detailed image analysis
 - Improves 3D and reformat visualization
- Neuro 3D Filter provides users the capability to filter head acquisition data using specially designed and optimized 3D

Item No. Qty

Description

filters.

Neuro 3D Filter is not available when ASiR is implemented.

Fast, User-Friendly Simultaneous Workflow:

- Advanced Workflow Platform, the next evolution of GE's workflow platform built to help you maximize productivity.
 - Delivers up to 16 images per second (lips) reconstruction
 - Up to 10 fps network transfer rates
 - Direct Multiplanar Reformats (DMPR) that enables the move from 2D review to prospective 3D review of sagittal, coronal and oblique planes automatically
 - Data Export and Interchange that allow you to easily share images with referring physicians and patients
- Includes reference protocols and the ability to customize your own for a total of 6,840 programmable protocols
- SmartPrep with Dynamic Transition allows low dose intermittent monitoring of intravenous contrast enhancement in a user-selected section of anatomy. With Dynamic Transition when the prescribed contrast enhancement is reached the system will automatically transition from the monitoring phase to the scan phase
- 10 Prospective Multiple Reconstructions: Up to 10 reconstructions can be pre-programmed as part of the scan protocol prior to acquisition. The operator can select different start/end location, slice thickness, interval, interval reconstruction algorithms and display fields of view for each reconstruction. Assisting to prospectively prescribing the image reconstructions needed, even for complex trauma exams and freeing the user up to focus on the patient
- Remote tilt from the operator console to increase exam speed
- Built-in breathing lights with a countdown timer, so the patient does not have to guess

how much longer to hold their breath

- New built-in 12-inch touch screen gantry display allows technologists to make personalized exams by displaying the patient's name on it. When not scanning, the video of relaxing scenes or cartoons may have a calming effect on children or patients of all ages
- By using the Default patient positioning on built-in 12-inch touch screen gantry display the bed provides automatic positioning according to the type of exam, reducing manual positioning and streamlining workflow
- In room start button mounted on gantry with countdown display, facilitates single technologist operation and improved departmental productivity
- GE software allows you to automate or build every task into

the protocols to increase throughput

- Has up to 250,000 uncompressed 512 x 2 image files storage capacity, and 3,520 scan rotations, or up to 1,500 scan data files, or up to 300 exams

Dose Management Leadership:

- OptiDose management features: new bowtie filters optimized for adult and pediatric body exams, full 3D dose modulation, color coding for kids, tracking collimator hardware and software for x-ray beam tracking to name a few of GE's dose optimization features, all based on the ALARA principle
- Dynamic Z-axis tracking provides automatic and continuous correction of the x-ray beam shape to block unused x-ray at the beginning and end of a helical scan to reduce unnecessary patient radiation
- 3D Dose modulation - Before the scan, clinicians must select the desired Noise Index as well as the minimum and maximum mA setting. The system automatically accounts for the changing dimensions of the patient's anatomy enabling patient to patient reproducibility in this aspect of image quality and real-time x-y-z during each scan
- Tracking collimator hardware and software for x-ray beam tracking to minimize patient dose
- Filtration of the x-ray beam is optimized independently for body and head applications
- DLP (dose length product), and dose efficiency display during scan prescription provides the patient's dose information to the operator
- Dose Reporting provides access to the CTDIvol and DLP with the patient record prior and post exam. DICOM Structured Dose Report is also supported.
- Dose Check provides the user with tools to help them manage CT dose in clinical practice and is based on the standard XR-25-2010 published by The Association of Electrical and Medical Imaging Equipment Manufacturers (NEMA). Dose Check provides the following:
 - Checking against a Notification Value if the estimated dose for the scan is above your site established value
 - Checking against an Alert Value where the user needs specific authority to continue the scan at the current estimated dose without changing the scan parameters if the estimated dose exceeds the alert value
 - The ability to define Alert Values for Adult and Pediatric with age threshold
 - Audit logging and review capabilities
 - Protocol Change Control capabilities

The Advanced Reconstruction breaks through existing limits on speed, image quality and flexibility to provide an optimized volumetric workflow solution from acquisition to final report and has the capability to deliver up to 16 full fidelity images per second (lips) reconstruction and 10 fps network transfer rates.

Clinical Benefits:

- CTA runoffs
- Thin slices fast; routine use of thin slices
- Organ coverage in arterial phase
- Long helical scans
- Multi-phase organ studies
- Improved multi-planar reformats with isotropic microvoxel imaging
- Fast scanning with outstanding image performance and GE's proprietary cross beam and hyperplane helical reconstruction algorithms
- System designed for optimization of z-axis resolution and dose with 0.625mm slice thickness

System Components:

Gantry:

- Advanced slip ring design continuously rotates the generator, Performix 40 X-ray tube, detector and Volara XT digital data acquisition system around the patient.
 - Aperture: 70 cm
 - Maximum SFOV: 50 cm
 - Rotational Speeds: 360 degrees in 0.4 0.5, 0.6, 0.7, 0.8, 0.9, 1.0 and 2.0 (axial) seconds
 - Tilt: +/- 30 degrees, speed 1 degree/sec
 - Remote tilt from operators console
 - Integrated breathing lights and countdown timer
 - Integrated 12-inch touch screen on gantry with workflow features
 - Integrated start scan button with countdown timer to indicate when x-ray will turn on
- Visual readout is easy to read from the tableside or from the operator console. Gantry tilt controls are located on the side of the gantry.

Laser Alignment Lights:

- Defined internal and external scan planes to +/- 1mm accuracy
- Operate over full range of gantry tilt
- Coronal light remains perpendicular to axial light as gantry tilts

Table:

- Cantilever design for easy access
- Vertical range: 43.0 cm to 99.1 cm
- Vertical scannable range: 79.1 cm to 99.1
- Horizontal range: 1,745 mm (VT1700 Table), or 2,045 mm (VT 2000 Table)
- Horizontal speed: up to 137.5 mm/sec
- Table load capacity: 227 kg (500 lb) +/- 0.25mm positional accuracy

X-ray Tube: Performix 40 metal-ceramic tube unit

- Performix 40 tube with 6.3 MHU of storage and capable of 72kW operation provides increased helical performance with greater patient throughput
- Wide range of technique (10 mA to 560 mA, in 5 mA increments) gives technologist and physician flexibility to tailor protocols to specific patient needs, while optimizing patient dose, and providing the power needed to perform a broad spectrum of examinations.
- Maximum anode heat storage capacity: 6.3 MHU
- Dual Focal Spots:
 - Small Focal Spot: 0.9 x 0.7 IEC60336:2005
 - Large Focal Spot: 1.2 x 1.1 IEC60336:2005
- Maximum power: 72 kW
- Beam collimated to 56 degree fan angle

High Voltage Generator: High Frequency on-board generator allows for continuous operation during scan.

- 72 kW Output Power
- kV: 80, 100, 120, 140 kV
- mA: 10 to 560 mA, 5 mA increments

Maximum mA for Each kV Selection (large focal spot):

- 400mA @ 80kV
- 480mA @ 100kV
- 560mA @ 120kV
- 515mA @ 140kV

V-Res Detector: The V-Res detector was designed for high performance imaging. V-Res detector benefits are:

- Solid 40mm coverage per rotation
- GE's exclusive patented detector material

Volara XT Digital DAS (Data Acquisition System): The Volara XT digital DAS dramatically reduces electrical noise for improved imaging performance.

- 2,460Hz maximum sample rate
- Effective analog to digital conversion

Optima CT660 Operator Console:

- 1,792GB of total system storage
- Up to 250,000 512 x 2 images and 3,520 scan rotations or up

| Item No. Qty. | Description |
|---------------|--|
| | <p>to 1,500 scan data files, or up to 300 exams</p> <ul style="list-style-type: none"> • 4.7 GB DVD-R/CD-R for DICOM interchange (not recommended as a long term archive) <p>Image Networking: Exams can be selected and moved between the Optima CT660 CT</p> <p>System and any imaging system supporting DICOM protocol for network send, receive and pull/inquiry.</p> <ul style="list-style-type: none"> • Standard Auto-configuring Ethernet • Direct Network Connection • Supports 1GB or 1000/100/10 BaseT <p>DICOM Conformance Standards</p> <ul style="list-style-type: none"> • DICOM Storage Service Class • Service Class User (SCU) for image send • Service Class Provider(SCP)for image receive • DICOM Query/Retrieve Service Class • DICOM Storage Commitment Class Push • DICOM Modality Worklist (incl. Performed Procedure Step) (through ConnectPro option) • DICOM Print <p>The Optima CT660 workflow platform is designed to deliver high performance in each of these tasks:</p> <ul style="list-style-type: none"> • SmartTools Simplifies Scan Setup and Includes All Reconstructions, Filming, Archiving, Transferring Prospectively • Workflow platform built on the LINUX operating system delivers up to 16 fps reconstruction and the fast network transfer rates of up to 10 fps • Data Export and Interchange allow you to easily share images with referring physicians and patients • Direct MPR that enables the move from 2D review to 3D image review of axial, sagittal, coronal and oblique planes automatically • Exam Split delivers the capability to split a series of patient images into separate groups for networking • Exam Rx desktop environment provides the clinical tools desired for fast, efficient control of patient studies. Exam Rx |

| Item No. Qty | Description |
|--------------|--|
| | <p>tools include patient scheduling and data entry, exam protocol selection, protocol viewing and editing, scan data acquisition, image display and routine analysis, AutoTransfer, AutoStore, and AutoFilm</p> <ul style="list-style-type: none"> ImageWorks is a desktop environment designed to take advantage of the Optima CT660 CT System advanced computer systems. Standard features include archive, network and manual film control, as well as some advanced image processing such as Direct multi-planar reformatting (DMPR), multi-projection volume rendering (MPVR) and display. The ImageWorks desktop also provides a gateway for DICOM 3.0 image transactions, either through a local area network, or via DICOM-formatted media Volume Viewer includes Volume Analysis, Volume Rendering and Navigator software. This combination allows the user to render volumetric data in three dimensions for use in analysis of patient condition, i.e. CT Angiography (CTA), gives more information on the spatial relationships of structures than standard 3D, allows the translucent visualization of structures for improved problem solving, can perform "virtual endoscopies" of air and contrast filled structures. Enables 3D reformats in any plane, ALL on the Xstream ready console. <p>Scan Modes: The Optima CT660 system can perform virtually any clinical application due to its wide variety of scan modes. Helical scan mode offers continuous 360 degree scanning with table incrementation and no interscan delay. Axial scan mode allows for up to 32 contiguous axial slices acquired simultaneously with each 360 degree rotation.</p> <ul style="list-style-type: none"> Helical scanning pitches: 0.516:1, 0.984:1, 1.375:1 Retrospective reconstruction image thicknesses: 32 x 0.625, 32 x 1.25, 64 x 0.625* <p>* Available only with Overlapped Reconstruction option (axial mode & 20mm coverage)</p> <p>Scan Enhancements:</p> <ul style="list-style-type: none"> Anatomical programmer: a ten region anatomical selector allows quick and easy access to user programmable protocols and a separate selector for adult and pediatric exams with greater than 6,840 protocol storage available |

| Item No. | Qty | Description |
|----------|-----|---|
| | | <ul style="list-style-type: none"> • Protocols include preset scan time, kV, mA, scan mode, image thickness and spacing, table speed, scan FOV, display FOV and center, recon algorithm, and special image acquisition and processing options like DMPR • Any scan parameters may be edited for each scan or all scans - either before or during an exam. The number of scans may also be easily changed • AutoScan: Automates longitudinal table movement and start of each scan • Auto-Voice: 3 preset (9 languages) and 17 user defined messages automatically deliver patient breathing instructions, especially useful for multiple helical scanning • Trauma Patient: Allows patient scans and image display/analysis without entering patient data before scanning • Reconstruction Algorithms: Soft Tissue, Standard, Detail, Chest, Bone, Bone Plus, Lung, and Edge <p>Warranty: The published Company warranty in effect on the date of shipment shall apply. The Company reserves the right to make changes. All specifications are subject to change. Regulatory compliance: This product is designed to comply with applicable standards under the radiation control for Health and Safety Act of 1968.</p> <p>Laser alignment devices contained within this product are appropriately labeled according to the requirements of the Center for Devices and Radiological Health.</p> <p>Siting Considerations: See the Pre-Installation manual for details of the siting requirements for the Optima CT660.</p> <p>This product is a CE-compliant device that satisfies IEC60601-1:1998 and applicable collateral and particular standards, including regulations regarding Electro-Magnetic Compatibility (EMC) and Electro-Magnetic Interference (EMI), pursuant to IEC-60601-1-2:2004.</p> |
| 2 | 1 | <p>English keyboard and Label Kit</p> <p>English Keyboard (Black) for CT systems and system labels</p> |

| Item No. Qty | | Description |
|--------------|---|--|
| 3 | 1 | <p>Optima CT660 Cable Set</p> <p>Optima CT660 Long Cable set</p> |
| 4 | 1 | <p>VT2000 Table</p> <p>The CT system 2000 table enables volume scanning. Key features of the VT 2000 table include: 500 lb weight capacity, 2000 mm scannable range, 175 mm/sec travel time, real-time position control to support advanced application such as SnapShot Pulse, VolumeShuttle, and Volume Helical Shuttle.</p> |
| 5 | 1 | <p>Freedom Workspace Desk - Adjustable</p> <p>Freedom Workspace is the desk designed with ergonomics. This table design enables the efficient use of space while enhancing clinical workflow and technologist comfort. Attributes: Fully adjustable monitor arms Adjustable height Flexible location of Console hardware</p> <p>Benefits: Improve patient visibility Clear path to patient More comfortable for technologist Improved ergonomics for technologist Sitting or Standing position Easy height adjustment</p> <p>Requirements: Tables are 3M apart from Operator console</p> |
| 6 | 1 | <p>ASIR - Adaptive Statistical Iterative Reconstruction for Optima CT660 system</p> <p>ASIR - Adaptive Statistical Iterative Reconstruction for Optima CT660 system</p> |
| 7 | 1 | <p>AW VolumeShare 5 with Two Flat Panel Monitors and 6GB of RAM</p> <p>AW VolumeShare 5 with Two Flat Panel Monitors and 6GB of RAM.</p> <p>AW VolumeShare 5 is a multi-modality image review, comparison and post processing workstation built with simplicity and power at its core. Powerful software is optimized to take advantage of state of the art 64 bit technology and multiple cores to ensure leading edge performance.</p> <p>AW VolumeShare 5 features include:</p> <p>Hardware:</p> <ul style="list-style-type: none"> • HP Z800 Workstation with Intel x5650 Six Core Xeon 2.66 GHz CPU with 8MB Shared L2 Cache / 1333 MHz Dual FSB |

| Item No.. | | Description |
|-----------|---|--|
| | | <ul style="list-style-type: none"> • 6GB DDR-3 1333 ECC DIMM • 300GB SAS 15,000rpm Hard Disk for OS and Apps. • 600GB SAS 15,000rpm Hard Disks for Image Data • 2 x 19" EIZO MX191 Monitors |
| | | <p>Software:</p> <ul style="list-style-type: none"> • Fast access to information you need through optional RIS integration & priors post-fetch • Efficient workflow through dynamic load, end review and Key Image Notes features • Optional productivity package to pre-process exams and allow up to 8 simultaneous sessions • Applications usage monitor to track usage of your system • Smart layouts with Volume Viewer General review protocol that optimizes comparison and single exam layouts • Enhanced multi-modality contouring tool with support for PET SUV's • Support for external DICOM USB media and preference management tool to exchange preferences across users • Support for optional, broad suite of multi-modality advanced applications |
| 8 | 1 | VS5 S/W ONLY UPGR.F/XW840 |
| 9 | 1 | <p>VolumeShare 5 Software Only Upgrade VESSELIQ SL AUTOBONE</p> <p>VesselIQ Xpress & AutoBone Xpress</p> <p>CT VesselIQ Xpress and Autobone Xpress is for AW VolumeShare5</p> <p>VesselIQ Xpress provides an optimized non-invasive application to analyze vascular anatomy and pathology and aid in determining treatment plans from a set of CTA images. This software supports the physician in:</p> <ul style="list-style-type: none"> • Assessment of aneurysms with or without thrombus (false lumen) for size and volume measurements with the capability to track the size and volume over time, stenosis analysis, pre/post stent and surgical planning and directional vessel tortuosity visualization. • Automatic tools for the segmentation of bony structures in the |

| Item No. | Description |
|----------|---|
| 101 | <p data-bbox="586 178 1291 279">brain and neck and other vascular areas for accurate identification of the vessels, single or double click vessel analysis.</p> <ul data-bbox="537 296 1312 751" style="list-style-type: none"> • Sizing the vessel, analyzing calcified and non-calcified plaque to determine the densities of plaque within a vessel, measure areas of abnormalities within a vessel (like stenosis, plaque, thrombus, dissection or leakage). • Semi-automated detection and segmentation of thrombus for subsequent measurements within the application. • Dedicated anatomy based protocols for improved workflow. • Compare a patients previous exam to their current exam in order to measure and track any changes over time of their vascular structures. • After review of the exams, there are multiple ways to film, archive and capture information for future review. <p data-bbox="518 781 771 808">System Requirements:</p> <ul data-bbox="537 812 841 837" style="list-style-type: none"> • AW VolumeShare5 <p data-bbox="518 861 1247 921">Note: All software are Non-Transferable to other hardware and are Non-Returnable.</p> <p data-bbox="532 932 1175 959">10 KVA Partial UPS for CT LightSpeed and LightSpeed PRO</p> <p data-bbox="522 978 1260 1005">2 Phase 10 KVA Partial UPS for CT Lightspeed and Lightspeed PRO</p> <p data-bbox="522 1026 1276 1453">The 2 Phase 10 KVA Partial System UPS kit has been specifically designed to coordinate with the BrightSpeed, LightSpeed and LightSpeed PRO 16 families of CT scanners. In the event of a power outage, a partial system UPS provides continuous back-up power to the scanner host and control computers, thus assuring no loss of usable scan data. In addition, critical circuits in the gantry and table remain powered which facilitate the safe removal of the patient from the scanner. If power is restored within the battery hold-up time, the operator can continue scanner operations without the need to reboot the system. When longer power outages are anticipated, the UPS provides time for the operator to complete an orderly shutdown of the system software.</p> <p data-bbox="522 1480 761 1507">FEATURES/BENEFITS</p> <ul data-bbox="537 1535 1252 1598" style="list-style-type: none"> • True double-conversion, online technology provides reliable operation and uninterrupted glitch free power. |

| Item No. | Qty | Description |
|----------|-----|---|
| | | <ul style="list-style-type: none"> Automatic voltage and frequency selection eases startup, i.e., 50 or 60 Hz compatible Integral Static Bypass switch means zero transfer time Integral Manual Bypass switch facilitates continued scanner operation while UPS is being serviced Single input connect utilized for both UPS input and static switch Maintains system electronics and allows critical scanner operations to continue for 10 minutes (typical) after loss of power Advanced Battery Management (ABM) software monitors / indicates battery health and doubles battery service life |
| | | SPECIFICATIONS |
| | | <ul style="list-style-type: none"> Dimensions (H x W x D): 32.7" x 12" x 32" Weight: 350 lbs. Rating: 10 kVA Input Voltage Range: 85-144V / ph; 2 Phase Output Frequency: 50 or 60 Hz, auto-sensing |
| | | COMPATIBILITY |
| | | <ul style="list-style-type: none"> HiSpeed Advantage-RP, CT/I, Lightspeed QXi, LightSpeed Plus, LightSpeed Ultra, LightSpeed 16, BrightSpeed Systems, LightSpeed Pro 16 and RT Systems, Discovery NM 670 (Nuc) |
| | | NOTES: |
| | | <ul style="list-style-type: none"> Customer is responsible for rigging and arranging for installation with a certified electrician ITEM IS NON-RETURNABLE AND NON-REFUNDABLE |
| 11 | 1 | <p>90 Amp Main Disconnect Panel for CT</p> <p>90 Amp Main Disconnect Panel for CT</p> <p>This 90 amp main disconnect panel for GEHC CT systems provides emergency shut down, undervoltage protection, overcurrent protection, local disconnect for the imaging system. It also reduces installation time and cost by providing a single-point power connection eliminating the need to mount and wire a number of individual components. The standardized design and testing</p> |

| Item No. | Qty | Description |
|----------|-----|--|
| | | <p>assures high product quality and system reliability, and it is UL and cUL listed for compliance with National Electric Code. Panel can be surface or semi-flush mounted and includes one remote emergency off push button. Customer is responsible for rigging and arranging for installation by a licensed electrician. ITEM IS NON-RETURNABLE and NON NON-REFUNDABLE Warranty Code: Y</p> |
| 12 | 1 | <p>CT Table Slicker with Cushion - VCT 2000 Systems (2-pc Set)</p> <p>Slicker - CT HD750 and VCT w/GT 2000 Table (2 Piece Set)</p> <p>FEATURES/BENEFITS</p> <ul style="list-style-type: none"> Two-piece, sealed slicker cushion set has comfort pads enclosed inside the slicker cover and extender cover Durable, clear PVC plastic cover facilitates faster, more thorough cleanup of blood and fluids Increase system uptime by protecting table from spills and particulate contaminants Thermo-sealed seams and flaps prevent contaminate buildup in hard to clean areas <p>COMPATIBILITY</p> <ul style="list-style-type: none"> VCT with GT 2000 Table, CT HD750 |
| 13 | 1 | <p>CT Footswitch Slicker - VCT 2000 & 1700 Systems</p> <p>Footswitch Slicker for CT HD750 and VCT Systems</p> <p>The footswitch slicker for CT VCT 2000 and 1700 systems is made of durable, clear PVC plastic that protects the footswitch and facilitates faster, more thorough cleanup of contamination caused by blood and other body fluids. Cover is held securely in place with Velcro...H</p> |
| 14 | 1 | <p>6 Day CT TiP Onsite System Training</p> <p>6 Day CT TIP Onsite System Training</p> <p>CT Onsite Training for a new CT system</p> <ul style="list-style-type: none"> One 4 day onsite visit to coincide with system start-up. One 2 day onsite follow-up visit 6-8 weeks post system start up. |

| Item No. | Qty | Description |
|----------|-----|---|
| 15 | 1 | <p data-bbox="526 174 1276 453">During the first visit, the applications specialist will work with the medical and technical staff on system operation and patient procedures. The training produces the best results when a dedicated core group of 2-4 CT technologists complete the session with a modified patient schedule. It is suggested that key physicians are available to participate in the protocol implementation and image quality review sessions. By the end of this visit, the core group should be able to perform the routine patient procedures.</p> <p data-bbox="526 478 1276 720">The 2 day revisit is suggested after the staff has run the system for 6-8 weeks, however this is flexible based on the site needs. The training will focus on the intermediate and advanced functions of the system or special needs of the customer. The training produces the best results when the same dedicated core group of 2-4 CT technologists from the initial visit complete the session with a modified patient schedule.</p> <p data-bbox="526 745 1230 810">This training program must be scheduled and completed within 12 months after the date of product delivery.</p> <p data-bbox="526 821 708 848">4 Days Ct Onsite</p> <p data-bbox="526 863 846 890">4 Days CT TiP Onsite Training</p> <p data-bbox="526 911 1211 1003">Four Days CT Onsite Training provided from 8AM to 5PM, Monday through Friday. Includes TELL expenses. Days provided consecutively.</p> <p data-bbox="526 1024 1224 1089">This training program must be scheduled and completed within 12 months after the date of product delivery.</p> |
| 16 | 1 | <p data-bbox="526 1115 976 1142">TiP HQ Class Optima CT660 - Full Service</p> <p data-bbox="526 1157 980 1184">TiP HQ Class Optima CT660 - Full Service</p> <p data-bbox="526 1205 1208 1270">3.5 day CT course held in the Milwaukee area. Includes travel and modest living expenses.</p> <p data-bbox="526 1291 1230 1356">This course is designed to introduce the technologist to the Optima CT660 system.</p> <p data-bbox="526 1371 1224 1436">This training program must be scheduled and completed within 12 months after the date of product delivery.</p> |
| 17 | 1 | <p data-bbox="526 1524 842 1551">CT LIGHTSPEED PRO ADV SER</p> <p data-bbox="526 1566 1224 1593">The LightSpeed Pro Advanced course is intended for engineers</p> |

| Item No. | Qty | Description |
|----------|-----|--|
| 18 | 1 | <p data-bbox="524 174 1247 275">servicing LightSpeed Pro 16, LightSpeed RT, and forward production LightSpeed 16/Ultra/Plus (starting in 2004) systems. This course must be taken within 2 years from the purchase date.</p> <p data-bbox="524 304 773 331">CT LS 7X St OPTIMA 660</p> <p data-bbox="524 361 1027 388">CT LightSpeed 7x and Optima 660 (Class/Lab)</p> <p data-bbox="524 411 1258 688">The CT LightSpeed 7x & Optima 660 course is a differences class and is intended for Engineers who have completed (R0026CT) LightSpeed Pro Training. It will equip the Engineer with system and subsystem theory and hands-on lab activities to address technical service issues for the 32/64-slice family of scanners (including LightSpeed VCT, LightSpeed VCT XT, LightSpeed VCT Select, and Optima 660. This training must be used within 2 years from the purchase date.</p> |
| 19 | 15 | <p data-bbox="524 726 824 753">Meals And Lodging Expense</p> <p data-bbox="524 777 1235 911">Meals and Lodging Expense has been developed to allow the customer the convenience of prepaying for their meals and lodging expenses when attending Technical Service Training at the GE Healthcare Institute located in Waukesha, WI.</p> <p data-bbox="524 940 1268 1113">The price of this convenience is based on a per day basis. Thus a quantity of 1 is equal to 1 day's meals and lodging expense. When purchasing the meals and lodging expense please be mindful of weekend days during the training stay and include 2 days to cover a weekend in the purchase quantity.</p> <p data-bbox="524 1136 1278 1379">Examples: A 5-day course needs a quantity of 5. Any course longer than 5 days should include 2 days to account for the weekend stay. Any course longer than 10 days will require an additional 4 days of the meals and lodging expense to cover the 2 weekends of the stay. Thus a 15-day course would have a quantity of 19 days to cover the 2 weekends of the stay. This expense must be used within 2 years from the purchase date.</p> <p data-bbox="524 1407 1235 1579">Three meals a day Monday thru Thursday, 2 meals on Friday, pluse breaks are provided in the onsite cafeteria. The GE Healthcare Institute cafeteria closes Friday after lunch and reopens Monday morning for breakfast. Weekend meals are the responsibility of the customer.</p> <p data-bbox="524 1606 1208 1633">Only for In-resident courses to be taken at the GE Healthcare</p> |

| Item No. | Qty | Description |
|----------|-----|--|
| | | Institute. |
| 20 | 2 | <p>Airfare Expense</p> <p>The AIRFARE EXPENSE has been developed to allow the customer the convenience to prepay their roundtrip Airfare expenses when attending Technical Service Training at the GE Healthcare Institute located in Waukesha, WI. To be used for engineers attending In-Resident Class/Lab courses for Diagnostic Imaging.</p> <p>Customer will make their Airfare arrangements thru the GE Travel Center. Specific directions will be provided to the customer upon confirmation of class. Please note that this expense must be used within 2 years of the purchase date</p> |
| 21 | 1 | <p>Lodging Weekend Expense</p> <p>Lodging Weekend Expense</p> <p>Weekend Lodging Expense is to cover Saturday and Sunday lodging expenses for those engineers who are staying at the Rivers Edge Condos while attending Diagnostic Imaging Biomed training at the Healthcare Institute. Please note that there are no meals included on the weekend. Must be used within 2 years from the purchase date.</p> |
| 22 | 1 | <p>CT Basic Physics/Instrumentation (web)</p> <p>CT Basic Physics/Instrumentation (Web)</p> <p>The CT Fundamentals Course is Designed for Service Engineers who have Little or No Familiarity with CT Systems. The Course Teaches General Processes, Concepts, and Equipment Used in CT Scanning. This Course is Delivered Via the internet as an online training course. This course must be taken within 2 years from the purchase date.</p> |
| 23 | 1 | <p>CT LIGHTSPEED PRO ADV SVC</p> <p>CT Lightspeed Pro Advanced Service (Web)</p> <p>Web course is 8 hours long</p> <p>Sales Description:</p> <p>Introduction to CT LightSpeed Pro system theory and subsystems</p> <p>Executive Summary:</p> |

| Item No. | Qty | Description |
|----------|-----|---|
| | | <p>This is a computer-based training course intended to prepare Service Engineers on basic system theory for the LightSpeed Pro product line.</p> <p>Course Competencies:</p> <p>The curriculum builds on concepts taught in CT Basic Physics and is a prerequisite for the CT LightSpeed Pro and Discovery ST in-resident training classes at the GE Healthcare Institute.</p> <p>Special Considerations:</p> <p>A functioning laptop computer with a CD-ROM reader, network card and a modem card is required for use during this course. The browser on the computer must be IE4 or Netscape 4.5 or higher. Minimum system requirements include 133 MHz Windows 95, NT 4.0 or higher 32 MB of RAM 16-bit color display adapter. Proof of completion of this eLearning course is necessary prior to attending any subsequent GE Healthcare In-Resident training. This course contains proprietary content. For customers attending this course, special paperwork is required to take this course. Please see the registration page for details on the enrollment process. This course must be taken within 2 years from the purchase date.</p> |
| 24 | 1 | <p>GLOBAL OPERATOR CONSOLE</p> <p>CT GLOBAL OPERATORS CONSOLE 3,4,& 5</p> <p>The Global Operators Console can be referred to as the Xtreme console as well. This is the current operator console for the CT LightSpeed and PET Discovery ST systems. This course must be taken within 2 years from the purchase date.</p> |
| 25 | 1 | <p>OPTIMA CT660 SERVICE (WEB)</p> <p>Optima CT660 Service (web)</p> <p>This upgrade course taken online is intended for Support Engineers who have previous LightSpeed VCT training. Topics covered include: New gantry display, new power saving mode, new gantry axial motor and control, new gantry covers removal and installation, safety awareness with gantry cover mounting hardware, new operators console (RIO), load from cold-Saturn detector. This course must be taken within 2 years from the purchase date or it expires without refund.</p> |

| Item No. | Qty. | Description |
|----------|------|--|
| 26 | 1 | <p data-bbox="524 254 935 279">Troubleshooting Basics Service (web)</p> <p data-bbox="524 296 935 321">Troubleshooting Basics Service (Web)</p> <p data-bbox="524 338 1274 657">This Course is Intended for Individuals Involved in Servicing Medical Equipment. By Taking This Course, You will Learn a Proven Process for Troubleshooting Problems with Medical Equipment. You will Also Learn How to Use Various Tools in a Troubleshooting Situation and How to Interpret Error Messages. This Course Does Not Address How to Troubleshoot Specific Products. It is Recommended That you Have Fundamental Training in a Modality Prior to Taking This Course. This course must be taken within 2 years from the purchase date.</p> |
| 27 | 1 | <p data-bbox="524 709 849 735">NETWORKING & DICOM BASIC</p> <p data-bbox="524 751 1065 777">Networking and Dicom Basic for DI Service (Web)</p> <p data-bbox="524 800 1247 903">Training will prepare engineers on configuring and troubleshooting networks, which use the DICOM protocol for transferring patient data and how to read and use DICOM Conformance Statements.</p> <p data-bbox="524 928 881 953">This course covers the following:</p> <ul data-bbox="544 982 1255 1255" style="list-style-type: none"> <li data-bbox="544 982 1255 1045">• Introduction to 7 layer OSI and 5 layer TCP/IP protocols (Basic model only) <li data-bbox="544 1062 1040 1087">• Identify hardware used in networking <li data-bbox="544 1104 1239 1167">• Review of the most used networking devices, cables, NIC, switch and routers <li data-bbox="544 1184 1138 1209">• Simple network connection with 2 to 5 devices <li data-bbox="544 1226 1114 1251">• Dicom definitions, theory and configuration <p data-bbox="524 1276 1255 1302">This course must be taken within 2 years from the purchase date.</p> |