

## *Accuray System Base Configuration and Additional Options*

---

### **A. TomoTherapy® TomoHDA™ Treatment System**

---

1. **TomoHDA™ Ring Gantry System**
  - 1.1. Treatment Delivery System
    - 1.1.1. 6MV 850MU/Minute Fixed Target Linear Accelerator (Linac)
    - 1.1.2. Primary Collimator System
    - 1.1.3. Fixed 1cm, 2.5cm and 5cm treatment beam widths
    - 1.1.4. Dynamic 2.5 cm and 5 cm treatment beam widths
    - 1.1.5. Binary Multi-Leaf Collimator System
    - 1.1.6. Dose Control System
  - 1.2. Image Guidance System
    - 1.2.1. High-Resolution MVCT Detector System
  - 1.3. Enclosure Assembly with 12" Positioning Control Panels
  - 1.4. Noise Eliminating Intercom System
  - 1.5. Laser Positioning System
  - 1.6. High-Performance Treatment Couch
2. **Operator Station**
  - 2.1. Operator Station and Status Console
    - 2.1.1. Operator Station Couch Control
3. **Standard Treatment Delivery Modes**
  - 3.1. TomoHelical™ treatment delivery mode
  - 3.2. TomoDirect™ treatment delivery mode
4. **System Interface**
  - 4.1. DICOM CT and RTStructure Set Import
  - 4.2. DICOM Export (CT, RTStructure, RTDose, RTPlan)
  - 4.3. OIS Connect™ (MOSAIQ®/ARIA® Interface)\*
5. **Treatment Planning System**
  - 5.1. Two (2) TomoTherapy Treatment Planning Workstations
  - 5.2. VoLO™ Optimization Cluster – High-Performance Option
  - 5.3. Planned Adaptive™ Software
6. **Remote Access**
  - 6.1. Remote Software Solutions Package†
7. **Data Management System**
  - 7.1. Data Server and SAN
    - 7.1.1. TomoGateway™ Node
    - 7.1.2. TomoLink™
  - 7.2. Data Management System (DMS) Software
    - 7.2.1. Patient Self Transfer
    - 7.2.2. Workstation Configuration Option for Patient Data Transfer
8. **Quality Assurance**
  - 8.1. TomoTherapy Quality Assurance (TQA™) Essentials
9. **User Documentation**
  - 9.1. TomoTherapy electronic User Documentation package

\* **OIS Connect.** Customer acknowledges that in order to take advantage of the OIS Connect interface functionality, Customer must separately purchase and install the corresponding version of the third party vendor software and any additional software/options required by that vendor.

† The Remote Software Solutions package includes two (2) Remote Planning user license and two (2) TomoPortal™ Remote Viewer user licenses. Additional user licenses may be purchased as Additional Options as described below.

---

**B. Additional Options Available for Purchase**

- The following Additional Options are available for purchase by Customers.
- By indicating a Quantity below, Customer may indicate any Additional Options Customer wishes to purchase, with the understanding that such selections will be subject to additional costs and are to be added to the “Total Due” amount indicated in Section A of the Agreement and Customer agrees to pay the applicable Unit Price(s), plus any taxes and other applicable fees.

| Quantity | Additional Options                                                                                                                                                                                                                                                             |  |  |   |
|----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|---|
| 0        | <b>1. Additional TomoTherapy Planning Workstation</b><br>Two (2) TomoTherapy treatment planning workstations are included with the base TomoHDA configuration; however, one (1) additional TomoTherapy workstation may be purchased for a maximum total of three workstations. |  |  |   |
| 0        | <b>2. Additional Planning Workstation Station Extended Network Option</b>                                                                                                                                                                                                      |  |  |   |
| 0        | <b>3. Network Data Storage</b>                                                                                                                                                                                                                                                 |  |  |   |
| 1        | <b>4. Power Conditioner System</b><br><i>Purchase of a Power Conditioner for the TomoHDA system is strongly recommended</i>                                                                                                                                                    |  |  |   |
| 0        | <b>5. HIS Workflow Connectivity</b>                                                                                                                                                                                                                                            |  |  |   |
| 0        | <b>6. DICOM Transit Dose Export</b>                                                                                                                                                                                                                                            |  |  |   |
| 0        | <b>7. StatRT™ Software</b>                                                                                                                                                                                                                                                     |  |  |   |
| 0        | <b>8. Data Services Package – SharePlan™</b>                                                                                                                                                                                                                                   |  |  |   |
| 0        | <b>9. Data Services – SharePlan™ One Year Maintenance Agreement</b>                                                                                                                                                                                                            |  |  |   |
| 0        | <b>10. Remote Software Solutions – Remote Planning- Third Client License Option</b>                                                                                                                                                                                            |  |  |   |
| 0        | <b>11. TomoPortal™ Client License Increase – from 2 to 4 users</b>                                                                                                                                                                                                             |  |  |   |
| 0        | <b>12. TomoPortal™ Client License Increase – from 2 to 8 users</b>                                                                                                                                                                                                             |  |  |   |
| 1        | <b>13. Standard TomoTherapy QA Package</b><br>Includes: ‘cheese’ phantom, solid water phantom and ion chambers                                                                                                                                                                 |  |  |   |
| 1        | <b>14. Beam Measurement and QA package</b><br>Includes: PC, software, TomoElectrometer and water tank                                                                                                                                                                          |  |  |   |
| 0        | <b>15. Beam Measurement and QA package (without water tank)</b><br>Includes: PC, software and TomoElectrometer                                                                                                                                                                 |  |  |   |
| 0        | <b>16. TomoElectrometer™</b>                                                                                                                                                                                                                                                   |  |  |   |
| 0        | <b>17. Calibrated Mini-Ion Chamber and Interconnect Cable</b>                                                                                                                                                                                                                  |  |  |   |
| 0        | <b>18. Calibrated CT Slice Ion Chamber, Buildup Cap and Jig, and Interconnect Cable</b>                                                                                                                                                                                        |  |  |   |
| 0        | <b>19. Mini-Ion Chamber Interconnect Cable - 30 meter</b>                                                                                                                                                                                                                      |  |  |   |
| 0        | <b>20. Film Digitizer Kit (with PC, Software and Vidar Scanner)</b>                                                                                                                                                                                                            |  |  | 0 |
| 0        | <b>21. Film Analysis Kit for Existing Vidar DosimetryPRO® Advantage or Flatbed Scanner</b>                                                                                                                                                                                     |  |  |   |
|          | <b>22. Tomo Quality Assurance (TQA™) Total Package</b>                                                                                                                                                                                                                         |  |  |   |
|          |                                                                                                                                                                                                                                                                                |  |  |   |