

VAMC BATTLE CREEK, MI
PO# 515-B30006

| Qty | Product Model Number and Description | | | | |
|-----|---|--|--|--|--|
| 1 | DISCOVERY W - BONE DENSITOMETER Discovery W One Pass Linear Fan-Beam Acquisition DXA Technology 128 Element High-Resolution Digital Detector Array Indexing Scan Table with Positioning Accessories Motorized Table and C-arm Internal Reference System Provides Automatic Continuous Calibration Computer Console Standard Computer Configuration (Minimum Specifications): Processor Intel Core 2 Duo (3.0GHz 6MB cache) 2GB memory 160GB hard drive DVD +/- RW drive 10 USB ports (4 front / 6 rear) Intel PRO 100 CT Gigabit Ethernet Microtower case High Resolution 17 inch Flat Panel Color Monitor Hewlett Packard Color Inkjet Printer Standard Software Configuration: Operating System Hologic APEX with Microsoft Windows 7 32 bit Professional Automatic PASS/FAIL Quality Control Express Exam Work-flow Management OneTime Auto Analysis ProTech with DXApro AccuView Automatic Hip Positioning QDR OnePage™ Report ISCD Compliant Reporting Colored Fracture Risk Indication FRAX WHO 10-Year FX Risk Patient/FRAX History Questionnaire Pediatric Reference Data DXA Practice Manager NHANES Body Composition Calibration NHANES Whole Body Bone and BCA Reference Data for Adults and Children Advanced Clinical Body Composition with InnerCore Visceral Fat Assessment and Rate-of-Change BCA Reporting and Color Whole Body Images Standard Software Configuration: Scan/Analysis Protocols: AP Lumbar Spine Proximal Femur Dual Hip Forearm Decubitus Lateral Hip Structure Analysis General Region of Interest | | | | |

| | | | | | |
|---|--|--|--|---|--|
| | Scoliotic Spine Automatic Low Density Spine Automatic Low Density Hip Pediatric Spine, Hip, and Forearm Infant spine Pediatric Whole Body Whole Body BMD Body Composition Analysis with Sub Region analysis Epress BMD Instant Vertebral Assessment (IVA-HD) 15-Second High Resolution Imaging Capability Quantitative Morphometry Physician's Viewer with IVA Toolbox MXApro – Computer Aided Fracture Assessment Additional Features Reposition/Rescan Feature Image Compare Mode for Serial Exams Single Energy Scan Switch Capability Window/Level Control for Image Optimization Patient Database Manager On-line Normative Reference Values and User-Definable Values Serial Trend Report with Automatic Rate of Change/Significance Calculation NHANES and Ethnic Reference Data Context Sensitive Help DICOM Connectivity Capability and DICOM Structured Reporting Integrated Physician's Viewer Reference Database Editor QDR Anthropomorphic Spine Phantom QDR for Windows User's Guide Two Days (12 Hours) of ECE Approved Application Training QDR Reference Manual Installation and a Twelve (12) Month Comprehensive Warranty | | | | |
| 1 | - 120 VOLT POWER SUPPLY | | | | |
| 1 | - IRIS Enterprise for systems on APEX 4.xx or higher IRIS Enterprise-APEX for systems on APEX 4.xx Includes DICOM Storage DICOM Modality Worklist Remote Physician's Viewer Physician's Report Writer HL7 Enterprise Data Management Installation and Twelve (12) Month Comprehensive Warranty on upgrade components only. | | | | |
| | | | | | |
| | | | | T | |
| | | | | | |
| | | | | T | |

| | | | | | |
|--|--|--|--|--|--|
| | | | | | |
|--|--|--|--|--|--|