

TRADE-IN : GE-OEC  
MODEL: UROVIEW 2800  
SERIAL: P6-0351-L  
ACQ DATE 11/19/2003  
EE#: 2572566

VAMC PINEVILLE, LA  
PO# 502-B35004

Qty	Description
-----	-------------

- 1      **OEC 9900 Elite Digital Mobile Standard C-arm VasMTS (Vascular MTS Platform with 30 fps digital disk and 12" I.I.)**

**High performance mobile C-arm for use in vascular and endovascular procedures. Provides flexible mechanical features, rotating anode x-ray tube, user-friendly touch screen user interface, and superb image quality.**

**The physical dimensions and the open design of the C-arm allow the system to be easily transported throughout the hospital and provide unobstructed imaging access around the patient and procedural table.**

**Digital Image Processing and Workstation**

- Dynamic Range Management (DRM) controls for image management
- Six pre-set imaging profiles for optimized anatomical capture
- Dual articulating, high resolution flat screen black & white monitors
- Integrated keyboard touchpad - cursor/tap controls
- DICOM 3.0 interface with send and query/retrieve
- CD/DVD on-board media storage, read/write compatible
- USB 2.0 mass storage device, write only, .jpg/.bmp/.avi file formats
- Automated sending of system log files to remote GE InSite 2.0 server

**Generator**

- 60kHz high frequency
- 15kW power
- Up to 120kVp
- Up to 75mA for radiographic film exposure
- Continuous high level fluoro (HLF) up to 20mA
- Digital spot up to 75mA
- Full power from standard wall outlet
- Patented battery buffered design

<b>Qty</b>	<b>Description</b>	<b>e</b>
	<b>X-ray Tube</b>	
	<ul style="list-style-type: none"> <li>- Rotating anode X-ray tube</li> <li>- 0.3mm and 0.6mm focal spots</li> <li>- Anode heat capacity: 300,000 HU (per IEC 60613)</li> <li>- Anode cooling rate: 85,000 HU/min.</li> <li>- Housing heat capacity: 1,600,000 HU</li> <li>- Standard housing cooling 22,500 HU/min.</li> </ul>	
	<b>Digital Image Rotation</b>	
	<ul style="list-style-type: none"> <li>- Digitally adjusts image display</li> <li>- Automatic image update</li> <li>- Image rotation</li> <li>- Image reversal (side-to-side)</li> <li>- Image invert (top-to-bottom)</li> <li>- Image positioning without additional exposure</li> </ul>	
	<b>PreView Collimator</b>	
	<ul style="list-style-type: none"> <li>- On-screen collimator position indication</li> <li>- PreView iris collimator</li> <li>- PreView Tungsten rotatable double leaf collimator</li> <li>- Adjusts collimators without X-ray exposure</li> </ul>	
	<b>Fluoro Mode</b>	
	<ul style="list-style-type: none"> <li>- kVp range: 40 - 120</li> <li>- mA range: 0.2 - 10</li> <li>- 1.0 - 20 HLF (high level fluoro)</li> <li>- Auto and manual fluoro modes</li> <li>- AutoTrak ABS varies mA, kVp, camera gain</li> </ul>	
	<b>Pulsed Fluoro Mode</b>	
	<ul style="list-style-type: none"> <li>- kVp range: 40 - 120</li> <li>- mA range: 0.2 - 10</li> <li>- Pulse rate: 1,2,4,8</li> <li>- Pulse width: 25 or 50ms</li> <li>- AutoTrak ABS, mA, kVp, camera gain</li> <li>- Reduces X-ray dose to patient and operator</li> </ul>	
	<b>High Level Pulsed Fluoro</b>	
	<ul style="list-style-type: none"> <li>- kVp range: 40 - 120</li> <li>- mA range: 1- 40</li> <li>- Pulse rate: 1, 2, 4, 8</li> <li>- Pulse width: 25 or 50ms</li> <li>- AutoTrak ABS, mA, kVp, and camera gain</li> </ul>	
	<b>Digital Cine Pulse Mode</b>	
	<ul style="list-style-type: none"> <li>- kVp range: 40 - 120</li> <li>- mA range: up to 150</li> <li>- Pulse rate: 15 or 30pps 60 Hz, 12 or 25pps 50Hz</li> <li>- Pulse width: 10ms</li> <li>- AutoTrak ABS, mA, kVp, camera gain</li> </ul>	
	<b>Digital Spot Mode</b>	
	<ul style="list-style-type: none"> <li>- kVp range: 40 - 120</li> <li>- mA range: Up to 75</li> <li>- Automatic exposure termination and automatic image save</li> </ul>	

Qty	Description
	<p><b>Radiographic Mode</b></p> <ul style="list-style-type: none"> <li>- mA range: Up to 75</li> <li>- mAs range: up to 300</li> <li>- Computer controlled exposure time</li> </ul> <p><b>12" Image Intensifier</b></p> <ul style="list-style-type: none"> <li>- Tri-mode 12"/9"/6" (31cm/23cm/15cm) image intensifier</li> <li>- Minimum central resolution (at the monitor): <ul style="list-style-type: none"> <li>- 12" (31cm): 1.6 lp/mm</li> <li>- 9" (23cm): 2.2 lp/mm</li> <li>- 6" (15cm): 2.6 lp/mm</li> </ul> </li> <li>- DQE: 65% (typical)</li> </ul> <p><b>Precision imaging with Dynamic Range Management (DRM) enhances features of interest while attenuating background noise.</b></p> <ul style="list-style-type: none"> <li>- Preset imaging Profiles <ul style="list-style-type: none"> <li>- 9800</li> <li>- General</li> <li>- Orthopedic</li> <li>- Spine</li> <li>- Vascular</li> <li>- Bolus Chase</li> </ul> </li> </ul> <p><b>AutoTrak(TM)</b></p> <ul style="list-style-type: none"> <li>- Automatic Brightness Stabilization (ABS)</li> <li>- Automatically seeks the subject anatomy anywhere within the imaging field and selects the optimum imaging technique</li> <li>- Automatically adjusts to anatomical size and location</li> <li>- Provides uniform image quality throughout entire image</li> <li>- Simplifies operation</li> </ul> <p><b>Image Quality</b></p> <ul style="list-style-type: none"> <li>- Smart Window</li> <li>- Dynamically senses the collimator position and automatically adjusts brightness and contrast to produce high image quality.</li> <li>- Smart Metal <ul style="list-style-type: none"> <li>- Allow user to adjust automatic brightness and contrast sensitivity levels to metal</li> <li>- Provides optimum image quality even when metal is introduced to the field</li> </ul> </li> <li>- Tungsten Collimator <ul style="list-style-type: none"> <li>- Denser collimator limits X-ray exposure area</li> <li>- Reduces scatter radiation</li> <li>- Improves image detail</li> </ul> </li> </ul> <p><b>Video Camera</b></p> <ul style="list-style-type: none"> <li>- High resolution 1k x 1k CCD camera</li> <li>- Full frame capture</li> <li>- Motorized rotation</li> <li>- On-screen orientation indicator (real-time feedback without fluoro)</li> <li>- Left-right image reversal</li> <li>- Top-bottom image invert</li> </ul> <p><b>Video Display</b></p> <ul style="list-style-type: none"> <li>- Dual display anti-glare, LCD flat panel monitors mounted on an articulating arm</li> </ul>

Qty	Description
-	22" horizontal travel 7° up/10° down Monitors viewable from all four sides of workstation Horizontal and vertical viewing angle 170° 1200 CD/M2 maximum brightness Touch screen system control 1280 x 1024 high resolution monitors
<b>Vascular MTS Platform</b>	
-	1k x 1k x 16 bit image processing Preset Imaging Profiles 9800 General Orthopedic Spine Vascular Bolus Chase Noise filter with on-screen indicator Minimal difference spatiotemporal noise filter (MDST) Real-time dynamic range management (DRM) Automatic digital brightness and contrast control Manual digital brightness and contrast control Negate mode Save and auto-save feature Swap and auto-swap feature Patient information Examination list Customized patient information Customize functions Workstation set-up Mainframe set-up Patient information set-up Date/time set-up DICOM interface set-up Last image hold 1000 image storage CD/DVD burner with DICOM viewer for displaying images on PC platforms -512 x 512 or 1k x 1k Integrated DICOM interface (store, print, worklist, and query/retrieve) HIPAA SecureView Password protection Blank screen function Delete all Zoom and roam function Image annotation Measurement software Real-time digital subtraction (DSA) Peak opacification Roadmapping Re-registration Variable landmarking Mask save/recall Motion tolerant subtraction (MTS) Digital cine pulse mode 30 pulses/sec 60Hz (25 pulses /sec 50Hz) Up to 150mA 10ms pulse width 30fps Dynamic digital disk 60 Hz (25 pulses/sec, 50Hz) Recording/playback rate: 1,2,4,8,15, 30fps 60 Hz (1, 2, 4, 8, 12, 25, 50 Hz)

<b>Qty</b>	<b>Description</b>
	- Recording time: 10 minutes @ 30fps 60 Hz (time depends on record frame rate) (25 fps, 50 Hz)

**User Interface**

- Entire system is computer controlled and software upgradeable
- Touchscreen control simplifies operation
- Automated system operation requires minimum operator interface
- Multi-functional controls
  - Footswitch
  - Hand-held control
- Simplified keyboard with integrated touchpad
- Multi-purpose image directory
  - Retrieve and review images
  - Compose hardcopy films
  - Copy images
- X-ray dose summary

**C-arm Mechanics**

- Counterbalanced, manual adjustment of orbital rotation, cranial-caudal rotation, wig-wag and horizontal motion
- Patented flip-flop C-arm reversal (SmartView)
- 31" free space, 28" depth in arc, 115° orbital rotation
- Dual, illuminated C-arm operator control panels

**OEC Clinical Excellence Onsite Training**

- Pre-training package with interactive CD-ROM
- Up to three days\* of in-service by our ARRT certified Clinical Imaging Specialists (CIS) during the one-year warranty period.
  - Post-training skills assessment & test
- Participants may be eligible for Continuing Education (CE) credits from the American Society of Radiologic Technologists\*\*

**Warranty**

- One Year Warranty

**Notes: WSP – Wireless Service Platform**

- The OEC 9900 comes with a device called the WSP which facilitates three capabilities:
  - 1) USB 2.0 mass storage , The USB functionality is write only, and saves images in .jpg/.bmp/.avi file formats. USB storage devices must be unencrypted and unprotected
  - 2) Automated sending of system log files to remote GE InSite 2.0 server: OEC joins GE Healthcare InSite 2.0 Remote Data Collection and Monitoring platform:
    - Also utilizes 128 bit Secure Sockets Layer (SSL) data encryption, HTTPS Port 443, and complies with your existing firewall rules and web proxies. Log files sent via the WSP do not contain patient information
  - 3) If purchased, Wireless DICOM. Wireless utilizes a Wireless LAN 802.11 b/g/n, 10/100MB device, runs a custom embedded Windows implementation, and currently supports WEP-64, WEP-128, WPA-PSK(TKIP), WPA2-PSK (AES-CCMP) security methods

**Notes: OEC Clinical Excellence Onsite Training**

- \* Onsite training provided from 8am to 5pm, Monday through Friday. Includes all CIS travel & living expenses.\*\* Training produces the best results when a dedicated core group of

Qty	Description
	<p>technologists complete the session. Those who complete the entire OEC Clinical Excellence curriculum should be competent to perform the tasks required for basic operation of the system. Competency will be measured through a skills assessment completed while the CIS is on-site.</p>
1	<b>30ft. Interconnect Cable</b>
1	<p><b>IDI Aspect 100UC Plus Fluoroscopic Procedure Table</b></p> <p><b>Multi-Purpose Mobile Urological Imaging Table with motorized Elevation, Trendelenburg Tilt Longitudinal and Transverse Tabletop Motions</b></p> <ul style="list-style-type: none"> <li>- Carbon fiber tabletop with low x-ray attenuation</li> <li>- 2" (5cm thick tabletop &amp; tabletop extension pads with memory foam construction and waterproof cover</li> <li>- Tabletop dimensions: 28"W x 46"L ( 71 x 117 cm)</li> <li>- Radiolucent tabletop extension: 24"W x 34"L (61 x 86 cm)</li> <li>- Imaging area, main tabletop: 25"W x 20"L (63.5 x 50.8 cm)</li> <li>- Imaging area, tabletop extension: 24"W x 34"L (61 x 86 cm)</li> <li>- Motorized Tabletop motions with hand control and foot control: <ul style="list-style-type: none"> <li>- Elevation 31.3" – 41.1" (79.4 – 104.4 cm)</li> <li>- Trendelenburg Tilt <math>\pm 15^\circ</math></li> <li>- Longitudinal tabletop travel = 10" (25.4cm)</li> <li>- Transverse tabletop travel = <math>\pm 3.5</math> (<math>\pm 8.9</math> cm, 17.8cm total)</li> </ul> </li> <li>- Built-in footswitch storage tray</li> <li>- Removable radiographic cassette tray, loads from either side of table</li> <li>- Standard US side-rails for accessory attachments on main tabletop</li> <li>- Patient weight capacity: <ul style="list-style-type: none"> <li>- With patient upper body on main tabletop: 525 lbs.(238 kg)</li> <li>- With patient upper body on extension: 400 lbs. (181 kg)</li> <li>- Load capacity on extension only: 250 lbs. (113.4 kg)</li> </ul> </li> <li>- Table weight: 575 lbs. (261 kg)</li> <li>- Locking swivel casters, 5" (12.7 cm), for easy transport</li> <li>- Electrical power required: 110 VAC, 60 Hz, or 230 VAC 50/60 Hz operation (must specify with order). 10 amp max. required at 110V.</li> <li>- Battery back-up operation</li> </ul> <p><b>Accessories included:</b></p> <ul style="list-style-type: none"> <li>- Tabletop Pad, 2" thick, waterproof</li> <li>- Tabletop Extension, with pad</li> <li>- Hand control pendant</li> <li>- Footswitch control</li> <li>- Drain Bag Support Hoop</li> <li>- Disposable drain bag and hose assembly, box of 20</li> <li>- Removable Radiation Shield Flap for perineal end of table</li> <li>- Patient Restraint Straps (2)</li> <li>- Armboard (1), rail mounted, with pad and strap</li> <li>- Disposable covers for footswitch, box of 50</li> </ul> <p><b>Warranty:</b></p> <ul style="list-style-type: none"> <li>- Image Diagnostics Inc. (IDI) products are warranted to be free from failures due to defects in materials or Workmanship for a</li> </ul>

Qty	Description
	<p>period of 3 years for parts and 1 year for labor.</p> <ul style="list-style-type: none"> <li>- During this period any warranty claims should be addressed directly with IDI.</li> </ul>
1	<p><b>NuBoom M2 1080P, High Definition System, All in one equipment management and visualization system.</b></p>
	<p><b>Includes:</b></p> <ul style="list-style-type: none"> <li>- Two (2) Overhead 24" HD Monitors</li> <li>- One 1-segment and One 2-segment booms</li> <li>- 7" Touch Screen Master Control</li> <li>- DOCS Software</li> <li>- Split-Screen compatibility on all monitors</li> <li>- Installation: Weekdays; 2 Days In-Service</li> <li>- Foot Pedal Basket Storage Kit - A 12"x8"x6" Utility Basket.</li> <li>- Suction Canister Kit - A 19" rail that can hold up to 4 Canisters.</li> <li>- CO2 Cylinder Basket (Qty 2) - Two CO2 cylinder baskets.</li> <li>- BOOM Mounted IV Hooks - A set of 4 hooks.</li> </ul>
	<p><b>Inputs:</b></p> <ul style="list-style-type: none"> <li>- 1 'S' Video</li> <li>- 1 DVI-I</li> <li>- 1 HDSDI 1080i</li> <li>- 1DVI-D 1080P</li> <li>- 1 VGA</li> <li>- 1 Composite</li> <li>- DVD-R (in/out. S-video)</li> </ul>
1	<p><b>12" I.I. laser aimer/localizer (with removable cross hairs) &amp; supplement</b></p>
1	<p><b>CI, REMOTE TOUCH PANEL CONFIGURATION, ACC</b></p>
1	<p><b>CI, I-POD SOUND UNIT</b></p>
1	<p><b>CI ITEM, SERV, MAINT, 12 INCH, TWELVE MONTH EXTENDED WARRANTY, 9900</b></p>
1	<p><b>CI ITEM, SERV, TRNG, 2 WEEK, 9900 TECHNICAL TRAINING</b></p>