

Item No. Qty

Description

1 1

Discovery XR656 Base System with TWO FlashPad Detectors

Discovery XR656 Base System with TWO FlashPad Detectors

The Discovery XR656 is an advanced digital radiographic imaging system that provides excellent image quality, a variety of image manipulation and post-processing tools as well as the option to utilize GE Healthcare's exclusive advanced clinical applications. Powered by FlashPad, GE's wireless digital detector, workflow is streamlined with TWO portable detectors that can be used in the table, wall-stand or freely to best accomodate most 2D exam requirements as well as advanced radiographic requirements.

At its core, the Discovery XR656 delivers improved productivity through an efficient and intuitive workflow. The base system comes equipped with a systems cabinet, acquisition review workstation, image processing tools, short-term storage and quick in-room viewing of images. Also included is a host of networking and connectivity options, dose reporting and system maintenance. These features are designed to make this system easy to use and reliable while providing high quality radiographic images in a digital environment.

The Discovery XR656's core feature set can help streamline even the toughest exam. Consistent image quality helps reduce retakes, and unique image processing helps save time. With DICOM connectivity, the Discovery XR656 connects easily to HIS, RIS or PACS for efficient data transfer.

The Discovery XR656 base system includes:

Two wireless Digital Detectors with batteries and optional tether.

- Single panel (non-tiled) amorphous silicon detector with a Cesium Iodide scintillator and two handgrips that facilitate easy positioning and a secure grip.
- Wireless connectivity through Ultrawide Band Technology for fast, efficient and secure communications between the detector and the base system.
- Image area 40.4 cm x 40.4 cm (15.9 in. x 15.9 in.)
- Active matrix 2022 x 2022 pixels.
- Pixel pitch 200 microns.

Item No. Qty.	Description
	<ul style="list-style-type: none"> • Typical upper dynamic range 7.8 mR @ RQA5. • Typical DQE 68 @ Olp @ RQA5. • Weighs 4.32 kg (9.52 lb.) without battery, 0.18 kg (0.40 lb.) battery weight. • Dimensions: L 580 mm, H 452 mm, T 24 mm. • Can support up to 160 kg (352 lb.) of distributed load. • Battery or tethered operated. Includes two rechargeable and exchangeable batteries and 7 m cable (4 m or 10 m optional) for optimal connectivity and power. • One 6:1 clip-on grid for FlashPad detectors. • Battery allows for 150 images to be taken in 3 hours. • QAP (Quality Assurance Procedure) • FlashPad detector can support Advanced Applications that are options available on the Discovery XR656

Acquisition Workstation

The Acquisition Workstation is the primary interface to the network and provides image post-processing capabilities. The System Controller Module provides single point control, directing and coordinating overall system operation, while monitoring all system modules automatically through software.

The Acquisition Workstation includes:

- Two 19-inch Flat Panel Monitors that help minimize desktop space required.
- Keyboard and Mouse
- CPU Tower with 6GB RAM, 320 GB of hard drive storage and capacity for over 22,000 images.
- 120/140 VAC, 50/60 Hz.
- Easy Image Manipulation and
- Image Display Tools which include:
 - Window width and level
 - Gray scale/contrast invert
 - Interpolated zoom and roam
 - Image flips (horizontal, vertical) with automatic indicator
 - Image Rotate - 90 degree increments
 - Free rotation - 360 degrees

Item No. Qty	Description
	<ul style="list-style-type: none"> – Image orientation management – Electronic Left/Right Markers – Free text annotation – Manual shuttering – Image Annotations and Measurement Tools – Multi-Resolution Post Processing – Customizable Image Processing to Match – User Preferences – CD-RW and DVD Drive for Image Archive – Image Viewer on Archive CD's and DVD's
	<p>Acquisition Control System</p> <ul style="list-style-type: none"> • Single Point System Control and Monitoring. • Auto-Protocol Programming: comes with default set of exam and view protocols and the ability to build an infinite number of exams or views through the editor function. User also has the ability to attach acquisition protocols. • System manager allows equipment error logging and provides resident power-up diagnostics. • kVp selections range from 40-150 in 1 kVp increments • Focal Spot Sizes of 0.6 and 1.2 mm - mA
	<p>Patient Management</p> <ul style="list-style-type: none"> • DICOM Modality Worklist for HIS/RIS - SCU (with programmable auto refresh) • Patient edit/auto-folding (Copy exam) • Patient Dose Reporting Calculated Dose Monitoring (mGy) • "Patient Directory" provides fast access to the image and exam database for case reviews and file management. • Detector Exposure Indicator: tool for detector dose feedback to ensure exposure was within normal limits. • Emergency Patient Exam Capability Emergency patient feature - allows user to open exam and acquire images without a worklist entry. • The Copy Patient function allows merging of the patient information with the exam images. • Fast Image Display - Average time for a partially processed image is approximately three seconds and less than eight

Item No. Qty	Description
	<p>seconds for a fully processed image. Times may vary based on how the detector is connected to the system (i.e. docked, tethered or wireless)</p> <ul style="list-style-type: none"> • Set of default adult and pediatric protocols allows quick selection of the appropriate techniques for common procedures/exams with the ability to define unlimited number of custom protocols. • Set of 4 Factory (GE pre-set) image processing selections (looks) optimized for each anatomical view with the ability to define multiple Custom looks for each anatomical view/patient size combination. • Automatic image storage and print with DICOM 3.0 and IHE Compliant networking, further increasing exam throughput and decreasing examination time for patients. <p>Image Quality and Dose</p> <ul style="list-style-type: none"> • Multi-resolution image processing capability. • Tissue Equalization used to correct over-penetrated and under-penetrated areas within the image. • Auto and manual image shuttering cropping tool. • Automated brightness/contrast setting (Smart Windowing) • Orthopedic Magnification/Print. • Detector Exposure Index (DEI) - dose tracking and QC metric. • Dose Area Product (DAP) - entrance dose metric. • Grid Line Reduction (selectable) • Intelligent Collimator Edge Detection automated, image based cropping/shuttering tool that relies solely on image information to locate the collimator edges present in the image. • Orthopedic Magnification/Print <p>DICOM and Standard Networking Capabilities</p> <p>Images may be transmitted manually or automatically through the DICOM interface to printers, archival devices, servers or review workstations.</p> <p>Please refer to the DICOM Conformance Statement for complete definition of</p> <p>supported DICOM connectivity services.</p>

Item No.	Qty	Description
		<p>DICOM and Standard Networking Capabilities include:</p> <ul style="list-style-type: none"> • Ethernet Network Link - DICOM 3.0. • DICOM Storage (with auto-send to different locations) • DICOM Storage Commitment (with programmable auto delete function) • DICOM Modality Worklist for HIS/RIS (with auto refresh) • DX/CR Worklist Filtering. • DICOM Media Interchange on DVD-R. • DICOM Modality Perform Procedure Step (MPPS) feedback to the HIS/RIS (SPS PPS) • DICOM Grayscale Print (with print layout at the console) • DICOM Query/Retrieve (retrieves images back from PACS) • System Access & Authorization Control to support HIPAA Compliance. • Full Range of Printing Options. • Numerous Layout and Format Options. <p>Systems Cabinet</p> <ul style="list-style-type: none"> • Built-in System Distribution Power Module and Circuit Breaker for Single Point Power Feed to Room Subsystems • Modular Designed X-ray Systems Based on a Digital Communications Network for Improved Reliability and Image Quality. <p>Quality Control/System Reliability Features</p> <p>Preventing customer experienced system failures and reducing unplanned system downtime are critical.</p> <p>The following features help to achieve these goals:</p> <ul style="list-style-type: none"> • Using the integrated system Quality Assurance Procedure (QAP), image quality checks can be easily performed by the customer. • The QAP includes a phantom, optimized for Digital Image Quality testing and is included with the system. • System changes are highlighted and can be corrected before they become a problem.
2	1	<p>Overhead Ceiling Tube Suspension</p> <p>Discovery XR656 Overhead Tube Suspension</p>

Item No. Qty	Description
	<p>The Overhead Tube Suspension (OTS) system with motorized movement delivers excellent levels of operational support for efficient operation and precise positioning.</p> <p>The Overhead Ceiling Tube Suspension Package includes:</p> <ul style="list-style-type: none"> • Auto-Positioning Package (included in base) <ul style="list-style-type: none"> – Auto-Positioning enables the users to select a predefined system position from the system console and automatically move the equipment by simply holding the "Auto Positioning" buttons. This feature reduces user fatigue and increases the productivity of the operator. – Auto-Positioning is controlled at the acquisition workstation or with the IR remote control, allowing the user to remain in the room while moving the system. – Pre-set positions at the table, wallstand and park position at various SIDs and vertical and horizontal orientations. – Auto-Positioning will incorporate angulation of the tube, longitudinal, lateral, rotational and vertical positioning of OTS, table detector longitudinal positioning, wallstand detector vertical positioning. – Auto-Positioning comes with 7 default positions and up to 10 additional user defined positions can be added to the system. • Tube and Detector Tracking Synchronized tube and detector tracking enables convenient workflow by keeping the tube correctly aligned with the wall-stand or table detector. • Motorized 5-axis movements Supported Positions: park, table 100 cm SID (head, center, foot), wallstand 100 cm and 180 cm CID, and horizontal wallstand 100 cm SID. • IR Remote - The infrared remote is an in-room control allowing the technologist greater flexibility & ability to pre-position the system automatically in preparation for the next clinical exam/view. • Auto or manual positioning with single lock release and auto detents for assisted manual positioning. • Touch-screen user interface with LCD screen display helps confirm patient data, review techniques, receptor selection

Item No.	Qty	Description
		and modify/confirm wireless detector association. The easy to Read, Auto Rotating user screen also includes a message readout line and easy to see light indicators.
		<ul style="list-style-type: none"> • Patient Side Touch Screen User provides the following functions to the user: <ul style="list-style-type: none"> – Lock, Detent Control – Field of View Image Size Selection – Collimator Field Light Selection – Technique Adjust (kVp, mAs) – Receptor Selection (table, wallstand, wireless or cassette) – Exam Inhibit Display – Collimator Manual Override – Position Display (Source-to-Image Distance, X-ray Tube Angle, Column Rotation Angle) – Display of Patient Name for In-Room Verification (this feature can be disabled)
3	1	2, 3, or 4 Meter Longitudinal Drive Belt Kit
		2, 3 or 4 Meter Longitudinal Drive Belt Kit
4	1	Inboard Rail Select
		2, 3 or 4 Meter Longitudinal Rail Select (Dependent on Room Size)
5	1	2, 3, or 4 Meter Bridge Select
		2, 3 or 4 Meter Bridge Select (Dependent on Room Size)
6	1	2, 3, or 4 Meter Bridge Cable Select
		2, 3 or 4 Meter Bridge Cable Select
7	1	2, 3, or 4 Meter Bridge Cable Drape Select
		2, 3 or 4 Meter Bridge Cable Drape Select
8	1	Required Language Labels
		Required System Language Labels - English
9	1	80kW High Frequency Generator
		80kW High Frequency Generator, 50 or 60Hz The high frequency 80kW power unit is designed for radiographic applications and

Item No.	Qty	Description
10	1	<p data-bbox="526 174 1292 201">utilizes microprocessor controlled power and parameter adjustment</p> <p data-bbox="526 222 686 249">Specifications:</p> <ul data-bbox="545 262 803 338" style="list-style-type: none"> <li data-bbox="545 262 803 289">• 1000mA at 80kVp <li data-bbox="545 306 803 333">• 800mA at 100kVp <p data-bbox="586 357 1045 384">Input Power: 360-480VAC, 3 Phase, 60Hz</p> <p data-bbox="586 407 1045 434">kVp Range: 40-150kVp,)kVp increments</p> <p data-bbox="526 472 954 499">kVp Accuracy: 3% +/-2kVp mA Range:</p> <ul data-bbox="545 512 954 588" style="list-style-type: none"> <li data-bbox="545 512 954 539">• Small Focal Spot: 10-320mA <li data-bbox="545 556 954 583">• Large Focal Spot: 160-1000mA <p data-bbox="526 609 1099 636">Discovery XR656 Table WITHOUT FlashPad Detector</p> <p data-bbox="526 657 1039 684">Discovery XR656 Table (No Detector Included)</p> <ul data-bbox="545 697 1187 976" style="list-style-type: none"> <li data-bbox="545 697 1076 724">• Table Base - Motor Driven Elevating Table <li data-bbox="545 737 1187 764">• Table Footprint - 1260 mm x 690 mm (50 in. x 27 in.) <li data-bbox="545 779 1089 806">• Elevating Range 57 cm-82 cm (20 in.-32 in.) <li data-bbox="545 821 1130 848">• Elevation from Min to Max Height <10 seconds <li data-bbox="545 863 1117 890">• Maximum Patient Weight 220kg (485 pounds) <li data-bbox="545 905 1040 932">• Foot Pedals on Front and Rear of Table <li data-bbox="545 947 935 974">• Two Safety Lockout Switches <p data-bbox="516 997 634 1024">Table Top</p> <ul data-bbox="545 1037 1235 1234" style="list-style-type: none"> <li data-bbox="545 1037 902 1064">• Eight Way Tabletop Motion <li data-bbox="545 1079 1182 1106">• Inherent Filtration <1.0 mm Al Equivalent at 100kVp <li data-bbox="545 1121 1235 1186">• Table Top Dimensions 88 cm Width x 229 cm Length (34.6 in. x 90.16 in.) <li data-bbox="545 1201 1089 1228">• Non-Protruding Edge Rails for Attachments <p data-bbox="516 1251 875 1278">Table Detector Servo Mechanics</p> <ul data-bbox="545 1291 1295 1604" style="list-style-type: none"> <li data-bbox="545 1291 1068 1318">• Longitudinal Tracking of Detector to Tube <li data-bbox="545 1333 954 1360">• Detector Travel 685 mm (27 in.) <li data-bbox="545 1375 1062 1402">• Lateral Patient Coverage 610 mm (24 in.) <li data-bbox="545 1417 1135 1444">• Longitudinal Patient Coverage 1830 mm (72 in.) <li data-bbox="545 1459 1295 1524">• Removable 100 cm Focus High Line Rate Stationary Grid 70 I/cm (198 I/in.) 12:1 Ratio <li data-bbox="545 1539 984 1566">• Automatic Exposure Control (AEC) <li data-bbox="545 1581 922 1608">• Cross Table Detector Holder

Item No.	Qty	Description
11	1	<p>Table Top FlashPad Lateral Detector Holder</p> <p>Table Top Lateral Detector Holder</p> <p>Wireless DR detector holder, designed specifically for GE, secures the detector in a vertical position on the tabletop for cross-table imaging.</p>
12	1	<p>Discovery XR656 Tilting Wallstand with Extended Length Arm</p> <p>Discovery XR656 Tilting Wallstand with Extended Length Arm</p> <p>Discovery XR656 Digital Tilting Wallstand is designed for radiography applications with the patient standing, sitting or lying on a gurney.</p> <p>The Discovery XR656 wallstand is designed for use with GE's exclusive FlashPad wireless digital detector, overhead tube suspension, ion chamber and removable non-reciprocating grid. The FlashPad wireless digital detector can be operated docked, tethered, or in a wireless digital cassette mode.</p> <ul style="list-style-type: none"> • The wallstand is motorized Electromagnetic braking secures vertical motion. • Motorized receptor tilting controlled with either IR remote control or hand switches located on the arm. 0 Vertical motorization of the wallstand with foot switch or IR remote control. • Auto-tracking from the overhead tube suspension. • Graphic outlines of image sizes and ion chamber scan areas on the front panel enhance accuracy and safety. Preparation is fast and simple for better patient throughput. • The wallstand tilts from -20 degrees to 90 degrees • The wallstand grids are removable from the side. • It is configurable for either left or right side insertion. • The wallstand comes with two removable hi-line rate fixed grids for optimum scatter cleanup and aluminum interspacing for image quality uniformity: <ul style="list-style-type: none"> - 100 cm (40 in.) focus grid with a SID range of 90 cm - 118 cm (70 lines/cm, 12:1 ratio) • 180 cm (72 in.) focus grid with a SID range of 145 cm - 245 cm (70 lines/cm, 13:1 ratio)

Item No.	Qty	Description
13	1	<ul style="list-style-type: none"> Automatic Exposure Control (AEC) utilizes three ion chamber sensors, which are mounted between the patient and digital detector. Automatic Exposure Control (AEC) on this extended arm utilizes three ion chamber sensors which are mounted between the patient and digital detector on the wallstand. The addition of the 4th ion chamber provides for supine chest and abdominal imaging when a stretcher is used (anytime the wallstand housing is horizontal). The standard right, left and center ion chambers are available when the wallstand is in the vertical position. Includes patient handgrips and a lateral support bar. <p>Mid-Range 130cm Focus Grid</p> <p>Mid-Range Wallstand Grid</p> <p>- 130 cm (51 inch) focus - 10:1 ratio, 70 lines/cm - Used for SID ranges from 90 to 190 cm</p>
14	1	<p>Dual Energy - Chest and Abdomen</p> <p>Dual Energy Subtraction Chest and Abdomen</p> <p>Dual Energy Subtraction is an Advanced Clinical Application.</p> <p>Provides three diagnostic images:</p> <ul style="list-style-type: none"> Standard Soft tissue only Bone only Technique is available for PA/AP chest and abdomen imaging.
15	1	<p>VolumeRAD Table and Wallstand - 4 Sweeps</p> <p>VolumeRAD at Wallstand and Table with Extended Arm (4 sweeps)</p> <p>VolumeRAD is an advanced application for radiographic X-ray imaging. It allows retrospective reconstruction of a number of planes (slices) from a series of low dose exposures acquired within a limited angular range. These slices show anatomical structures at different depths</p> <p>This application supports the following four sweeps at wallstand:</p> <ul style="list-style-type: none"> Table Horizontal sweep of the tube over the table for prone or

Item No.	Qty	Description
16	1	<p data-bbox="574 178 756 205">supine imaging</p> <ul data-bbox="537 216 1289 443" style="list-style-type: none"> <li data-bbox="537 216 1230 243">• Vertical tube sweep at the wallstand for upright imaging <li data-bbox="537 254 1289 354">• Horizontal sweep of the tube over the stretcher for prone or supine imaging (requires the selection of the extended arm wallstand) <li data-bbox="537 380 1230 443">• Lateral tube sweep at the wallstand for cross table imaging over the stretcher for prone or supine imaging <p data-bbox="509 468 1089 495">Repeat and Reject Analysis for the Discovery XR656</p> <p data-bbox="509 510 1089 537">Repeat and Reject Analysis for the Discovery XR656</p> <p data-bbox="509 562 1256 735">The Repeat and Reject Analysis (RRA) software package is a quality control (QC) application available for the Discovery XR656 that allows for repeat or reject images to be captured and categorized by technologist. It is designed to help track and analyze the X-ray repeat rate</p> <ul data-bbox="529 751 1284 1129" style="list-style-type: none"> <li data-bbox="529 751 1284 852">• The easy-to-use operator interface helps technologists classify each image they reject and to select a specific reason for the repeat /image rejection. <li data-bbox="529 869 1284 970">• RRA can be a helpful teaching tool because it includes links to actual JPEGs of the rejected images to help the user analyze why the image was rejected. <li data-bbox="529 987 1243 1050">• The RRA application tracks the rejected image data by operator, exam type, date and reason code. <li data-bbox="529 1066 1260 1129">• Reports can be exported in DVD, CD or USB format for ease of use.
17	1	<p data-bbox="509 1163 987 1190">Auto Field of View for the Discovery XR656</p> <p data-bbox="509 1205 987 1232">Auto Field of View for the Discovery XR656</p> <p data-bbox="509 1257 1268 1392">Auto Field of View enables the user to pre define the collimation size on an individual view basis. Helps improve workflow by reducing the time it takes to make the final collimation adjustment for a particular exam.</p>
18	1	<p data-bbox="509 1428 1011 1455">Auto Protocol Assist for the Discovery XR656</p> <p data-bbox="509 1480 1011 1507">Auto Protocol Assist for the Discovery XR656</p> <p data-bbox="509 1524 1224 1587">Auto-Protocol Assist - the Discovery XR656 system will automatically transition directly to the Acquire screen when the</p>

Item No.	Qty	Description
19	1	<p data-bbox="526 170 1279 310">protocol code downloaded from the HIS/RIS (automatically performed with worklist refresh) matches the exam code contained in the protocol database. This tool eliminates the user steps required to select patient exam types and initiate an exam.</p> <p data-bbox="526 338 1224 365">Auto Image Paste for the Discovery XR656 Wallstand and Table</p> <p data-bbox="526 392 1219 420">Wallstand Auto Image Paste for Spine and Long Bone Imaging</p> <ul data-bbox="537 434 1295 1052" style="list-style-type: none"> <li data-bbox="537 434 1295 533">• Fully automated acquisition and processing of a series of images with user defined start and stop locations on the anatomical regions of interest. <li data-bbox="537 552 1295 688">• Average acquisition time for a 3-image exam (90 cm coverage) is less than 10 seconds. Image pasting and processing time for a 3-image exam is less than 15 seconds from the last exposure. <li data-bbox="537 707 1247 772">• Allows 2 to 5 images to be pasted together with a maximum range of 150 cm. <li data-bbox="537 791 1240 856">• Includes imaging of the spine for scoliosis evaluation and imaging of the legs for orthopedic evaluations. <li data-bbox="537 875 1287 974">• Supports anatomies/view combinations of Spine Antero-posterior, Spine Postero anterior, Spine Lateral, Leg Antero-posterior, Leg Postero-anterior. <li data-bbox="537 993 1247 1052">• Includes a patient stand with screen to keep the patient comfortable during acquisition. <p data-bbox="526 1073 1162 1100">Table Auto Image Paste for Spine and Long Bone Imaging</p> <ul data-bbox="537 1115 1295 1644" style="list-style-type: none"> <li data-bbox="537 1115 1295 1213">• Fully automated acquisition and processing of a series of images with user defined start and stop locations on the anatomical regions of interest. <li data-bbox="537 1232 1295 1369">• Average acquisition time for a 3-image exam (90 cm coverage) is less than 10 seconds. Image pasting and processing time for a 3-image exam is less than 15 seconds from last exposure. <li data-bbox="537 1388 1247 1453">• Allows 2 to 3 images to be pasted together with a maximum range of 100 cm. <li data-bbox="537 1472 1240 1537">• Includes imaging of the spine for scoliosis evaluation and imaging of the legs for orthopedic evaluations. <li data-bbox="537 1556 1287 1644">• Supports anatomies/view combinations of Spine Antero-posterior, Spine Postero anterior, Spine Lateral, Leg Antero-posterior.

Item No.	Qty	Description
20	1	<p data-bbox="516 216 1182 243">25 KAIC X-Ray Main Disconnect Panel 80 Amp, 480 V / 208 V</p> <p data-bbox="516 268 1198 296">25 KAIC X-Ray Main Disconnect Panel 80 Amp, 480 V / 208 V</p> <p data-bbox="516 321 756 348">FEATURES/BENEFITS</p> <ul data-bbox="529 373 1287 1003" style="list-style-type: none"> Serves as the main power disconnect between the X-Ray system and the facility 480V or 208V power source Provides emergency shut down, undervoltage protection and overcurrent protection for the X-Ray power distribution cabinet Standardized design provides a platform for future upgrades of the system Offers a number of advantages by combining a variety of individual components into a single pre-engineered and factory tested panel UL and cUL listed for compliance with NEC Article 100 and Article 110-3 Remote emergency off pushbutton located by X-Ray control provides immediate shut down of the entire system to comply with NEC required disconnecting means Surface or semi-flush mounting <p data-bbox="508 1020 703 1047">SPECIFICATIONS</p> <ul data-bbox="529 1066 1247 1213" style="list-style-type: none"> Dimensions (H x W x D): 48" x 20" x 6.68" Weight: 80 lbs. Mounting: via keyhole slots; Width is 16" on centers, Height is 45.5" on centers <p data-bbox="508 1236 691 1264">COMPATIBILITY</p> <ul data-bbox="529 1285 927 1312" style="list-style-type: none"> GE Three Phase X-Ray generators <p data-bbox="516 1337 609 1365">NOTES:</p> <ul data-bbox="529 1411 1230 1516" style="list-style-type: none"> Customer is responsible for rigging and arranging for installation with a certified electrician ITEM IS NON-RETURNABLE AND NON-REFUNDABLE
21	1	<p data-bbox="516 1549 829 1577">6 Day X-ray System Training</p> <p data-bbox="508 1591 797 1619">6 Day XR System Training</p>

Item No.	Qty	Description
22	1	<p data-bbox="524 201 1192 264">One 4 day and one 2 day TIP Onsite Training visits for the X-ray system.</p> <p data-bbox="524 289 1114 317">Includes TELL expenses. Days provided consecutively.</p> <p data-bbox="524 342 1230 405">This training program must be scheduled and completed within 12 months after the date of product delivery.</p> <p data-bbox="524 436 1153 464">VolumeRAD TIP Training: 4 Days Onsite (3 Days + 1 Day)</p> <p data-bbox="524 489 1149 516">VolumeRAD TiP Training: 4 Days Onsite (3 Days + 1 Day)</p> <p data-bbox="524 541 1170 604">One 3 day and one 1 day TIP onsite applications training visit covering the VolumeRAD software functionality.</p> <p data-bbox="524 630 1271 657">Onsite training provided from 8AM to 5PM, Monday through Friday.</p> <p data-bbox="524 672 786 699">Includes TELL expenses.</p> <p data-bbox="524 716 1230 779">This training program must be scheduled and completed within 12 months after the date of product delivery.</p>
23	1	<p data-bbox="524 814 1170 842">Auto Image Paste TiP Training: 2 Consecutive Days Onsite</p> <p data-bbox="524 867 1170 894">Auto Image Paste TiP Training: 2 Consecutive Days Onsite</p> <p data-bbox="524 919 1175 982">Two consecutive days of TiP Applications Onsite X-ray training covering the Auto Image Paste function.</p> <p data-bbox="524 1008 1271 1035">Onsite training provided from 8AM to 5PM, Monday through Friday.</p> <p data-bbox="524 1050 786 1077">Includes TELL expenses.</p> <p data-bbox="524 1094 1230 1157">This training program must be scheduled and completed within 12 months after the date of product delivery.</p>
24	1	<p data-bbox="524 1192 834 1220">WALLSTAND CABLE SELECT</p> <p data-bbox="524 1245 919 1272">Connection cables for the wallstand</p>
25	1	<p data-bbox="524 1308 760 1335">TABLE CABLE SELECT</p> <p data-bbox="524 1350 870 1377">Connection cables for the table</p>