

VAMC SALISBURY, NC
PO# 659-B40007

Qty	Product Model Number and Description			
1	<p>SELENIA DIMENSIONS TOMO / AWS 5000 - Selenia Dimensions Tomosynthesis with AWS 5000</p> <p>X-Ray Gantry: Constant Potential, High Frequency, Inverter Type Generator kV Range – 20 kV to 49 kV in 1 kV increments mAs Range – 3.0 mAs to 500 mAs mA Range – 10 mA to 200 mA, Lg Focal Spot; 10 mA to 50 mA, Sm Focal Spot</p> <p>X-Ray Tube: Tungsten, Bi-Angular, High Speed, High Heat Capacity Filtration – 0.050 mm Rhodium (Rh); 0.050 mm Silver (Ag); 0.70 mm Aluminum (Al) Focal Spot Size – 0.3 mm (Lg); 0.1 mm (Sm)</p> <p>Digital Image Receptor: Amorphous Selenium, TFT HTC Anti-Scatter Grid w/Auto-Retraction for Geometric Magnification and Tomosynthesis Views Image Receptor Size – Single Plate, 24 cm x 29 cm Image Matrix Sizes – 2560 x 3328 (18 cm x 24 cm); 3328 x 4096 (24 cm x 29 cm) Pixel Size – 0.070 mm Limiting Spatial Resolution – 7.1 lp/mm Ambient Operating Temperature – 20°C to 30°C</p> <p>Standard Acquisition Workstation (AWS 5000): CPU – Multi-Core Intel Based CPU Memory – 16 GB RAM Hard Drive – 1TB min Storage Media – DVD +/- R/W Image Display Technology – Standard, 2 MP Color LCD Display User Interface Display – 1.2 MP Color LCD Display Operating System – Operating System – Microsoft Windows 7 / 64 bit Local Image Store Capacity – Approximately 9,000 Screening Mammography Studies</p> <p>Standard Compression Paddles: 24 cm x 29 cm Screening Paddle 18 cm x 24 cm Screening Paddle 10 cm Contact Paddle 10 cm Magnification Paddle</p> <p>Accessories: Face Shields (2) Magnification Platform Dual Function Footswitches (2)</p> <p>DICOM Services: Worklist, Storage, Storage Commitment, Query/Retrieve, Print</p> <p>IHE Profiles: Scheduled Workflow, Patient Information Reconciliation, Mammography Image</p> <p>Licenses: Diagnostic Imaging License Dynamic Tube Head Motion License</p>			

	<p>Tomosynthesis Modality AWS License (1)</p> <p>If softcopy review of tomosynthesis images is via SecurView, a SecurView Tomosynthesis license for each system must be included in your order and the following SecurView recommendations apply:</p> <p>Software: SecurView v7.3.1 or higher Standalone or Client Computer (applies to DX or RT): Dell Precision T7500 Manager Computer (applies to DX or RT): SecurView Dell PowerEdge T610 Client Display: Barco Nio, Coronis or High Brightness LCD Displays (CRTs not supported) Keypad: Tomo capable keypad</p> <p>Please consult your Hologic Technical Sales Specialist for more details on the above (US only).</p> <p>Training Requirements:</p> <p>A pre-requisite for the Tomosynthesis Modality is ACR or State FFDM accreditation. If your site has not yet been FFDM accredited, you will need to apply to the ACR or State for FFDM certification. Once that is received, you may contact the MSQA for the FFDM certification extension program for tomosynthesis at 301-796-5710.</p> <p>Manuals: User Manual, Service and Maintenance Manual, Quality Control Manual</p> <p>Warranty: Standard One-Year Parts and Labor Warranty; Two-Year Prorated Manufacturer's Warranty on X-Ray Tube</p> <p>NOTE: BEFORE WE CAN ACCEPT YOUR ORDER FOR A SELENIA DIMENSIONS TOMOSYNTHESIS SYSTEM ("SYSTEM") OR TOMOSYNTHESIS OPTION FOR AN EXISTING SELENIA DIMENSIONS SYSTEM, A TOMOSYNTHESIS ENTERPRISE CUSTOMER SITE ASSESSMENT SURVEY ("SURVEY") MUST BE COMPLETED WITH YOUR HOLOGIC REPRESENTATIVE TO ASSURE YOUR CONFIGURATION NEEDS ARE ACCURATELY ADDRESSED. THE ADDITIONAL INFORMATION OBTAINED FROM THE SURVEY MAY MODIFY THE REQUIRED COMPONENTS OF YOUR SYSTEM AND THEREFORE YOUR TOTAL PRICE MAY BE ADJUSTED IN A FINAL QUOTE.</p>
1	<p>Initial Applications - Dimensions 2D - Initial Applications per site - Dimensions 2D Three (3) days of Applications Training for up to 5 technologists per site.</p> <p>Applications must be completed within 12 months of equipment shipment.</p> <p>Please note: Cancellations must be made 48 hours prior to the end of the business week before your scheduled applications to avoid cancellation fees.</p> <p>Initial Applications Added Value: \$5,100.00</p>
1	<p>Initial Applications – Dimensions Tomo - Tech - Initial Applications – Dimensions Tomo Modality - Technologist Two (2) days of Applications Training for up to 5 technologists per site. Applications must be completed within 12 months of equipment shipment. This training cannot be performed until FFDM certification extension is received for tomosynthesis. Please note: Cancellations must be made 48 hours prior to the end of the business week before your scheduled applications to avoid cancellation fees. Initial Applications Added Value: \$3,800</p>
1	<p>Initial Education – Dimensions Tomo - Physicist - Initial Education – Dimensions Tomo Modality - Physicist Included in the price are 8 hours of Hologic tomosynthesis educational training for up to 2 physicists. This consists of 3 hours of Online training; and, 5 hours of QC training with a Hologic Field Engineer. Education must be completed within 12 months of equipment shipment. Visit www.hologic.com/medical-professionals for a list of Hologic educational opportunities. Added Value: \$1,500.00</p>
1	<p>Initial Education – Dimensions Tomo - Radiologist - Initial Education – Dimensions Tomo Modality - Radiologist (Requires SecurViewDX) Included in the price of your first system are 8 hours of Hologic Tomosynthesis educational training for up to 5 radiologists per site. This consists of 8 hours of live, online image review training. This training requires a SecurViewDX workstation (licensed for Tomosynthesis) and a co-located computer with internet connection; ONLY two radiologists to a SecurViewDX per Web Ex Training. Please refer to SecurView hardware/software minimums.</p> <p>Education must be completed within 12 months of equipment shipment. Excludes travel and expenses. Visit www.hologic.com/medical-professionals for a list of Hologic educational opportunities. Added Value: \$7,500.00</p>
1	<p>Standard 2 MP Color Display - Standard 2 MP Color Display for AWS 5000 Standard 2 MP Color Display for AWS 5000</p>
1	<p>ASY-04165 - Bar Code Scanner for Selenia Dimensions with AWS 5000 The barcode scanner option offers an enhancement to the workflow of the AWS 5000. It promotes ease and accuracy in the selection of patients from a worklist and can expedite the imaging procedure.</p>
1	<p>ASY-04194 - Expanded Paddle Kit for Selenia Dimensions (AWS 5000 Packages)</p>

	<p>The optional expanded compression paddle kit includes the following extra paddles:</p> <p>(a) Small breast screening paddle (b) Frameless spot paddle (c) 7.5 cm spot contact paddle (d) 7.5 cm spot magnification paddle</p> <p>This compression paddle kit is designed to enhance both screening and diagnostic patient imaging procedures for the Selenia Dimensions with Acquisition Workstation 5000 or Avia packages.</p>	
1	<p>PRF Field Guide - Pink Ribbon Facility Field Guide</p> <p>A comprehensive business development guide, provided by Hologic, full of proven marketing materials and communication strategies to support your individual business-building efforts.</p> <p>Value Added List Price: \$5,000.00</p>	
1	<p>MammoPad Introductory Package - MammoPad Introductory Package</p> <p>MammoPad Introductory Package Contain:</p> <p>Starter supply of twenty-five (25) 29 x 30 cm MammoPad Breast Cushions MammoPad Initial Order Form</p> <p>Educational Information: "Technical Analysis of Radiolucent Breast Cushions" "Marketing Digital Mammography and a Mammography Breast Cushion to Increase Capacity Volume"</p>	
1	<p>Trade In Program - Selenia Base (3+ Yrs) - Used Equipment Trade-in Discount - Selenia Base (aged 3+ yrs)</p> <p>Used Equipment Trade-in Discount (based on FMV) applicable to Selenia or Selenia Dimensions orders. Indicate the trade-in system serial number with the order.</p>	
1	<p>R2UP-2002 - Cenova Digital Port License Transfer</p> <p>Provides license transfer from one configuration to another. Allows an existing ImageChecker CAD port or Quantra port to be migrated from one digital manufacturer FFDM model to another manufacturer's FFDM model. Also allows an existing ImageChecker CAD port or Quantra port to be migrated from one Cenova server to another, preserving the investment in R2 ImageChecker CAD.</p> <p>Includes:</p> <ul style="list-style-type: none"> • Support for transferring one ImageChecker CAD port or Quantra port: <ul style="list-style-type: none"> • Between manufacturers (e.g. GE to Hologic) • Between models (e.g. Hologic Selenia to Hologic Selenia Dimensions) • From one Cenova server to another <p>NOTE: Cenova serial number must be provided.</p>	
1	<p>WS-LIC-1003-NC - SecurView (DX/RT) Tomosynthesis License</p> <p>This license enables the display of tomosynthesis images on a SecurView DX or RT workstation that is in use at a customer site.</p> <p>It is provided at no charge as part of a promotional program for new Dimensions 3D customers or for customers purchasing the Tomosynthesis Option on an existing Dimensions system (2D to 3D).</p> <p>NOTE: the system serial number of the workstation must be supplied. If the workstation is part of a multi-workstation cluster, the serial number of the Manager must also be supplied.</p>	
2	<p>WS-HW-1002-SC - SV DX or RT Workstation Computer Upgrade with Windows 7.0 (Under Service Contract)</p> <p>The SecurView DX or RT workstation computer upgrade kit allows an existing SecurView DX or RT standalone or client workstation to be upgraded to the latest computer hardware (Dell T7600). This kit is for customers who are ARE covered under a Service Contract, and requires trade-in of existing computer. This product is discounted due to the existing service contract that is expected to be renewed; therefore it is sold as is with NO warranty.</p> <p>SPECIFICATIONS:</p> <p>High End, Windows Based System equipped with a Quad-Core Processor (2.0GHz or higher) 16GB High Speed RAM (3) 2TB SAS drives in RAID 5 array 10/100/1000 Base T Ethernet Network Interface</p> <p>REQUIREMENTS:</p> <p>Replaces any SecurView DX or RT standalone or client. SecurView DX systems equipped with older CRT displays and video cards will have to be upgraded to the newer LCD displays and video cards by ordering the appropriate Product ID.</p> <p>NOTE: Price includes trade-in of existing computer. The warranty information in this product description shall supersede any warranty language in the attached terms.</p>	
1	<p>ASY-02452 - Dimensions Needle Localization Kit (incl. crosshairs and open paddles for both contact and mag)</p> <p>Needle Localization Kit for Selenia Dimensions</p> <p>Includes:</p> <p>Contact Localization Crosshairs (for both contact and magnification procedures) 10 cm Open Localization Paddle 10 cm Open Localization Paddle for Magnification (ASY-02385,ASY-02258,ASY-01993 & ASY-02163)</p>	
3	<p>ASY-04662 - Dimensions Paddle Rack</p>	

Paddle storage rack designed for use with all Selenia Dimensions configurations. Rack holds approximately 3-4 Dimensions paddles. Hologic recommends a minimum of 3 racks for Dimensions with AWS 8000 and a minimum of 2 racks for Dimensions with AWS 5000.

Approximate paddle rack dimensions:

Length: 36 1/4"
Depth: 4" (from the wall out)
Height: 7"

Pitch of 12" to 14" when multiple racks are installed over each other.

Price includes wall mounting bracket component. Does not include installation.