

VA INDIANAPOLIS, IN
PO# 583-B40088

Item No.	Qty	Description
		Optima XR220amx - 30Kw
1	1	<p>Optima XR220amx 30kW</p> <p>Optima XR220amx Digital Mobile Radiographic system - with 30kW generator</p> <p>The Optima XR220amx is a self-contained battery operated mobile radiographic digital X-Ray imaging system designed for performing radiographic exams at the point of care</p> <p>Key Features</p> <ul style="list-style-type: none">o 30kW generatoro Wireless Digital Detector with 6:1 removable grid, back-up tether, QAP (Quality Assurance Proceedure)o Dose Area Product Meter (DAP)o Capable of 100-240V nominal, 50/60 Hz operationo Stand-by mode to eliminate boot up cycles and allow exposure within 25 secondso Exposures can be taken and processed while the unit is chargingo Detector battery charges automatically while while the detector is in the bino Optimized GUI - Technique, image acquisition and display tools in a single integrated user interfaceo The detector can be used in additional wireless enabled GE radiographic systems: please refer to the current literature for system compatibility <p>Productivity</p> <ul style="list-style-type: none">o Up to 1,200 w of power available to minimize

Item No. Qty	Description
	<p>charge time</p> <ul style="list-style-type: none"> o System can be driven within 4 seconds of activation o Pre-programmed techniques per anatomy and patient size o Systems can be used without the detector o Modality Perform Procedure Step (MPPS; SPS/PPS configurable) o Automated and customizable image transfer to PACS and printers o Can reprocess images post acquisition and during an exam o Usage reporting tools by individuals and user groups o System Health dashboard for system status o Bin stores detector and grid o Built-in storage for cleaning wipes, gloves and lead apron o Self-propelled single drive handle control with variable speed of up to 5 km/h (3.1 mph on flat surfaces) forward and reverse to automatically adjusts to the operators pace <p>Wireless Digital Detector Specifications</p> <ul style="list-style-type: none"> o Detector battery can take up to 45 exposures per hour and provide enough power for 3 hours of use on a single charge o Single panel (non-tiled) amorphous silicon detector with a Cesium Iodide scintillator o Image area 40.4cm x 40.4cm (15.9in x 15.9in) o Active matrix 2022 x 2022 pixels o 8mb raw image file size o Pixel Pitch 200 microns o Typical upper dynamic range 7.8mR

Item No. Qty	Description
	<ul style="list-style-type: none"> o Typical DQE @ Olp/mm: (68%) o Two handgrips o Dimensions: L 23.1in., H 17.8in., o T 0.94in. (L 580mm, H 452mm, T 24mm) o Wireless point-to-point network between the system and detector for transferring image data <ul style="list-style-type: none"> – Communication over wide 500MHz channels to achieve very high data rates – Designed to co-exist with 802.11 networks without interference – Frequency: 3.1-10.6 GHz Max Power Output: -41.3 dBm – Max PHY Data rate: 480 Mbps – Effective Throughput: 30-70 Mbps <p>Worklist can be retrieved from HIS/RIS systems and images can be transmitted through the DICOM interface to printers, archival devices (PACS) servers or review workstations</p> o RJ45 10/100/1000 Base T Ethernet port <p>Please refer to the DICOM conformance statement for complete definition of supported DICOM services.</p> <p>Generator</p> <ul style="list-style-type: none"> o 300 mA maximum o kVp and mAs controls o Less than 2% low frequency ripple o Frequency: greater than 100 kHz, super resonant inverter with varying frequency <p>X-ray Source</p> <ul style="list-style-type: none"> o Nominal Tube Voltage (radiographic) - 150kV o Nominal Focal Spot size (IEC 60336) <ul style="list-style-type: none"> – Large Focus - 1.3 mm

Item No.	Qty	Description
2	1	<ul style="list-style-type: none"> – Small Focus - 0.6 mm o Anode Rotation Speed (minimal): 3200 min o Permanent Filtration: 0.9 mm A1/75 kV IEC60522: 1999 o Maximum X-ray Tube Current – Large Focus: 500 mA – Small Focus: 200 mA o Maximum Continuous Heat Dissipation: Without Air-circulator: 170W (238 HU/s) Collimator A pair of independent collimator blades control the X-ray field o 180 lux (1000 Lumen/mt²) light field lamp o The collimator rotates plus and minus 180 degrees with detents at -180, -90, 0, +90 and +180 degrees o Full 43cm x 43cm (17 in.) coverage at a 100cm SID The column may be rotated up to plus or minus 270 degrees from the park position o Drive Inhibit keypad access o Password protected access to patient information for compliance with confidentiality regulations o Automatic safety brake: Operator must hold drive handle to allow system movement o Integrated front bumper stops unit and activates brakes when activated
		<p>Wireless Connectivity</p> <p>Wireless Connectivity for Optima XR220amx and Optima XR200amx</p> <p>802.11 a/b/g n-compatible wireless connectivity</p>

Item No.	Qty	Description
		to hospital network Wi-Fi Certified Compatible with: o 802.11i, Wi-Fi Protected Access 2 (WPA2), WPA 802.1X o AES - TKIP o 64-, 128-WEP o VPN: IPSec - IKE o Management Frame Protection (MFP) EAP Types: - LEAP - LEAP + 128-WEP - LEAP + WPA - EAP - TLS - EAP-TTLS/MSCHAPv2 - EAP-FAST - PEAP-GTC - PEAP/MSCHAPV2
3	1	Auto Protocol Assist Auto Protocol Assist for Optima XR200amx and Optima XR220amx
4	1	Repeat/Reject Analysis Repeat/Reject Analysis for Optima XR220amx/Upgraded Optima XR200amx
5	1	Sterile Protective Drapes - Detector Drape Sterile Protective Drapes - Detector Drape
6	1	Optima XR220amx Training: 3 Days Onsite Optima XR220amx Training: 3 Days Onsite (2 Days + 1 Day) One 2 day and one 1 day TIP onsite training

Item No.	Qty	Description
7	1	<p data-bbox="518 234 821 263">visit for Optima XR220amx.</p> <p data-bbox="518 283 1008 343">Includes TELL expenses. Days provided in two customer visits.</p> <p data-bbox="518 375 1045 482">This training program must be scheduled and completed within 12 months after the date of product delivery.</p> <p data-bbox="518 507 1029 536">X-ray TiP Training, 2 Consecutive Days Onsite</p> <p data-bbox="518 555 1024 584">X-ray TiP Training, 2 Consecutive Days Onsite</p> <p data-bbox="518 607 1076 671">Two Day X-ray Onsite Training provided from 8AM to 5PM, Monday through Friday.</p> <p data-bbox="518 696 928 760">Includes TELL expenses. Days provided consecutively.</p> <p data-bbox="518 787 899 936">This training program must be scheduled and completed within 12 months after the date of product delivery.</p>
8	1	<p data-bbox="518 963 743 992">OPTIMA XR220 AMX</p> <p data-bbox="518 1017 1073 1081">Optima XR220AMX, XR200AMX, & Brivo XR285AMX Full Service Training (Class/Lab) (3 Days)</p> <p data-bbox="518 1102 1104 1253">This service training class covers all three mobile X-ray systems Optima XR220AMX, XR200AMX, and Brivo XR285AMX. Engineer must have completed XR Basic Service R0182RY & R0181RY</p> <p data-bbox="518 1280 1053 1386">or equivalent experience before attending this course. This course must be taken within 2 years from the purchase date.</p>
9	3	<p data-bbox="524 1485 833 1514">Meals And Lodging Expense</p> <p data-bbox="524 1531 1437 1630">Meals and Lodging Expense has been developed to allow the customer the convenience of prepaying for their meals and lodging expenses when attending Technical Service Training at the GE Healthcare Institute located in Waukesha, WI.</p> <p data-bbox="524 1653 1495 1682">The price of this convenience is based on a per day basis. Thus a quantity of 1 is equal</p>

Item No.	Qty	Description
		<p>to 1 day's meals and lodging expense. When purchasing the meals and lodging expense please be mindful of weekend days during the training stay and include 2 days to cover a weekend in the purchase quantity.</p> <p>Examples: A 5-day course needs a quantity of 5. Any course longer than 5 days should include 2 days to account for the weekend stay. Any course longer than 10 days will require an additional 4 days of the meals and lodging expense to cover the 2 weekends of the stay. Thus a 15-day course would have a quantity of 19 days to cover the 2 weekends of the stay. This expense must be used within 2 years from the purchase date.</p> <p>Three meals a day Monday thru Thursday, 2 meals on Friday, plus breaks are provided in the onsite cafeteria. The GE Healthcare Institute cafeteria closes Friday after lunch and reopens Monday morning for breakfast. Weekend meals are the responsibility of the customer.</p> <p>Only for In-resident courses to be taken at the GE Healthcare Institute.</p>
10	1	<p>Airfare Expense</p> <p>The AIRFARE EXPENSE has been developed to allow the customer the convenience to prepay their roundtrip Airfare expenses when attending Technical Service Training at the GE Healthcare Institute located in Waukesha, WI. To be used for engineers attending In-Resident Class/Lab courses for Diagnostic Imaging.</p> <p>Customer will make their Airfare arrangements thru the GE Travel Center. Specific directions will be provided to the customer upon confirmation of class. Please note that this expense must be used within 2 years of the purchase date</p>
11	1	<p>OPTIMAXRAMX WIRELSSSVTRNG</p> <p>OPTIMAXRAMX WIRELSSSVTRNG</p>
12	1	<p>XR CHALLENGE EXAM</p> <p>X-ray Basic Service Challenge Exam</p>

Item No. Qty	Description
	<p data-bbox="516 171 1078 623">This exam is for engineers with considerable X-ray service experience and wish to test out of of the X-ray Basic Service courses(R0181RY and R0182RY). This exam will cover the competencies of the in-resident course. If you fail the challenge exam, you must attend the in-resident course and pass the end of course test to receive credit. Note: Cost of the challenge exam will apply as a credit towards the X-ray Basic Service course.</p> <p data-bbox="516 658 1049 733">This exam must be taken within 2 years from the purchase date.</p>