

XR-U/S, VAMC AUSTIN, TX

PO# 674-B40011

Line #	Part #	Description	Qty
1		<b>CX50 CompactXtreme General Imaging</b>	1
		<b>Interface:</b> 15.0 inch high resolution display with wide viewing angle Quick Keys and Active Mode Laptop style Alphanumeric QWERTY keyboard 8 TGCs and 2 LGCs Ergonomic carrying handle Includes AC adapter , power cord and system battery pack 2 USB flash drives on system 80 GB hard drive Internal DVD RW drive	
		<b>Architecture:</b> All-digital compact broadband beamformer, Microfine 2D focusing with Dynamic Focal Tuning that includes Advanced X-Res signal processing, 170 dB full time input dynamic range 18,432 digitally-processed channels, Continuously variable steering in 2D, color and Doppler modes 2D Opt signal processing with 4X multi-line parallel processing and frequency compounding.	
		<b>Intelligent Controls:</b> The CX50 has been designed to make portable exams easy and efficient. With a single button, iSCAN technology automatically samples data for a new level of 2D and Doppler optimization iSCAN one-touch Intelligent Optimization, iSCAN one-touch Intelligent Color Optimization, iSCAN Doppler one-touch optimization.	
		<b>Transducers:</b> Supports Compact family of transducers featuring PureWave imaging technology in the S5-1, CX7-2t, C5-1, D5CWC and C9-3V. Also supports the high resolution S12-4, S8-3, C8-5 and L12-3 transducers. All transducers provide breakthrough frequency bandwidths and array configurations. These transducers also have ergonomically designed lightweight flexible cables and compact connectors.	
		<b>Modes:</b> 2D M-mode Anatomical M-mode Color M-mode Pulsed Wave Doppler Color Power Angio (CPA) Continuous Wave Doppler Invert and Color Invert Color compare mode Dual mode Duplex for simultaneous 2D and Doppler 2D Optimization Signal Processing Live Compare Tissue Harmonic Imaging (THI) Reconstructed zoom with pan (read zoom)	

Write zoom  
 Pulse Inversion Harmonic imaging  
 Adaptive Doppler  
 Adaptive Color Doppler  
 Color Tissue Doppler imaging  
 Pulsed Wave Tissue Doppler imaging  
 Active Native Data - manipulation of image data  
 Cineloop review  
 Acquisition, storage, and display in real-time and duplex modes of up to 500 frames  
 On-board workstation-class data management with thumbnail previews and storage of images, loops, and reports. Retrospective and prospective clip capture to internal drive or removable media  
 Integrated DVD/CD burning capability for storage of images or export in DICOM, JPEG and .avi for PC compatibility. Philips DICOM viewer option to imbed in media transfer for easy viewing of study on most PCs.  
 Maintenance and Serviceability  
 Remote Access for Expedient Clinical and Technical Support  
 Flexible Service Agreements  
 Clinical Application and Educational Support  
 Scheduled Preventative Maintenance and System Optimization

**2**
**DICOM Package**
**1**

**Dicom**

Provides DICOM 3.0 network Print and Store, Performed Procedure Step (PPS), and Modality Worklist functionality. Networking capability supportable in both wired and wireless environments

**Dicom Structured Reporting**

Provides adult echo and vascular DICOM Structured Reporting. (Requires appropriate adult echo or vascular clinical options.

**3**
**Radiology Clinical Option Package**
**1**

This clinical option package includes the Abdominal, Vascular, OB/GYN, Musculoskeletal and Small Parts clinical options.

**Abdominal Clinical Option**

This clinical option includes abdominal Tissue Specific Imaging software, SonoCT and Freehand 3D for abdominal applications. This clinical option also includes analysis and reporting packages for abdominal applications. Allows operation for abdominal applications of the C5-1, S5-1 and L12-3 transducers

**Vascular Clinical Option**

This clinical option includes Tissue Specific Imaging software and SonoCT for Cerebrovascular, Peripheral vascular, abdominal vascular and Transcranial applications. This clinical option also includes in depth analysis and reporting packages for vascular applications. Freehand 3D is also provide within this clinical option. Allows operation for vascular applications of the C5-1, S5-1, L12-3 and D5cwc transducers

**OB/GYN Clinical Option**

This clinical option includes OB/GYN Tissue Specific Imaging software, SonoCT and Freehand 3D for OB/GYN applications. This clinical option also includes in depth analysis and reporting packages for OB/GYN applications. Allows operation for OB/GYN applications of the C5-1 and C9-3V transducers.

**Musculoskeletal Clinical Option**

This clinical option includes MSK Tissue Specific Imaging software, SonoCT and Freehand 3D for MSK applications. This clinical option also includes analysis and reporting packages for MSK

applications. Allows operation for MSK applications of the L12-3.

#### Small Parts Clinical Option

This clinical option includes small parts Tissue Specific Imaging software, SOnoCT and Freehand 3D for a wide range of small parts applications (eg. Breast, thyroid, testical) . This clinical option also includes analysis and reporting packages for small parts applications. Allows operation for small parts applications of the L12-3 and C5-1.

4	<b>Cart with Multi-port Adapter</b>	<b>1</b>
	Highly mobile cart that features hardware module to support transducer switching among up to three imaging transducers at the touch of a button. Includes: 4 swivel wheels with 2 locking casters, rear handle, micro-positioning grips, quick-connect tray, storage shelf, footrest, internal isolation transformer, B&W printer brackets, integrated transducer connector holder, gel holders and cable management. Includes USB hub for additional connectivity.	
5	<b>USA Power Cord</b>	<b>1</b>
6	<b>C5-1 Broadband Curved Array Transducer</b>	<b>1</b>
	C5-1 broadband, Curved Array PureWave crystal transducer with 5 to 1 Mhz extended operating frequency range for abdominal, obstetrical, gynecological, interventional, vascular, contrast acute care and regional anesthesia applications.	
7	<b>Transducer L12-5, 50</b>	<b>1</b>
	High frequency linear array transducer with 12 to 5 MHz extended operating frequency range for small parts, vascular, abdominal, and MSK applications. (Compact connector)	
8	<b>Transducer, L15-7io</b>	<b>1</b>
	High frequency compact linear array transducer with 15 to 7 MHz extended operating frequency range for cardiac and vascular surgery, superficial vascular, and MSK applications. (Compact connector)	
9	<b>English Manual</b>	<b>1</b>
	Operation Manual	
10	<b>1 Day PAS Onsite</b>	<b>1</b>
	<b>1 Day PAS Onsite</b> - Ultrasound system or upgrade onsite training provided by a PAS (Product Applications Specialist) for specific system applications or upgrades; not per modality. <i>Education is provided Monday - Friday during normal business hours.</i> Note: Philips Healthcare personnel are not responsible for actual patient contact or operation of equipment during education sessions except to demonstrate proper equipment operation. The training sessions should be attended by the appropriate healthcare professional as identified by the department director. <i>Repeat training for staff non-attendance will not be accepted.</i> Site must be patient-ready to meet training expectations. All onsite training day expires within 90 days from system or upgrade installation date. Exceptions are for 3D Stress onsite training (which expires 9 months from system or upgrade installation date) and Fusion & Needle Navigation onsite training (which expires 180 days from system or upgrade installation date).	
11	<b>1st SVC Manual for Gov</b>	<b>1</b>