

675-B47013  
ORLANDO, FL.

Line #	Part #	Description	Qty	Each	Price
--------	--------	-------------	-----	------	-------

1	**NNAP204	CX50 2D xMATRIX	1		
---	-----------	-----------------	---	--	--

**Interface:**

- 15.0 inch high resolution display with wide viewing angle
- Quick Keys and Active Mode
- Laptop style Alphanumeric QWERTY keyboard
- 8 TGCs and 2 LGCs
- Ergonomic carrying handle
- Includes AC adapter , power cord and system battery pack
- 2 USB flash drives on system
- 80 GB hard drive
- Internal DVD RW drive

**Architecture:**

All-digital compact broadband beamformer, Microfine 2D focusing with Dynamic Focal Tuning that includes Advanced X-Res signal processing, 170 dB full time input dynamic range 18,432 digitally-processed channels, Continuously variable steering in 2D, color and Doppler modes 2D Opt signal processing with 4X multi-line parallel processing and frequency compounding.

**Intelligent Controls:**

The CX50 has been designed to make portable exams easy and efficient. With a single button, iSCAN technology automatically samples data for a new level of 2D and Doppler optimization iSCAN one-touch Intelligent Optimization, iSCAN one-touch Intelligent Color Optimization, iSCAN Doppler one-touch optimization.

**Transducers:**

Supports Compact family of transducers featuring PureWave imaging technology in the S5-1, CX7-2t, C5-1, D5CWC. Also supports the high resolution S12-4, S8-3, C8-5 and L12-3 transducers. All transducers provide breakthrough frequency bandwidths and array configurations. These transducers also have ergonomically designed lightweight flexible cables and compact connectors.

**Modes:**

- 2D
- M-mode
- Anatomical M-mode
- Color M-mode
- Pulsed Wave Doppler
- Color Power Angio (CPA)
- Continuous Wave Doppler
- Invert and Color Invert
- Color compare mode

## 100662 CX50 CompactXtreme

Line #	Part #	Description	Qty	Each	Price
		<ul style="list-style-type: none"> <li>• Dual mode</li> <li>• Duplex for simultaneous 2D and Doppler</li> <li>• 2D Optimization Signal Processing</li> <li>• Live Compare</li> <li>• Tissue Harmonic Imaging (THI)</li> <li>• Reconstructed zoom with pan (read zoom)</li> <li>• Write zoom</li> <li>• Pulse Inversion Harmonic imaging</li> <li>• Adaptive Doppler</li> <li>• Adaptive Color Doppler</li> <li>• Color Tissue Doppler imaging</li> <li>• Pulsed Wave Tissue Doppler imaging</li> <li>• Active Native Data - manipulation of image data</li> <li>• DICOM Networking <ul style="list-style-type: none"> <li>• Ethernet @100Mb/second Includes DICOM wired and wireless "G" and "N".</li> <li>• Provides DICOM 3.0 network print and store and storage commitment</li> <li>• Performed Procedure Step (pps)</li> <li>• Modality Worklist</li> </ul> </li> <li>• DICOM Structured Reporting</li> </ul>			

### Cineloop review

- On-board workstation-class data management with thumbnail previews and storage of images, loops, and reports. Retrospective and prospective clip capture to internal drive or removable media  
Integrated DVD/CD burning capability for storage of images or export in DICOM, JPEG and .avi for PC compatibility. Philips DICOM viewer option to imbed in media transfer for easy viewing of study on most PCs.
- Maintenance and Serviceability
- Remote Access for Expedient Clinical and Technical Support
- Flexible Service Agreements
- Clinical Application and Educational Support

### Cinical Education

**1 Day PAS Onsite** - Ultrasound system or upgrade onsite training provided by a PAS (Product Applications Specialist) for specific system applications or upgrades; not per modality. *Education is provided Monday - Friday during normal business hours.* Note: Philips Healthcare personnel are not responsible for actual patient contact or operation of equipment during education sessions except to demonstrate proper equipment operation. The training sessions should be attended by the appropriate healthcare professional as identified by the department director. *Repeat training for staff non-attendance will not be accepted.* Site must be patient-ready to meet training expectations. All onsite training day expires within 90 days from system or upgrade installation date. Exceptions are for 3D Stress onsite training (which expires 9 months from system or upgrade installation date) and Fusion & Needle Navigation onsite training (which expires 180 days from system or upgrade installation date).

## 100662 CX50 CompactXtreme

Line #	Part #	Description	Qty	Each	Price
<p>***THE NUMBER OF ONSITE TUITIONS YOU RECIEVE MAY VARY BASED ON PURCHASED OPTIONS. PLEASE CONSULT YOUR SALES REPRESENTATIVE FOR FURTHER DETAILS***</p> <p><i>All Tuitions must be registered prior to the expiration date. The course chosen must be taken within 90 days of expiration.</i></p>					
2	**NUSN033	<b>Adult Echo Clinical Option</b>	1		
<p>Tissue Specific Imaging software for adult cardiac ultrasound applications. Display optimization software with Tissue Specific presets for adult cardiac imaging and Doppler applications. Analysis software package includes cardiac imaging protocol measurements and configurable reports and finding codes. Active native data for post-process optimization and advanced XRES adaptive image processing for improved tissue conspicuity. iSCAN intelligent one-button optimization for adaptive gain compensation in 2D, Doppler, Tissue Doppler Imaging and LVO contrast functions. Includes Live compare mode, cardiac High-Q Automatic Doppler Analysis and respiration waveform from chest impedance. Allows operation of S5-1, CX X7-2t and D2cwc transducers.</p>					
3	**NUSM019	<b>Exam Protocols/Smart Exam</b>	1		
<p>Fast Exam provides easy to use, customizable guides that help the clinician complete studies on every patient. An on screen menu, guides the clinician through the required views for a specific exam type, automatically enters annotation and build the patient report.</p>					
4	**FUS5042	<b>Travel Case</b>	1		
<p>Padded wheeled travel case with front zipper pockets &amp; telescoping handle. Includes a retractable stacking extention for transporting additional supplies and a customized transducer &amp; gel carrying bag for up to 4 transducers. Dimensions: Accommodates CX50 imaging system, AC adapter, separate transducer carry bag for up to 4 transducers and gel,</p>					
5	**FUS5047	<b>Cart</b>	1		
<p>Highly mobile cart that includes: 4 swivel wheels with 4 locking casters, rear handle, micropositioning grips, quick-connect tray, utility drawer, storage shelf, footrest, internal isolation transformer, B&amp;W printer brackets, integrated transducer connector holder, gel holders and cable management. Includes USB hub for additional connectivity.</p>					
6	**FUS5101	<b>USA Power Cord</b>	1		
7	**FUS5083	<b>D2cwc Static Transducer</b>	1		
<p>Non-imaging 2 MHz PW/CW Doppler transducer for cardiac applications</p>					
8	**FUS5120	<b>S5-1 Transducer</b>	1		
<p>Features compact connector designed for reliability and improved ergonomics. Compatible with both EPIQ and CX50 systems. Manufactured in accordance with the European Union's Restriction of Hazardous Substances (RoHS) directive.</p> <p>PureWave crystal Sector array transducer with 5 to 1 MHz extended operating frequency range for adult cardiology, abdominal, vascular, TCD and Acute Care.</p>					
9	**FUS7000	<b>English Manual</b>	1		
<p>Operation Manual</p>					
10	**FNA8204	<b>1st SVC Manual for Gov</b>	1		