

Exhibit 6 – Manufacturing Equipment

The PCC has the following manufacturing equipment:

Blenders:

Blender Brand	Blender Type	Size
Globe Pharma Mini-Blender -	V-Shell	1 and 2 quart
Keith	V-Shell	8 quart
Keith	V-Shell	16 quart
2 Cubic Foot	V-Shell	2 Cubic Foot
MVP V-Shell	V-Shell	3 Cubic Foot
AIM Ribbon Blender	Ribbon blender	5 Cubic Foot

Mill

Machine	Function
Fitzpatrick D6A Comminutor	grinds tablets or course raw materials into powder

Sieve Sizes

8" Mesh	12" Mesh
8	12
10	18
11	20
14	30
20	35
30	40
40	80
60	N/A
140	N/A
200	N/A

Capsule Fillers Powder:

Machine	Capsule Tooling Sizes	Maximum Line Speed
Bosch 701 GKF	00, 0, 1, 2, 3 and 4	700/min
Bosch 702 GKF	00, 0, 1, 2, 3 and 4	700/min
Bosch 702 GKF "mini-bowl"	00, 0, 1, 2, 3 and 4	Small capacity bowl for pilot fills. Bowl fill volume is approximately 350 mL \pm depending on product density
Ultra 8 II	00, 0, 1, 2, 3 and 4, DBAA	2,000 to 5,000/hour depending on product. Semi-automated capsule filler typically used when running matching placebo for over-encapsulated products.
Sitco Minicap	DBAA, 00, 0, and 1	Manual filling. Maybe 300/hour

Capsule Filler Pellets

Machine	Capsule Tooling Sizes	Maximum Line Speed
Bosch 702 GKF	00	500/min

Capsule Filler, hard gelatin capsules, non aqueous Liquids:

Machine	Capsule Tooling Sizes	Maximum Line Speed
Shionogi	Size 0	400/min

Over-encapsulation of Tablets/Capsules

Machine	Capsule Tooling Sizes	Maximum Line Speed
Ultra 8 II	00, 0, 1, 2, 3 and 4, DBAA	2,000 to 5,000/hour depending on product. Semi-automated capsule filler typically used when running matching placebo for over-encapsulated products or for overencapsulating an oblong tablet or a capsule in another capsule.
Minicap	DBAA	Essentially manual filling. Maybe 300/hour

Tablet Press

Machine	Capsule Tooling Sizes	Maximum Line Speed
Korsch Tablet Press	Tooling specific to product	≤ 800/minute

Coaters

Machine	Capsule Tooling Sizes	Maximum Capacity
Vector VP-1355 Hi-Coater System	VPC-1355 Perforated	13 Liters
O'Hara	8.5" pan	750 g max tabs
	10 " pan	1 Kg max tabs
	15" pan	3 Kgs max tabs

Spray Guns for Coaters

Machine	Spray Gun
Vector VP-1355 Hi-Coater System	
O'Hara	

Sotax AT4 In-Line, In-Process Tablet and Capsule Tester

Machine	Capsule Tooling Sizes	Maximum Line Speed
Sotax AT4	Adjusts to product size	N/A. Will sample at intervals designated by operator/batch record.