

NM BONE DEN, VAMC DENVER, CO

PO# 554-B50008

Qty	Product Description
1	<p>HORIZON-W Advanced Fan-Beam DXA Technology One Pass™ Acquisition Technique; Multi-Detector Array Scanning Method New 128 high-resolution multi element ceramic digital detector array New High Frequency X-ray Generator Indexing Scan Table with Positioning Accessories Motorized Table and C-arm Dynamic Internal Reference System for Continuous Calibration Computer Console QDR Anthropomorphic Spine Phantom</p> <p>Standard Computer Package Computer Workstation with Dual Core 3 GHz Windows® 7 Professional 250 GB hard drive 2 GB RAM 19" Widescreen LCD Monitor HP Professional Series Color DeskJet® printer DVD RAM drive</p> <p>Standard Configuration: Hologic APEX 5™ Operating System Automatic PASS/FAIL Quality Control Express BMD™ 10 Second Acquisition Single Energy Scan Display Capability Window/Level Control for Image Optimization</p> <p>APEX Productivity Tools Express Exam™ Workflow Management OneTime™ Auto Analysis with Histogram ProTech with DXApro™</p>

	<p>Auto Hip Positioning Reposition/Rescan Feature Automatic Scan Comparison for Serial Exams Least Significant Change Configuration</p> <p>Horizon Advance Reporting Solutions QDR OnePage™ Report with Rate of Change Assessment FRAX® 10 Year Fracture Assessment New Dual Hip™ Report New integrated Report Writer DX™</p> <p>Horizon Scan and Analysis Protocols AP Lumbar Spine with Automatic Low Density Analysis and Scoliosis Analysis Proximal Femur, Automatic Low Density Analysis and Hip Structure Analysis (HSA) Dual Hip™ Forearm Whole Body BMD Advanced Body Composition™ Analysis with InnerCore™ Visceral Fat Assessment IVA HD with Image Pro High Resolution Imaging Capability Quantitative Morphometry Integrated Physician's Viewer with MXApro™ New Atypical Femur Fracture Assessment (AFF) High Resolution Imaging Capability Pediatric Analysis for Spine, Femur, Forearm, Infant Spine, and Whole Body</p> <p>DXA Practice Manager Patient Database Manager On-line Normative Reference Values and User-Definable Values NHANES and Ethnic and Pediatric Reference Data Context Sensitive Help General Region of Interest DAP Reporting</p> <p>QDR Reference Manual CD Two Days (12 Hours) of ECE Approved Application Training Installation and a Twelve (12) Month Comprehensive Warranty</p>
1	18.5" FLAT PANEL MONITOR
1	ASSY, P4 COMPUTER, PCIE, W7
1	FRAX OPTION FOR APEX
1	ATYPICAL FEMUR FRACTURE OPTION
1	APEX VISCERAL FAT OPTION
1	<p>DXA TRAINING, INITIAL, 1.5 EIGHT-HOUR DAYS, 1SITE, MAX 4 STAFF</p> <ul style="list-style-type: none"> • One and a half (1.5) consecutive days (total of 12 hours) of Applications Training for up to (4) Operators per site. • Applications must be completed within 12 months of equipment shipment. • Please note: Cancellation must be made 48 hours prior to your scheduled applications to avoid cancellation fees
1	<p>INKJET-COLOR HP Office Inkjet-Color Printer (CMP-00993) included in all APEX 4.XX Windows 7 Upgrades</p>
1	<p>KIT, PCIE CABLE, INTERFACE 30' 30 FT CABLE INTERFACE Kit</p>
1	POWER KIT FOR 120VAC
1	<p>APEX IRIS CONNECTIVITY SUITE W/ENTERPRISE DATA MANAGEMENT IRIS Enterprise-APEX for systems on APEX 4.xx Includes DICOM Storage DICOM Modality Worklist Remote Physician's Viewer Physician's Report Writer HL7 Enterprise Data Management</p> <p>Installation and Twelve (12) Month Comprehensive Warranty for IRIS Enterprise</p>

Warranty Information:
12 Months Parts and Labor
