

CHIEF, A & MM B53035
 V.A. Medical Center
 800 ZORN AVE
 LOUISVILLE, KY 40206

P.O.# 603-B53035

RADIOLOGY/ROOM B-140

Qty	Product Number	Product Description
1	HORIZON-A	<p>HORIZON-A Advanced Fan-Beam DXA Technology One Pass™ Acquisition Technique; Multi-Detector Array Scanning Method New 216 high-resolution multi element ceramic digital detector array New High Frequency X-ray Generator Indexing Scan Table with Positioning Accessories Motorized Table and Rotating C-arm Dynamic Internal Reference System for Continuous Calibration Computer Console QDR Anthropomorphic Spine Phantom</p> <p>Standard Computer Package Computer Workstation with Dual Core 3 GHz Windows® 7 Professional 250 GB hard drive 2 GB RAM 19" Widescreen LCD Monitor HP Professional Series Color DeskJet® printer DVD RAM drive</p> <p>Standard Configuration: Hologic APEX 5™ Operating System Automatic PASS/FAIL Quality Control Express BMD™ 10 Second Acquisition Single Energy Scan Display Capability Window/Level Control for Image Optimization</p> <p>APEX Productivity Tools Express Exam™ Workflow Management OneTime™ Auto Analysis with Histogram ProTech with DXApro™ Auto Hip Positioning</p>

	<p>Reposition/Rescan Feature Automatic Scan Comparison for Serial Exams Least Significant Change Configuration</p> <p>Horizon Advance Reporting Solutions QDR OnePage™ Report with Rate of Change Assessment FRAX® 10 Year Fracture Assessment New Dual Hip™ Report New integrated Report Writer DX™</p> <p>Horizon Scan and Analysis Protocols AP Lumbar Spine with Automatic Low Density Analysis and Scoliosis Analysis Supine Lateral Spine with Baseline Compensation Proximal Femur, Automatic Low Density Analysis and Hip Structure Analysis (HSA) Dual Hip™ Forearm Whole Body BMD Advanced Body Composition™ Analysis with InnerCore™ Visceral Fat Assessment IVA HD with Image Pro High Resolution Imaging Capability Quantitative Morphometry Integrated Physician's Viewer with MXApro™ New Atypical Femur Fracture Assessment (AFF) High Resolution Imaging Capability Pediatric Analysis for Spine, Femur, Forearm, Infant Spine and Whole Body</p> <p>DXA Practice Manager Patient Database Manager On-line Normative Reference Values and User-Definable Values NHANES and Ethnic and Pediatric Reference Data Context Sensitive Help General Region of Interest DAP Reporting</p> <p>QDR Reference Manual CD Three Days (16 Hours) of ECE Approved Application Training Installation and a Twelve (12) Month Comprehensive Warranty</p>
1	18.5" FLAT PANEL MONITOR
1	ASSY, P4 COMPUTER, PCIE, W7
1	FRAX OPTION FOR APEX
1	ATYPICAL FEMUR FRACTURE OPTION
1	APEX VISCERAL FAT OPTION
1	<p>DXA A TRAINING, INITIAL, 2 EIGHT-HOUR DAYS , 1 SITE, MAX 4 STAFF</p> <ul style="list-style-type: none"> • Two (2) consecutive days or can be split over three (3) consecutive days (16 hours total) of Applications Training for up to (4) Operators per site • Applications must be completed within 12 months of equipment shipment. • Please note: Cancellation must be made 48 hours prior to your scheduled applications to avoid cancellation fees
1	<p>INKJET-COLOR HP Office Inkjet-Color Printer (CMP-00993) included in all APEX 4.XX Windows 7 Upgrades</p>
1	<p>KIT, PCIE CABLE, INTERFACE 30' 30 FT CABLE INTERFACE Kit</p>
1	POWER KIT FOR 120VAC
1	<p>APEX IRIS CONNECTIVITY SUITE W/ENTERPRISE DATA MANAGEMENT IRIS Enterprise-APEX for systems on APEX 4.xx Includes DICOM Storage DICOM Modality Worklist Remote Physician's Viewer Physician's Report Writer HL7 Enterprise Data Management</p> <p>Installation and Twelve (12) Month Comprehensive Warranty for IRIS Enterprise</p>
1	<p>APEX DATABASE SYNCHRONIZER Requires APEX 3.2 and higher Stand Alone Application Enables synchronization of 2 or more QDR databases on a network</p>

1	<p>Allows sharing of patient information between Hologic QDR scanners</p> <p>RESEARCH OPTION KIT APEX</p> <p>Prosthetic Hip option allows the scanning of hips that include metal, for example a femur with a prosthesis. The metal components of the implant are removed from the bone mass measurement.</p> <p>Small Animal option contains rat whole body, body composition and region high resolution applications. Includes tissue calibration bar for body composition measurements</p> <p>Infant Whole body and Body Composition option is designed to be used on children weighing approximately 6.6 lbs-22 lbs. Whole Body Systems only.</p>
1	<p>DXA BIOMED TRAINING, 5 DAYS, 1 ENGINEER</p> <p>This 5 day program course provides the information necessary to maintain, troubleshoot, and repair Hologic's current bone densitometry systems. The course is conducted at Hologic's corporate headquarters Massachusetts and is a combination of lecture/lab work with approximately sixty percent hands-on instruction.</p>
2	<p>DXA TRAINING, ADDITIONAL, 1 SIX-HOUR DAY, 1 SITE</p> <ul style="list-style-type: none"> • Additional Apps Training – Extra (6) hours

Warranty Information:
12 Months Parts and Labor