

WHSE/PAD B54011
V.A. Medical Center
3801 Miranda Ave
Building 50
Palo Alto, CA 94304

P.O.# 640-B54011

Qty	Product Description
1	<p>SECURVIEW-DX 400 (STD/CLNT) SecurView DX is cleared under a 510(k) submission for display, manipulation, and interpretation of FFDM images (from any vendor), as well as other DICOM multi-modality breast images.</p> <p>System Specifications:</p> <p>High End Workstation High End, six core, Xeon Processor Windows Based System (2.0GHz or higher) 16GB High Speed RAM (3) 2TB SAS drives in RAID 5 array 10/100/1000 Base T Ethernet Network Interface</p> <p>Display Hardware Barco High Brightness Dual 5MP Display (Barco model MDMG-5221) with "blue base". Luminance - 1000 cd/m2. Adjustable viewing angle, I-Luminate feature. Active screen size (diagonal) 541 mm each. Dimensions with stand: 408x582-642 x 250 mm each display (this model has two displays): High End 10-bit Medical Grade Display Card for high brightness</p> <p>Attached Devices Multi-function Mouse Backlit Keyboard Dedicated Workflow Keypad Designed for Mammography Applications Power UPS</p> <p>Software Windows 7 Operating System Dedicated Mammography Based Reviewing Software (User Specific Hanging Protocols, Workflow, etc) Integrated CAD Display (CAD Software and Systems sold separately) Automated Display QC Software</p>

		<p>DICOM Connectivity - Print Store, Query/Retrieve, CAD SR Objects, Grayscale Software Presentation State Objects</p> <p>User Manual Installation Manual</p> <p>Options: Advanced Multimodality Module, Application Synchronization Module, MPE Software, SecurView DX Manager, Barcode Scanner</p> <p>The SecurView software may include certain third party software products and libraries. In purchasing this product, the customer agrees and accepts that such software is licensed to the customer under the existing terms and conditions of the applicable license.</p>
1		SVDX LICENSE & ACCESSORY KIT - US
1		UPS FOR SVDX, SVRT, OR MANAGERS - US
1		<p>SECURVIEW DX TRAINING, INITIAL, 3 DAYS, 1 SITE, MAX 5 RADS</p> <p>* Three (3) Consecutive days of applications Training for up to 5 Radiologist per site.</p> <p>* Applications must be completed within 12 months of equipment shipment.</p> <p>* Please note: Cancellation must be made 48 hours prior to the end of the business week before your scheduled applications to avoid cancellation fees</p> <p>* Initial Applications Added Value: \$5,100.00</p>

Warranty Information:
12 Months Parts and Labor