

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0001	Surgical Drill	Ergonomic hand piece - powered by electrical power console - Accepts Wires from 0.7mm – 2.0mm.	MicroAire/5995U	Installation/assembly required	Non-Disposable	EA	2
E0-0002	Surgical Drill	Electric Drill - ability to 2mm reciprocating stroke at least 4000 cycle per minute - ability to adjust speed throttle on hand piece by thumb - powered by electric console - ability to sterilize instrument and electric cable	MicroAire/PAL-650	Installation/assembly required	Non-Disposable	EA	2
E0-0003	Surgical Drill	Electric hand piece - used for Small bone surgery - ability of Variable-speed trigger- ability for precise control- ability for Forward/reverse control on instrument base	MicroAire/1641	Installation/assembly required	Non-Disposable	EA	2
E0-0004	Surgical Drill	Pneumatic Drill- speed ability At least 750 rpm Saw - torque ability 50 in. lbs. (5.6 Nm) - Used for Large bone surgeries -ability to connect to hall style air hose connector -ability to use Large bone high torques reaming couplers- drill couplers- Wire and pin driver - capable of Reamer speed of At least 250 rpm - capable of reamer torque: 150 in. lbs. (16.9 Nm)	MicroAire/7105-200	Installation/assembly required	Non-Disposable	EA	2
E0-0005	Surgical Drill	Pneumatic drill - ability to connect with Hall air hose connector- capable of 12 position rotating head -ability to blade locking hub. ability to use Large power oscillating saw blades - capable of Saw speed of 13-000 cpm- capable arc of 4.5°- ability to latch-style blade locking system	MicroAire/7205-200	Installation/assembly required	Non-Disposable	EA	2
E0-0006	Surgical Drill	Pneumatic drill - capable of using Large power reciprocating saw blades - ability to connect with Hall air hose connector- capable of twist lock blade hub -a ability of 360° blade rotation - capable of Saw speed of 14-000 cpm- capable of Blade stroke of 3.2 mm	MicroAire/7405-200	Installation/assembly required	Non-Disposable	EA	2
E0-0007	Surgical Drill	Cordless Hand piece - capable of Dual-trigger that provides forward/reverse/oscillate operation - ability for variable-speed triggers - ability to use in Single-hand operations- ability to be sterilize- used for Small bone surgeries	MicroAire/6643	Installation/assembly required	Non-Disposable	EA	2
E0-0008	Surgical Drill	Pneumatic Hand piece - ability to connect to MicroAire Style Air Hose Connector ability to use couplers for drilling- reaming- pinning- and osteotomies accessories	MicroAire/6640-100	Installation/assembly required	Non-Disposable	EA	2
E0-0009	Surgical Drill	Pneumatic Hand piece - ability to connect to Hall Style Air Hose Connector - ability to use couplers for drilling- reaming- pinning- and osteotomies	MicroAire/6640-200	Installation/assembly required	Non-Disposable	EA	2

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0010	Surgical Cart	Ability for for tumescent liposuction- Pumps cable of infiltration and Aspiration - Ability for touchscreen interface- capable of an integrated control system -ability of lockable rolling casters -capable of Adjustable 50 – 700ml/min infiltration pump - ability to operate manually or with Footswitch - capable of adjustable IV pole - capable of Aspiration pumpb flow rate of 4.2 SCFM- Capable of adjusting vaccum pressure up to 29 inHg - ability to use with any handpiece - ability to control Aspirations on/off by handpiece- footswitch or console - capable of dual-1/3HP pumps - capable of dry rotary-vane components - ability for digital LCD touchscreen controls - capable of multiple canister-mounting positions - capble of adjustable equipment storage to include front and rear baskets	MicroAire/ASP-1021	Installation/assembly required	Non-Disposable	EA	2
E0-0011	ADAPTOR	Ability to connect Toomey to Luer Lock Syringe- compatable surgical cart and lyposuction pumps	MicroAire/ASP-ADP-2		Disposable	EA	5
E0-0012	ADAPTOR	Ability to connect Toomey to Leur Lock Needle - Compatible surgical cart and lyposuction pumps	MicroAire/ASP-ADP-3		Disposable	EA	5
E0-0013	WASTE CANISTER	Ability to hold 1200 CC - liner-free - compatible with surgical cart and lyposuction pumps	MicroAire/ASP-CAN-1200	Installation/assembly required	Disposable	10PK	4
E0-0014	WASTE CANNISTER	Ability to hold 1200 CC of waste -ability to be liner-free - compatible with surgical cart and lyposuction pumps	MicroAire/ASP-CAN-2000	Installation/assembly required	Disposable	10PK	4
E0-0015	CLAMP	Stainless steel - C-Clamp- attaches filter stand to table or stand	MicroAire/ASP-CAN-2C	Installation/assembly required	Disposable	EA	2
E0-0016	CANISTER STAND	Ability to hold 4 canister up to 2000 CC - compatible with surgical cart and lyposuction pumps - ability to be reused - ability to be autoclavable - ability to clamp to mayo stand or surgical table	MicroAire/ASP-CAN-2R	Installation/assembly required	Disposable	EA	2
E0-0017	Filter	3 liter filter- dual chamber-ability to fat harvest and separate - ability to maintain high cell viability - disposable - ability of a buoyance filtration	MicroAire/ASP-CAN-2S	Installation/assembly required	Disposable	1PK	10
E0-0018	CABLE	Ability to connect surgical cart and lyposuction pump to electric power source	MicroAire/ASP-CBL-1020	Installation/assembly required	Disposable	EA	2

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0019	CABLE	Ability to connect surgical cart and lyposuction pump to electric power source	MicroAire/ASP-CBL-5020	Installation/assembly required	Disposable	EA	2
E0-0020	BIOFILTER ASSEMBLY	Compatible with surgical cart and lyposuction pumps	MicroAire/ASP-FLTR	Installation/assembly required	Disposable	EA	4
E0-0021	FOOT PEDAL	Ability to operate lyposuction machine - Compatible with surgical cart and lyposuction pumps	MicroAire/ASP-FOOT-1	Installation/assembly required	Disposable	EA	2
E0-0022	CANISTER RING	Ability to hold the 1200 CC waste canister to the cart	MicroAire/ASP-RNG-1200	Installation/assembly required	Disposable	EA	8
E0-0023	CANISTER RING	Ability to hold the 2000 CC waste canister to the cart	MicroAire/ASP-RNG-2000	Installation/assembly required	Disposable	EA	8
E0-0024	SHELF	Compatible with surgical cart and ability to hold power console	MicroAire/ASP-SHLF-1	Installation/assembly required	Disposable	EA	2
E0-0025	TUBING	Aspiration - Sterile - 10 foot Long	MicroAire/ASP-TB-ASP		Disposable	10PK	5
E0-0026	TUBING	Non Sterile - Vacuum - 18inches - connects the biofilter to waste cannistors	MicroAire/ASP-TB-VAC		Disposable	4PK	5
E0-0027	Infiltration Pump	Ability to inject local anesthesia for lyposuction techniques - Compatible with surgical cart - Compatible with lyposuction pumps - Compatible with electric power console	MicroAire/ASP-TUM	Installation/assembly required	Non-Disposable	EA	2
E0-0028	Implant	Bioabsorbable - capable of use in plastic surgery - Three 3.0mm tines - ultra thin platform - Capable of absorption rate from 9-12 months	MicroAire/CFD-010-0167		Disposable	EA	5
E0-0029	Implant	Capable of plastic surgery in Midface - 4.5mm bone - bioabsorbable - capable absorption rate of approximately 9-12 months	MicroAire/CFD-020-0197		Disposable	EA	5
E0-0030	Implant	Capable of Midface plastic surgery - 4.5mm - Absorption rate ability from 9-12 months	MicroAire/CFD-060-0197		Disposable	EA	5
E0-0031	Implant	Used for brow bone surgery - 3.5mm - bioabsorbable - includes disposable drill bit and insertion tool	MicroAire/CFD-080-0167		Disposable	EA	5
E0-0032	Implant	Capable of plastic surgery for the brow - 3.0mm - bioabsorbable - includes disposable drill bit and insertion tool	MicroAire/MicroAire/CFD-080-0177		Disposable	EA	5
E0-0033	Ribbon	Bioabsorbable fixation device - 2.5mm tine height - Capable of 90 day absorption rate	MicroAire/CFD-13005		Disposable	5PK	5
E0-0034	Surgical Wire	Consists of two needles - connected by a three-foot section of fine-gauge- surgical-grade steel -ability to twist	MicroAire/CFD-19005		Disposable	5PK	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0035	Implant	Bioabsorbable brow fixation device - 2 Pack - 3.5mm tine height	MicroAire/CFD-22102		Disposable	2PK	5
E0-0036	Implant	Bioabsorbable brow fixation device - 2 Pack - 3.0mm tine height	MicroAire/CFD-22202		Disposable	2PK	5
E0-0037	Implant	Bioabsorbable implants - Three 3.5mm tines and an ultra thin platform - Absorption rate ability approximately 90 days	MicroAire/CFD-23102		Disposable	2PK	5
E0-0038	Implant	Bioabsorbable- Three 3.0mm tines - ultra thin platform - absorption rate ability approximately 90 days	MicroAire/CFD-23202		Disposable	2PK	5
E0-0039	INSTRUMENT	insertion instrument - used in brow lift procedures	MicroAire/CFD-010-4305		Disposable	EA	5
E0-0040	Drill bit	drill bit with stopper	MicroAire/CFD-010-4313		Disposable	EA	5
E0-0041	INSTRUMENT KIT	Includes insertion tool- drill bits and sterilization tray with lid for use with bioabsorbable forehead implants	MicroAire/CFD-010-4320		Disposable	EA	5
E0-0042	INSTRUMENT TRAY	Tray for forehead bioabsorbable implants	MicroAire/CFD-010-4325		Disposable	EA	5
E0-0043	- Tap TOOL	Midface Bone - Tap Tool for use with bioabsorbable implants	MicroAire/CFD-020-4200		Disposable	EA	5
E0-0044	Installation tool	Midface bone installation Tool for use with bioabsorbable implants	MicroAire/CFD-020-4250		Disposable	EA	5
E0-0045	DRILLBIT	1.25mm	MicroAire/CFD-020-4302		Disposable	EA	5
E0-0046	Instrument kit	Contains Forehead insertion tool- drill bits- sterilization tray with lid and compatible with bioabsorbable implants	MicroAire/CFD-020-4320		Disposable	EA	5
E0-0047	Installation tool	Midface bone installation Tool for use with bioabsorbable implants	MicroAire/CFD-020-4325		Disposable	EA	5
E0-0048	Clipper	Screw Clipper used in midface plastic surgery	MicroAire/CFD-020-4400		Disposable	EA	5
E0-0049	Hand drill	manual- ability to implant bioabsorbable implants	MicroAire/CFD-899-4000		Non-Disposable	EA	2
E0-0050	Blade	ergonomic handle- 35 degree v-shaped blade - ability to remove strips of epithelial tissue - Ability to control the depth of the strips	MicroAire/PL-4006		Disposable	6PK	5
E0-0051	Blade	ergonomic handle- 55 degree v-shaped blade - Ability to remove strips of epithelial tissue - Ability control the depth of the strips	MicroAire/PL-4556		Disposable	6PK	5
E0-0052	ELECTRIC CONSOLE	Power source for handpiece - ability to use dual handpiece cable - 110V	MicroAire/5020	Installation/assembly required	Non-Disposable	EA	2

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0053	Cable	Electric cable - ability to connect power source to Small power handpiece power equipment	MicroAire/1006-PALE	Installation/assembly required	Non-Disposable	EA	2
E0-0054	Cable	Electric cable - ability to connect to power counsel and Large hand pieces	MicroAire/5006-PAL	Installation/assembly required	Non-Disposable	EA	2
E0-0055	Cap	Ability to cover connections on instrument cable during cleaning	MicroAire/CAP-600E		Non-Disposable	EA	20
E0-0056	TUBING	Aspiration tubing- 12' -Sterile	MicroAire/PAL-1200		Disposable	5PK	5
E0-0057	Cannula	2.4mm tube 7.5mm Port - 15CM Long	MicroAire/PAL-240LM		Disposable	EA	5
E0-0058	Cannula	2.4mm tube 7.5mm Port - 22CM Long	MicroAire/PAL-240LS		Disposable	EA	5
E0-0059	Cannula	2.4mm tube 5mm Port - 22CM Long	MicroAire/PAL-240SS		Disposable	EA	5
E0-0060	Cannula	2.4mm tube 7.5mm Port - 15CM Long	MicroAire/PAL-244LM		Disposable	EA	5
E0-0061	Cannula	2.4mm tube 3mm Port - 15CM Long	MicroAire/PAL-24S1		Disposable	EA	5
E0-0062	Cannula	3.0mm tube 8mm Port - 30CM Long	MicroAire/PAL-300LL		Disposable	EA	5
E0-0063	CANNULA	3.0mm- 8mm Port - 22CM Long	MicroAire/PAL-300LS		Disposable	EA	5
E0-0064	Cannula	3.0mm- 6mm Port - 22CM Long	MicroAire/PAL-300SS		Disposable	EA	5
E0-0065	Cannula	3.0mm- double Mercedes- 30CM Long	MicroAire/PAL-3010LL		Disposable	EA	5
E0-0066	Cannula	3.0mm- double Mercedes- 22CM Long	MicroAire/PAL-3010LS		Disposable	EA	5
E0-0067	Cannula	3.0mm X 22CM Long 4.1mm Port	MicroAire/PAL-3011LS		Disposable	EA	5
E0-0068	Cannula	3.0mm X 22CM Long 5.8mm Port	MicroAire/PAL-3012LS		Disposable	EA	5
E0-0069	Cannula	3.0mm tube 10mm Port - 30CM Long	MicroAire/PAL-302LL		Disposable	EA	5
E0-0070	Cannula	3.0mm tube 8mm Port - 22CM Long	MicroAire/PAL-302LS		Disposable	EA	5
E0-0071	Cannula	3.0mm tube 8mm Port - 30CM Long	MicroAire/PAL-303LL		Disposable	EA	5
E0-0072	Cannula	3.0mm tube 8mm Port - 22CM Long	MicroAire/PAL-303LS		Disposable	EA	5
E0-0073	Cannula	3.0mm tube 8mm Port - 30CM Long	MicroAire/PAL-304LL		Disposable	EA	5
E0-0074	Cannula	3.0mm bent tube	MicroAire/PAL-304LLB		Disposable	EA	5
E0-0075	Cannula	3.0mm tube 8mm Port - 15CM Long	MicroAire/PAL-304LM		Disposable	EA	5
E0-0076	Cannula	3.0mm tube 8mm Port - 22CM Long	MicroAire/PAL-304LS		Disposable	EA	5
E0-0077	Cannula	3.0mm flared tube 30CM Long	MicroAire/PAL-308LL		Disposable	EA	5
E0-0078	Cannula	3.0mm flared tube 22CM Long	MicroAire/PAL-308LS		Disposable	EA	5
E0-0079	Cannula	3.0mm TriPort tube 15CM Long	MicroAire/PAL-309LM		Disposable	EA	5
E0-0080	Cannula	3.0mm TriPort tube 22CM Long	MicroAire/PAL-309LS		Disposable	EA	5
E0-0081	Cannula	3.0mm multiPort tube 30CM Long	MicroAire/PAL-30MLL		Disposable	EA	5
E0-0082	Cannula	3.0mm multiPort tube 15CM Long	MicroAire/PAL-30MLM		Disposable	EA	5
E0-0083	Cannula	3.0mm multiPort tube 22CM Long	MicroAire/PAL-30MLS		Disposable	EA	5
E0-0084	Cannula	3.0mm tube 1.5mm Port -15CM Long	MicroAire/PAL-30S2		Disposable	EA	5
E0-0085	Cannula	3.0mm tube 1.5mm Port -15CM Long	MicroAire/PAL-30S3		Disposable	EA	5
E0-0086	Cannula	3.6mm offset double tube 30CM Long	MicroAire/PAL-3610SL		Disposable	EA	5
E0-0087	Cannula	3.5mm TriPort tube 30CM Long	MicroAire/PAL-369LL		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0088	Cannula	3.5mm TriPort tube 15CM Long	MicroAire/PAL-369LM		Disposable	EA	5
E0-0089	Cannula	3.5mm TriPort tube 22CM Long	MicroAire/PAL-369LS		Disposable	EA	5
E0-0090	Cannula	4.0mm offset double tube 30CM Long	MicroAire/PAL-4010SL		Disposable	EA	5
E0-0091	Cannula	4mm X 30CM- 5.6mm Port tube	MicroAire/PAL-4011LL		Disposable	EA	5
E0-0092	Cannula	4mm X 30CM- 6.4mm Port tube	MicroAire/PAL-4012LL		Disposable	EA	5
E0-0093	Cannula	4.0mm tube 10mm Port - 30CM Long	MicroAire/PAL-402LL		Disposable	EA	5
E0-0094	Cannula	4.0mm tube 10mm Port - 22CM Long	MicroAire/PAL-402LS		Disposable	EA	5
E0-0095	Cannula	4.0mm tube 10mm Port - 30CM Long	MicroAire/PAL-403LL		Disposable	EA	5
E0-0096	Cannula	4.0mm tube 10mm Port - 22CM Long	MicroAire/PAL-403LS		Disposable	EA	5
E0-0097	Cannula	4.0mm tube 10mm Port - 30CM Long	MicroAire/PAL-404LL		Disposable	EA	5
E0-0098	Cannula	4.0mm bent tube	MicroAire/PAL-404LLB		Disposable	EA	5
E0-0099	Cannula	4.0mm turbo tube - 22CM Long Mercedes tip	MicroAire/PAL-404LLT		Disposable	EA	5
E0-0100	Cannula	4.0mm tube 10mm Port - 22CM Long	MicroAire/PAL-404LS		Disposable	EA	5
E0-0101	Cannula	4.0mm turbo tube - 22CM Long Mercedes tip	MicroAire/PAL-404LST		Disposable	EA	5
E0-0102	Cannula	4.0mm tube 10mm Port - 30CM Long	MicroAire/PAL-405LL		Disposable	EA	5
E0-0103	Cannula	4.0mm turbo tube - 30CM Long double Mercedes tip	MicroAire/PAL-405LLT		Disposable	EA	5
E0-0104	Cannula	4.0mm tube 10mm Port - 22CM Long	MicroAire/PAL-405LS		Disposable	EA	5
E0-0105	Cannula	4.0mm tube 10mm Port - 30CM Long	MicroAire/PAL-406LL		Disposable	EA	5
E0-0106	Cannula	4.0mm tube 10mm Port - 22CM Long	MicroAire/PAL-406LS		Disposable	EA	5
E0-0107	Cannula	4.0mm tube 10mm Port - 30CM Long	MicroAire/PAL-407LL		Disposable	EA	5
E0-0108	Cannula	4.0mm helixed tube	MicroAire/PAL-407LS		Disposable	EA	5
E0-0109	Cannula	4.0mm tube 30CM Long	MicroAire/PAL-408LL		Disposable	EA	5
E0-0110	Cannula	4.0mm bent tube	MicroAire/PAL-408LLB		Disposable	EA	5
E0-0111	Cannula	4.0mm flared tube 22CM Long	MicroAire/PAL-408LS		Disposable	EA	5
E0-0112	STERILIZATION CASE	forlyposuction hand pieces -ability to holds cable -ability to hold up to 10 cannulas - ability for sterilization	MicroAire/PAL-500		Non-Disposable	EA	5
E0-0113	CANNULA	5.0mm tube 10mm Port - 30CM Long	MicroAire/PAL-504LL		Disposable	EA	5
E0-0114	CANNULA	5.0mm tube 10mm Port - 22CM Long	MicroAire/PAL-504LS		Disposable	EA	5
E0-0115	CANNULA	5.0mm tube 10mm slot- 30CM Long	MicroAire/PAL-505LL		Disposable	EA	5
E0-0116	CANNULA	5.0mm turbo tube 30CM Long double Mercedes tip	MicroAire/PAL-505LLT		Disposable	EA	5
E0-0117	CANNULA	5.0mm tube 10mm Port - 22CM Long	MicroAire/PAL-505LS		Disposable	EA	5
E0-0118	CANNULA	5.0mm tube 10mm Port -30CM Long	MicroAire/PAL-506LL		Disposable	EA	5
E0-0119	CANNULA	5.0mm tube 10mm Port - 22CM Long	MicroAire/PAL-506LS		Disposable	EA	5
E0-0120	CANNULA	5.0mm tube 10mm Slot- 22CM Long	MicroAire/PAL-507LL		Disposable	EA	5
E0-0121	CANNULA	5.0mm helixed tri-Port tube 30CM Long	MicroAire/PAL-507LLT		Disposable	EA	5
E0-0122	CANNULA	5.0mm tube 10mm Port - 22CM Long	MicroAire/PAL-507LS		Disposable	EA	5
E0-0123	CANNULA	5.0mm diameterMETER- 7.0mm diameterMETER AT FLARE- 30CM USABLE LENGTH- 10mm Port LENGTH tube	MicroAire/PAL-508LL		Disposable	EA	5
E0-0124	CANNULA	5.0mm bent tube	MicroAire/PAL-508LLB		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0125	CANNULA	5.0mm diameterMETER- 7.0mm diameterMETER AT FLARE- 22CM USABLE LENGTH- 10mm Port LENGTH tube	MicroAire/PAL-508LS		Disposable	EA	5
E0-0126	LOCK HUB	Ability to use cannulae tubing from multiple manufacturers	MicroAire/PAL-700		Non-Disposable	EA	5
E0-0127	RELEASE KEY	Abiliyt to release tube from hand piece device	MicroAire/PAL-710		Non-Disposable	EA	5
E0-0128	Wand	Manual hand piece - ability to use in delicate plastic procedures	MicroAire/PAL-730		Non-Disposable	EA	5
E0-0129	TUBING	Aspiration Tubing - 9' - Sterile	MicroAire/PAL-900		Disposable	6PK	5
E0-0130	CANNULA	2.4mm diameterMETER- 15CM USABLE LENGTH- 7.5mm Port LENGTH tube	MicroAire/PAL-R240LM		Disposable	EA	5
E0-0131	CANNULA	2.4mm diameterMETER- 22CM USABLE LENGTH- 7.5mm Port LENGTH tube	MicroAire/PAL-R240LS		Disposable	EA	5
E0-0132	CANNULA	2.4mm diameterMETER- 22CM USABLE LENGTH- 5.0mm Port LENGTH tube	MicroAire/PAL-R240SS		Disposable	EA	5
E0-0133	CANNULA	2.4mm diameterMETER- 15CM USABLE LENGTH- 7.5mm Port LENGTH tube	MicroAire/PAL-R244LM		Disposable	EA	5
E0-0134	CANNULA	2.4mm diameterMETER- 15CM USABLE LENGTH- 1 Port 3.0mm Port LENGTH tube	MicroAire/PAL-R24S1		Disposable	EA	5
E0-0135	CANNULA	3.0mm diametermeter 30CM usable length- 6.0mm Port length	MicroAire/PAL-R300LL		Disposable	EA	5
E0-0136	CANNULA	3.0mm diametermeter 22CM usable length- 6.0mm Port length	MicroAire/PAL-R300LS		Disposable	EA	5
E0-0137	CANNULA	3.0mm diametermeter 22CM usable length- 6.0mm Port length	MicroAire/PAL-R300SS		Disposable	EA	5
E0-0138	CANNULA	3.0mm diametermeter 30CM usable length- 10mm Port length	MicroAire/PAL-R3010LL		Disposable	EA	5
E0-0139	CANNULA	3.0mm diametermeter 22CM usable length- 10mm Port length	MicroAire/PAL-R3010LS		Disposable	EA	5
E0-0140	CANNULA	3.0mm diametermeter 30CM usable length- 8.0mm Port length	MicroAire/PAL-R302LL		Disposable	EA	5
E0-0141	CANNULA	3.0mm diametermeter 22CM usable length- 8.0mm Port length	MicroAire/PAL-R302LS		Disposable	EA	5
E0-0142	CANNULA	3.0mm diametermeter 30CM usable length- 8.0mm Port length	MicroAire/PAL-R303LL		Disposable	EA	5
E0-0143	CANNULA	3.0mm diametermeter 22CM usable length- 8.0mm Port length	MicroAire/PAL-R303LS		Disposable	EA	5
E0-0144	CANNULA	3.0mm diametermeter 22CM usable length- 8.0mm Port length	MicroAire/PAL-R304LL		Disposable	EA	5
E0-0145	CANNULA	3.0mm bent tube	MicroAire/PAL-R304LLB		Disposable	EA	5
E0-0146	CANNULA	3.0mm diametermeter 15CM usable length- 8.0mm Port length	MicroAire/PAL-R304LM		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0147	CANNULA	3.0mm diametermeter 22CM usable length- 8.0mm Port length	MicroAire/PAL-R304LS		Disposable	EA	5
E0-0148	CANNULA	3.0mm diametermeter 4.2 diametermeter at flare 30CM usable length- 9.7mm Port length	MicroAire/PAL-R308LL		Disposable	EA	5
E0-0149	CANNULA	3.0mm diametermeter 4.2 diametermeter at flare 22CM usable length- 9.7mm Port length	MicroAire/PAL-R308LS		Disposable	EA	5
E0-0150	CANNULA	3.0mm diametermeter 22CM usable length- 1.6mm Port length	MicroAire/PAL-R309LM		Disposable	EA	5
E0-0151	CANNULA	3.0mm diametermeter 22CM usable length- 1.6mm Port length	MicroAire/PAL-R309LS		Disposable	EA	5
E0-0152	CANNULA	3.0mm diametermeter 30CM usable length- 1.9mm Port length	MicroAire/PAL-R30MLL		Disposable	EA	5
E0-0153	CANNULA	3.0mm diametermeter 15CM usable length- 1.9mm Port length	MicroAire/PAL-R30MLM		Disposable	EA	5
E0-0154	CANNULA	3.0mm diametermeter 22CM usable length- 1.9mm Port length	MicroAire/PAL-R30MLS		Disposable	EA	5
E0-0155	CANNULA	3.0mm diametermeter 15CM usable length- 2 Port s- 1.5mm Port length	MicroAire/PAL-R30S2		Disposable	EA	5
E0-0156	CANNULA	3.0mm diametermeter 15CM usable length- 3 Port s- 1.5mm Port length	MicroAire/PAL-R30S3		Disposable	EA	5
E0-0157	CANNULA	3.6mm diametermeter 30CM usable length- 4.5mm and 3.5mm Port lengths	MicroAire/PAL-R3610SL		Disposable	EA	5
E0-0158	Cannula	3.6mm diametermeter 30CM usable length- 4.1mm and 2.0mm Port lengths	MicroAire/PAL-R369LL		Disposable	EA	5
E0-0159	Cannula	3.6mm diametermeter 15CM usable length- 4.1mm and 2.0mm Port lengths	MicroAire/PAL-R369LM		Disposable	EA	5
E0-0160	Cannula	3.6mm diametermeter 22CM usable length- 4.1mm and 2.0mm Port lengths	MicroAire/PAL-R369LS		Disposable	EA	5
E0-0161	Cannula	4.0mm diametermeter 30CM usable length- 5.1mm and 4.2mm Port lengths	MicroAire/PAL-R4010SL		Disposable	EA	5
E0-0162	Cannula	4.0mm X 40CM- 5.6mm Port Tube	MicroAire/PAL-R4011XL		Disposable	EA	5
E0-0163	Cannula	4.0mm X 40CM- 6.4mm Port Tube	MicroAire/PAL-R4012XL		Disposable	EA	5
E0-0164	Cannula	4.0mm diametermeter 30CM usable length- 10mm Port tube lengths	MicroAire/PAL-R402LL		Disposable	EA	5
E0-0165	Cannula	4.0mm diametermeter 22CM usable length- 10mm Port tube lengths	MicroAire/PAL-R402LS		Disposable	EA	5
E0-0166	Cannula	4.0mm diametermeter 30CM usable length- 10mm Port tube lengths	MicroAire/PAL-R403LL		Disposable	EA	5
E0-0167	Cannula	4.0mm diametermeter 22CM usable length- 10mm Port tube lengths	MicroAire/PAL-R403LS		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0168	Cannula	4.0mm diametermeter 30CM usable length- 10mm Port tube lengths	MicroAire/PAL-R404LL		Disposable	EA	5
E0-0169	Cannula	4.0mm BENT tube lengths	MicroAire/PAL-R404LLB		Disposable	EA	5
E0-0170	Cannula	4.0mm diametermeter 22CM usable length- 10mm Port tube lengths	MicroAire/PAL-R404LS		Disposable	EA	5
E0-0171	Cannula	4.0mm diametermeter 40CM usable length- 10mm Port tube lengths	MicroAire/PAL-R404XL		Disposable	EA	5
E0-0172	Cannula	4.0mm diametermeter 30CM usable length- 10mm Port tube lengths	MicroAire/PAL-R405LL		Disposable	EA	5
E0-0173	Cannula	4.0mm diametermeter 22CM usable length- 10mm Port tube lengths	MicroAire/PAL-R405LS		Disposable	EA	5
E0-0174	Cannula	4.0mm diametermeter 30CM usable length- 10mm Port tube lengths	MicroAire/PAL-R406LL		Disposable	EA	5
E0-0175	Cannula	4.0mm diametermeter 22CM usable length- 10mm Port tube lengths	MicroAire/PAL-R406LS		Disposable	EA	5
E0-0176	Cannula	4.0mm diametermeter 30CM usable length- 10mm Port tube lengths	MicroAire/PAL-R407LL		Disposable	EA	5
E0-0177	Cannula	4.0mm diametermeter 22CM usable length- 10mm Port tube lengths	MicroAire/PAL-R407LS		Disposable	EA	5
E0-0178	Cannula	4.0mm diameterMETER- 5.6mm diameterMETER AT FLARE- 30CM USABLE LENGTH- 9.7mm Port LENGTH tube lengths	MicroAire/PAL-R408LL		Disposable	EA	5
E0-0179	Cannula	4.0mm BENT tube lengths	MicroAire/PAL-R408LLB		Disposable	EA	5
E0-0180	Cannula	4.0mm diameterMETER- 22CM USABLE LENGTH- 9.7mm Port LENGTH tube lengths	MicroAire/PAL-R408LS		Disposable	EA	5
E0-0181	Cannula	5.0mm diameterMETER- 30CM USABLE LENGTH- 10mm Port LENGTH tube lengths	MicroAire/PAL-R504LL		Disposable	EA	5
E0-0182	Cannula	5.0mm diameterMETER- 22CM USABLE LENGTH- 10mm Port LENGTH tube lengths	MicroAire/PAL-R504LS		Disposable	EA	5
E0-0183	Cannula	5.0mm diameterMETER- 30CM USABLE LENGTH- 10mm Port LENGTH tube lengths	MicroAire/PAL-R505LL		Disposable	EA	5
E0-0184	Cannula	5.0mm diameterMETER- 22CM USABLE LENGTH- 10mm Port LENGTH tube lengths	MicroAire/PAL-R505LS		Disposable	EA	5
E0-0185	Cannula	5.0mm diameterMETER- 30CM USABLE LENGTH- 10mm Port LENGTH tube lengths	MicroAire/PAL-R506LL		Disposable	EA	5
E0-0186	Cannula	5.0mm diameterMETER- 22CM USABLE LENGTH- 10mm Port LENGTH tube lengths	MicroAire/PAL-R506LS		Disposable	EA	5
E0-0187	Cannula	5.0mm diameterMETER- 30CM USABLE LENGTH- 10mm Port LENGTH tube lengths	MicroAire/PAL-R507LL		Disposable	EA	5
E0-0188	Cannula	5.0mm diameterMETER- 22CM USABLE LENGTH- 10mm Port LENGTH tube lengths	MicroAire/PAL-R507LS		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0189	Cannula	5.0mm diameterMETER- 7.0mm diameterMETER AT FLARE- 30CM USABLE LENGTH- 9.7mm Port LENGTH tube lengths	MicroAire/PAL-R508LL		Disposable	EA	5
E0-0190	Cannula	5.0mm BENT tube lengths	MicroAire/PAL-R508LLB		Disposable	EA	5
E0-0191	Cannula	5.0mm diameterMETER- 7.0mm diameterMETER AT FLARE- 22CM USABLE LENGTH- 9.7mm Port LENGTH tube lengths	MicroAire/PAL-R508LS		Disposable	EA	5
E0-0192	Arthroscope handpiece	2.9mm 0 degrees - lightweight aluminum handpiece - pistol-grip style Compatable with endoscopic carpal tunnel release system - Ability to slides into position from the back - Ability to click-locks into place - Ability for the dual-function latch is lifted to insert and hold the blade assembly - Ability to twist lock the rotational aspect of the blade and scope - Ability to rotate left or right the blade assembly within the handpiece - Ability to connect a 3.0mm eyepiece endoscope connects to a Standard camera coupler.	MicroAire/81014		Non-Disposable	EA	2
E0-0193	ENDOSCOPE	3.0mm Eyepiece - Ability to fit into the arthroscope handpiece - autoclavable	MicroAire/81025		Non-Disposable	EA	2
E0-0194	Dilator	Size Small- Capable of dilating the tunnel for carpal tunnel surgery	MicroAire/81026		Non-Disposable	EA	5
E0-0195	Dilator	Size Medium -Capable dilating the tunnel for carpal tunnel surgery	MicroAire/81027		Non-Disposable	EA	5
E0-0196	ELEVATOR	Ability to separate synovium from deep side of ligament	MicroAire/81029		Non-Disposable	EA	5
E0-0197	STERILIZATION TRAY	Made of medical-grade polymer - Includes silicone insert and stainless steel latches - Ability multiple sterilization cycles - Ability to hold handpieces- elvator- dilator- endoscope and adapters instruments - Ability to include space for surgical instrumentation and autoclavable camera & light cables - Tray dimensions: 24.89cm x 52.32cm x 5.08cm	MicroAire/81040		Non-Disposable	EA	15
E0-0198	SLEEVE	Scope Protective Sleeve (Non-Autoclavable)	MicroAire/81057		Non-Disposable	EA	5
E0-0199	Dilator	Coequal to Disposable Blade Assembly	MicroAire/81061		Non-Disposable	EA	5
E0-0200	Adapter	Compatable with Endoscope Light Post Adapter - Ability to coneect with Storz / Olympus style fitting for lenses	MicroAire/81151		Non-Disposable	EA	5
E0-0201	Adapter	Compatible with Endoscope Light Post Adapter abiity to connect with Wolf/ Dyonics style fitting for lenses	MicroAire/81152		Non-Disposable	EA	5
E0-0202	Blade	Compatible with Arthroscope handpiece - Ability to retract the blade - Include a window at the tip of the disposable blade - Ability to provide a direct view of the transverse carpal ligament - Ability to adjust for right- or left-handed surgeons - Single Pack	MicroAire/81010-1		Disposable	PG	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0203	Blade	Compatible with Arthroscope handpiece - Ability to retract the blade - Include a window at the tip of the disposable blade - Ability for a direct view of the transverse carpal ligament - The ability to adjust for right- or left-handed surgeons - 6 per Pack	MicroAire/81010-6		Disposable	6PK	5
E0-0204	Scope	Ability for minimally-invasive incision - Ability for direct connection between the endoscope and camera	MicroAire/81020N		Non-Disposable	PG	2
E0-0205	Handpiece	Drill/Reamer - Compatible with 14.4 volt batteries - Ability for high torque - Ability for quick connec reamers - Compatiable with drilling couplers- Wire and pin Driver - Ability to use 110V Electric Power	MicroAire/7505	Installation/assembly required	Non-Disposable	EA	8
E0-0206	Handpiece	Oscillating Saw - Ability for Large Power - Ability for high and low speed settings - Ability for 8 position rotating head - Ability for secure blade locking hub - compatible with 14.4 volt batteries - Ability to use 110V Electric Power	MicroAire/7506	Installation/assembly required	Non-Disposable	EA	8
E0-0207	Handpiece	Reciprocating Saw - Ability for twist lock blade mount - Ability for secure blade locking hub - compatible with 14.4 volt batteries - electric powered - Compatible with reciprocating saw blades - Ability to use 110V Electric Power	MicroAire/7507	Installation/assembly required	Non-Disposable	EA	8
E0-0208	Handpiece	Sternum Saw Compatible with 14.4v batteries - Ability of keyless blade hub - Ability to attach sternum guard - Ability to use for mediametern sternotomies	MicroAire/7800	Installation/assembly required	Non-Disposable	EA	8
E0-0209	KEY	T handle design- 1/4 inch	MicroAire/4100-030		Non-Disposable	EA	5
E0-0210	Air Hose	Lenth:10' -Ability to vent at the diffuser end - Double-layer vulcanized rubber supPort ed with nylon mesh - Abilty to reduce instrument noise up to 40% - Connection is compatible with a scharder style connector	MicroAire/7000-000	Installation/assembly required	Non-Disposable	EA	2
E0-0211	Air Hose	Lenth: 20' - ability to vent at the diffuser end - Double-layer vulcanized rubber supPort ed with nylon mesh - Ability to reduce instrument noise up to 40%- Connection compatible with a scharder style connector	MicroAire/7000-020	Installation/assembly required	Non-Disposable	EA	2
E0-0212	Power Cable	Electric - Compatible with electric powered handpieces	MicroAire/7007-000	Installation/assembly required	Non-Disposable	EA	2
E0-0213	Air Hose	10' pneumatic hose - Connector compatible with Hall style connector to a scharder style connector	MicroAire/7013-000	Installation/assembly required	Non-Disposable	EA	2
E0-0214	Air Hose	20' pneumatic hose - Connectin compatible withh Hall style connector	MicroAire/7013-020	Installation/assembly required	Non-Disposable	EA	2

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0215	Coupler	At least At least 250 rpm - Compatible with Large bone drill reamers - Ability to connect with Hall style connector	MicroAire/7100-001		Disposable	EA	5
E0-0216	Coupler	At least At least 750 rpm - Compatible with Large bone Drill Reamers	MicroAire/7100-002		Disposable	EA	5
E0-0217	Coupler	1/4" Coupler w/ key - At least 750rpm - Compatible with Large bone Drill Reamers	MicroAire/7100-003		Disposable	EA	5
E0-0218	Coupler	1/4" Keyless Coupler- At least 750 rpm- use with Large bone drill Reamers. NOT recommended for high torque applications running in reverse	MicroAire/7100-004		Disposable	EA	5
E0-0219	Coupler	1/8" Keyless Coupler- At least At least 750 rpm - Compatible use with Large bone handpieces	MicroAire/7100-005		Disposable	EA	5
E0-0220	Coupler	Ability for Quick Connection - at least 750rpm - use with Large bone handpieces - Compatable with synthesis style twist drills and adapters	MicroAire/7100-006		Disposable	EA	5
E0-0221	Coupler	At least At least 250 rpm - Compatible with Large bone handpieces device	MicroAire/7100-008		Disposable	EA	5
E0-0222	Coupler	1/4" Coupler w/ key - At least At least 250 rpm - Compatable with Large bone handpieces	MicroAire/7100-010		Disposable	EA	5
E0-0223	Coupler	At least 250 rpm - Compatible with Large bone handpiece instruments	MicroAire/7100-011		Disposable	EA	5
E0-0224	PIN DRIVER	Compatible pins from 1.6mm to 2.8mm.	MicroAire/7100-045		Disposable	EA	5
E0-0225	PIN DRIVER	Compatible pins from 3.2mm to 4.0mm	MicroAire/7100-050		Disposable	EA	5
E0-0226	Coupler	At least At least 750 rpm - Compatable with Large bone handpiece reamers	MicroAire/7100-051		Disposable	EA	5
E0-0227	Coupler	At least 750 rpm- Compatible with Large bone Drill Reamers	MicroAire/7100-058		Disposable	EA	5
E0-0228	Coupler	At least 250 rpm- Compatible with Large bone Drill Reamers	MicroAire/7100-111		Disposable	EA	5
E0-0229	Coupler	Non-Reversing - Ability for spring style IM flexible reamers - Ability to prevent binding - Compatable with Large bone hand pieces.	MicroAire/7100-112		Disposable	EA	5
E0-0230	Rechargeable battery housing	Contains rechargeable batteries - Sterile housing - auoclavable	MicroAire/7500-615	Installation/a ssembly required	Disposable	EA	8
E0-0231	Rechargeable Battery	14.4 volt - Non-Sterile - Compatible with battery housing	MicroAire/7500-620	Installation/a ssembly required	Non-Disposable	EA	8
E0-0232	Charger	single unit - Compatible with aseptic battery pack	MicroAire/7500-625	Installation/a ssembly required	Non-Disposable	EA	8

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0233	Shield	Sterile - Compatible with aseptic battery charger	MicroAire/7500-630	installation/assembly required	Non-Disposable	EA	8
E0-0234	Battery Charger	Ability to charge 14.4-volt Nickel Metal Hydride (NiMH) rechargeable-autoclavable- or aseptic-transfer batteries - Ability to charge up to four batteries at one time - Ability to charge different size batteries at the same time - Ability for no memory-based charging -Ability to charge in 90 minutes or less	MicroAire/7500-700	Installation/assembly required	Non-Disposable	EA	8
E0-0235	STERILIZATION TRAY	Compatible with automated washer disinfectors systems - Aluminum outer case and lid - Nylon coated brackets - Inner tray is made from medical-grade stainless steel -Ability to use Inner tray in the Sterile field - Inner tray is Capable of holding three hand pieces- fourteen drive couplers- and two instruments - Dimensions: 4.28cm x 38.50cm x 19.86cm	MicroAire/7505-040	Installation/assembly required	Non-Disposable	EA	5
E0-0236	Coupler	Ability to accept Wires/pins 1.6mm - 2.8mm.	MicroAire/7505-045		Disposable	EA	5
E0-0237	Coupler	Ability to accept pins from 3.0mm - 4.0mm	MicroAire/7505-050		Disposable	EA	5
E0-0238	Battery Charger	14.4 volts - Compatible with Large bone handpiece instruments	MicroAire/7505-710	Installation/assembly required	Non-Disposable	EA	8
E0-0239	Guard	Stainless Steel - Compatible with Sternum saw blade	MicroAire/7800-002		Disposable	EA	5
E0-0240	Nitrogen Regulator	Single Out-put- Ability to connect with schrader style connection	MicroAire/9500-000	Installation/assembly required	Non-Disposable	EA	2
E0-0241	Adapter	Provides the ability for a Stryker/Zimmer-3m- MicroAire style connector to connect to a Stryker style regulator	MicroAire/9500-001		Disposable	EA	5
E0-0242	Connector	High pressure - Ability to handle air pressures up to 15-000 psi - Ability to connect to Hall- 3M- or MicroAire style connections	MicroAire/9500-002		Disposable	EA	5
E0-0243	Adapter	Compatible with Air hose -Ability to connect with Hall/Stryker style connections	MicroAire/9500-003		Disposable	EA	5
E0-0244	Adapter	Compatible with Air Hose - Ability to connect Zimmer style connections to Ohio style connections	MicroAire/9500-004		Disposable	EA	5
E0-0245	Manifold	Compatible with Air Hose -Ability to connect 3 hoses - Ability to connect to Schrader style connectors (Hall-3M- or MicroAire) styles	MicroAire/9500-011	Installation/assembly required	Non-Disposable	EA	2
E0-0246	Manifold	Compatible with Air Hose - Ability to connect 3 hoses - ability to connect with Stryker style connections	MicroAire/9500-012	Installation/assembly required	Non-Disposable	EA	2

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0247	Manifold	Compatible with AIR HOSE - Ability to connect 3 hoses - Ability to connect 2 SCHRADER style CONNECTORS and 1 STRYKER style CONNECTOR	MicroAire/9500-013	Installation/assembly required	Non-Disposable	EA	2
E0-0248	Manifold	Compatible with AIR HOSE - Ability to connect to 3 hoses - Ability to connect 1 SCHRADER style CONNECTOR and 2 STRYKER style CONNECTORS	MicroAire/9500-014	Installation/assembly required	Non-Disposable	EA	2
E0-0249	Nitrogen Regulator	Dual Out put - Ability to connect to MicroAire- Hall- and 3M connection style air hoses -Ability to show Main tank and pressure output readings	MicroAire/9502-000	Installation/assembly required	Non-Disposable	EA	2
E0-0250	Brush	drill bit with brush for drill/reamer systems	MicroAire/9600-061		Disposable	EA	5
E0-0251	Brush	drill bit with burs and blades	MicroAire/9600-062		Disposable	EA	5
E0-0252	Brush	Drill bit with brush for Wire driver	MicroAire/9600-063		Disposable	EA	5
E0-0253	Brush	Brush For Cleaning Drill Collets	MicroAire/9600-064		Disposable	EA	5
E0-0254	BRUSH	FOR CLEANING - FLEXIBLE - REAMER CANNULATIONS -19" Long	MicroAire/9600-161		Disposable	EA	5
E0-0255	Cabinet	16 Drawers - Include divider s- Ability to hold burs and blade attachments	MicroAire/9700-050	Installation/assembly required	Non-Disposable	EA	5
E0-0256	STERILIZATION CASE	Ability to clean Large power tool battery systems - Ability to hold 11 drive couplers and pin driver and one power cable - autoclavable	MicroAire/9720-000	Installation/assembly required	Non-Disposable	EA	5
E0-0257	Oil	Capable of lubricating instruments	MicroAire/9916-000		Disposable	EA	5
E0-0258	Cap	WASHER DISINFECTOR CAP - Ability to connect to MICROAIRE STYLE CONNECTOR	MicroAire/CAP-100		Disposable	EA	5
E0-0259	Cap	WASHER DISINFECTOR CAP - Ability to connect to Hall STYLE CONNECTOR	MicroAire/CAP-200		Disposable	EA	5
E0-0260	Saw module	At least 30-000 CPM -6° arc - Compatible with Small bone power motor modules.	MicroAire/1922		Non-Disposable	EA	5
E0-0261	Drill module	At least 30-000 RPM - Ability at least 6.0 in./oz. of torque -Includes built-in Medium bur guard - Compatible with For use with Small bone motor modules	MicroAire/1930		Non-Disposable	EA	5
E0-0262	Drill module	20 Degree angled - Capable of at least 30-000 RPM - Capable of 6.0 in./oz. of torque - Includes built-in Medium bur guard - Compatible with Small bone motor modules-	MicroAire/1932		Non-Disposable	EA	5
E0-0263	Saw module	Reciprocating - Ability at least 18-000 CPM - Compatible with Small bone motor modules	MicroAire/1945		Non-Disposable	EA	5
E0-0264	Saw module	Wrenchless - Capable of at least 18-000 CPM - Ability to connect with hall style connections - Compatible with Small bone motor modules	MicroAire/1950		Non-Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0265	Saw module	Ability to use hex locking key - Ability at least 18-000 CPM - Ability to connect with Hall style connector- Compatible with Small bone 1000 motor modules	MicroAire/1955		Non-Disposable	EA	5
E0-0266	Saw module	Oscillating-Ability at least 30-000 CP - Compatible with Small bone motor modules	MicroAire/1970		Non-Disposable	EA	5
E0-0267	Saw module	Oscillating -Compatible for Foot Surgery - Ability at least 30-000 CPM - Compatible with Small bone motor modules	MicroAire/1972		Non-Disposable	EA	5
E0-0268	Saw module	Wrenchless - Ability at least 24-000 CPM - Compatible to connect with hall style connections - Compatible with Small bone motor modules	MicroAire/1976		Non-Disposable	EA	5
E0-0269	Drill module	Compatible with Jacob style connector with key - Capable of at least 1-670 RPM - Compatible with Small bone motor modules - Ability to accepts up to 3.5mm diameter twist drills	MicroAire/1980		Disposable	EA	5
E0-0270	Drill module	Compatible with Synthes Style connection - Ability at least 1-670 RPM - Compatible with Small bone motor modules	MicroAire/1990+F254		Disposable	EA	5
E0-0271	Wire driver	Ability to operate at speed of at least of 1-400 RPM - Automatic grip mechanism - Ability to hold Wire while being driven - Ability to automatically advance Wire when the instrument collar is pulled back- Ability for locking grip mechanism for Wire removal - Compatible with Non-- Threaded K-Wires with guide sized between .028mm and .062mm in diametermeter	MicroAire/1995		Disposable	EA	5
E0-0272	Power Console	Ability to operate two hand pieces simultaneously - Ability for user interface instrument control - Includdes intergrated irrigation pump - Compatible with Small power instruments- micro power instruments- Driver and Solo handpiece	MicroAire/5025	Installation/a ssembly required	Disposable	EA	2
E0-0273	FOOT PEDAL	Electric - Compatible with Small power Electric Instrument Control Console - Two pedals- Ability to provide variable speed control for forward or reverse rotation - Ability to provide irrigation control from the foot pedal	MicroAire/5401	Installation/a ssembly required	Non-Disposable	EA	2
E0-0274	Surgical Drill	Dual Trigger - Pistol-Grip handpiece - Compatible with couplers for drilling- pinning- cutting- and debridement - Compatible with Small power electric console	MicroAire/5641	Installation/a ssembly required	Non-Disposable	EA	2
E0-0275	Saw module	Sagittal- Ability at least 30-000 CPM - Ability at least 6° arc - Compatable with Small power electric console motor module	MicroAire/5922		Disposable	EA	5
E0-0276	Drill module	Capable of at least 30-000 RPM - Capable of at least 6.0 in./oz. of torque - Includes Built-in Medium bur guard	MicroAire/5930		Disposable	EA	5
E0-0277	Saw module	Reciprocating - Capable of at least 16-000 CPM	MicroAire/5945		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0278	Saw module	Sagittal - Keyless - Ability to connect with hall style connections- Capable of at least 18-000 CPM	MicroAire/5950		Disposable	EA	5
E0-0279	Saw module	Sagittal - Hex locking key - Compatible with hall style connection- Capable of at least 18-000 CPM	MicroAire/5955		Disposable	EA	5
E0-0280	Saw module	Oscillating -Ability for use in Foot Surgery - Capable of at least 30-000 CPM -Compatible with Small power handpiece	MicroAire/5972		Disposable	EA	5
E0-0281	Drill module	Jacob style connection with key -Capable of at least 1-670 RPM - Compatible with Small power hand piece	MicroAire/5980		Disposable	EA	5
E0-0282	Drill module	Compatible with Synthes Style connection- Capable of at least 1-670 RPM - Compatible with Small power hand piece	MicroAire/5990		Disposable	EA	5
E0-0283	Motor Module	Ability to switch heads during a procedure without disconnecting the hose or cable - Includes hand throttle - Pneumatic powered- Ability to connect to MicroAire style connectors	MicroAire/1000-100		Disposable	EA	4
E0-0284	Motor Module	Ability to switch heads during a procedure without disconnecting the hose or cable- Include hand throttle - Pneumatic powered -Ability to connect to Linvatec and Hall style connectors	MicroAire/1000-200		Disposable	EA	4
E0-0285	MOTOR MODULE	Ability to switch heads during a procedure without disconnecting the hose or cable - Includes hand throttle -Electric powered -Ability to connect to MicroAire style connectors	MicroAire/1000ET		Non-Disposable	EA	4
E0-0286	Cable	Connecting Cable from Instruments to power console	MicroAire/1006-1000		Non-Disposable	EA	2
E0-0287	Cable	Electric Connecting cable from instruments to power console	MicroAire/1006-6000		Non-Disposable	EA	2
E0-0288	Bur Guard	Small surgical drill - Ability to attache to drill module - Size Long	MicroAire/1100-005		Disposable	EA	5
E0-0289	Bur Guard	Small surgical drill - Ability to attache to drill module- Size extra Long	MicroAire/1100-006		Disposable	EA	5
E0-0290	Nose Shield	Ability to prevent blood and liquids from entering drill - Ability to protect against heat	MicroAire/1100-008		Disposable	10PK	5
E0-0291	Wire Driver	Pencil style -Ability to accepts K-Wires .028" - .062" with Wire guard - Cannulated - Ability to use with multiple Wire lengths up to 9" - Compatible to connect with MicroAire style Air Hose Connector	MicroAire/1625-100		Non-Disposable	EA	5
E0-0292	Wire driver	Pencil style -Ability to use K-Wires .028" - .062" with Wire Guide - Cannulated - Ability to use with multiple Wire lengths up to 9" - Compatible to connect with Hall style Connector	MicroAire/1625-200		Non-Disposable	EA	5
E0-0293	Jacob KEY	Jacob Key Style- Size 5/32 -Compatible with 1645-003/2180/6680 Couplers	MicroAire/1645-004		Non-Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0294	LOCKING TOOL	Wrench - Ability to lock saw attachments Oscillating Saw attachments	MicroAire/1745-001		Non-Disposable	EA	5
E0-0295	Drill module	Capable of at least 90-000 RPM - Included built-in Medium bur guard	MicroAire/2130-000		Non-Disposable	EA	5
E0-0296	LOCKING TOOL	Hex wrench - Ability to Blade Lock Compatible with HALL-STYLE connections - Compatible with SAGITTAL SAW	MicroAire/2250-001		Non-Disposable	EA	5
E0-0297	Impactor	Pneumatic- Ability to use for General Surgery - Compatabile with Micro Hose Connector	MicroAire/2500-100		Non-Disposable	EA	5
E0-0298	Impactor	Pneumatic - Ability to use in General Surgery - Compatible with Hall Style Hose Connector	MicroAire/2500-200		Non-Disposable	EA	5
E0-0299	Blade	Chisel-like knife- Ability to cut bone - Straight- 50 X10mm Ability to use with Impactor	MicroAire/2500-203		Disposable	EA	5
E0-0300	Blade	Chisel-like knife - Ability to cut bone - Straight - 80 X6mm - Compatible with Impactor	MicroAire/2500-205		Disposable	EA	5
E0-0301	Blade	Chisel-like knife - Ability to cut bone - Straight - 80 X10mm-Compatible with Impactor	MicroAire/2500-206		Disposable	EA	5
E0-0302	Blade	Chisel-like knife - Ability to cut bone - Straight - 70 X15mm - Compatible with Impactor	MicroAire/2500-209		Disposable	EA	5
E0-0303	Blade	Chisel-like knife - Ability to cut bone - Curved - 80 X10mm- Compatible with Impactor	MicroAire/2500-216		Disposable	EA	5
E0-0304	Blade	Chisel with a concave blade- Curved - 80 X 10mm- Compatible with Impact	MicroAire/2500-230		Disposable	EA	5
E0-0305	Blade	Chisel-like knife - Ability to cut bone - Straight - 80 X6mm - Compatible with Impactor -Ability to use in the Notch Plasty procedure	MicroAire/2500-316		Disposable	EA	5
E0-0306	STERILIZATION CASE	Molded Polymer - For modular instruments	MicroAire/2651-000		Non-Disposable	EA	5
E0-0307	Bur Guard	Extended surgical drill -Compatible with high speed drill module	MicroAire/2910-011		Disposable	EA	5
E0-0308	Bur Guard	Medium surgical drill- Compatible with high speed drill module	MicroAire/2910-012		Disposable	EA	5
E0-0309	Motor Module	Pneumatic - Compatible with MicroAire style air hose connector -	MicroAire/5000-100		Non-Disposable	EA	5
E0-0310	Motor Module	Pneumatic - Ability to connect with Hall style air hose connector	MicroAire/5000-200		Non-Disposable	EA	5
E0-0311	Motor Module	Electric -Without hand throttle - compatable with Electric power Console- foot pedal and foot pedal cable - Compatible with Small bone instruments/accessories	MicroAire/5000E		Non-Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0312	Motor Module	Electric - Includes variable speedhand throttle - Compatible with Electric power Console - Compatible with all Small bone modular attachments. Ability to use with foot pedal	MicroAire/5000ET		Non-Disposable	EA	5
E0-0313	STERILIZATION CASE	case for Small bone modular pieces and hand piece with couplers - autoclavable	MicroAire/5000SC		Non-Disposable	EA	5
E0-0314	Cable	Electric cable connecting motor module and Electric Power Control Console.	MicroAire/5006-1000	Installation/assembly required	Non-Disposable	EA	2
E0-0315	Cable	Electric cable for connecting motor module with Electric Power Control Console.	MicroAire/5006-5000	Installation/assembly required	Non-Disposable	EA	2
E0-0316	CABLE	Provides power between foot pedal and power console	MicroAire/5025-5401	Installation/assembly required	Non-Disposable	EA	2
E0-0317	Tubing	Irrigation tubing - Compatible with electric/pneumatic motors system	MicroAire/5601-IT		Disposable	10PK	5
E0-0318	Coupler	K-Wire driver coupler -Compatible with Wires from .028" to .062" (0.6mm - 1.6mm) - Capable of at least 1600 rpm - Capable of at least 135 in. oz. torque	MicroAire/5641-045		Disposable	EA	5
E0-0319	Coupler	Pin driver - Compatible with pins from .078" to .125" (2.0mm - 3.2mm)	MicroAire/5641-050		Disposable	EA	5
E0-0320	STERILIZATION CASE	Medical grade polymer - Ability to hold handpiece with cable in base of case - Ability for include insert tray to hold 10 couplers and Wire/pin collets	MicroAire/5641SC		Non-Disposable	EA	5
E0-0321	FOOT CONTROL SWITCH	compatible with Electric Power Console	MicroAire/6401-000	Installation/assembly required	Non-Disposable	EA	2
E0-0322	Air Hose	10' pneumatic - Ability to connect with MicroAire style connector to a scharder style connector	MicroAire/9000-000	Installation/assembly required	Non-Disposable	EA	4
E0-0323	Air Hose	20' pneumatic hose - Ability to connect with MicroAire style connector to a scharder style connector	MicroAire/9000-020	Installation/assembly required	Non-Disposable	EA	4
E0-0324	THROTTLE RING	Attachment to finger throttle on hand pieces - Ability to control foot pad with hand piece -package of 10	MicroAire/9010-000		Disposable	10PK	5
E0-0325	Air Hose	10' - pneumatic - Compatible with Hall style connection	MicroAire/9013-000	Installation/assembly required	Non-Disposable	EA	4

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0326	AIR PRESSURE TESTER	Ability to test air pressor for Zimmer/hall style Fittings/connectors	MicroAire/9013-005		Non-Disposable	EA	
E0-0327	Air Hose	20' - pneumatic - Compatible with Hall style connection .	MicroAire/9013-020	Installation/assembly required	Non-Disposable	EA	4
E0-0328	Control Valve	Ability to limit pressure to 100PSI -Compatible with Schrader style Connector	MicroAire/9050-000		Non-Disposable	EA	4
E0-0329	Foot Pedal	Pneumatic - Includes 10 Ft. Hose - Compatible with schrader style connector	MicroAire/9100-000		Non-Disposable	EA	4
E0-0330	Air Hose	10' Extension hose - Compatible with Schrader style connection to Schrader style connector	MicroAire/9300-010		Non-Disposable	EA	4
E0-0331	Air Hose	15' Extension hose - Compatible with Schrader style connecting to Schrader style connector	MicroAire/9300-015		Non-Disposable	EA	4
E0-0332	Air Hose	20' Extension hose - Compatible with Schrader style connector to Schrader style connector	MicroAire/9300-020		Non-Disposable	EA	4
E0-0333	Air Hose	10' Extension hose - Compatible with Schrader style connector to Ohio style connector	MicroAire/9300-030		Non-Disposable	EA	4
E0-0334	Air Hose	15' Extension hose - Compatible with Stryker style connector to Stryker style connector	MicroAire/9300-045		Non-Disposable	EA	4
E0-0335	Sterilization Case	Holds one handpiece and one air hose - Stainless steel -18cm x 25cm	MicroAire/9650-000		Non-Disposable	EA	5
E0-0336	Sterilization Case	Medical-grade polymer /aluminum - Ability for the Base to hold 3 hand pieces and hose - Ability for Insert tray to have compartments for 11 couplers and pin driver - Ability for additional space for other accessories - 25cm x 50cm x 10cm	MicroAire/9700-000		Non-Disposable	EA	5
E0-0337	Coupler	Reciprocating saw -Capable of at least 16-000 CPM - Compatible with Small bone reciprocating Blades and Coarse file s	MicroAire/6400		Disposable	EA	5
E0-0338	Wire Driver	Capable of at least 1600 RPM - size: 0.7mm – 1.6mm -Compatible with Wires from .028" to .062" (0.6mm - 1.6mm)	MicroAire/6645		Disposable	EA	5
E0-0339	Pin Driver	Capable of at least 1600 RPM - Size: 1.7mm – 3.2mm - Compatible for pins from .078" to .125" (2.0mm - 3.2mm)	MicroAire/6650		Disposable	EA	5
E0-0340	Coupler	Capable of at least 1600 rpm - Capable of at least 135 in. oz. torque - Compatible with Trinkle style twist drills- screwdrivers- and adapters	MicroAire/6660		Disposable	EA	5
E0-0341	Coupler	Capable of at least 500 rpm - Capable of at least 31 in/lb. torque	MicroAire/6665		Disposable	EA	5
E0-0342	Coupler	Sagittal saw style - Includes key	MicroAire/6670		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0343	Coupler	Keyless - Sagittal saw module - Ability for a simple push button hub - Capable of at least 18-000 CPM	MicroAire/6672		Disposable	EA	5
E0-0344	Coupler	Heavy Duty - Sagittal Saw Module - Capable of at least Operating Speed of 18-000 CPM	MicroAire/6673		Disposable	EA	5
E0-0345	Coupler	5/32" - Jacob Drill style- Capable of at least operating speed of 1600 rpm - Capable of at least 135 in oz. torque - compabtible with pins and dirll bits up to 1/8"	MicroAire/6680		Disposable	EA	5
E0-0346	Coupler	1/4" - High Torque - Jacob Style chuck Drill & Reamer -Capable of at least 500 rpm -Capable of at least 31 in/lb. torque	MicroAire/6685		Disposable	EA	5
E0-0347	Coupler	1/4" (6.4mm) Jacob style Drill - Includes Key -Capable of at least 1600 rpm -Capable of at least 135 in.oz. torque	MicroAire/6686		Disposable	EA	5
E0-0348	Coupler	Reamer - Ability to connect to Zimmer style connection coupler - Capable of at least 500 rpm -Capable of at least 31 in/lb. torque	MicroAire/6687		Disposable	EA	5
E0-0349	Coupler	1/4" (6.4mm) Keyless Drill - Capable of operating speed of 1600 rpm	MicroAire/6688		Disposable	EA	5
E0-0350	Coupler	high torque - reamer style couple - Capable of at least operating speed of 500 rpm - Capable of at least 31 in/lb. torque	MicroAire/6689		Disposable	EA	5
E0-0351	Coupler	Capable of using Synthes style connection - Capable of at least operating speed 1600 rpm - Capable of at least 135 in. oz. torque	MicroAire/6690		Disposable	EA	5
E0-0352	Wire Guard	Sterilized cover	MicroAire/6640-002		Disposable	EA	5
E0-0353	Rechargeable Battery	Large Battery - 14.4 volts - Compatible with saw handpieces	MicroAire/6640-710		Non-Disposable	EA	8
E0-0354	Sterilization Case	Aluminum outer case - Polymer lid - Stainless steel inner tray - Ability to autoclave in automated washer disinfectors systems - Nylon coated brackets- Ability to use Inner tray in the Sterile field - Tray is capable of holding two hand pieces- eight drive couplers- saws- and Wire and pin driver - 23cm x 36cm x 17cm	MicroAire/6643-040		Non-Disposable	EA	5
E0-0355	Coupler	K- Wire driver - Compatible with Wires from .028" to .062 (.6mm - 1.6mm)	MicroAire/6643-045		Non-Disposable	EA	2

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0356	Coupler	Pin driver - Compatible with Wires from .078" to .125" (0.2mm - 3.2mm)	MicroAire/6643-050		Disposable	EA	5
E0-0357	Coupler	Ability to connect with Hudson style connection -Capable of at least 1600 RPM Operating Speed	MicroAire/6665-M1		Disposable	EA	5
E0-0358	Saw Blade	Small - Sagittal - 12.7mm X 19.8mm X 0.5mm	MicroAire/1200-001		Disposable	EA	5
E0-0359	Saw Blade	Small - Sagittal - 12.7mm X 24.3mm X 0.5mm	MicroAire/1200-003		Disposable	EA	5
E0-0360	Saw Blade	Small- Sagittal - 6.9mm X 19.8mm X 0.5mm	MicroAire/1200-005		Disposable	EA	5
E0-0361	Saw blade	Small - Sagittal - 4.1mm X 17mm X 0.5mm	MicroAire/1200-006		Disposable	EA	5
E0-0362	saw Blade	Small - Sagittal - 4.3mm X 15mm X 0.5mm	MicroAire/1200-007		Disposable	EA	5
E0-0363	Saw Blade	Small - Sagittal - 4.3mm X 15mm X 0.5mm	MicroAire/1200-008		Disposable	EA	5
E0-0364	saw Blade	Small - Sagittal - 10.2mm X 17.8mm X 0.5mm	MicroAire/1200-009		Disposable	EA	5
E0-0365	Saw Blade	Small - Sagittal - 10.2mm X 13.7mm X 0.5mm	MicroAire/1200-010		Disposable	EA	5
E0-0366	Saw Blade	Small - Sagittal - 8.1mm X 13.7mm X 0.5mm	MicroAire/1200-011		Disposable	EA	5
E0-0367	Saw Blade	Small - Sagittal - 9.1mm X 13.7mm X 0.5mm	MicroAire/1200-012		Disposable	EA	5
E0-0368	Saw Blade	Small - Sagittal - 12.7mm X 29mm X 0.5mm	MicroAire/1200-013		Disposable	EA	5
E0-0369	Saw Blade	Small - Sagittal - 7.5mm X 15mm X 0.5mm	MicroAire/1200-014		Disposable	EA	5
E0-0370	Saw Blade	Small - Sagittal - 5.8mm X 18mm X 0.5mm	MicroAire/1200-017		Disposable	EA	5
E0-0371	Saw Blade	Small - Sagittal - 10.2mm X 18.3mm X 0.5mm	MicroAire/1200-018		Disposable	EA	5
E0-0372	Saw Blade	Small - Sagittal - 6.9mm X 28.7mm X 0.5mm	MicroAire/1200-019		Disposable	EA	5
E0-0373	Saw Blade	Small - 6.9mm X 13.5mm X 0.5mm	MicroAire/1200-019A		Disposable	EA	5
E0-0374	Saw Blade	Micro - Sagittal - 9.4mm X 25.5mm X 0.5mm	MicroAire/1200-038		Disposable	EA	5
E0-0375	Saw Blade	Small - Reciprical - 8mm X 25.4mm X 0.7mm	MicroAire/1400-001		Disposable	EA	5
E0-0376	Saw Blade	Small - Reciprical - 8mm X 38.1mm X 0.7mm	MicroAire/1400-002		Disposable	EA	5
E0-0377	Saw Blade	Small - Reciprical - 4.3mm X 12.7mm X 0.9mm	MicroAire/1400-003		Disposable	EA	5
E0-0378	saw Blade	Small - Reciprical - 4.3mm X 12.7mm X 0.9mm	MicroAire/1400-004		Disposable	EA	5
E0-0379	saw Blade	Small - Reciprical - 5mm X 22.9mm X 0.6mm	MicroAire/1400-005		Disposable	EA	5
E0-0380	saw Blade	Small - Reciprical - 9mm X 40mm X 0.8mm	MicroAire/1400-006		Disposable	EA	5
E0-0381	saw Blade	Micro - Oscillating - 10.2mm X 11mm X 0.5mm	MicroAire/1700-001		Disposable	EA	5
E0-0382	Saw Blade	Micro - Oscillating - 10.2mm X 12mm X 0.5mm	MicroAire/1700-002		Disposable	EA	5
E0-0383	Saw Blade	Micro - Oscillating - 13mm X 16.5mm X 0.5mm	MicroAire/1700-003		Disposable	EA	5
E0-0384	Saw Blade	Micro - Oscillating - 10.2mm X 11mm X 0.5mm	MicroAire/1700-007		Disposable	EA	5
E0-0385	Saw Blade	Micro - Oscillating - 10.2mm X 12mm X 0.5mm	MicroAire/1700-008		Disposable	EA	5
E0-0386	Saw Blade	Micro - Oscillating - 10.2mm X 7.7mm X 0.5mm	MicroAire/1700-009		Disposable	EA	5
E0-0387	Saw Blade	Micro - Oscillating - 10.2mm X 9.3mm X 0.5mm	MicroAire/1700-010		Disposable	EA	5
E0-0388	Saw Blade	Micro - Oscillating - 10.2mm X 10.2mm X 0.5mm	MicroAire/1700-012		Disposable	EA	5
E0-0389	saw Blade	Micro - Oscillating - 9.1mm X 7mm X 0.5mm	MicroAire/1700-013		Disposable	EA	5
E0-0390	Saw Blade	Compatible wit Synthes style - Sagittal - 19mm X 95mm X 0.9mm	MicroAire/AO-107		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0391	Saw Blade	Compatible wit Synthes style - Sagittal - 25mm X 95mm X 0.9mm	MicroAire/AO-108		Disposable	EA	5
E0-0392	Saw Blade	Compatible wit Synthes style - Sagittal - 19mm X 95mm X1.27mm	MicroAire/AO-114		Disposable	EA	5
E0-0393	saw blade	Compatible wit Synthes style - Sagittal -25mm X 95mm X1.27mm	MicroAire/AO-115		Disposable	EA	5
E0-0394	saw Blade	Compatible wit Synthes style - Sagittal - 19mm X 95mm X1.35mm	MicroAire/AO-116		Disposable	EA	5
E0-0395	saw Blade	Compatible wit Synthes style - Sagittal - 25mm X 95mm X1.27mm	MicroAire/AO-117		Disposable	EA	5
E0-0396	saw Blade	Compatible wit Synthes style - Sagittal - 12.5mm X 90mm X.9mm	MicroAire/AO-118		Disposable	EA	5
E0-0397	Saw Blade	Cast Cutter BLADE - 19mm X 63.5mm X1mm -Non-Sterile	MicroAire/CC-1000NS		Disposable	EA	5
E0-0398	Saw Blade	Large - Oscillating - 31.2mm X 60mm X 1mm	MicroAire/DO-701		Disposable	EA	5
E0-0399	Saw Blade	Large - Oscillating - 15.3mm X 30mm X 1mm	MicroAire/DO-703		Disposable	EA	5
E0-0400	Saw Blade	Large - Oscillating - 9.2mm X 30mm X 1mm	MicroAire/DO-704		Disposable	EA	5
E0-0401	Saw Blade	Large - Oscillating - 13.2mm X 60mm X 1mm	MicroAire/DO-705		Disposable	EA	5
E0-0402	Saw Blade	Large - Oscillating - 15.2mm X 50mm X 1.6mm	MicroAire/MO-121		Disposable	EA	5
E0-0403	Saw Blade	Large - Oscillating - 31.5mm X 50mm X 1.6mm	MicroAire/MO-122		Disposable	EA	5
E0-0404	Saw Blade	Large - Oscillating - 15.2mm X 50mm X 1.1mm	MicroAire/MO-123		Disposable	EA	5
E0-0405	Saw Blade	Large - Oscillating - 31.5mm X 50mm X 1.1mm	MicroAire/MO-124		Disposable	EA	5
E0-0406	Saw Blade	Large - Oscillating - 6.6mm X 20mm X 1mm	MicroAire/MO-125		Disposable	EA	5
E0-0407	Saw Blade	Large - Oscillating - 6.6mm X 20mm X 1mm	MicroAire/MO-126		Disposable	EA	5
E0-0408	Saw Blade	Large - Oscillating - 31.5mm X 50mm X 1.6mm	MicroAire/MO-127		Disposable	EA	5
E0-0409	Saw Blade	Large - Oscillating - 31.5mm X 50mm X 0.9mm	MicroAire/MO-128		Disposable	EA	5
E0-0410	Saw Blade	Large - Oscillating - 15.5mm X 50mm X 0.9mm	MicroAire/MO-129		Disposable	EA	5
E0-0411	Saw blade	Mini - Sagittal - 6.4mm X 15.7mm X 1mm	MicroAire/MS-121		Disposable	EA	5
E0-0412	saw Blade	Mini - Sagittal - 6.4mm X 31.7mm X 1mm	MicroAire/MS-122		Disposable	EA	5
E0-0413	Saw Blade	Mini - Sagittal - 10.4mm X 15.7mm X 1mm	MicroAire/MS-123		Disposable	EA	5
E0-0414	Saw Blade	Mini - Sagittal - 10.4mm X 31.7mm X 1mm	MicroAire/MS-124		Disposable	EA	5
E0-0415	Saw Blade	Mini - Sagittal - 15.5mm X 15.7mm X 1mm	MicroAire/MS-125		Disposable	EA	5
E0-0416	Saw Blade	Mini - Sagittal - 6.4mm X 15.7mm X 0.7mm	MicroAire/MS-131		Disposable	EA	5
E0-0417	Saw Blade	Mini - Sagittal - 6.4mm X 31.7mm X 0.7mm	MicroAire/MS-132		Disposable	EA	5
E0-0418	Saw Blade	Mini - Sagittal - 10.4mm X 31.7mm X 0.7mm	MicroAire/MS-134		Disposable	EA	5
E0-0419	Saw Blade	Small - Oscillating - 10.7mm X 20.3mm X 1mm	MicroAire/SO-008		Disposable	EA	5
E0-0420	Saw Blade	Ability to use in Autopsy - 61mm X 20mm X 0.9mm	MicroAire/SO-010		Disposable	EA	5
E0-0421	Saw Blade	Ability to use in Autopsy - 61mm X 20mm X 0.9mm -Non-Sterile	MicroAire/SO-010NS		Disposable	EA	5
E0-0422	Saw Blade	Large - Oscillating - 21.8mm X 43mm X 1mm	MicroAire/SO-016		Disposable	EA	5
E0-0423	Saw Blade	Ability to use in Autopsy - 96.5mm X 41.3mm X 1mm	MicroAire/SO-020		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0424	Saw Blade	Small - Oscillating - 9.7mm X 25.4mm X 0.7mm	MicroAire/SO-103		Disposable	EA	5
E0-0425	Saw Blade	Large - Oscillating -Dual Tooth - 11mm X 90mm X 0.89mm	MicroAire/SO11-089-090		Disposable	EA	5
E0-0426	Saw Blade	Large - Oscillating -Dual Tooth - 1mm X 90mm X 1.19mm	MicroAire/SO11-119-090		Disposable	EA	5
E0-0427	Saw Blade	Large - Oscillating -Dual Tooth - 11mm X 90mm X 1.27mm	MicroAire/SO11-127-090		Disposable	EA	5
E0-0428	Saw Blade	Large - Oscillating -Dual Tooth - 11mm X 90mm X 1.35mm	MicroAire/SO11-135-090		Disposable	EA	5
E0-0429	Saw Blade	Large - Oscillating - 20.6mm X 31.2mm X 1mm	MicroAire/SO-116		Disposable	EA	5
E0-0430	Saw Blade	Ability to use in Sternum - 42mm X 45mm X 1.1mm	MicroAire/SO-125		Disposable	EA	5
E0-0431	Saw Blade	Large - Oscillating -Dual Tooth - 13mm X 90mm X 1.19mm	MicroAire/SO13-119-090		Disposable	EA	5
E0-0432	Saw Blade	Large - Oscillating -Dual Tooth - 18mm X 75mm X 0.89mm	MicroAire/SO18-089-075		Disposable	EA	5
E0-0433	Saw Blade	Large - Oscillating -Dual Tooth - 18mm X 90mm X 0.89mm	MicroAire/SO18-089-090		Disposable	EA	5
E0-0434	Saw Blade	Large - Oscillating -Dual Tooth - 18mm X 90mm X 1.19mm	MicroAire/SO18-119-090		Disposable	EA	5
E0-0435	Saw Blade	Large - Oscillating -Dual Tooth - 18mm X 75mm X 1.27mm	MicroAire/SO18-127-075		Disposable	EA	5
E0-0436	Saw Blade	Large - Oscillating -Dual Tooth - 18mm X 90mm X 1.27mm	MicroAire/SO18-127-090		Disposable	EA	5
E0-0437	Saw Blade	Large - Oscillating -Dual Tooth - 18mm X 100mm X 1.27mm	MicroAire/SO18-127-100		Disposable	EA	5
E0-0438	Saw Blade	Large - Oscillating -Dual Tooth - 18mm X 90mm X 1.35mm	MicroAire/SO18-135-090		Disposable	EA	5
E0-0439	Saw Blade	Large - Oscillating -Dual Tooth - 25mm X 75mm X 0.89mm	MicroAire/SO25-089-075		Disposable	EA	5
E0-0440	Saw Blade	Large - Oscillating -Dual Tooth - 25mm X 90mm X 0.89mm	MicroAire/SO25-089-090		Disposable	EA	5
E0-0441	Saw Blade	Large - Oscillating -Dual Tooth - 25mm X 90mm X 1.07mm	MicroAire/SO25-107-090		Disposable	EA	5
E0-0442	Saw Blade	Large - Oscillating -Dual Tooth - 25mm X 90mm X 1.19mm	MicroAire/SO25-119-090		Disposable	EA	5
E0-0443	Saw Blade	Large - Oscillating -Dual Tooth - 25mm X 90mm X 1.27mm	MicroAire/SO25-127-090		Disposable	EA	5
E0-0444	Saw Blade	Large - Oscillating -Dual Tooth - 25mm X 100mm X 1.27mm	MicroAire/SO25-127-100		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0445	Saw Blade	Large - Oscillating -Dual Tooth - 25mm X 90mm X 1.35mm	MicroAire/SO25-135-090		Disposable	EA	5
E0-0446	Saw Blade	Large - Oscillating -Dual Tooth - 25mm X 90mm X 1.47mm	MicroAire/SO25-147-090		Disposable	EA	5
E0-0447	Saw Blade	Small - Oscillating - 5.5mm X 12mm X 0.7mm	MicroAire/SO-410		Disposable	EA	5
E0-0448	Saw Blade	Small - Oscillating - 9.7mm X 12mm X 0.7mm	MicroAire/SO-411		Disposable	EA	5
E0-0449	Saw Blade	Small - Oscillating - 5.5mm X 20mm X 0.7mm	MicroAire/SO-412		Disposable	EA	5
E0-0450	Saw Blade	Small - Oscillating - 9.6mm X 20mm X 0.7mm	MicroAire/SO-413		Disposable	EA	5
E0-0451	Saw Blade	Small - Oscillating - 5.5mm X 25mm X 0.7mm	MicroAire/SO-414		Disposable	EA	5
E0-0452	Saw Blade	Crescentic Style Blade - 9.1mm X 16.5mm X 0.6mm	MicroAire/SO-415		Disposable	EA	5
E0-0453	Saw Blade	Crescentic Style Blade - 16.5mm X 24.1mm X 0.6mm	MicroAire/SO-416		Disposable	EA	5
E0-0454	Saw Blade	Crescentic Style Blade - 20.3mm X 24.1mm X 0.6mm	MicroAire/SO-417		Disposable	EA	5
E0-0455	Saw Blade	Mandibular Style Blade - Oscillating - 11.2mm X 12mm X 0.8mm	MicroAire/SO-434		Disposable	EA	5
E0-0456	Saw Blade	Small - Oscillating - 5.5mm X 12mm X 0.6mm	MicroAire/SO-510		Disposable	EA	5
E0-0457	Saw Blade	Small - Oscillating - 5.5mm X 18mm X 0.6mm	MicroAire/SO-512		Disposable	EA	5
E0-0458	Saw Blade	Small - Oscillating - 9.5mm X 20mm X 0.6mm	MicroAire/SO-513		Disposable	EA	5
E0-0459	Saw Blade	Small - Oscillating - 9.5mm X 31.2mm X 0.6mm	MicroAire/SO-515		Disposable	EA	5
E0-0460	Saw Blade	Large - Oscillating - 25.4mm X 70mm X 1mm	MicroAire/SO-600		Disposable	EA	5
E0-0461	Saw Blade	Large - Oscillating - 19.5mm X 90mm X 1.27mm	MicroAire/SO-601		Disposable	EA	5
E0-0462	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1mm	MicroAire/SO-602		Disposable	EA	5
E0-0463	Saw Blade	Large - Oscillating - 25.4mm X 70mm X 1mm	MicroAire/SO-605		Disposable	EA	5
E0-0464	Saw Blade	Large - Oscillating - 19.1mm X 85mm X 1mm	MicroAire/SO-609		Disposable	EA	5
E0-0465	Saw Blade	Large - Oscillating - 25.4mm X 70mm X 1.4mm	MicroAire/SO-610		Disposable	EA	5
E0-0466	Saw Blade	SYSTEM 6Oscillating - 25mm X 73mm X 1mm	MicroAire/SO6-100		Disposable	EA	5
E0-0467	Saw Blade	SYSTEM 6Oscillating - 25mm X 89.5mm X 1mm	MicroAire/SO6-102		Disposable	EA	5
E0-0468	Saw Blade	SYSTEM 6Oscillating - 12.5mm X 89.5mm X 1mm	MicroAire/SO6-103		Disposable	EA	5
E0-0469	Saw Blade	SYSTEM 6Oscillating - 25mm X 73mm X 1mm	MicroAire/SO6-105		Disposable	EA	5
E0-0470	Saw Blade	SYSTEM 6Oscillating - 19mm X 84.5mm X 1mm	MicroAire/SO6-109		Disposable	EA	5
E0-0471	Saw Blade	Large - Oscillating - 25.4mm X 85mm X 1.47mm	MicroAire/SO-611		Disposable	EA	5
E0-0472	Saw Blade	Oscillating -25mm X 73mm X 1.47mm	MicroAire/SO6-110		Disposable	EA	5
E0-0473	Saw Blade	Large - Oscillating - Dual tooth - 11mm X 75mm X 0.89mm	MicroAire/SO611-089-075		Disposable	EA	5
E0-0474	Saw Blade	Large - Oscillating - Dual tooth - 11mm X 90mm X 0.89mm	MicroAire/SO611-089-090		Disposable	EA	5
E0-0475	Saw Blade	Oscillating - 25mm X 79.5mm X 1.47mm	MicroAire/SO6-111		Disposable	EA	5
E0-0476	Saw Blade	Large - Oscillating - Dual tooth - 11mm X 90mm X 1.07mm	MicroAire/SO611-107-090		Disposable	EA	5
E0-0477	Saw Blade	Large - Oscillating - Dual tooth - 11mm X 90mm X 1.19mm	MicroAire/SO611-119-090		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0478	Saw Blade	Large - Oscillating - Dual tooth - 11mm X 90mm X 1.27mm	MicroAire/SO611-127-090		Disposable	EA	5
E0-0479	Saw Blade	Large - Oscillating - Dual tooth - 11mm X 90mm X 1.35mm	MicroAire/SO611-135-090		Disposable	EA	5
E0-0480	Saw Blade	Large - Oscillating - Dual tooth - 11mm X 75mm X 1.47mm	MicroAire/SO611-147-075		Disposable	EA	5
E0-0481	Saw Blade	Oscillating -16.5mm X 34mm X 1.73mm	MicroAire/SO6-113		Disposable	EA	5
E0-0482	Saw Blade	Oscillating - 18.5mm X 73mm X 0.89mm	MicroAire/SO6-115		Disposable	EA	5
E0-0483	Saw Blade	Oscillating -25mm X 79.5mm X 1.24mm	MicroAire/SO6-118		Disposable	EA	5
E0-0484	Saw Blade	Oscillating -18.5mm X 60.5mm X 0.99mm	MicroAire/SO6-120		Disposable	EA	5
E0-0485	Saw Blade	Oscillating -40.5mm X 47.5mm X 1.1mm	MicroAire/SO6-125		Disposable	EA	5
E0-0486	Saw Blade	Oscillating -17.5mm X 47.5mm X 1mm	MicroAire/SO6-131		Disposable	EA	5
E0-0487	Saw Blade	Oscillating - Dual tooth -13mm X 90mm X 1.19mm	MicroAire/SO613-119-090		Disposable	EA	5
E0-0488	Saw Blade	Oscillating -28mm X 80mm X 1mm	MicroAire/SO6-137		Disposable	EA	5
E0-0489	Saw Blade	Oscillating -29mm X 58mm X 1mm	MicroAire/SO6-140		Disposable	EA	5
E0-0490	Saw Blade	Oscillating -5.8mm X 24mm X 0.8mm	MicroAire/SO6-144		Disposable	EA	5
E0-0491	Saw Blade	Oscillating - 9mm X 24mm X 0.8mm	MicroAire/SO6-145		Disposable	EA	5
E0-0492	Saw Blade	Oscillating -9mm X 34.5mm X 0.8mm	MicroAire/SO6-148		Disposable	EA	5
E0-0493	Saw Blade	Large - Oscillating - 19.1mm X 70mm X 0.9mm	MicroAire/SO-615		Disposable	EA	5
E0-0494	Saw Blade	Oscillating -12.5mm X 70mm X 1mm	MicroAire/SO6-150		Disposable	EA	5
E0-0495	Saw Blade	SAW BLADE- 25mm X 85mm X 1.19mm	MicroAire/SO6-151		Disposable	EA	5
E0-0496	Saw Blade	Oscillating -12.5mm X 81.5mm X 1.19mm	MicroAire/SO6-152		Disposable	EA	5
E0-0497	Saw Blade	Oscillating - 12.5mm X 73mm X 1.47mm	MicroAire/SO6-155		Disposable	EA	5
E0-0498	Saw Blade	Oscillating - 12.5mm X 79.5mm X 1.47mm	MicroAire/SO6-156		Disposable	EA	5
E0-0499	Saw Blade	Oscillating -12.5mm X 81.5mm X 1.27mm	MicroAire/SO6-158		Disposable	EA	5
E0-0500	Saw Blade	Oscillating -25mm X 73mm X 1.3mm	MicroAire/SO6-160		Disposable	EA	5
E0-0501	Saw Blade	Oscillating -9.5mm X 60.5mm X 1mm	MicroAire/SO6-175		Disposable	EA	5
E0-0502	Saw Blade	Oscillating - 25mm X 79.5mm X 1.5mm	MicroAire/SO6-176		Disposable	EA	5
E0-0503	Saw Blade	Oscillating -25mm X 86.5mm X 1.32mm	MicroAire/SO6-177		Disposable	EA	5
E0-0504	Saw Blade	Large - Oscillating - 25.4mm X 85mm X 1.27mm	MicroAire/SO-618		Disposable	EA	5
E0-0505	Saw Blade	Oscillating - Dual tooth - 18mm X 75mm X 0.89mm	MicroAire/SO618-089-075		Disposable	EA	5
E0-0506	Saw Blade	Large - Oscillating - Dual tooth - 8mm X 90mm X 0.97mm	MicroAire/SO618-097-090		Disposable	EA	5
E0-0507	Saw Blade	Large - Oscillating - Dual tooth - 18mm X 90mm X 1.19mm	MicroAire/SO618-119-090		Disposable	EA	5
E0-0508	Saw Blade	Large - Oscillating - Dual tooth - 18mm X 75mm X 1.27mm	MicroAire/SO618-127-075		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0509	Saw Blade	Large - Oscillating - Dual tooth - 18mm X 90mm X 1.27mm	MicroAire/SO618-127-090		Disposable	EA	5
E0-0510	Saw Blade	Large - Oscillating - Dual tooth - 18mm X 100mm X 1.27mm	MicroAire/SO618-127-100		Disposable	EA	5
E0-0511	Saw Blade	Large - Oscillating - Dual tooth - 18mm X 90mm X 1.35mm	MicroAire/SO618-135-090		Disposable	EA	5
E0-0512	Saw Blade	Large - Oscillating - Dual tooth - 18mm X 105mm X 1.35mm	MicroAire/SO618-135-105		Disposable	EA	5
E0-0513	Saw Blade	Large - Oscillating - Dual tooth - 18mm X 105mm X 1.47mm	MicroAire/SO618-147-105		Disposable	EA	5
E0-0514	Saw Blade	Oscillating - 18.5mm X 83.5mm X 1.24mm	MicroAire/SO6-182		Disposable	EA	5
E0-0515	Saw Blade	Oscillating -18.5mm X 100mm X 1.24mm	MicroAire/SO6-182-10		Disposable	EA	5
E0-0516	Saw Blade	Oscillating -29mm X 83.5mm X 1.24mm	MicroAire/SO6-183		Disposable	EA	5
E0-0517	Saw Blade	Oscillating - 18.5mm X 104mm X 1.24mm	MicroAire/SO6-185		Disposable	EA	5
E0-0518	Saw Blade	Oscillating -25mm X 81.5mm X 1.27mm	MicroAire/SO6-189		Disposable	EA	5
E0-0519	Saw Blade	Oscillating -25mm X 81.5mm X 1.04mm	MicroAire/SO6-193		Disposable	EA	5
E0-0520	Saw Blade	Oscillating -48.5mm X 90mm X 0.96mm	MicroAire/SO6-195		Disposable	EA	5
E0-0521	Saw Blade	Oscillating -48.5mm X 61mm X 1.08mm	MicroAire/SO6-197		Disposable	EA	5
E0-0522	Saw Blade	Large - Oscillating - 19.1mm X 58mm X 0.9mm	MicroAire/SO-620		Disposable	EA	5
E0-0523	Saw Blade	Oscillating - 25mm X 98.5mm X 1.24mm	MicroAire/SO6-218		Disposable	EA	5
E0-0524	Saw Blade	Large - Oscillating - 19.5mm X 95mm X 1.27mm	MicroAire/SO-623		Disposable	EA	5
E0-0525	Saw Blade	Ability to use in Sternum - 41.9mm X 45mm X 1.1mm	MicroAire/SO-625		Disposable	EA	5
E0-0526	Saw Blade	Large - Oscillating - Dual tooth - 25mm X 75mm X 0.64mm	MicroAire/SO625-064-075		Disposable	EA	5
E0-0527	Saw Blade	Large - Oscillating - Dual tooth - 25mm X 75mm X 0.89mm	MicroAire/SO625-089-075		Disposable	EA	5
E0-0528	Saw Blade	Large - Oscillating - Dual tooth - 25mm X 90mm X 0.89mm	MicroAire/SO625-089-090		Disposable	EA	5
E0-0529	Saw Blade	Large - Oscillating - Dual tooth - 25mm X 90mm X 0.97mm	MicroAire/SO625-097-090		Disposable	EA	5
E0-0530	Saw Blade	Large - Oscillating - Dual tooth - 25mm X 90mm X 1.07mm	MicroAire/SO625-107-090		Disposable	EA	5
E0-0531	Saw Blade	Large - Oscillating - Dual tooth - 25mm X 90mm X 1.19mm	MicroAire/SO625-119-090		Disposable	EA	5
E0-0532	Saw Blade	Large - Oscillating - Dual tooth - 25mm X 90mm X 1.27mm	MicroAire/SO625-127-090		Disposable	EA	5
E0-0533	Saw Blade	Large - Oscillating - Dual tooth - 25mm X 100mm X 1.27mm	MicroAire/SO625-127-100		Disposable	EA	5
E0-0534	Saw Blade	Large - Oscillating - Dual tooth - 25mm X 90mm X 1.35mm	MicroAire/SO625-135-090		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0535	Saw Blade	Large - Oscillating - Dual tooth - 25mm X 90mm X 1.47mm	MicroAire/SO625-147-090		Disposable	EA	5
E0-0536	Saw Blade	Oscillating - 25mm X 89.5mm X 0.89mm	MicroAire/SO6-302		Disposable	EA	5
E0-0537	Saw Blade	Tooth Blade - 25.4mm X 90mm X 0.89mm	MicroAire/SO-6302		Disposable	EA	5
E0-0538	Saw Blade	Oscillating -25mm X 73mm X 0.89mm	MicroAire/SO6-305		Disposable	EA	5
E0-0539	Saw Blade	Tooth Blade - 25.4mm X 75mm X 0.89mm	MicroAire/SO-6305		Disposable	EA	5
E0-0540	Saw Blade	Oscillating - 19mm X 84.5mm X 0.89mm	MicroAire/SO6-309		Disposable	EA	5
E0-0541	Saw Blade	Tooth Blade - 19.1mm X 90mm X 0.89mm	MicroAire/SO-6309		Disposable	EA	5
E0-0542	Saw Blade	Oscillating -25mm X 73mm X 1.24mm	MicroAire/SO6-310		Disposable	EA	5
E0-0543	Saw Blade	Tooth Blade - 25.4mm X 80mm X 1.24mm	MicroAire/SO-6318		Disposable	EA	5
E0-0544	Saw Blade	Tooth Blade - 12.7mm X 85mm X 1.19mm	MicroAire/SO-6352		Disposable	EA	5
E0-0545	Saw Blade	Oscillating -12.5mm X 79mm X 1.24mm	MicroAire/SO6-356		Disposable	EA	5
E0-0546	Saw Blade	Tooth Blade - 12.7mm X 80mm X 1.24mm	MicroAire/SO-6356		Disposable	EA	5
E0-0547	Saw Blade	Ability to use in Sternum - 80mm X 28mm X 1mm	MicroAire/SO-637		Disposable	EA	5
E0-0548	Saw Blade	Tooth Blade - 19mm X 70mm X 1.37mm	MicroAire/SO-6380		Disposable	EA	5
E0-0549	Saw Blade	Sagittal -19mm X 90mm X 1.37mm	MicroAire/SO6-390		Disposable	EA	5
E0-0550	Saw Blade	Large - Oscillating - 19mm X 90mm X 1.35mm	MicroAire/SO-6390		Disposable	EA	5
E0-0551	Saw Blade	Large - Oscillating - 30.5mm X 60mm X 1mm	MicroAire/SO-640		Disposable	EA	5
E0-0552	Saw Blade	Large - Oscillating - 9.7mm X 25mm X 0.8mm	MicroAire/SO-645		Disposable	EA	5
E0-0553	Saw Blade	Large - Oscillating - 9mm X 34mm X 0.8mm	MicroAire/SO-648		Disposable	EA	5
E0-0554	Saw Blade	Large - Oscillating - 12.7mm X 70mm X 1mm	MicroAire/SO-650		Disposable	EA	5
E0-0555	Saw Blade	Sagittal - 25.4mm X 85mm X 1.19mm	MicroAire/SO-651		Disposable	EA	5
E0-0556	Saw Blade	Sagittal - 25.4mm X 95mm X 1.19mm	MicroAire/SO-651L		Disposable	EA	5
E0-0557	Saw Blade	Large - Oscillating - 12.7mm X 70mm X 1.47mm	MicroAire/SO-655		Disposable	EA	5
E0-0558	Saw Blade	Large - Oscillating - 12.7mm X 80mm X 1.47mm	MicroAire/SO-656		Disposable	EA	5
E0-0559	Saw Blade	Oscillating - 12.7mm X 94mm X 1.35mm	MicroAire/SO-658		Disposable	EA	5
E0-0560	Saw Blade	Tooth Blade - 12.7mm X 90mm X 1.37mm	MicroAire/SO-6580		Disposable	EA	5
E0-0561	Saw Blade	Large - Oscillating - 25.4mm X 73mm X 1.3mm	MicroAire/SO-660		Disposable	EA	5
E0-0562	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.27mm	MicroAire/SO-661		Disposable	EA	5
E0-0563	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.37mm	MicroAire/SO-663		Disposable	EA	5
E0-0564	Saw Blade	Large - Oscillating - 9.5mm X 60.5mm X 1mm	MicroAire/SO-675		Disposable	EA	5
E0-0565	Saw Blade	SAGITTAL SAW BLADE- 25.4mm X 79mm X 1.63mm	MicroAire/SO-676		Disposable	EA	5
E0-0566	Saw Blade	Large - Oscillating - 25.4mm X 85mm X 1.35mm	MicroAire/SO-677		Disposable	EA	5
E0-0567	Saw Blade	Large - Oscillating - 19.1mm X 85mm X 1.24mm	MicroAire/SO-682		Disposable	EA	5
E0-0568	Saw Blade	SAGITTAL SAW BLADE- 30.5mm X 83mm X 1.24mm	MicroAire/SO-683		Disposable	EA	5
E0-0569	Saw Blade	Large - Oscillating - 19mm X 105mm X 1.24mm	MicroAire/SO-685		Disposable	EA	5
E0-0570	Saw Blade	SAGITTAL SAW BLADE- 25.4mm X 85mm X 1.37mm	MicroAire/SO-689		Disposable	EA	5
E0-0571	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 0.9mm	MicroAire/SO-696		Disposable	EA	5
E0-0572	Saw Blade	SAGITTAL SAW BLADE- 25.4mm X 85mm X 1.19mm	MicroAire/SO-699		Disposable	EA	5
E0-0573	Saw Blade	Large - Oscillating - 25mm X 98mm X 1.24mm	MicroAire/SO-718		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0574	Saw Blade	Large - Oscillating - 12.7mm X 85mm X 1mm	MicroAire/SO-750		Disposable	EA	5
E0-0575	Saw Blade	Large - Oscillating - 25.4mm X 85mm X 1.2mm	MicroAire/SO-805		Disposable	EA	5
E0-0576	Saw Blade	Small - Oscillating - 6.6mm X 18.5mm X 0.38mm	MicroAire/SP-051		Disposable	EA	5
E0-0577	Saw Blade	Small - Oscillating - 9mm X 10mm X 0.38mm	MicroAire/SP-101		Disposable	EA	5
E0-0578	Saw Blade	Oscillating - 9mm X 10mm X 0.51mm	MicroAire/SP-101A		Disposable	EA	5
E0-0579	Saw Blade	Small - Oscillating - 16.8mm X 18.5mm X 0.38mm	MicroAire/SP-102		Disposable	EA	5
E0-0580	Saw Blade	Oscillating - 16.8mm X 18.5mm X 0.51mm	MicroAire/SP-102A		Disposable	EA	5
E0-0581	Saw Blade	Small - Oscillating - 7mm X 15mm X 0.38mm	MicroAire/SP-103		Disposable	EA	5
E0-0582	Saw Blade	Oscillating - 7mm X 15mm X 0.51mm	MicroAire/SP-103A		Disposable	EA	5
E0-0583	Saw Blade	Small - Oscillating - 16.5mm X 34.5mm X 0.38mm	MicroAire/SP-104		Disposable	EA	5
E0-0584	Saw Blade	Oscillating - 16.5mm X 34.5mm X 0.51mm	MicroAire/SP-104A		Disposable	EA	5
E0-0585	Saw Blade	Small - Oscillating - 9mm X 18.5mm X 0.38mm	MicroAire/SP-105		Disposable	EA	5
E0-0586	Saw Blade	Oscillating - 9mm X 18.5mm X 0.51mm	MicroAire/SP-105A		Disposable	EA	5
E0-0587	Saw Blade	Small - Oscillating - 13mm X 34.5mm X 0.38mm	MicroAire/SP-106		Disposable	EA	5
E0-0588	Saw Blade	Oscillating - 13mm X 34.5mm X 0.51mm	MicroAire/SP-106A		Disposable	EA	5
E0-0589	Saw Blade	Small - Oscillating - 13mm X 39mm X 0.38mm	MicroAire/SP-108		Disposable	EA	5
E0-0590	Saw Blade	Oscillating - 13mm X 39mm X 0.51mm	MicroAire/SP-108A		Disposable	EA	5
E0-0591	Saw Blade	Small - Oscillating - 9mm X 24.6mm X 0.38mm	MicroAire/SP-111		Disposable	EA	5
E0-0592	Saw Blade	Small - Oscillating - 9mm X 24.6mm X 0.51mm	MicroAire/SP-111A		Disposable	EA	5
E0-0593	Saw Blade	Small - Oscillating - 7mm X 18.5mm X 0.38mm	MicroAire/SP-114		Disposable	EA	5
E0-0594	Saw Blade	Oscillating - 7mm X 18.5mm X 0.51mm	MicroAire/SP-114A		Disposable	EA	5
E0-0595	Saw Blade	Small - Oscillating - 7mm X 29.5mm X 0.38mm	MicroAire/SP-115		Disposable	EA	5
E0-0596	Saw Blade	Oscillating - 7mm X 29.5mm X 0.51mm	MicroAire/SP-115A		Disposable	EA	5
E0-0597	Saw Blade	Oscillating - 9mm X 12mm X 0.38mm	MicroAire/SP-121		Disposable	EA	5
E0-0598	Saw Blade	Small - Oscillating - 9mm X 31mm X 0.38mm	MicroAire/SP-125		Disposable	EA	5
E0-0599	Saw Blade	Oscillating - 9mm X 31mm X 0.51mm	MicroAire/SP-125A		Disposable	EA	5
E0-0600	Saw Blade	Oscillating - 9mm X 20mm X 0.64mm	MicroAire/SP-150		Disposable	EA	5
E0-0601	Saw Blade	Oscillating - 9.5mm X 11mm X 0.64mm	MicroAire/SP-151		Disposable	EA	5
E0-0602	Saw Blade	Oscillating - 9mm X 16mm X 0.64mm	MicroAire/SP-152		Disposable	EA	5
E0-0603	Saw Blade	Oscillating - 9mm X 18.5mm X 0.38mm	MicroAire/SP-155		Disposable	EA	5
E0-0604	Saw Blade	Oscillating - 7mm X 20mm X 0.38mm	MicroAire/SP-194		Disposable	EA	5
E0-0605	Saw Blade	Small - Oscillating - 12mm X 34.5mm X 0.38mm	MicroAire/SP-206		Disposable	EA	5
E0-0606	Saw Blade	Oscillating - 12mm X 34.5mm X 0.51mm	MicroAire/SP-206A		Disposable	EA	5
E0-0607	Saw Blade	Small - Oscillating - 5.5mm X 18.5mm X 0.38mm	MicroAire/SP-212		Disposable	EA	5
E0-0608	Saw Blade	Oscillating - 5.5mm X 18.5mm X 0.5mm	MicroAire/SP-212A		Disposable	EA	5
E0-0609	Saw Blade	Oscillating - 19mm X 41mm X 0.38mm	MicroAire/SP-214		Disposable	EA	5
E0-0610	Saw Blade	Oscillating - 9mm X 10mm X 0.51mm	MicroAire/SP-220A		Disposable	EA	5
E0-0611	Saw Blade	Oscillating - 9mm X 12mm X 0.38mm	MicroAire/SP-221		Disposable	EA	5
E0-0612	Saw Blade	Oscillating - 9mm X 12mm X 0.51mm	MicroAire/SP-221A		Disposable	EA	5
E0-0613	Saw Blade	Oscillating - 17mm X 34.5mm X 0.51mm	MicroAire/SP-224A		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0614	Saw Blade	Small - Oscillating - 9mm X 31mm X 0.38mm	MicroAire/SP-225		Disposable	EA	5
E0-0615	Saw Blade	Oscillating - 9mm X 31mm X 0.51mm	MicroAire/SP-225A		Disposable	EA	5
E0-0616	Saw Blade	Oscillating - 9mm X 18.5mm X 0.25mm	MicroAire/SP-23-105		Disposable	EA	5
E0-0617	Saw Blade	Oscillating - 9mm X 25mm X 0.25mm	MicroAire/SP-23-111		Disposable	EA	5
E0-0618	Saw Blade	Oscillating - 5.5mm X 18.5mm X 0.25mm	MicroAire/SP-23-412		Disposable	EA	5
E0-0619	Saw Blade	Oscillating - 5.5mm X 25mm X 0.25mm	MicroAire/SP-23-414		Disposable	EA	5
E0-0620	Saw Blade	Oscillating - 5.5mm X 9mm X 0.9mm	MicroAire/SP-234A		Disposable	EA	5
E0-0621	Saw Blade	Oscillating - 9mm X 18.5mm X 0.38mm	MicroAire/SP-255		Disposable	EA	5
E0-0622	Saw Blade	Small - Oscillating - 5.5mm X 11.5mm X 0.38mm	MicroAire/SP-410		Disposable	EA	5
E0-0623	Saw Blade	Oscillating - 5.5mm X 11.5mm X 0.51mm	MicroAire/SP-410A		Disposable	EA	5
E0-0624	Saw Blade	SAW BLADE - 19mm X 0.64mm X 41.6mm- DEPTH MARKS AT 23mm & 38mm	MicroAire/SP-410R		Disposable	EA	5
E0-0625	Saw Blade	Small - Oscillating - 9mm X 11.5mm X 0.38mm	MicroAire/SP-411		Disposable	EA	5
E0-0626	Saw Blade	Oscillating - 9mm X 11.5mm X 0.51mm	MicroAire/SP-411A		Disposable	EA	5
E0-0627	Saw Blade	Small - Oscillating - 5.5mm X 18.5mm X 0.38mm	MicroAire/SP-412		Disposable	EA	5
E0-0628	Saw Blade	Oscillating - 5.5mm X 18.5mm X 0.51mm	MicroAire/SP-412A		Disposable	EA	5
E0-0629	Saw Blade	Oscillating - 2.5mm X 18.5mm X 0.38mm	MicroAire/SP-413		Disposable	EA	5
E0-0630	Saw Blade	Oscillating - 2.5mm X 18.5mm X 0.51mm	MicroAire/SP-413A		Disposable	EA	5
E0-0631	Saw Blade	Small - Oscillating - 5.5mm X 25mm X 0.38mm	MicroAire/SP-414		Disposable	EA	5
E0-0632	Saw Blade	Oscillating - 5.5mm X 25mm X 0.51mm	MicroAire/SP-414A		Disposable	EA	5
E0-0633	Saw Blade	Oscillating - 16.5mm X 34.5mm X 0.64mm	MicroAire/SP-504		Disposable	EA	5
E0-0634	Saw Blade	Oscillating - 13mm X 34.5mm X 0.64mm	MicroAire/SP-506		Disposable	EA	5
E0-0635	Saw Blade	Oscillating - 13mm X 39mm X 0.64mm	MicroAire/SP-508		Disposable	EA	5
E0-0636	Saw Blade	Oscillating - 9mm X 24.6mm X 0.64mm	MicroAire/SP-511		Disposable	EA	5
E0-0637	Saw Blade	Oscillating - 13.4mm X 39mm X 0.64mm	MicroAire/SP-515		Disposable	EA	5
E0-0638	Saw Blade	Oscillating - 9mm X 31mm X 0.51mm	MicroAire/SP-519A		Disposable	EA	5
E0-0639	Saw Blade	Oscillating - 6.5mm X 18.5mm X 0.51mm	MicroAire/SP-51A		Disposable	EA	5
E0-0640	Saw Blade	Oscillating - 9mm X 35mm X 0.64mm	MicroAire/SP-520A		Disposable	EA	5
E0-0641	Saw Blade	Oscillating - 9mm X 12mm X 0.64mm	MicroAire/SP-521		Disposable	EA	5
E0-0642	Saw Blade	Oscillating - 9mm X 25mm X 0.51mm	MicroAire/SP-521A		Disposable	EA	5
E0-0643	Saw Blade	Oscillating - 9mm X 25mm X 0.64mm	MicroAire/SP-522A		Disposable	EA	5
E0-0644	Saw Blade	Oscillating - 5.5mm X 25mm X 0.64mm	MicroAire/SP-524		Disposable	EA	5
E0-0645	Saw Blade	Oscillating - 9mm X 31mm X 0.64mm	MicroAire/SP-525		Disposable	EA	5
E0-0646	Saw Blade	Oscillating - 13mm X 39mm X 0.51mm	MicroAire/SP-526A		Disposable	EA	5
E0-0647	Saw Blade	Oscillating - 9mm X 8mm X 0.64mm	MicroAire/SP-528		Disposable	EA	5
E0-0648	Saw Blade	Oscillating -9mm X 10mm X 0.64mm	MicroAire/SP-529		Disposable	EA	5
E0-0649	Saw Blade	Oscillating -9mm X 18mm X 0.64mm	MicroAire/SP-530		Disposable	EA	5
E0-0650	Saw Blade	Ability to use in Sternum - 6.4mm X 32.4mm X 1mm	MicroAire/SR-027		Disposable	EA	5
E0-0651	Saw Blade	Ability to use in Sternum - 6.1mm X 32.3mm X 1mm	MicroAire/SR-100		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0652	Saw Blade	Ability to use in Sternum - 6.1mm X 32.3mm X 1mm- Non-Sterile	MicroAire/SR-100NS		Disposable	EA	5
E0-0653	Saw Blade	Ability to use in Sternum - 6.1mm X 32mm X 1mm	MicroAire/SR-101		Disposable	EA	5
E0-0654	Saw Blade	Ability to use in Sternum - 6.1mm X 46.7mm X 1mm	MicroAire/SR-101L		Disposable	EA	5
E0-0655	Saw Blade	Large - Reciprocating - 6.2mm X 60.2mm X 1mm	MicroAire/SR-250		Disposable	EA	5
E0-0656	Saw Blade	Large - Reciprocating - Offset 6.4mm X 47.5mm X 1mm	MicroAire/SR-251		Disposable	EA	5
E0-0657	Saw Blade	Large - Reciprocating - 12.7mm X 78.7mm X 1mm	MicroAire/SR-275		Disposable	EA	5
E0-0658	Saw Blade	Large - Reciprocating - Offset 12.7mm X 68mm X 1.15mm	MicroAire/SR-277		Disposable	EA	5
E0-0659	Saw Blade	Large - Reciprocating - Offset 12.7mm X 68mm X 1.37mm	MicroAire/SR-281		Disposable	EA	5
E0-0660	Saw Blade	Large - Reciprocating - 6.4mm X 78.7mm X 1mm	MicroAire/SR-300		Disposable	EA	5
E0-0661	Saw Blade	Large - Reciprocating - 11.2mm X 77.7mm X 1.07mm	MicroAire/SR-325		Disposable	EA	5
E0-0662	Saw Blade	Large - Reciprocating - Offset 11mm X 77.7mm X 1.6mm	MicroAire/SR-326		Disposable	EA	5
E0-0663	Saw Blade	Reciprocating - 9.4mm X 87.8mm X 1mm	MicroAire/SR-350		Disposable	EA	5
E0-0664	Saw Blade	Small RECIP BLADE- 6.3mm X 27mm X 0.5mm	MicroAire/SR-51-007		Disposable	EA	5
E0-0665	Saw Blade	SM. EXTENDED RECIP BLADE- 6.3mm X 27mm X 0.5mm	MicroAire/SR-51-007L		Disposable	EA	5
E0-0666	Saw Blade	Small RECIP BLADE- 6.3mm X 27mm X 1mm	MicroAire/SR-51-007T		Disposable	EA	5
E0-0667	Saw Blade	Small RECIP BLADE-Short- 6.3mm X 33.5mm X 0.5mm	MicroAire/SR-51-137-103		Disposable	EA	5
E0-0668	Saw Blade	BLADE-Short- 6.3mm X 25mm X 0.5mm	MicroAire/SR-51-137-106		Disposable	EA	5
E0-0669	Saw Blade	Small RECIP BLADE-Short- 6.3mm X 14.5mm X 0.5mm	MicroAire/SR-51-137-11		Disposable	EA	5
E0-0670	Saw Blade	Small RECIP BLADE- Long - 6.3mm X 33.5mm X 0.5mm	MicroAire/SR-51-137-12		Disposable	EA	5
E0-0671	Saw Blade	Small RECIP BLADE-Short- 6.3mm X 7mm X 0.5mm	MicroAire/SR-51-137-120		Disposable	EA	5
E0-0672	Saw Blade	Small RECIP BLADE-Short- 6.3mm X 14.5mm X 0.5mm	MicroAire/SR-51-137-121		Disposable	EA	5
E0-0673	Saw Blade	VERT. OSTEO BLADE-Short- 6.3mm X 34mm X 0.5mm	MicroAire/SR-51-137-20		Disposable	EA	5
E0-0674	Saw Blade	Small RECIP BLADE-Short- 6.3mm X 16mm X 0.5mm	MicroAire/SR-51-137-44		Disposable	EA	5
E0-0675	Saw Blade	Small RECIP BLADE- Long - 6.3mm X 23mm X 0.5mm	MicroAire/SR-51-337-107		Disposable	EA	5
E0-0676	Saw Blade	Small RECIP BLADE-Short- 6.3mm X 16mm X 0.5mm	MicroAire/SR-51-337-111		Disposable	EA	5
E0-0677	Saw Blade	Small RECIP BLADE- Long - 6.3mm X 33.5mm X 0.5mm	MicroAire/SR-51-337-12		Disposable	EA	5
E0-0678	Saw Blade	Small RECIP BLADE- Long - 6.3mm X 14.5mm X 0.5mm	MicroAire/SR-51-337-121		Disposable	EA	5
E0-0679	Saw Blade	Small RECIP BLADE- Long - 6.3mm X 26mm X 0.5mm	MicroAire/SR-51-337-123		Disposable	EA	5
E0-0680	Saw Blade	VERT. OSTEO BLADE- Long - 6.3mm X 34mm X 0.5mm	MicroAire/SR-51-337-20		Disposable	EA	5
E0-0681	Saw Blade	SM. ROUNDED RECIP BLADE- 6mm X 23mm X 0.5mm	MicroAire/SR-51-707		Disposable	EA	5
E0-0682	Saw Blade	Small RECIP BLADE- 6.3mm X 26mm X 0.5mm	MicroAire/SR-51-723		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0683	Saw Blade	Small RECIP BLADE- 6.3mm X 22.5mm X 0.5mm	MicroAire/SR-51-733		Disposable	EA	5
E0-0684	Saw Blade	SM. EXTENDED RECIP BLADE- 6.3mm X 22.5mm X 0.5mm	MicroAire/SR-51-733L		Disposable	EA	5
E0-0685	Saw Blade	Ability to use in Sternum - 6mm X 31.5mm X 0.8mm	MicroAire/TR-032		Disposable	EA	5
E0-0686	Saw Blade	Ability to use in Sternum - 6mm X 31.5mm X 0.8mm- Non-Sterile	MicroAire/TR-032NS		Disposable	EA	5
E0-0687	Saw Blade	Micro - Sagittal - 9.5mm X 25.5mm X 1mm	MicroAire/XZS-038		Disposable	EA	5
E0-0688	Saw Blade	Micro - Sagittal - 9.5mm X 25.5mm X 1mm	MicroAire/XZS-338		Disposable	EA	5
E0-0689	Saw Blade	Small - Oscillating - 5.8mm X 12mm X 0.6mm	MicroAire/ZO-060		Disposable	EA	5
E0-0690	Saw Blade	Small - Oscillating - 9.7mm X 12mm X 0.6mm	MicroAire/ZO-061		Disposable	EA	5
E0-0691	Saw Blade	Small - Oscillating - 5.8mm X 18.5mm X 0.6mm	MicroAire/ZO-062		Disposable	EA	5
E0-0692	Saw Blade	Small - Oscillating - 9.7mm X 18.5mm X 0.6mm	MicroAire/ZO-063		Disposable	EA	5
E0-0693	Saw Blade	Small - Oscillating - 5.8mm X 28.5mm X 0.6mm	MicroAire/ZO-064		Disposable	EA	5
E0-0694	Saw Blade	Small - Oscillating - 9.7mm X 28.5mm X 0.6mm	MicroAire/ZO-065		Disposable	EA	5
E0-0695	Saw Blade	Small - Oscillating - 18.5mm X 28.5mm X 0.6mm	MicroAire/ZO-066		Disposable	EA	5
E0-0696	Saw Blade	Small - Oscillating - 5.8mm X 13.5mm X 0.6mm	MicroAire/ZO-067		Disposable	EA	5
E0-0697	Saw Blade	Small - Oscillating - 9.7mm X 13.5mm X 0.6mm	MicroAire/ZO-068		Disposable	EA	5
E0-0698	Saw Blade	Small - Oscillating - 9.7mm X 13.5mm X 0.6mm	MicroAire/ZO-069		Disposable	EA	5
E0-0699	Saw Blade	CRESCENTIC SAW BLADE- 9.1mm X 11.9mm X 0.6mm	MicroAire/ZO-070		Disposable	EA	5
E0-0700	Saw Blade	CRESCENTIC SAW BLADE- 14mm X 21.6mm X 0.6mm	MicroAire/ZO-071		Disposable	EA	5
E0-0701	Saw Blade	CRESCENTIC SAW BLADE- 22.9mm X 27.7mm X 0.6mm	MicroAire/ZO-072		Disposable	EA	5
E0-0702	Saw Blade	Small - Oscillating - 9.7mm X 31.2mm X 0.6mm	MicroAire/ZO-073		Disposable	EA	5
E0-0703	Saw Blade	Ability to use in Sternum - 45mm X 16mm X 0.7mm	MicroAire/ZO-076		Disposable	EA	5
E0-0704	Saw Blade	Small - Oscillating - 9.4mm X 40.1mm X 0.6mm	MicroAire/ZO-077		Disposable	EA	5
E0-0705	Saw Blade	Ability to use in Sternum - 53.3mm X 30mm X 0.6mm	MicroAire/ZO-079		Disposable	EA	5
E0-0706	Saw Blade	DUALCUT OSCILLATING BLADE- 11mm X 90mm X 0.89mm	MicroAire/ZO11-089-090		Disposable	EA	5
E0-0707	Saw Blade	DUALCUT OSCILLATING BLADE- 11mm X 90mm X 1.19mm	MicroAire/ZO11-119-090		Disposable	EA	5
E0-0708	Saw Blade	DUALCUT OSCILLATING BLADE- 11mm X 90mm X 1.27mm	MicroAire/ZO11-127-090		Disposable	EA	5
E0-0709	Saw Blade	Small - Oscillating - 9.7mm X 27.6mm X 0.6mm	MicroAire/ZO-165		Disposable	EA	5
E0-0710	Saw Blade	DUALCUT OSCILLATING BLADE- 18mm X 90mm X 0.89mm	MicroAire/ZO18-089-090		Disposable	EA	5
E0-0711	Saw Blade	DUALCUT OSCILLATING BLADE- 18mm X 90mm X 1.27mm	MicroAire/ZO18-127-090		Disposable	EA	5
E0-0712	Saw Blade	DUALCUT OSCILLATING BLADE- 18mm X 100mm X 1.27mm	MicroAire/ZO18-127-100		Disposable	EA	5
E0-0713	Saw Blade	DUALCUT OSCILLATING BLADE- 18mm X 90mm X 1.35mm	MicroAire/ZO18-135-090		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0714	Saw Blade	DUALCUT OSCILLATING BLADE- 25mm X 75mm X 0.89mm	MicroAire/ZO25-089-075		Disposable	EA	5
E0-0715	Saw Blade	DUALCUT OSCILLATING BLADE- 25mm X 90mm X 0.89mm	MicroAire/ZO25-089-090		Disposable	EA	5
E0-0716	Saw Blade	DUALCUT OSCILLATING BLADE- 25mm X 90mm X 1.07mm	MicroAire/ZO25-107-090		Disposable	EA	5
E0-0717	Saw Blade	DUALCUT OSCILLATING BLADE- 25mm X 90mm X 1.19mm	MicroAire/ZO25-119-090		Disposable	EA	5
E0-0718	Saw Blade	DUALCUT OSCILLATING BLADE- 25mm X 90mm X 1.27mm	MicroAire/ZO25-127-090		Disposable	EA	5
E0-0719	Saw Blade	DUALCUT OSCILLATING BLADE- 25mm X 100mm X 1.27mm	MicroAire/ZO25-127-100		Disposable	EA	5
E0-0720	Saw Blade	DUALCUT OSCILLATING BLADE- 25mm X 90mm X 1.35mm	MicroAire/ZO25-135-090		Disposable	EA	5
E0-0721	Saw Blade	DUALCUT OSCILLATING BLADE- 25mm X 90mm X 1.47mm	MicroAire/ZO25-147-090		Disposable	EA	5
E0-0722	Saw Blade	SAGITTAL SAW BLADE- 30.5mm X 50mm X 1mm	MicroAire/ZO-640		Disposable	EA	5
E0-0723	Saw Blade	Large - Oscillating - 48mm X 63mm X 0.9mm	MicroAire/ZO-7050		Disposable	EA	5
E0-0724	Saw Blade	Large - Oscillating - 13mm X 64mm X 0.9mm	MicroAire/ZO-7052		Disposable	EA	5
E0-0725	Saw Blade	Large - Oscillating - 19.5mm X 63mm X 0.9mm	MicroAire/ZO-7053		Disposable	EA	5
E0-0726	Saw Blade	Large - Oscillating - 32mm X 68mm X 1.5mm	MicroAire/ZO-7055		Disposable	EA	5
E0-0727	Saw Blade	Ability to use in Sternum - 53mm X 35mm X 0.6mm	MicroAire/ZO-7079		Disposable	EA	5
E0-0728	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.4mm	MicroAire/ZO-7111		Disposable	EA	5
E0-0729	Saw Blade	Large - Oscillating - 19.5mm X 95mm X 1.27mm	MicroAire/ZO-7123		Disposable	EA	5
E0-0730	Saw Blade	Ability to use in Sternum - 42mm X 50mm X 1.1mm	MicroAire/ZO-7125		Disposable	EA	5
E0-0731	Saw Blade	Ability to use in Sternum - 83mm X 30mm X 1mm	MicroAire/ZO-7137		Disposable	EA	5
E0-0732	Saw Blade	Large - Oscillating - 11mm X 30mm X 1.1mm	MicroAire/ZO-7145		Disposable	EA	5
E0-0733	Saw Blade	Large - Oscillating - 33.3mm X 63mm X 0.9mm	MicroAire/ZO-7151		Disposable	EA	5
E0-0734	Saw Blade	Large - Oscillating - 19.5mm X 71mm X 1.35mm	MicroAire/ZO-7154		Disposable	EA	5
E0-0735	Saw Blade	Large - Oscillating - 33.5mm X 64mm X 0.9mm	MicroAire/ZO-7156		Disposable	EA	5
E0-0736	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.27mm	MicroAire/ZO-7161		Disposable	EA	5
E0-0737	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.37mm	MicroAire/ZO-7163		Disposable	EA	5
E0-0738	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.47mm	MicroAire/ZO-7165		Disposable	EA	5
E0-0739	Saw Blade	Large - Oscillating - 32mm X 63mm X 0.9mm	MicroAire/ZO-7166		Disposable	EA	5
E0-0740	Saw Blade	Large - Oscillating - 19.5mm X 90mm X 1.2mm	MicroAire/ZO-7167		Disposable	EA	5
E0-0741	Saw Blade	Large - Oscillating - 19.5mm X 71mm X 0.9mm	MicroAire/ZO-7172		Disposable	EA	5
E0-0742	Saw Blade	Large - Oscillating - 10mm X 85mm X 1.27mm	MicroAire/ZO-7180		Disposable	EA	5
E0-0743	Saw Blade	Large - Oscillating - 19.5mm X 90mm X 1.27mm	MicroAire/ZO-7201		Disposable	EA	5
E0-0744	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1mm	MicroAire/ZO-7205		Disposable	EA	5
E0-0745	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.35mm	MicroAire/ZO-7216		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0746	Saw Blade	STABLECUT(R) RIBBED BLADE- 25.6mm X 90mm X 1.3mm	MicroAire/ZO-72222		Disposable	EA	5
E0-0747	Saw Blade	Ability to use in Sternum - 51mm X 19mm X 0.7mm	MicroAire/ZO-7225		Disposable	EA	5
E0-0748	Saw Blade	Ability to use in Sternum - 63.5mm X 19mm X 0.7mm	MicroAire/ZO-7226		Disposable	EA	5
E0-0749	Saw Blade	Large - Oscillating - 25.4mm X 80mm X 1.4mm	MicroAire/ZO-7240		Disposable	EA	5
E0-0750	Saw Blade	Large - Oscillating - 19.5mm X 80mm X 1.4mm	MicroAire/ZO-7241		Disposable	EA	5
E0-0751	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.4mm	MicroAire/ZO-7242		Disposable	EA	5
E0-0752	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 0.9mm	MicroAire/ZO-7243		Disposable	EA	5
E0-0753	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.5mm	MicroAire/ZO-7244		Disposable	EA	5
E0-0754	Saw Blade	Large - Oscillating - 13mm X 70mm X 1.5mm	MicroAire/ZO-7245		Disposable	EA	5
E0-0755	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 0.9mm	MicroAire/ZO-7246		Disposable	EA	5
E0-0756	Saw Blade	Large - Oscillating - 13mm X 90mm X 1mm	MicroAire/ZO-7250		Disposable	EA	5
E0-0757	Saw Blade	Large - Oscillating - 12mm X 89mm X 1.35mm	MicroAire/ZO-7251		Disposable	EA	5
E0-0758	Saw Blade	OSC. BLADE- 25mm X 90mm X 1.19mm	MicroAire/ZO-72622		Disposable	EA	5
E0-0759	Saw Blade	OSC. BLADE- 25.4mm X 90mm X 1.27mm	MicroAire/ZO-72623		Disposable	EA	5
E0-0760	Saw Blade	OSC. BLADE- 25mm X 90mm X 1.35mm	MicroAire/ZO-72625		Disposable	EA	5
E0-0761	Saw Blade	OSC. BLADE- 25.4mm X 90mm X 0.89mm	MicroAire/ZO-72689		Disposable	EA	5
E0-0762	Saw Blade	OSC. BLADE- 21mm X 90mm X 1mm	MicroAire/ZO-72709		Disposable	EA	5
E0-0763	Saw Blade	OSC. BLADE- 13mm X 90mm X 1.19mm	MicroAire/ZO-72719		Disposable	EA	5
E0-0764	Saw Blade	OSC. BLADE- 19.1mm X 105mm X 1.27mm	MicroAire/ZO-72730		Disposable	EA	5
E0-0765	Saw Blade	OSC. BLADE- 21.3mm X 90mm X 1.35mm	MicroAire/ZO-72737		Disposable	EA	5
E0-0766	Saw Blade	Large - Oscillating - 19.5mm X 86mm X 0.9mm	MicroAire/ZO-7274		Disposable	EA	5
E0-0767	Saw Blade	Large - Oscillating - 19.5mm X 71mm X 0.9mm	MicroAire/ZO-7277		Disposable	EA	5
E0-0768	Saw Blade	Large OSC. SAW BLADE	MicroAire/ZO-7277T		Disposable	EA	5
E0-0769	Saw Blade	Large - Oscillating - 19.5mm X 80mm X 0.9mm	MicroAire/ZO-7278		Disposable	EA	5
E0-0770	Saw Blade	OSC. BLADE- 25.4mm X 95mm X 0.89mm	MicroAire/ZO-7302		Disposable	EA	5
E0-0771	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.2mm	MicroAire/ZO-7305		Disposable	EA	5
E0-0772	Saw Blade	OSC. BLADE- 19.1mm X 90mm X 0.89mm	MicroAire/ZO-7309		Disposable	EA	5
E0-0773	Saw Blade	OSC. BLADE- 25.4mm X 90mm X 1.19mm	MicroAire/ZO-73416		Disposable	EA	5
E0-0774	Saw Blade	RIBBED BLADE- 25.4mm X 90mm X 1.47mm	MicroAire/ZO-73417		Disposable	EA	5
E0-0775	Saw Blade	RIBBED BLADE- 25.4mm X 90mm X 1.37mm	MicroAire/ZO-73418		Disposable	EA	5
E0-0776	Saw Blade	RIBBED BLADE- 25.4mm X 90mm X 1.19mm	MicroAire/ZO-73419		Disposable	EA	5
E0-0777	Saw Blade	OSC. BLADE- 12.7mm X 90mm X 1.27mm	MicroAire/ZO-73420		Disposable	EA	5
E0-0778	Saw Blade	OSC. BLADE- 12.7mm X 90mm X 1mm	MicroAire/ZO-73421		Disposable	EA	5
E0-0779	Saw Blade	OSC. BLADE- 12.7mm X 90mm X 1.37mm	MicroAire/ZO-73422		Disposable	EA	5
E0-0780	Saw Blade	OSC. BLADE- 12.7mm X 90mm X 0.89mm	MicroAire/ZO-73489		Disposable	EA	5
E0-0781	Saw Blade	OSC. BLADE- 21.3mm X 90mm X 1.19mm	MicroAire/ZO-73613		Disposable	EA	5
E0-0782	Saw Blade	OSC. BLADE- 21.3mm X 90mm X 1.27mm	MicroAire/ZO-73616		Disposable	EA	5
E0-0783	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.5mm	MicroAire/ZO-7377		Disposable	EA	5
E0-0784	Saw Blade	Large - Oscillating - 19mm X 85mm X 1.35mm	MicroAire/ZO-7380		Disposable	EA	5
E0-0785	Saw Blade	Large - Oscillating - 19mm X 90mm X 1.35mm	MicroAire/ZO-7390		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0786	Saw Blade	Large - Oscillating - 19.5mm X 90mm X 0.9mm	MicroAire/ZO-7595		Disposable	EA	5
E0-0787	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 0.9mm	MicroAire/ZO-7596		Disposable	EA	5
E0-0788	Saw Blade	Large - Oscillating - 25.4mm X 90mm X 1.19mm	MicroAire/ZO-7651		Disposable	EA	5
E0-0789	Saw Blade	Large - Oscillating - 13mm X 90mm X 1.19mm	MicroAire/ZO-7654		Disposable	EA	5
E0-0790	Saw Blade	OSC. BLADE- 19.1mm X 85mm X 1.27mm	MicroAire/ZO-7682		Disposable	EA	5
E0-0791	Saw Blade	Large - Oscillating - 19.1mm X 105mm X 1.19mm	MicroAire/ZO-7685		Disposable	EA	5
E0-0792	Saw Blade	OSC. BLADE- 12.7mm X 94mm X 0.89mm	MicroAire/ZO-7689		Disposable	EA	5
E0-0793	Saw Blade	OSC. BLADE- 25mm X 105mm X 1.24mm	MicroAire/ZO-7718		Disposable	EA	5
E0-0794	Saw Blade	Small - Reciprical - 7.6mm X 25.4mm X 0.6mm	MicroAire/ZR-020		Disposable	EA	5
E0-0795	Saw Blade	Small - Reciprical - 9.4mm X 39.1mm X 0.6mm	MicroAire/ZR-021		Disposable	EA	5
E0-0796	Saw Blade	Ability to use in Sternum - 9.3mm X 53.6mm X 1mm	MicroAire/ZR-031		Disposable	EA	5
E0-0797	Saw Blade	Ability to use in Sternum - 9.4mm X 34mm X 1mm	MicroAire/ZR-032		Disposable	EA	5
E0-0798	Saw Blade	Ability to use in Sternum - 9.4mm X 34mm X 1mm	MicroAire/ZR-032M		Disposable	EA	5
E0-0799	Saw Blade	Ability to use in Sternum - 9.4mm X 34mm X 1mm- Non-Sterile	MicroAire/ZR-032MNS		Disposable	EA	5
E0-0800	Saw Blade	Ability to use in Sternum - 9.4mm X 34mm X 1mm- Non-Sterile	MicroAire/ZR-032NS		Disposable	EA	5
E0-0801	Saw Blade	STERNUM BLADE- 9.4mm X 34mm X 1mm	MicroAire/ZR-033M		Disposable	EA	5
E0-0802	Saw Blade	Small - Reciprical - 5.8mm X 6.8mm X 0.8mm	MicroAire/ZR-050		Disposable	EA	5
E0-0803	Saw Blade	Small - Reciprical - 5mm X 17.2mm X 0.8mm	MicroAire/ZR-051		Disposable	EA	5
E0-0804	Saw Blade	Small - Reciprical - 6.1mm X 23.1mm X 0.8mm	MicroAire/ZR-052		Disposable	EA	5
E0-0805	Saw Blade	Small - Reciprical - 8.1mm X 25.1mm X 0.8mm	MicroAire/ZR-053		Disposable	EA	5
E0-0806	Saw Blade	Large - Reciprocating - 12.2mm X 87.9mm X 1mm	MicroAire/ZR-058		Disposable	EA	5
E0-0807	Saw Blade	Large - Reciprocating - 9.4mm X 52.8mm X 1mm	MicroAire/ZR-059		Disposable	EA	5
E0-0808	Saw Blade	Large - Reciprocating - 9.4mm X 38.1mm X 0.64mm	MicroAire/ZR-060		Disposable	EA	5
E0-0809	Saw Blade	Large - Reciprocating - 12.2mm X 80mm X 1.78mm	MicroAire/ZR-061		Disposable	EA	5
E0-0810	Saw Blade	Large - Reciprocating - 12.7mm X 73.5mm X 0.9mm	MicroAire/ZR-073		Disposable	EA	5
E0-0811	Saw Blade	RECIP. SAW BLADE- RIBBED- 12.7mm X 76mm X 0.9mm	MicroAire/ZR-079		Disposable	EA	5
E0-0812	Saw Blade	Small - Reciprical - 8.1mm X 25.4mm X 0.5mm	MicroAire/ZR-120		Disposable	EA	5
E0-0813	Saw Blade	Small - Reciprical - 7.4mm X 23.6mm X 0.6mm	MicroAire/ZR-154		Disposable	EA	5
E0-0814	Saw Blade	Reciprocating - 9mm X 52.8mm X 1mm	MicroAire/ZR-159		Disposable	EA	5
E0-0815	Saw Blade	Large - Reciprocating - 12.8mm X 38.1mm X 0.5mm	MicroAire/ZR-160		Disposable	EA	5
E0-0816	Saw Blade	Reciprocating - 6.4mm X 22.9mm X 0.6mm	MicroAire/ZR-733		Disposable	EA	5
E0-0817	Saw Blade	Micro - Sagittal - 5.8mm X 15mm X 0.6mm	MicroAire/ZS-015		Disposable	EA	5
E0-0818	Saw Blade	Micro - Sagittal - 9.7mm X 24.4mm X 0.6mm	MicroAire/ZS-018		Disposable	EA	5
E0-0819	Saw Blade	Micro - Sagittal - 19.1mm X 41mm X 0.6mm	MicroAire/ZS-027		Disposable	EA	5
E0-0820	Saw Blade	Micro - Sagittal - 9.4mm X 7.1mm X 0.6mm	MicroAire/ZS-028		Disposable	EA	5
E0-0821	Saw Blade	Micro - Sagittal - 9.4mm X 13.5mm X 0.6mm	MicroAire/ZS-029		Disposable	EA	5
E0-0822	Saw Blade	Micro - Sagittal - 19.1mm X 25.4mm X 0.6mm	MicroAire/ZS-030		Disposable	EA	5
E0-0823	Saw Blade	Micro - Sagittal - 19.1mm X 41mm X 0.6mm	MicroAire/ZS-031		Disposable	EA	5
E0-0824	Saw Blade	Micro - Sagittal - 5.8mm X 25.4mm X 0.6mm	MicroAire/ZS-032		Disposable	EA	5
E0-0825	Saw Blade	Micro - Sagittal - 14mm X 41mm X 0.6mm	MicroAire/ZS-033		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0826	Saw Blade	Micro - Sagittal - 5.8mm X 8.9mm X 1mm	MicroAire/ZS-034		Disposable	EA	5
E0-0827	Saw Blade	Micro - Sagittal - 5.8mm X 17.3mm X 1mm	MicroAire/ZS-035		Disposable	EA	5
E0-0828	Saw Blade	Micro - Sagittal - 14mm X 16mm X 1mm	MicroAire/ZS-036		Disposable	EA	5
E0-0829	Saw Blade	Micro - Sagittal - 14mm X 25.4mm X 0.6mm	MicroAire/ZS-037		Disposable	EA	5
E0-0830	Saw Blade	Micro - Sagittal - 9.4mm X 25.4mm X 0.6mm	MicroAire/ZS-038		Disposable	EA	5
E0-0831	Saw Blade	Micro - Sagittal - 9.4mm X 25.4mm X 0.6mm	MicroAire/ZS-038A		Disposable	EA	5
E0-0832	Saw Blade	Micro - Sagittal - 9.4mm X 25.4mm X 0.6mm- Non-Sterile	MicroAire/ZS-038NS		Disposable	EA	5
E0-0833	Saw Blade	Micro - Sagittal - 4.5mm X 25.4mm X 0.6mm	MicroAire/ZS-039		Disposable	EA	5
E0-0834	Saw Blade	Micro - Sagittal - 9.4mm X 41mm X 0.6mm	MicroAire/ZS-040		Disposable	EA	5
E0-0835	Saw Blade	Micro - Sagittal - 5.8mm X 33.5mm X 0.6mm	MicroAire/ZS-041		Disposable	EA	5
E0-0836	Saw Blade	Micro - Sagittal - 4.3mm X 35.6mm X 0.6mm	MicroAire/ZS-042		Disposable	EA	5
E0-0837	Saw Blade	Micro - Sagittal - 9.5mm X 25.4mm X 0.6mm	MicroAire/ZS-043		Disposable	EA	5
E0-0838	Saw Blade	Micro - Sagittal - 7mm X 13mm X 0.4mm	MicroAire/ZS-103		Disposable	EA	5
E0-0839	Saw Blade	Micro - Sagittal - 9.4mm X 15mm X 0.6mm	MicroAire/ZS-104		Disposable	EA	5
E0-0840	Saw Blade	Micro - Sagittal - 9mm X 17mm X 0.4mm	MicroAire/ZS-105		Disposable	EA	5
E0-0841	Saw Blade	19.1mm X 40mm X 0.6mm	MicroAire/ZS-110		Disposable	EA	5
E0-0842	Saw Blade	Micro - Sagittal - 7mm X 16mm X 0.4mm	MicroAire/ZS-114		Disposable	EA	5
E0-0843	Saw Blade	Micro - Sagittal - 7mm X 29.5mm X 0.4mm	MicroAire/ZS-115		Disposable	EA	5
E0-0844	Saw Blade	Micro - Sagittal - 9mm X 31mm X 0.4mm	MicroAire/ZS-125		Disposable	EA	5
E0-0845	Saw Blade	Micro - Sagittal - 9mm X 31mm X 0.4mm- Non-Sterile	MicroAire/ZS-125NS		Disposable	EA	5
E0-0846	Saw Blade	Micro - Sagittal	MicroAire/ZS-125NS25		Disposable	25PK	5
E0-0847	Saw Blade	Micro - Sagittal - 10.4mm X 8.1mm X 0.5mm	MicroAire/ZS-128		Disposable	EA	5
E0-0848	Saw Blade	Micro - Sagittal - 10.4mm X 13.5mm X 0.5mm	MicroAire/ZS-129		Disposable	EA	5
E0-0849	Saw Blade	Micro - Sagittal - 9mm X 31mm X 0.4mm	MicroAire/ZS-225		Disposable	EA	5
E0-0850	Saw Blade	Micro - Sagittal - 9.4mm X 15mm X 0.6mm	MicroAire/ZS-304		Disposable	EA	5
E0-0851	Saw Blade	Small - Sagittal - 9.5mm X 25.4mm X 0.5mm	MicroAire/ZS-318		Disposable	EA	5
E0-0852	Saw Blade	Micro - Sagittal - 19.1mm X 41mm X 0.6mm	MicroAire/ZS-327		Disposable	EA	5
E0-0853	Saw Blade	Small - Sagittal - 9.7mm X 14mm X 0.56mm -70 Degree angled radiometer	MicroAire/ZS-329		Disposable	EA	5
E0-0854	Saw Blade	Small - Sagittal - 19mm X 25.5mm X 0.6mm	MicroAire/ZS-330		Disposable	EA	5
E0-0855	Saw Blade	Small - Sagittal - 19.5mm X 41mm X 0.6mm	MicroAire/ZS-331		Disposable	EA	5
E0-0856	Saw Blade	Small - Sagittal - 5.6mm X 25.4mm X 0.6mm	MicroAire/ZS-332		Disposable	EA	5
E0-0857	Saw Blade	Small - Sagittal - 14mm X 41mm X 0.6mm	MicroAire/ZS-333		Disposable	EA	5
E0-0858	Saw Blade	Micro - Sagittal - 5.8mm X 9mm X 1mm	MicroAire/ZS-334		Disposable	EA	5
E0-0859	Saw Blade	Small - Sagittal - 14mm X 16mm X 1mm	MicroAire/ZS-336		Disposable	EA	5
E0-0860	Saw Blade	Small - Sagittal - 14mm X 25.4mm X 0.6mm	MicroAire/ZS-337		Disposable	EA	5
E0-0861	Saw Blade	Small - Sagittal - 9.4mm X 25.4mm X 0.5mm	MicroAire/ZS-338		Disposable	EA	5
E0-0862	Saw Blade	Small - Sagittal - 4.5mm X 25.4mm X 0.6mm	MicroAire/ZS-339		Disposable	EA	5
E0-0863	Saw Blade	Small - Sagittal - 9.5mm X 41mm X 0.5mm	MicroAire/ZS-340		Disposable	EA	5
E0-0864	Saw Blade	Small - Sagittal - 6mm X 34mm X 0.5mm	MicroAire/ZS-341		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0865	Saw Blade	Small - Sagittal - 4mm X 34mm X 0.5mm	MicroAire/ZS-342		Disposable	EA	5
E0-0866	Saw Blade	Small - Sagittal - 9.5mm X 25.4mm X 0.6mm	MicroAire/ZS-343		Disposable	EA	5
E0-0867	Saw Blade	Small - Oscillating - 5.5mm X 11mm X 0.5mm	MicroAire/ZS-360		Disposable	EA	5
E0-0868	Saw Blade	Small - Oscillating - 10mm X 11mm X 0.5mm	MicroAire/ZS-361		Disposable	EA	5
E0-0869	Saw Blade	Small - Oscillating - 5.5mm X 17.5mm X 0.5mm	MicroAire/ZS-362		Disposable	EA	5
E0-0870	Saw Blade	Small - Oscillating - 10mm X 17.5mm X 0.5mm	MicroAire/ZS-363		Disposable	EA	5
E0-0871	Saw Blade	Small - Oscillating - 5.5mm X 27.5mm X 0.5mm	MicroAire/ZS-364		Disposable	EA	5
E0-0872	Saw Blade	Small - Oscillating - 10mm X 27.5mm X 0.5mm	MicroAire/ZS-365		Disposable	EA	5
E0-0873	Saw Blade	Small - Oscillating - 18.5mm X 27mm X 0.5mm	MicroAire/ZS-366		Disposable	EA	5
E0-0874	Saw Blade	Small - Oscillating - 5.8mm X 13.5mm X 0.6mm	MicroAire/ZS-367		Disposable	EA	5
E0-0875	Saw Blade	Small - Oscillating - 9.6mm X 14mm X 0.6mm	MicroAire/ZS-368		Disposable	EA	5
E0-0876	Saw Blade	Small - Oscillating - 10mm X 33.5mm X 0.5mm	MicroAire/ZS-369		Disposable	EA	5
E0-0877	Saw Blade	9mm X 13.2mm X 0.6mm	MicroAire/ZS-370		Disposable	EA	5
E0-0878	Saw Blade	13.5mm X 23.1mm X 0.6mm	MicroAire/ZS-371		Disposable	EA	5
E0-0879	Saw Blade	18mm X 29mm X 0.6mm	MicroAire/ZS-372		Disposable	EA	5
E0-0880	Saw Blade	13.5mm X 29mm X 0.6mm	MicroAire/ZS-373		Disposable	EA	5
E0-0881	Saw Blade	Small -8mm X 12mm X 0.6mm	MicroAire/ZS-382		Disposable	EA	5
E0-0882	Saw Blade	Small -9.5mm X 10mm X 0.6mm	MicroAire/ZS-390		Disposable	EA	5
E0-0883	Saw Blade	Small - 9.5mm X 15mm X 0.6mm	MicroAire/ZS-395		Disposable	EA	5
E0-0884	Saw Blade	Micro - Sagittal - 5.5mm X 9.5mm X 0.4mm	MicroAire/ZS-410		Disposable	EA	5
E0-0885	Saw Blade	Micro - Sagittal - 5.5mm X 18mm X 0.4mm	MicroAire/ZS-412		Disposable	EA	5
E0-0886	Saw Blade	Micro - Sagittal - 5.5mm X 22mm X 0.4mm	MicroAire/ZS-414		Disposable	EA	5
E0-0887	Saw Blade	14mm X 40mm X 0.6mm- RadiameterI Distraction	MicroAire/ZS-510		Disposable	EA	5
E0-0888	Bur	Barrel - 14 FlutesX 6mm diameter	MicroAire/K-01		Disposable	EA	5
E0-0889	Bur	Egg Style -10 FlutesX 4mm diameter	MicroAire/K-02		Disposable	EA	5
E0-0890	Bur	Egg Style -14 FlutesX 6mm diameter	MicroAire/K-03		Disposable	EA	5
E0-0891	Bur	Pear Style - 14 FlutesX 6mm diameter	MicroAire/K-04		Disposable	EA	5
E0-0892	Bur	Round Cutting Style -16 FlutesX 8mm diameter	MicroAire/K-05		Disposable	EA	5
E0-0893	Bur	Flame - 14 FlutesX 6mm diameter	MicroAire/K-08		Disposable	EA	5
E0-0894	Bur	Round Cutting Style - 6 FlutesX 1.6mm diameter	MicroAire/K-09		Disposable	EA	5
E0-0895	Bur	Round Cutting Style - 8 FlutesX 2.3mm diameter	MicroAire/K-10		Disposable	EA	5
E0-0896	Bur	Round Cutting Style - 10 FlutesX 3.1mm diameter	MicroAire/K-11		Disposable	EA	5
E0-0897	Bur	Round Cutting Style - 12 FlutesX 4mm diameter	MicroAire/K-12		Disposable	EA	5
E0-0898	Bur	Round Cutting Style - 12 FlutesX 4.7mm diameter	MicroAire/K-13		Disposable	EA	5
E0-0899	Bur	Round Cutting Style - 8 FlutesX 4mm diameter	MicroAire/K-137		Disposable	EA	5
E0-0900	Bur	Round Cutting Style - 16 FlutesX 6mm diameter	MicroAire/K-14		Disposable	EA	5
E0-0901	Bur	Round Cutting Style - 6 FlutesX 2.4mm diameter	MicroAire/K-147		Disposable	EA	5
E0-0902	Bur	Round Cutting Style - 14 FlutesX 7mm diameter	MicroAire/K-15		Disposable	EA	5
E0-0903	Bur	Wire Pass Drill - 2 FlutesX 1.1mm diameter	MicroAire/K-16		Disposable	EA	5
E0-0904	Bur	Wire Pass Drill - 2 FlutesX 1.3mm diameter	MicroAire/K-17		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0905	Bur	Wire Pass Drill - 2 FlutesX 1.6mm diameter	MicroAire/K-18		Disposable	EA	5
E0-0906	Bur	Wire Pass Drill - 2 FlutesX 2mm diameter	MicroAire/K-19		Disposable	EA	5
E0-0907	Bur	BUR- FISSURE - 8 FlutesX 3.2mm diameter	MicroAire/K-20		Disposable	EA	5
E0-0908	Bur	Serrated - 4 FlutesX 2.3mm diameter	MicroAire/K-21		Disposable	EA	5
E0-0909	Bur	4 FlutesX 2.2mm diameter -used for harvesting bone cores	MicroAire/K-22		Disposable	EA	5
E0-0910	Bur	used for harvesting bone cores - 4 FlutesX 2.2mm diameter- Non-Sterile	MicroAire/K-22NS		Disposable	EA	5
E0-0911	Drill Bit	2 Flutes X 1.5mm diameter	MicroAire/K-23		Disposable	EA	5
E0-0912	Drill Bit	- Tapered Side - 3 FlutesX 2.1mm diameter	MicroAire/K-25		Disposable	EA	5
E0-0913	Bur	4 FlutesX 2mm diameter	MicroAire/K-26		Disposable	EA	5
E0-0914	Bur	6 FlutesX 4mm diameter	MicroAire/K-27		Disposable	EA	5
E0-0915	Bur	3 FlutesX 2.1mm diameter	MicroAire/K-28		Disposable	EA	5
E0-0916	Bur	Ability to side cut	MicroAire/K-32		Disposable	EA	5
E0-0917	Bur	Barrel style14 Flutes X 6mm diameter	MicroAire/K-33		Disposable	EA	5
E0-0918	Bur	Egg Style -14 Flutes X 6mm diameter	MicroAire/K-34		Disposable	EA	5
E0-0919	Bur	Pear Style - 14 Flutes X 6mm diameter	MicroAire/K-35		Disposable	EA	5
E0-0920	Bur	Round Cutting Style - 15 Flutes X 6mm diameter	MicroAire/K-36		Disposable	EA	5
E0-0921	Bur	Egg Style -10 Flutes X 4mm diameter	MicroAire/K-38		Disposable	EA	5
E0-0922	Drill Bit	Round - Diamond - Otology - 1mm diameter	MicroAire/K-40		Disposable	EA	5
E0-0923	Drill Bit	Round - Diamond - Otology - 1.5mm diameter	MicroAire/K-41		Disposable	EA	5
E0-0924	Drill Bit	Round - Diamond - Otology - 2mm diameter	MicroAire/K-42		Disposable	EA	5
E0-0925	Drill Bit	Round - Diamond - Otology - 2.5mm diameter	MicroAire/K-43		Disposable	EA	5
E0-0926	Drill Bit	Round - Diamond - Otology - 3mm diameter	MicroAire/K-44		Disposable	EA	5
E0-0927	Drill Bit	Round - Diamond - Otology - 3.5mm diameter	MicroAire/K-45		Disposable	EA	5
E0-0928	Drill Bit	Round - Diamond - Otology - 4mm diameter	MicroAire/K-46		Disposable	EA	5
E0-0929	Drill Bit	Round - Diamond - Otology - 4.5mm diameter	MicroAire/K-47		Disposable	EA	5
E0-0930	Drill Bit	Round - Diamond - Otology - 5mm diameter	MicroAire/K-48		Disposable	EA	5
E0-0931	Bur	Round - Steel - Otology - 6 Flutes X 2.4mm diameter	MicroAire/K-50		Disposable	EA	5
E0-0932	Bur	Round - Steel - Otology - 6 Flutes X 3.2mm diameter	MicroAire/K-51		Disposable	EA	5
E0-0933	Bur	Round - Steel - Otology - 6 Flutes X 4mm diameter	MicroAire/K-52		Disposable	EA	5
E0-0934	Bur	Round - Steel - Otology - 8 Flutes X 4mm diameter	MicroAire/K-53		Disposable	EA	5
E0-0935	Bur	Round - Steel - Otology - 6 Flutes X 4.8mm diameter	MicroAire/K-54		Disposable	EA	5
E0-0936	Bur	Round - Steel - Otology - 8 Flutes X 4.8mm diameter	MicroAire/K-55		Disposable	EA	5
E0-0937	Bur	Round - Steel - Otology - 8 Flutes X 5.5mm diameter	MicroAire/K-56		Disposable	EA	5
E0-0938	Bur	Round - Steel - Otology - 8 Flutes X 6.3mm diameter	MicroAire/K-57		Disposable	EA	5
E0-0939	Bur	Round - Steel - Otology - 8 Flutes X 6.3mm diameter	MicroAire/K-58		Disposable	EA	5
E0-0940	Bur	Round - Steel - Otology - 12 Flutes X 1mm diameter	MicroAire/K-60		Disposable	EA	5
E0-0941	Bur	Round - Steel - Otology - 18 Flutes X 2mm diameter	MicroAire/K-62		Disposable	EA	5
E0-0942	Bur	Round - Steel - Otology - 20 Flutes X 3mm diameter	MicroAire/K-64		Disposable	EA	5
E0-0943	Bur	Round - Steel - Otology - 25 Flutes X 4mm diameter	MicroAire/K-65		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0944	Bur	Carbide - Ball Cutter - 6 Flutes X 2.3mm diameter	MicroAire/K-70		Disposable	EA	5
E0-0945	Bur	Carbide - Side Cutter - 6 Flutes X 1.6mm diameter	MicroAire/K-71		Disposable	EA	5
E0-0946	Bur	Carbide - Side Cutter - 6 Flutes X 1.2mm diameter	MicroAire/K-72		Disposable	EA	5
E0-0947	Bur	Carbide - Side Cutter - 6 Flutes X 1.6mm diameter	MicroAire/K-73		Disposable	EA	5
E0-0948	Bur	Carbide - Side Cutter - 6 Flutes X 2.1mm diameter	MicroAire/K-74		Disposable	EA	5
E0-0949	Bur	Carbide - Ball Cutter - 6 Flutes X 1.8mm diameter	MicroAire/K-75		Disposable	EA	5
E0-0950	Bur	Carbide - Side Cutter - 6 Flutes X 1.2mm diameter	MicroAire/K-76		Disposable	EA	5
E0-0951	Bur	Carbide - Side Cutter - 6 Flutes X 1.8mm diameter	MicroAire/K-77		Disposable	EA	5
E0-0952	Bur	Carbide - Cross-Cutter - 8 Flutes X 2.1mm diameter	MicroAire/K-78		Disposable	EA	5
E0-0953	Tool	2 Flutes147mm 3.2mm diameter	MicroAire/MRAM-8		Disposable	EA	5
E0-0954	Tool	Ability to use in Disecting	MicroAire/MRG8-330		Disposable	EA	5
E0-0955	Bur	Carbide - End Cutter - 6 Flutes- Medium - 0.8mm diameter	MicroAire/ZB-100		Disposable	EA	5
E0-0956	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 1mm diameter	MicroAire/ZB-101		Disposable	EA	5
E0-0957	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 1.5mm diameter	MicroAire/ZB-102		Disposable	EA	5
E0-0958	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 1mm diameter	MicroAire/ZB-103		Disposable	EA	5
E0-0959	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 1.5mm diameter	MicroAire/ZB-104		Disposable	EA	5
E0-0960	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 2mm diameter	MicroAire/ZB-106		Disposable	EA	5
E0-0961	Bur	Used for harvesting bone cores - 4 Flutes X 2.3mm diameter	MicroAire/ZB-112		Disposable	EA	5
E0-0962	Bur	Carbide - Side Cutter - 6 Flutes- Medium -2mm diameter	MicroAire/ZB-113		Disposable	EA	5
E0-0963	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 2mm diameter	MicroAire/ZB-114		Disposable	EA	5
E0-0964	Bur	Carbide - Round Cutting - 6 Flutes- Medium - 1mm diameter	MicroAire/ZB-121		Disposable	EA	5
E0-0965	Bur	Carbide - Round Cutting - 6 Flutes- Medium - 1.5mm diameter	MicroAire/ZB-122		Disposable	EA	5
E0-0966	Bur	Carbide - Round Cutting - 6 Flutes- Medium - 2mm diameter	MicroAire/ZB-124		Disposable	EA	5
E0-0967	Bur	Carbide - Round Cutting - 6 Flutes- Medium - 2.4mm diameter	MicroAire/ZB-125		Disposable	EA	5
E0-0968	Bur	Carbide - Round Cutting - 8 Flutes- Medium - 3mm diameter	MicroAire/ZB-126		Disposable	EA	5
E0-0969	Bur	Carbide - Round Cutting - 8 Flutes- Medium - 4mm diameter	MicroAire/ZB-128		Disposable	EA	5
E0-0970	Bur	Carbide - Round Cutting - 10 Flutes- Medium - 5mm diameter	MicroAire/ZB-130		Disposable	EA	5
E0-0971	Bur	Carbide - Round Cutting - 10 Flutes- Medium - 6.5mm diameter	MicroAire/ZB-132		Disposable	EA	5
E0-0972	Bur	Carbide - Oval Cutting - 8 Flutes- Medium - 4mm diameter	MicroAire/ZB-136		Disposable	EA	5
E0-0973	Bur	Carbide - Oval Cutting - 10 Flutes- Medium - 5.5mm diameter	MicroAire/ZB-138		Disposable	EA	5
E0-0974	Bur	Reamer - 4 Flutes- Medium - 2mm diameter	MicroAire/ZB-142		Disposable	EA	5
E0-0975	Bur	Reamer - 5 Flutes- Medium - 3mm diameter	MicroAire/ZB-143		Disposable	EA	5
E0-0976	Bur	Reamer - 6 Flutes- Medium - 4mm diameter	MicroAire/ZB-144		Disposable	EA	5
E0-0977	Bur	Reamer - 8 Flutes- Medium - 5mm diameter	MicroAire/ZB-145		Disposable	EA	5
E0-0978	Bur	Wire Pass Drill - 2 Flutes- Medium - 1.1mm diameter	MicroAire/ZB-147		Disposable	EA	5
E0-0979	Bur	Wire Pass Drill - 2 Flutes- Medium - 1.5mm diameter	MicroAire/ZB-148		Disposable	EA	5
E0-0980	Bur	Wire Pass Drill - 2 Flutes- Medium - 2mm diameter	MicroAire/ZB-149		Disposable	EA	5
E0-0981	Bur	Carbide - Barrel -8 Flutes- Medium - 2mm diameter	MicroAire/ZB-152		Disposable	EA	5
E0-0982	Bur	Carbide - Barrel -12 Flutes- Medium - 4mm diameter	MicroAire/ZB-154		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-0983	Drill Bit	Round - Ball Cutter - Diamond - 2mm diameter	MicroAire/ZB-160		Disposable	EA	5
E0-0984	Drill Bit	Round - Ball Cutter - Diamond - 3mm diameter	MicroAire/ZB-162		Disposable	EA	5
E0-0985	Drill Bit	Round - Ball Cutter - Diamond - 4mm diameter	MicroAire/ZB-164		Disposable	EA	5
E0-0986	Drill Bit	Round - Ball Cutter - Diamond - 6mm diameter	MicroAire/ZB-168		Disposable	EA	5
E0-0987	Drill Bit	Oval Cutter - Diamond - 4mm diameter	MicroAire/ZB-170		Disposable	EA	5
E0-0988	Drill Bit	Ability to remove or rount out scars - Diamond - 4mm diameter	MicroAire/ZB-180		Disposable	EA	5
E0-0989	Drill Bit	Ability to remove or rount out scars - Diamond - 7mm diameter	MicroAire/ZB-181		Disposable	EA	5
E0-0990	Bur	Ability to remove or rount out scars - Cylinder - Medium - 10.5mm diameter	MicroAire/ZB-182		Disposable	EA	5
E0-0991	Bur	Ability to remove or rount out scars -- Taper - 10.2mm diameter	MicroAire/ZB-183		Disposable	EA	5
E0-0992	Bur	Ability to remove or rount out scars -- Taper - 5mm diameter	MicroAire/ZB-184		Disposable	EA	5
E0-0993	Bur	Carbide - 4 Flutes- Medium - 1.9mm diameter	MicroAire/ZB-190		Disposable	EA	5
E0-0994	Bur	Shannon 44 - 3 Flutes - Medium - 2.3mm diameter	MicroAire/ZB-192		Disposable	EA	5
E0-0995	Bur	Carbide -4 Flutes- Long - 1.9mm diameter	MicroAire/ZB-200		Disposable	EA	5
E0-0996	Bur	Carbide - Side Cutter - 6 Flutes - Long - 1.7mm diameter	MicroAire/ZB-203		Disposable	EA	5
E0-0997	Bur	Carbide - Side Cutting - 6 Flutes - Long - 1.5mm diameter	MicroAire/ZB-204		Disposable	EA	5
E0-0998	Bur	Carbide - Side Cutting - 6 Flutes - Long - 2mm diameter	MicroAire/ZB-206		Disposable	EA	5
E0-0999	Bur	Carbide - Side Cutting - 6 Flutes - Long - 2mm diameter	MicroAire/ZB-214		Disposable	EA	5
E0-1000	Bur	Carbide - Round Cutting -6 Flutes - Long - 1mm diameter	MicroAire/ZB-220		Disposable	EA	5
E0-1001	Bur	Carbide - Round Cutting -6 Flutes - Long - 1.5mm diameter	MicroAire/ZB-222		Disposable	EA	5
E0-1002	Bur	Carbide - Round Cutting -6 Flutes - Long - 2mm diameter	MicroAire/ZB-224		Disposable	EA	5
E0-1003	Bur	Carbide - Round Cutting - 8 Flutes - Long - 3mm diameter	MicroAire/ZB-226		Disposable	EA	5
E0-1004	Bur	Carbide - Round Cutting - 8 Flutes - Long - 4mm diameter	MicroAire/ZB-228		Disposable	EA	5
E0-1005	Bur	Carbide - Round Cutting - 10 Flutes- Long - 5mm diameter	MicroAire/ZB-230		Disposable	EA	5
E0-1006	Bur	Carbide - Round Cutting - 10 Flutes- Long - 6.5mm diameter	MicroAire/ZB-232		Disposable	EA	5
E0-1007	Bur	Carbide - Oval Cutting - 8 Flutes - Long - 4mm diameter	MicroAire/ZB-236		Disposable	EA	5
E0-1008	Bur	Carbide - Oval Cutting	MicroAire/ZB-236NS25		Disposable	25PK	5
E0-1009	Bur	Carbide - Oval Cutting - 10 Flutes- Long - 5.5mm diameter	MicroAire/ZB-238		Disposable	EA	5
E0-1010	Bur	Carbide - Oval Cutting	MicroAire/ZB-238NS25		Disposable	25PK	5
E0-1011	Bur	Carbide - Barrel - 8 Flutes - Long - 2mm diameter	MicroAire/ZB-252		Disposable	EA	5
E0-1012	Bur	Carbide - Barrel -12 Flutes- Long - 4mm diameter	MicroAire/ZB-254		Disposable	EA	5
E0-1013	Bur	Round - Ball Cutter - Diamond - Long - 2mm diameter	MicroAire/ZB-260		Disposable	EA	5
E0-1014	Bur	Round - Ball Cutter - Diamond - Long - 3mm diameter	MicroAire/ZB-262		Disposable	EA	5
E0-1015	Bur	Round - Ball Cutter - Diamond - Long - 4mm diameter	MicroAire/ZB-264		Disposable	EA	5
E0-1016	Bur	Round - Ball Cutter - Diamond - Long - 5mm diameter	MicroAire/ZB-266		Disposable	EA	5
E0-1017	Bur	Round - Ball Cutter - Diamond - Long - 6mm diameter	MicroAire/ZB-268		Disposable	EA	5
E0-1018	Bur	Oval Cutter - Diamond - Long - 4mm diameter	MicroAire/ZB-270		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1019	Bur	Shannon 44 - 3 Flutes- Long - 2.3mm diameter	MicroAire/ZB-292		Disposable	EA	5
E0-1020	Bur	Carbide - Side Cutter - 6 Flutes- X-Long - 2mm diameter	MicroAire/ZB-306		Disposable	EA	5
E0-1021	Bur	Carbide - Round Cutting - 6 Flutes- X-Long - 1.5mm diameter	MicroAire/ZB-322		Disposable	EA	5
E0-1022	Bur	Carbide - Round Cutting - 6 Flutes- X-Long - 2mm diameter	MicroAire/ZB-324		Disposable	EA	5
E0-1023	Bur	Carbide - Round Cutting - 8 Flutes- X-Long - 3mm diameter	MicroAire/ZB-326		Disposable	EA	5
E0-1024	Bur	Carbide - Round Cutting - 8 Flutes- X-Long - 4mm diameter	MicroAire/ZB-328		Disposable	EA	5
E0-1025	Bur	Carbide - Round Cutting - 10 Flutes- X-Long - 5mm diameter	MicroAire/ZB-330		Disposable	EA	5
E0-1026	Bur	Carbide - Round Cutting	MicroAire/ZB-330MNS		Disposable	EA	5
E0-1027	Bur	Carbide - Oval Cutting - 8 Flutes- X-Long - 4mm diameter	MicroAire/ZB-336		Disposable	EA	5
E0-1028	Bur	Carbide - Oval Cutting	MicroAire/ZB-336MNS		Disposable	EA	5
E0-1029	Bur	Carbide - Oval Cutting - 10 Flutes- X-Long - 5.5mm diameter	MicroAire/ZB-338		Disposable	EA	5
E0-1030	Drill Bit	Round - Ball Cutter - Diamond - X-Long - 2mm diameter	MicroAire/ZB-360		Disposable	EA	5
E0-1031	Drill Bit	Round Ball Cutter - Diamond -	MicroAire/ZB-360S		Disposable	EA	5
E0-1032	Drill Bit	Round - Ball Cutter - Diamond - X-Long - 3mm diameter	MicroAire/ZB-362		Disposable	EA	5
E0-1033	Drill Bit	Round - Ball Cutter - Diamond - X-Long - 4mm diameter	MicroAire/ZB-364		Disposable	EA	5
E0-1034	Drill Bit	Round - Ball Cutter - Diamond - X-Long - 5mm diameter	MicroAire/ZB-366		Disposable	EA	5
E0-1035	Bur	Wire Pass Drill - 2 Flutes -Short- 1mm diameter	MicroAire/ZB-447		Disposable	EA	5
E0-1036	Bur	Wire Pass Drill - 2 FlutesSHORT 1.5mm diameter	MicroAire/ZB-448		Disposable	EA	5
E0-1037	Bur	Wire Pass Drill - 2 Flutes-Short- 2mm diameter	MicroAire/ZB-449		Disposable	EA	5
E0-1038	Bur	Cebotome - 2 Flutes- 6mm diameter	MicroAire/ZB-582		Disposable	EA	5
E0-1039	Bur	Cebotome - 10 Flutes- 7mm diameter	MicroAire/ZB-584		Disposable	EA	5
E0-1040	Bur	Cebotome - 10 Flutes- 8mm diameter	MicroAire/ZB-585		Disposable	EA	5
E0-1041	Bur	Cebotome - 10 Flutes- 8mm diameter	MicroAire/ZB-586		Disposable	EA	5
E0-1042	Bur	Carbide - Otology - Short - 10 Flutes - 30mm - 6.5mm diameter	MicroAire/ZB-617		Disposable	EA	5
E0-1043	Bur	Carbide - Otology - 10 Flutes- 30mm- 5mm diameter	MicroAire/ZB-618		Disposable	EA	5
E0-1044	Bur	Carbide - Otology - 8 Flutes- 30mm- 4mm diameter	MicroAire/ZB-619		Disposable	EA	5
E0-1045	Bur	Carbide - Otology - 10 Flutes38mm 6.5mm diameter	MicroAire/ZB-625		Disposable	EA	5
E0-1046	Bur	Carbide - Otology - 10 Flutes- 38mm- 5mm diameter	MicroAire/ZB-626		Disposable	EA	5
E0-1047	Bur	Carbide - Otology - 8 Flutes- 38mm- 4mm diameter	MicroAire/ZB-627		Disposable	EA	5
E0-1048	Bur	Carbide - Otology - 8 Flutes- 38mm- 3mm diameter	MicroAire/ZB-628		Disposable	EA	5
E0-1049	Bur	Otology - 6 Flutes- 48mm- 1.5mm diameter	MicroAire/ZB-640		Disposable	EA	5
E0-1050	Bur	Otology - 6 Flutes- 48mm- 2mm diameter	MicroAire/ZB-642		Disposable	EA	5
E0-1051	Bur	Carbide - Otology - 10 Flutes- 48mm- 5mm diameter	MicroAire/ZB-644		Disposable	EA	5
E0-1052	Bur	Carbide - Otology - 8 Flutes- 48mm- 4mm diameter	MicroAire/ZB-645		Disposable	EA	5
E0-1053	Bur	Carbide - Otology - 8 Flutes- 48mm- 3mm diameter	MicroAire/ZB-646		Disposable	EA	5
E0-1054	Bur	Diamond - Otology - 38mm- 5mm diameter	MicroAire/ZB-647		Disposable	EA	5
E0-1055	Bur	Diamond - Otology - 38mm- 4mm diameter	MicroAire/ZB-648		Disposable	EA	5
E0-1056	Bur	Diamond - Otology - 48mm- 5mm diameter	MicroAire/ZB-650		Disposable	EA	5
E0-1057	Bur	Diamond - Otology - 48mm- 4mm diameter	MicroAire/ZB-652		Disposable	EA	5
E0-1058	Bur	Diamond - Otology - 48mm- 1mm diameter	MicroAire/ZB-658		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1059	Bur	Diamond - Otology - 48mm- 3mm diameter	MicroAire/ZB-659		Disposable	EA	5
E0-1060	Bur	Diamond - Otology - 48mm- 0.5mm diameter	MicroAire/ZB-661		Disposable	EA	5
E0-1061	Bur	Diamond - Otology - 48mm- 2mm diameter	MicroAire/ZB-662		Disposable	EA	5
E0-1062	Bur	Diamond - Otology - 48mm- 1.5mm diameter	MicroAire/ZB-663		Disposable	EA	5
E0-1063	Bur	Otology - 6 Flutes- 48mm- 1mm diameter	MicroAire/ZB-664		Disposable	EA	5
E0-1064	Bur	Otology - 4 Flutes- 48mm- 0.5mm diameter	MicroAire/ZB-666		Disposable	EA	5
E0-1065	Bur	Ability to remove or round out scars - 30mm- 4mm diameter	MicroAire/ZB-680		Disposable	EA	5
E0-1066	Bur	Ability to remove or round out scars - 30mm- 10mm diameter	MicroAire/ZB-683		Disposable	EA	5
E0-1067	Bur	Carbide - Round Cutting - Medium- 2.3mm diameter	MicroAire/ZB-70		Disposable	EA	5
E0-1068	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 1.6mm diameter	MicroAire/ZB-71		Disposable	EA	5
E0-1069	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 1.2mm diameter	MicroAire/ZB-72		Disposable	EA	5
E0-1070	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 1.6mm diameter	MicroAire/ZB-73		Disposable	EA	5
E0-1071	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 2.1mm diameter	MicroAire/ZB-74		Disposable	EA	5
E0-1072	Bur	Carbide - Round Cutting - 6 Flutes- Medium - 1.8mm diameter	MicroAire/ZB-75		Disposable	EA	5
E0-1073	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 1.2mm diameter	MicroAire/ZB-76		Disposable	EA	5
E0-1074	Bur	Carbide - Side Cutter - 6 Flutes- Medium - 1.2mm diameter	MicroAire/ZB-77		Disposable	EA	5
E0-1075	Bur	Carbide - Cross-Cutter - 8 Flutes- Medium - 2.1mm diameter	MicroAire/ZB-78		Disposable	EA	5
E0-1076	Cutter	Femoral Mill	MicroAire/2135-020		Disposable	1PK	5
E0-1077	Tip	Shower Spray - Includes splash shield	MicroAire/201-00-37		Disposable	6PK	5
E0-1078	Tip	Fan Spray - Includes splash shield	MicroAire/201-00-44		Disposable	6PK	5
E0-1079	Tip	Femoral Spray	MicroAire/201-00-46		Disposable	6PK	5
E0-1080	Tip	Shower Spary	MicroAire/201-00-50		Disposable	6PK	5
E0-1081	Tip	Fan Spay	MicroAire/201-00-51		Disposable	6PK	5
E0-1082	Tip	Femoral Brush	MicroAire/201-00-52		Disposable	6PK	5
E0-1083	Shield	Splash Shield - 6 INCH (15.2 CM)	MicroAire/201-00-80		Disposable	6PK	5
E0-1084	COMPONENT KIT	Includes Fully Disposable Pulse Lavage- Includes handpiece and battery	MicroAire/201-02-00		Disposable	6PK	5
E0-1085	KNEE KIT	Includes Fully Disposable Pulse Lavage- Includes handpiece and battery - Includes Shower Pray Tip with splash shield	MicroAire/201-03-00		Disposable	6PK	5
E0-1086	HIP KIT	includes Fully Disposable Pulse Lavage- handpiece and battery- FEMORAL SPRAY TIP SHOWER SPRAY TIP WITH SPLASH SHIELD	MicroAire/201-04-00		Disposable	6PK	5
E0-1087	FAN KIT	Includes Pulse Lavage - Includes handpiece - Includes battery - Fan Spray with with splach shield	MicroAire/201-05-00		Disposable	6PK	5
E0-1088	Splash Shield	6" - Ability to use for irrigation	MicroAire/4700-015		Disposable	6PK	5
E0-1089	Splash Shield	4" - Ability to use for irrigation	MicroAire/4700-016		Disposable	6PK	5
E0-1090	Splash Shield	1.6" - Cone- Ability to use for irrigation	MicroAire/4740-013		Disposable	6PK	5
E0-1091	Nozzle	4" - Single Orifice Nozzle	MicroAire/5740-004		Disposable	6PK	5
E0-1092	Nozzle	4" - Single Orifice nozzle with suction capability	MicroAire/5740-004S		Disposable	6PK	5
E0-1093	Nozzle	Shower-Head	MicroAire/5740-005		Disposable	6PK	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1094	Nozzle	Shower Head - Ability to suction	MicroAire/5740-005S		Disposable	6PK	5
E0-1095	Nozzle	Femoral Spray	MicroAire/5740-006		Disposable	6PK	5
E0-1096	Nozzle	Femoral Spray - Ability to suction	MicroAire/5740-006S		Disposable	6PK	5
E0-1097	Nozzle	Angled - Shower Head	MicroAire/5740-008		Disposable	6PK	5
E0-1098	Nozzle	Extended - Femoral Spray	MicroAire/5740-009		Disposable	6PK	5
E0-1099	Nozzle	Extended - Femoral Spray - Ability to provide suction	MicroAire/5740-009S		Disposable	6PK	5
E0-1100	Nozzle	Femoral Brush - Ability to provide Suction	MicroAire/5740-010S		Disposable	6PK	5
E0-1101	Tubing	Compatible for irrigation -Compatible with Hand pieces	MicroAire/5740-020		Disposable	6PK	5
E0-1102	Tubing	Ability to use with suction - Ability to use with for irrigation	MicroAire/5740-020S		Disposable	6PK	5
E0-1103	Tubing	Includes nozzle set - Ability to use with irrigation	MicroAire/5740-060		Disposable	6PK	5
E0-1104	Tubing	Ability to use with irrigation - Includes Standard Tubing- Include Showerhead Nozzle- Includes Femoral Spray Nozzle	MicroAire/5740-061		Disposable	6PK	5
E0-1105	Tubing	Ability to use with irrigation - Includes Standard Tubing - Includes Showerhead Nozzle - Includes Femoral Spray Nozzle.	MicroAire/5740-061S		Disposable	6PK	5
E0-1106	Tubing	Abiliyt to use with irrigation - Includes Tubing - Includes Angled Showerhead Nozzle.	MicroAire/5740-062		Disposable	6PK	5
E0-1107	Tubing	Ability to use with irrigation - Includes Tubing - Included Angled Showerhead Nozzle.	MicroAire/5740-062S		Disposable	6PK	5
E0-1108	Handpiece	Pneumatic - Abililty to use for irrigation - Compatible with MicroAire style hose connector	MicroAire/5740-100		Non-Disposable	EA	4
E0-1109	Handpiece	Pneumatic -Ability to use with irrigation - Compatible with Hall style hose connector	MicroAire/5740-200		Non-Disposable	EA	4
E0-1110	Handpiece	Pneumatic - Ability to use with irrigation -Compatible with MicroAire style hose connector	MicroAire/5740-600		Non-Disposable	EA	4
E0-1111	Modular Connector	Ability to debrid bone surface - Ability to clean trauma - Ability to clean soft tissue wounds	MicroAire/6640-474		Non-Disposable	EA	4
E0-1112	Coupler	Compatible with pulse lavage- Compatible with Large bone Drill Reamers- Compatible with irrigation pulse lavage tubing sets.	MicroAire/7100-474		Disposable	EA	5
E0-1113	Rasp	Coarse file - Cross CutTear - 5.3mm X 11mm	MicroAire/1400-011		Disposable	EA	5
E0-1114	Rasp	Coarse file - Straight Cottle - 5.6mm X 15.5mm	MicroAire/1400-013		Disposable	EA	5
E0-1115	Rasp	Coarse file - Small - Straight Cut - 5.5mm X 9.5mm	MicroAire/SR-080		Disposable	EA	5
E0-1116	Rasp	Coarse file - Large - Straight Cut - 6.5mm X 12.5mm	MicroAire/SR-081		Disposable	EA	5
E0-1117	Rasp	Coarse file - Small - Cross Cut- 5.5mm X 9.5mm	MicroAire/SR-082		Disposable	EA	5
E0-1118	Rasp	Coarse file - Large - Cross Cut- 6.5mm X 12.5mm	MicroAire/SR-083		Disposable	EA	5
E0-1119	Rasp	Small Tear - Coarse file - 5mm X 10.7mm	MicroAire/SR-51-080		Disposable	EA	5
E0-1120	Rasp	Large Tear - Coarse file - 7mm X 12.7mm	MicroAire/SR-51-081		Disposable	EA	5
E0-1121	Rasp	Small Tear - Cross Cut -Coarse file - 5mm X 11mm	MicroAire/SR-51-082		Disposable	EA	5
E0-1122	Rasp	Large Tear - Cross Cut - Coarse file - 7mm X 14mm	MicroAire/SR-51-083		Disposable	EA	5
E0-1123	Rasp	Coarse file - Small- Straight Cut- 5.5mm X 9.5mm	MicroAire/ZR-080		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1124	Rasp	Coarse file - Large- Straight Cut- 6.5mm X 12.5mm	MicroAire/ZR-081		Disposable	EA	5
E0-1125	Rasp	Coarse file - Small- Cross Cut- 5.5mm X 9.5mm	MicroAire/ZR-082		Disposable	EA	5
E0-1126	Rasp	Coarse file - Large- Cross Cut- 6.5mm X 12.5mm	MicroAire/ZR-083		Disposable	EA	5
E0-1127	Rasp	Coarse file - Large- Cross Cut- 6.5mm X 12.5mm- Non-Sterile	MicroAire/ZR-083NS		Disposable	EA	5
E0-1128	Rasp	Coarse file -Small-Intramedullary- 5.3mm X 24.5mm	MicroAire/ZR-084		Disposable	EA	5
E0-1129	Rasp	Coarse file -Large-Intramedullary- 9.5mm X 30.5mm	MicroAire/ZR-085		Disposable	EA	5
E0-1130	Drill Bit	Stainless Stell - Twist - Dental- 1.1mm diameter X 49mm	MicroAire/8051-011		Disposable	EA	5
E0-1131	Drill Bit	Stainless Stell - Twist - Dental- 1.3mm diameter X 49mm	MicroAire/8051-013		Disposable	EA	5
E0-1132	Drill Bit	Stainless Stell - Twist - Dental- 1.5mm diameter X 49mm	MicroAire/8051-015		Disposable	EA	5
E0-1133	Drill Bit	Stainless Stell - Twist - Dental- 1.8mm diameter X 49mm	MicroAire/8051-018		Disposable	EA	5
E0-1134	Drill Bit	Stainless Stell - Twist - Dental- 2mm diameter X 65mm	MicroAire/8051-020		Disposable	EA	5
E0-1135	Drill Bit	Stainless Stell - Twist - Dental- 2mm diameter X 65mm- Non-Sterile	MicroAire/8051-020NS		Disposable	EA	5
E0-1136	Drill Bit	Stainless Stell - Twist - Dental- 2.2mm diameter X 49mm	MicroAire/8051-021		Disposable	EA	5
E0-1137	Drill Bit	Stainless Stell - Twist - Dental- 2.7mm diameter X 79mm	MicroAire/8051-027		Disposable	EA	5
E0-1138	Drill Bit	Stainless Stell - Twist - Dental- 3.2mm diameter X 79mm	MicroAire/8051-032		Disposable	EA	5
E0-1139	Drill Bit	Stainless Stell - Twist - Dental- 3.5mm diameter X 79mm	MicroAire/8051-035		Disposable	EA	5
E0-1140	Drill Bit	Stainless Stell - Twist - Small- 2mm diameter X 65mm	MicroAire/8051-120		Disposable	EA	5
E0-1141	Drill Bit	Stainless Stell - Twist - Small- 2.7mm diameter X 79mm	MicroAire/8051-127		Disposable	EA	5
E0-1142	Drill Bit	Stainless Stell - Twist - AO Type- 1.1mm diameter X 35mm	MicroAire/8053-011		Disposable	EA	5
E0-1143	Drill Bit	Stainless Stell - Twist - AO Type- 1.5mm diameter X 60mm	MicroAire/8053-015		Disposable	EA	5
E0-1144	Drill Bit	Stainless Stell - Twist - AO Type- 1.5mm diameter X 85mm	MicroAire/8053-016		Disposable	EA	5
E0-1145	Drill Bit	Stainless Stell - Twist - AO Type- 1.5mm diameter X 85mm- Non-Sterile	MicroAire/8053-016NS		Disposable	EA	5
E0-1146	Drill Bit	Stainless Stell - Twist - AO Type- 2mm diameter X 75mm	MicroAire/8053-019		Disposable	EA	5
E0-1147	Drill Bit	Stainless Stell - Twist - AO Type- 2mm diameter X 75mm- Non-Sterile	MicroAire/8053-019NS		Disposable	EA	5
E0-1148	Drill Bit	Stainless Stell - Twist - Jacob- 2mm diameter X 27mm	MicroAire/8053-020		Disposable	EA	5
E0-1149	Drill Bit	Stainless Stell - Twist - AO Type- 2mm diameter X 100mm	MicroAire/8053-021		Disposable	EA	5
E0-1150	Drill Bit	Stainless Stell - Twist - AO Type- 2mm diameter X 100mm- Non-Sterile	MicroAire/8053-021NS		Disposable	EA	5
E0-1151	Drill Bit	Stainless Stell - Twist - AO Type- 2.5mm diameter X 155mm	MicroAire/8053-023		Disposable	EA	5
E0-1152	Drill Bit	Stainless Stell - Twist - AO Type- 2.5mm diameter X 155mm- Non-Sterile	MicroAire/8053-023NS		Disposable	EA	5
E0-1153	Drill Bit	Stainless Stell - Twist - Jacob- 2.5mm diameter X 60mm	MicroAire/8053-024		Disposable	EA	5
E0-1154	Drill Bit	Stainless Stell - Twist - AO Type- 2.5mm diameter X 85mm	MicroAire/8053-025		Disposable	EA	5
E0-1155	Drill Bit	Stainless Stell - Twist - AO Type- 2.5mm diameter X 85mm- Non-Sterile	MicroAire/8053-025NS		Disposable	EA	5
E0-1156	Drill Bit	Stainless Stell - Twist - AO Type- 2.7mm diameter X 75mm	MicroAire/8053-026		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1157	Drill Bit	Stainless Stell - Twist - AO Type- 2.7mm diameter X 75mm- Non-Sterile	MicroAire/8053-026NS		Disposable	EA	5
E0-1158	Drill Bit	Stainless Stell - Twist - AO Type- 2.7mm diameter X 100mm	MicroAire/8053-028		Disposable	EA	5
E0-1159	Drill Bit	Stainless Stell - Twist - AO Type- 3.2mm diameter X 170mm	MicroAire/8053-029		Disposable	EA	5
E0-1160	Drill Bit	Stainless Stell - Twist - AO Type- 3.2mm diameter X 120mm	MicroAire/8053-031		Disposable	EA	5
E0-1161	Drill Bit	Stainless Stell - Twist - AO Type- 3.2mm diameter X 120mm- Non-Sterile	MicroAire/8053-031NS		Disposable	EA	5
E0-1162	Drill Bit	Stainless Stell - Twist - Jacob- 3.2mm diameter X 85mm	MicroAire/8053-032		Disposable	EA	5
E0-1163	Drill Bit	Stainless Stell - Twist - AO Type- 3.5mm diameter X 85mm	MicroAire/8053-035		Disposable	EA	5
E0-1164	Drill Bit	Stainless Stell - Twist - AO Type- 3.5mm diameter X 85mm- Non-Sterile	MicroAire/8053-035NS		Disposable	EA	5
E0-1165	Drill Bit	Stainless Stell - Twist - AO Type	MicroAire/8053-035NS25		Disposable	25PK	5
E0-1166	Drill Bit	Stainless Stell - Twist - Jacob- 3.5mm diameter X 60mm	MicroAire/8053-036		Disposable	EA	5
E0-1167	Drill Bit	Stainless Stell - Twist - AO Type- 3.5mm diameter X 170mm	MicroAire/8053-037		Disposable	EA	5
E0-1168	Drill Bit	Stainless Stell - Twist - AO Type- 4.5mm diameter X 120mm	MicroAire/8053-044		Disposable	EA	5
E0-1169	Drill Bit	Stainless Stell - Twist - AO Type- 4.5mm diameter X 120mm- Non-Sterile	MicroAire/8053-044NS		Disposable	EA	5
E0-1170	Drill Bit	Stainless Stell - Twist - AO Type- 4.5mm diameter X 170mm	MicroAire/8053-048		Disposable	EA	5
E0-1171	Drill Bit	Stainless Stell - Twist - AO Type- 6mm diameter X 170mm	MicroAire/8053-060		Disposable	EA	5
E0-1172	Drill Bit	Stainless Stell - Twist - AO Type- 6mm diameter X 170mm- Non-Sterile	MicroAire/8053-060NS		Disposable	EA	5
E0-1173	Drill Bit	Stainless Steel - Tap- 1.5mm diameter- 1.1mm DRILL- 1.5mm SCREW	MicroAire/8053-115		Disposable	EA	5
E0-1174	Drill Bit	Stainless Steel - Tap- 2.0mm diameter- 1.5mm DRILL- 2.0mm SCREW	MicroAire/8053-119		Disposable	EA	5
E0-1175	Drill Bit	Stainless Steel - Tap- 2.7mm diameter- 2.0mm DRILL- 2.7mm SCREW	MicroAire/8053-126		Disposable	EA	5
E0-1176	Drill Bit	Stainless Steel - Tap- 3.5mm diameter- 2.5mm DRILL- 3.5mm SCREW	MicroAire/8053-132		Disposable	EA	5
E0-1177	Drill Bit	Stainless Steel - Tap- 3.5mm diameter- 2.5mm DRILL- 3.5mm SCREW- Non-Sterile	MicroAire/8053-132NS		Disposable	EA	5
E0-1178	Drill Bit	Stainless Steel - Tap- 3.5mm diameter-X-Long - 2.5mm DRILL- 3.5mm SCREW	MicroAire/8053-133		Disposable	EA	5
E0-1179	Drill Bit	Stainless Steel - Tap- 3.5mm diameter-X-Long - 2.5mm DRILL- 3.5mm SCREW- Non-Sterile	MicroAire/8053-133NS		Disposable	EA	5
E0-1180	Drill Bit	Stainless Steel - Tap- 4.0mm diameter- 2.5mm DRILL- 4.0mm SCREW	MicroAire/8053-134		Disposable	EA	5
E0-1181	Drill Bit	Stainless Steel - Tap- 4.0mm diameter- 2.5mm DRILL- 4.0mm SCREW- Non-Sterile	MicroAire/8053-134NS		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1182	Drill Bit	Stainless Steel - Tap- 4.5mm diameter- 3.2mm DRILL- 4.5mm SCREW	MicroAire/8053-146		Disposable	EA	5
E0-1183	Drill Bit	Stainless Steel - Tap- 4.5mm diameter- 3.2mm DRILL- 4.5mm SCREW- Non-Sterile	MicroAire/8053-146NS		Disposable	EA	5
E0-1184	Drill Bit	Stainless Steel - Tap- 4.5mm diameter-X-Long - 3.2mm DRILL- 4.5mm SCREW	MicroAire/8053-148		Disposable	EA	5
E0-1185	Drill Bit	Stainless Steel - Tap- 4.5mm diameter-X-Long - 3.2mm DRILL- 4.5mm SCREW- Non-Sterile	MicroAire/8053-148NS		Disposable	EA	5
E0-1186	Drill Bit	Stainless Stell - Twist - Standard - 1.6mm diameter X 127mm	MicroAire/8054-010		Disposable	EA	5
E0-1187	Drill Bit	Stainless Stell - Twist - Standard - 1.6mm diameter X 127mm- Non-Sterile	MicroAire/8054-010NS		Disposable	EA	5
E0-1188	Drill Bit	Stainless Stell - Twist - Standard - 2mm diameter X 127mm	MicroAire/8054-011		Disposable	EA	5
E0-1189	Drill Bit	Stainless Stell - Twist - Standard - 2mm diameter X 127mm- Non-Sterile	MicroAire/8054-011NS		Disposable	EA	5
E0-1190	Drill Bit	Stainless Stell - Twist - Standard - 2.4mm diameter X 127mm	MicroAire/8054-012		Disposable	EA	5
E0-1191	Drill Bit	Stainless Stell - Twist - Standard - 2.4mm diameter X 127mm- Non-Sterile	MicroAire/8054-012NS		Disposable	EA	5
E0-1192	Drill Bit	Stainless Stell - Twist - Standard - 2.8mm diameter X 127mm	MicroAire/8054-013		Disposable	EA	5
E0-1193	Drill Bit	Stainless Stell - Twist - Standard - 2.8mm diameter X 127mm- Non-Sterile	MicroAire/8054-013NS		Disposable	EA	5
E0-1194	Drill Bit	Stainless Stell - Twist - Standard - 3.2mm diameter X 127mm	MicroAire/8054-014		Disposable	EA	5
E0-1195	Drill Bit	Stainless Stell - Twist - Standard - 3.2mm diameter X 127mm- Non-Sterile	MicroAire/8054-014NS		Disposable	EA	5
E0-1196	Drill Bit	Stainless Stell - Twist - Standard - 3.6mm diameter X 127mm	MicroAire/8054-015		Disposable	EA	5
E0-1197	Drill Bit	Stainless Stell - Twist - Standard - 3.6mm diameter X 127mm- Non-Sterile	MicroAire/8054-015NS		Disposable	EA	5
E0-1198	Drill Bit	Stainless Stell - Twist - Standard - 4mm diameter X 127mm	MicroAire/8054-016		Disposable	EA	5
E0-1199	Drill Bit	Stainless Stell - Twist - Standard - 4mm diameter X 127mm- Non-Sterile	MicroAire/8054-016NS		Disposable	EA	5
E0-1200	Drill Bit	Stainless Stell - Twist - Standard - 4.5mm diameter X 127mm	MicroAire/8054-017		Disposable	EA	5
E0-1201	Drill Bit	Stainless Stell - Twist - Standard - 4.5mm diameter X 127mm- Non-Sterile	MicroAire/8054-017NS		Disposable	EA	5
E0-1202	Drill Bit	Stainless Stell - Twist - Standard - 4.8mm diameter X 127mm	MicroAire/8054-018		Disposable	EA	5
E0-1203	Drill Bit	Stainless Stell - Twist - Standard - 4.8mm diameter X 127mm- Non-Sterile	MicroAire/8054-018NS		Disposable	EA	5
E0-1204	Drill Bit	Stainless Stell - Twist - Standard - 6.4mm diameter X 127mm	MicroAire/8054-019		Disposable	EA	5
E0-1205	Drill Bit	Stainless Stell - Twist - Standard - 6.4mm diameter X 127mm- Non-Sterile	MicroAire/8054-019NS		Disposable	EA	5
E0-1206	Drill Bit	Stainless Stell - Twist - Standard - 2mm diameter X 127mm	MicroAire/8054-020		Disposable	EA	5
E0-1207	Drill Bit	Stainless Stell - Twist - Standard - 2.8mm diameter X 127mm	MicroAire/8054-021		Disposable	EA	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1208	Drill Bit	Stainless Stell - Twist - Standard - 3.2mm diameter X 127mm	MicroAire/8054-022		Disposable	EA	5
E0-1209	Drill Bit	Stainless Stell - Twist - Standard - 3.6mm diameter X 127mm	MicroAire/8054-023		Disposable	EA	5
E0-1210	Drill Bit	Stainless Stell - Twist - Standard - 4.5mm diameter X 127mm	MicroAire/8054-024		Disposable	EA	5
E0-1211	Wire Drill	K-Wire- 0.7mm(.028")X 140mm(5.5")- Non-Sterile	MicroAire/1600-001NS		Disposable	6PK	5
E0-1212	Wire Drill	K-Wire- 0.9mm(.035")X 140mm(5.5")- Non-Sterile	MicroAire/1600-002NS		Disposable	6PK	5
E0-1213	Wire Drill	K-Wire- 1.1mm(.045")X 140mm(5.5")- Non-Sterile	MicroAire/1600-003NS		Disposable	6PK	5
E0-1214	Wire Drill	K-Wire- 1.6mm(.062")X 140mm(5.5")- Non-Sterile	MicroAire/1600-004NS		Disposable	6PK	5
E0-1215	Wire Guard	Wire Guide- Red - Compatible with 0.7mm(0.028") Wire- Non-Sterile	MicroAire/1600-009NS		Disposable	6PK	5
E0-1216	Wire Guard	Wire Guide- Blue- Compatible with 0.9mm(0.035") Wire- Non-Sterile	MicroAire/1600-010NS		Disposable	6PK	5
E0-1217	Wire Guard	Wire Guide- Yellow- Compatible with 1.1mm(0.045") Wire- Non-Sterile	MicroAire/1600-011NS		Disposable	6PK	5
E0-1218	Wire Guard	Wire Guide - Black - Compatible with 1.6mm(0.062") Wire- Non-Sterile	MicroAire/1600-012NS		Disposable	6PK	5
E0-1219	Wire Guard	Wire Guide - White- Compatible with 1.4mm(0.054") Wire- Non-Sterile	MicroAire/1600-013NS		Disposable	6PK	5
E0-1220	Wire Drill	K-Wire - Includes Wire Guide- 0.7mm(0.028") X 140mm(5.5)	MicroAire/1600-020		Disposable	6PK	5
E0-1221	Wire Drill	K-Wire - Includes Wire Guide- 0.9mm(0.035") X 140mm(5.5)	MicroAire/1600-021		Disposable	6PK	5
E0-1222	Wire Drill	K-Wire - Includes Wire Guide- 1.1mm(0.045") X 140mm(5.5)	MicroAire/1600-022		Disposable	6PK	5
E0-1223	Wire Drill	K-Wire - Includes Wire Guide- 1.6mm(0.062") X 140mm(5.5)	MicroAire/1600-023		Disposable	6PK	5
E0-1224	Wire Drill	K-Wire .028" X 4"- Trocar points at both ends	MicroAire/1600-428		Disposable	6PK	5
E0-1225	Wire Drill	K-Wire .028" X 4"- Trocar points at both ends- Non-Sterile	MicroAire/1600-428NS		Disposable	6PK	5
E0-1226	Wire Drill	K-Wire .035" X 4"- Trocar points at both ends	MicroAire/1600-435		Disposable	6PK	5
E0-1227	Wire Drill	K-Wire .035" X 4"- Trocar points at both ends- Non-Sterile	MicroAire/1600-435NS		Disposable	6PK	5
E0-1228	Wire Drill	K-Wire .035" X 4" - Threaded - Trocar points at both ends	MicroAire/1600-435T		Disposable	6PK	5
E0-1229	Wire Drill	K-Wire .035" X 4" - Threaded - Trocar points at both ends- Non-Sterile	MicroAire/1600-435TNS		Disposable	6PK	5
E0-1230	Wire Drill	K-Wire .045" X 4"- Trocar points at both ends	MicroAire/1600-445		Disposable	6PK	5
E0-1231	Wire Drill	K-Wire .045" X 4"- Trocar points at both ends- Non-Sterile	MicroAire/1600-445NS		Disposable	6PK	5
E0-1232	Wire Drill	K-Wire .045" X 4" - Threaded- Trocar points at both ends	MicroAire/1600-445T		Disposable	6PK	5
E0-1233	Wire Drill	K-Wire .045" X 4" - Threaded- Trocar points at both ends- Non-Sterile	MicroAire/1600-445TNS		Disposable	6PK	5
E0-1234	Wire Drill	K-Wire .054" X 4"- Trocar points at both ends	MicroAire/1600-454		Disposable	6PK	5
E0-1235	Wire Drill	K-Wire .054" X 4"- Trocar points at both ends- Non-Sterile	MicroAire/1600-454NS		Disposable	6PK	5
E0-1236	Wire Drill	K-Wire .054" X 4" - Threaded- Trocar points at both ends	MicroAire/1600-454T		Disposable	6PK	5
E0-1237	Wire Drill	K-Wire .062" X 4"- Trocar points at both ends	MicroAire/1600-462		Disposable	6PK	5
E0-1238	Wire Drill	K-Wire .062" X 4"- Trocar points at both ends- Non-Sterile	MicroAire/1600-462NS		Disposable	6PK	5
E0-1239	Wire Drill	K-Wire .062" X 4" - Threaded- Trocar points at both ends	MicroAire/1600-462T		Disposable	6PK	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1240	Wire Drill	K-Wire .062" X 4" - Threaded- Trocar points at both ends- Non-Sterile	MicroAire/1600-462TNS		Disposable	6PK	5
E0-1241	Wire Drill	K-Wire .028" X 6"- Trocar points at both ends	MicroAire/1600-628		Disposable	6PK	5
E0-1242	Wire Drill	K-Wire .028" X 6"- Trocar points at both ends- Non-Sterile	MicroAire/1600-628NS		Disposable	6PK	5
E0-1243	Wire Drill	K-Wire .035" X 6"- Trocar points at both ends	MicroAire/1600-635		Disposable	6PK	5
E0-1244	Wire Drill	K-Wire .035" X 6"- Trocar points at both ends- Non-Sterile	MicroAire/1600-635NS		Disposable	6PK	5
E0-1245	Wire Drill	K-Wire .035" X 6" - Threaded- Trocar points at both ends	MicroAire/1600-635T		Disposable	6PK	5
E0-1246	Wire Drill	K-Wire .045" X 6"- Trocar points at both ends.	MicroAire/1600-645		Disposable	6PK	5
E0-1247	Wire Drill	K-Wire .045" X 6"- Trocar points at both ends- Non-Sterile	MicroAire/1600-645NS		Disposable	6PK	5
E0-1248	Wire Drill	K-Wire .045"X 6" - Threaded- Trocar points at both ends	MicroAire/1600-645T		Disposable	6PK	5
E0-1249	Wire Drill	K-Wire .045"X 6" - Threaded- Trocar points at both ends- Non-Sterile	MicroAire/1600-645TNS		Disposable	6PK	5
E0-1250	Wire Drill	K-Wire .054" X 6"- Trocar points at both ends	MicroAire/1600-654		Disposable	6PK	5
E0-1251	Wire Drill	K-Wire .054" X 6"- Trocar points at both ends- Non-Sterile	MicroAire/1600-654NS		Disposable	6PK	5
E0-1252	Wire Drill	K-Wire .054"X 6" - Threaded- Trocar points at both ends	MicroAire/1600-654T		Disposable	6PK	5
E0-1253	Wire Drill	K-Wire .054"X 6" - Threaded- Trocar points at both ends- Non-Sterile	MicroAire/1600-654TNS		Disposable	6PK	5
E0-1254	Wire Drill	K-Wire .062" X 6"- Trocar points at both ends	MicroAire/1600-662		Disposable	6PK	5
E0-1255	Wire Drill	K-Wire .062" X 6"- Trocar points at both ends	MicroAire/1600-662NS		Disposable	6PK	5
E0-1256	Wire Drill	K-Wire .062"X 6" - Threaded- Trocar points at both ends	MicroAire/1600-662T		Disposable	6PK	5
E0-1257	Wire Drill	K-Wire .062"X 6" - Threaded- Trocar points at both ends- Non-Sterile	MicroAire/1600-662TNS		Disposable	6PK	5
E0-1258	Wire Drill	K-Wire .035" X 9"- Trocar points at both ends	MicroAire/1600-935		Disposable	6PK	5
E0-1259	Wire Drill	K-Wire- .035" X 9"-Single Trocar Point	MicroAire/1600-9355		Disposable	6PK	5
E0-1260	Wire Drill	K-Wire- .035" X 9"-Single Trocar Point- Non-Sterile	MicroAire/1600-9355NS		Disposable	6PK	5
E0-1261	Wire Drill	K-Wire .035" X 9"- Trocar points at both ends- Non-Sterile	MicroAire/1600-935NS		Disposable	6PK	5
E0-1262	Wire Drill	K-Wire .035" X 9" - Threaded- Trocar points at both ends	MicroAire/1600-935T		Disposable	6PK	5
E0-1263	Wire Drill	K-Wire .035" X 9" - Threaded- Trocar points at both ends- Non-Sterile	MicroAire/1600-935TNS		Disposable	6PK	5
E0-1264	Wire Drill	K-Wire .045" X 9"- Trocar points at both ends	MicroAire/1600-945		Disposable	6PK	5
E0-1265	Wire Drill	K-Wire- .045" X 9"-Single Trocar Point	MicroAire/1600-9455		Disposable	6PK	5
E0-1266	Wire Drill	K-Wire- .045" X 9"-Single Trocar Point- Non-Sterile	MicroAire/1600-9455NS		Disposable	6PK	5
E0-1267	Wire Drill	K-Wire .045" X 9"- Trocar points at both ends- Non-Sterile	MicroAire/1600-945NS		Disposable	6PK	5
E0-1268	Wire Drill	K-Wire .045" X 9" - Threaded- Trocar points at both ends	MicroAire/1600-945T		Disposable	6PK	5
E0-1269	Wire Drill	K-Wire .045" X 9" - Threaded- Trocar points at both ends- Non-Sterile	MicroAire/1600-945TNS		Disposable	6PK	5
E0-1270	Wire Drill	K-Wire .054" X 9"- Trocar points at both ends	MicroAire/1600-954		Disposable	6PK	5
E0-1271	Wire Drill	K-Wire- .054" X 9"-Single Trocar Point	MicroAire/1600-9545		Disposable	6PK	5
E0-1272	Wire Drill	K-Wire .054" X 9"- Trocar points at both ends- Non-Sterile	MicroAire/1600-954NS		Disposable	6PK	5
E0-1273	Wire Drill	K-Wire .054"X 9" - Threaded- Trocar points at both ends	MicroAire/1600-954T		Disposable	6PK	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1274	Wire Drill	K-Wire .054" X 9" - Threaded- Trocar points at both ends- Non-Sterile	MicroAire/1600-954TNS		Disposable	6PK	5
E0-1275	Wire Drill	K-Wire .062" X 9"- Trocar points at both ends	MicroAire/1600-962		Disposable	6PK	5
E0-1276	Wire Drill	K-Wire- .062" X 9"-Single Trocar Point	MicroAire/1600-9625		Disposable	6PK	5
E0-1277	Wire Drill	K-Wire- .062" X 9"-Single Trocar Point- Non-Sterile	MicroAire/1600-9625NS		Disposable	6PK	5
E0-1278	Wire Drill	K-Wire .062" X 9" - Trocar points at both ends- Non-Sterile	MicroAire/1600-962NS		Disposable	6PK	5
E0-1279	Wire Drill	K-Wire .062" X 9" - Threaded- Trocar points at both ends	MicroAire/1600-962T		Disposable	6PK	5
E0-1280	Wire Drill	K-Wire .062" X 9" - Threaded- Trocar points at both ends- Non-Sterile	MicroAire/1600-962TNS		Disposable	6PK	5
E0-1281	Protector	Pin- .028" (.7mm)- Yellow	MicroAire/1601-028		Disposable	4PK	5
E0-1282	Protector	Pin- .035" (.9mm)- Blue	MicroAire/1601-035		Disposable	4PK	5
E0-1283	Head	Head - Light Blue	MicroAire/1601-035NS		Disposable	EA	5
E0-1284	Protector	Pin- .045" (1.1mm)- White	MicroAire/1601-045		Disposable	4PK	5
E0-1285	Pin	Head - White	MicroAire/1601-045NS		Disposable	EA	5
E0-1286	Protector	Pin- .062" (1.6mm)- Green	MicroAire/1601-062		Disposable	4PK	5
E0-1287	Head	Head- GREEN	MicroAire/1601-062NS		Disposable	EA	5
E0-1288	Protector	Pin- 1/8" (3.2mm)- Lavendar	MicroAire/1601-180		Disposable	4PK	5
E0-1289	Protector	Pin- 3/16" (4.8mm)- Orange	MicroAire/1601-316		Disposable	4PK	5
E0-1290	Protector	Pin- 3/32" (2.4mm)- Grey	MicroAire/1601-332		Disposable	4PK	5
E0-1291	Protector	Pin- 5/32" (4.0mm) - Navy Blue	MicroAire/1601-532		Disposable	4PK	5
E0-1292	Protector	Pin- 5/64" (2.0mm)- Pink	MicroAire/1601-564		Disposable	4PK	5
E0-1293	Protector	Pin- 7/64" (2.8mm)- Tan	MicroAire/1601-764		Disposable	4PK	5
E0-1294	Protector	Pin- 9/64" (3.6mm)-Brown	MicroAire/1601-964		Disposable	4PK	5
E0-1295	Wire Guide	K-Wire- 0.7mm(0.028") X 102mm(4")	MicroAire/1604-128		Disposable	6PK	5
E0-1296	Wire Guide	K-Wire- 0.7mm(0.028") X 102mm(4")- Non-Sterile	MicroAire/1604-128NS		Disposable	6PK	5
E0-1297	Wire Guide	K-Wire- 0.9mm(0.035") X 102mm(4")	MicroAire/1604-135		Disposable	6PK	5
E0-1298	Wire Guide	K-Wire- 0.9mm(0.035") X 102mm(4")- Non-Sterile	MicroAire/1604-135NS		Disposable	6PK	5
E0-1299	Wire Guide	K-Wire- 1.1mm(0.045") X 102mm(4")	MicroAire/1604-145		Disposable	6PK	5
E0-1300	Wire Guide	K-Wire- 1.1mm(0.045") X 102mm(4")- Non-Sterile	MicroAire/1604-145NS		Disposable	6PK	5
E0-1301	Wire Guide	K-Wire- 1.4mm(0.054") X 102mm(4")	MicroAire/1604-154		Disposable	6PK	5
E0-1302	Wire Guide	K-Wire- 1.4mm(0.054") X 102mm(4")- Non-Sterile	MicroAire/1604-154NS		Disposable	6PK	5
E0-1303	Wire Guide	K-Wire- 1.6mm(0.062") X 102mm(4")	MicroAire/1604-162		Disposable	6PK	5
E0-1304	Wire Guide	K-Wire- 1.6mm(0.062") X 102mm(4")- Non-Sterile	MicroAire/1604-162NS		Disposable	6PK	5
E0-1305	Pin	Pin	MicroAire/1605-020NS		Disposable	10PK	5
E0-1306	Pin	Fine - Surgical- 5/64" X 9"- TrocarPoint on both ends	MicroAire/1620-109		Disposable	6PK	5
E0-1307	Pin	Fine - Surgical- 5/64" X 9"- TrocarPoint on both ends- Non-Sterile	MicroAire/1620-109NS		Disposable	6PK	5
E0-1308	Pin	Fine - Surgical- 5/64" X 9"- Threaded- TrocarPoint on both ends	MicroAire/1620-109T		Disposable	6PK	5
E0-1309	Pin	Fine - Surgical- 5/64" X 9"- Threaded- TrocarPoint on both ends- Non-Sterile	MicroAire/1620-109TNS		Disposable	6PK	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1310	Pin	Fine - Surgical- 5/64" X 9"- Trocar Points on one end	MicroAire/1620-509		Disposable	6PK	5
E0-1311	Pin	Fine - Surgical- 5/64" X 9"- Trocar Points on one end- Non-Sterile	MicroAire/1620-509NS		Disposable	6PK	5
E0-1312	Pin	Fine - Surgical- 5/64" X 9"- Threaded- Trocar Points on one end	MicroAire/1620-509T		Disposable	6PK	5
E0-1313	Pin	Fine - Surgical- 5/64" X 9"- Threaded- Trocar Points on one end- Non-Sterile	MicroAire/1620-509TNS		Disposable	6PK	5
E0-1314	Pin	Fine - Surgical- 5/64" X 9"-Diamond - Point on one end	MicroAire/1620-609		Disposable	6PK	5
E0-1315	Pin	Fine - Surgical- 5/64" X 9"-Diamond - Point on one end- Non-Sterile	MicroAire/1620-609NS		Disposable	6PK	5
E0-1316	Pin	Fine - Surgical- 5/64" X 9"- Threaded-Diamond - Point on one end	MicroAire/1620-609T		Disposable	6PK	5
E0-1317	Pin	Fine - Surgical- 5/64" X 9"- Threaded-Diamond - Point on one end- Non-Sterile	MicroAire/1620-609TNS		Disposable	6PK	5
E0-1318	Pin	Fine - Surgical- 5/64" X 9"-Diamond -Point on both ends	MicroAire/1620-709		Disposable	6PK	5
E0-1319	Pin	Fine - Surgical- 5/64" X 9"-Diamond -Point on both ends- Non-Sterile	MicroAire/1620-709NS		Disposable	6PK	5
E0-1320	Pin	Fine - Surgical- 5/64" X 9"- Threaded-Diamond -Point on both ends	MicroAire/1620-709T		Disposable	6PK	5
E0-1321	Pin	Fine - Surgical- 5/64" X 9"- Threaded-Diamond -Point on both ends- Non-Sterile	MicroAire/1620-709TNS		Disposable	6PK	5
E0-1322	Pin	Fine - Surgical- 3/32" X 9"- TrocarPoint on both ends	MicroAire/1624-109		Disposable	6PK	5
E0-1323	Pin	Fine - Surgical- 3/32" X 9"- TrocarPoint on both ends- Non-Sterile	MicroAire/1624-109NS		Disposable	6PK	5
E0-1324	Pin	Fine - Surgical- 3/32" X 9"- Threaded- TrocarPoint on both ends	MicroAire/1624-109T		Disposable	6PK	5
E0-1325	Pin	Fine - Surgical- 3/32" X 9"- Threaded- TrocarPoint on both ends- Non-Sterile	MicroAire/1624-109TNS		Disposable	6PK	5
E0-1326	Pin	Fine - Surgical- 3/32" X 9"- Trocar Points on one end	MicroAire/1624-509		Disposable	6PK	5
E0-1327	Pin	Fine - Surgical- 3/32" X 9"- Trocar Points on one end- Non-Sterile	MicroAire/1624-509NS		Disposable	6PK	5
E0-1328	Pin	Fine - Surgical 3/32" X 9"- Trocar Points on one end- Non-Sterile	MicroAire/1624-509NS25		Disposable	25PK	5
E0-1329	Pin	Fine - Surgical- 3/32" X 9"- Threaded- Trocar Points on one end	MicroAire/1624-509T		Disposable	6PK	5
E0-1330	Pin	Fine - Surgical- 3/32" X 9"- Threaded- Trocar Points on one end- Non-Sterile	MicroAire/1624-509TNS		Disposable	6PK	5
E0-1331	Pin	Fine - Surgical- 3/32" X 9"-Diamond - Point on one end	MicroAire/1624-609		Disposable	6PK	5
E0-1332	Pin	Fine - Surgical- 3/32" X 9"-Diamond - Point on one end- Non-Sterile	MicroAire/1624-609NS		Disposable	6PK	5
E0-1333	Pin	Fine - Surgical- 3/32" X 9"- Threaded-Diamond - Point on one end	MicroAire/1624-609T		Disposable	6PK	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1334	Pin	Fine - Surgical- 3/32" X 9"- Threaded-Diamond - Point on one end- Non-Sterile	MicroAire/1624-609TNS		Disposable	6PK	5
E0-1335	Pin	Fine - Surgical- 3/32" X 9"-Diamond -Point on both ends	MicroAire/1624-709		Disposable	6PK	5
E0-1336	Pin	Fine - Surgical- 3/32" X 9"-Diamond -Point on both ends- Non-Sterile	MicroAire/1624-709NS		Disposable	6PK	5
E0-1337	Pin	Fine - Surgical- 3/32" X 9"- Threaded-Diamond -Point on both ends	MicroAire/1624-709T		Disposable	6PK	5
E0-1338	Pin	Fine - Surgical- 3/32" X 9"- Threaded-Diamond -Point on both ends- Non-Sterile	MicroAire/1624-709TNS		Disposable	6PK	5
E0-1339	Pin	Fine - Surgical- 7/64" X 9"- TrocarPoint on both ends	MicroAire/1628-109		Disposable	6PK	5
E0-1340	Pin	Fine - Surgical- 7/64" X 9"- TrocarPoint on both ends- Non-Sterile	MicroAire/1628-109NS		Disposable	6PK	5
E0-1341	Pin	Fine - Surgical- 7/64" X 9"- Threaded- TrocarPoint on both ends	MicroAire/1628-109T		Disposable	6PK	5
E0-1342	Pin	Fine - Surgical- 7/64" X 9"- Threaded- TrocarPoint on both ends- Non-Sterile	MicroAire/1628-109TNS		Disposable	6PK	5
E0-1343	Pin	Fine - Surgical- 7/64" X 9"- Trocar Points on one end	MicroAire/1628-509		Disposable	6PK	5
E0-1344	Pin	Fine - Surgical- 7/64" X 9"- Trocar Points on one end- Non-Sterile	MicroAire/1628-509NS		Disposable	6PK	5
E0-1345	Pin	Fine - Surgical- 7/64" X 9"- Threaded- Trocar Points on one end	MicroAire/1628-509T		Disposable	6PK	5
E0-1346	Pin	Fine - Surgical- 7/64" X 9"- Threaded- Trocar Points on one end- Non-Sterile	MicroAire/1628-509TNS		Disposable	6PK	5
E0-1347	Pin	Fine - Surgical- 7/64" X 9"-Diamond - Point on one end	MicroAire/1628-609		Disposable	6PK	5
E0-1348	Pin	Fine - Surgical- 7/64" X 9"-Diamond - Point on one end- Non-Sterile	MicroAire/1628-609NS		Disposable	6PK	5
E0-1349	Pin	Fine - Surgical- 7/64" X 9"- Threaded-Diamond - Point on one end	MicroAire/1628-609T		Disposable	6PK	5
E0-1350	Pin	Fine - Surgical- 7/64" X 9"- Threaded-Diamond - Point on one end- Non-Sterile	MicroAire/1628-609TNS		Disposable	6PK	5
E0-1351	Pin	Fine - Surgical- 7/64" X 9"-Diamond -Point on both ends	MicroAire/1628-709		Disposable	6PK	5
E0-1352	Pin	Fine - Surgical- 7/64" X 9"-Diamond -Point on both ends- Non-Sterile	MicroAire/1628-709NS		Disposable	6PK	5
E0-1353	Pin	Fine - Surgical- 7/64" X 9"- Threaded-Diamond -Point on both ends	MicroAire/1628-709T		Disposable	6PK	5
E0-1354	Pin	Fine - Surgical- 7/64" X 9"- Threaded-Diamond -Point on both ends- Non-Sterile	MicroAire/1628-709TNS		Disposable	6PK	5
E0-1355	Pin	Fine - Surgical- 1/8" X 9"- TrocarPoint on both ends	MicroAire/1632-109		Disposable	6PK	5
E0-1356	Pin	Fine - Surgical- 1/8" X 9"- TrocarPoint on both ends- Non-Sterile	MicroAire/1632-109NS		Disposable	6PK	5
E0-1357	Pin	Fine - Surgical- 1/8" X 9"- Threaded- TrocarPoint on both ends	MicroAire/1632-109T		Disposable	6PK	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1358	Pin	Fine - Surgical- 1/8" X 9"- Threaded- TrocarPoint on both ends- Non-Sterile	MicroAire/1632-109TNS		Disposable	6PK	5
E0-1359	Pin	Fine - Surgical- 1/8" X 9"- Trocar Points on one end	MicroAire/1632-509		Disposable	6PK	5
E0-1360	Pin	Fine - Surgical- 1/8" X 9"- Trocar Points on one end- Non-Sterile	MicroAire/1632-509NS		Disposable	6PK	5
E0-1361	Pin	Fine - Surgical- 1/8" X 9"- Threaded- Trocar Points on one end	MicroAire/1632-509T		Disposable	6PK	5
E0-1362	Pin	Fine - Surgical- 1/8" X 9"- Threaded- Trocar Points on one end- Non-Sterile	MicroAire/1632-509TNS		Disposable	6PK	5
E0-1363	Pin	Fine - Surgical- 1/8" X 9"-Diamond - Point on one end	MicroAire/1632-609		Disposable	6PK	5
E0-1364	Pin	Fine - Surgical- 1/8" X 9"-Diamond - Point on one end- Non-Sterile	MicroAire/1632-609NS		Disposable	6PK	5
E0-1365	Pin	Fine - Surgical- 1/8" X 9"- Threaded-Diamond - Point on one end	MicroAire/1632-609T		Disposable	6PK	5
E0-1366	Pin	Fine - Surgical- 1/8" X 9"- Threaded-Diamond - Point on one end- Non-Sterile	MicroAire/1632-609TNS		Disposable	6PK	5
E0-1367	Pin	Fine - Surgical- 1/8" X 9"-Diamond -Point on both ends	MicroAire/1632-709		Disposable	6PK	5
E0-1368	Pin	Fine - Surgical- 1/8" X 9"-Diamond -Point on both ends- Non-Sterile	MicroAire/1632-709NS		Disposable	6PK	5
E0-1369	Pin	Fine - Surgical- 1/8" X 9"- Threaded-Diamond -Point on both ends	MicroAire/1632-709T		Disposable	6PK	5
E0-1370	Pin	Fine - Surgical- 1/8" X 9"- Threaded-Diamond -Point on both ends- Non-Sterile	MicroAire/1632-709TNS		Disposable	6PK	5
E0-1371	Pin	Fine - Surgical- 9/64" X 9"- TrocarPoint on both ends	MicroAire/1636-109		Disposable	6PK	5
E0-1372	Pin	Fine - Surgical- 9/64" X 9"- TrocarPoint on both ends- Non-Sterile	MicroAire/1636-109NS		Disposable	6PK	5
E0-1373	Pin	Fine - Surgical- 9/64" X 9"- Threaded- TrocarPoint on both ends	MicroAire/1636-109T		Disposable	6PK	5
E0-1374	Pin	Fine - Surgical- 9/64" X 9"- Threaded- TrocarPoint on both ends- Non-Sterile	MicroAire/1636-109TNS		Disposable	6PK	5
E0-1375	Pin	Fine - Surgical- 9/64" X 9"- Trocar Points on one end	MicroAire/1636-509		Disposable	6PK	5
E0-1376	Pin	Fine - Surgical- 9/64" X 9"- Trocar Points on one end- Non-Sterile	MicroAire/1636-509NS		Disposable	6PK	5
E0-1377	Pin	Fine - Surgical- 9/64" X 9"- Threaded- Trocar Points on one end	MicroAire/1636-509T		Disposable	6PK	5
E0-1378	Pin	Fine - Surgical- 9/64" X 9"- Threaded- Trocar Points on one end- Non-Sterile	MicroAire/1636-509TNS		Disposable	6PK	5
E0-1379	Pin	Fine - Surgical- 9/64" X 9"-Diamond - Point on one end	MicroAire/1636-609		Disposable	6PK	5
E0-1380	Pin	Fine - Surgical- 9/64" X 9"-Diamond - Point on one end- Non-Sterile	MicroAire/1636-609NS		Disposable	6PK	5
E0-1381	Pin	Fine - Surgical- 9/64" X 9"- Threaded-Diamond - Point on one end	MicroAire/1636-609T		Disposable	6PK	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1382	Pin	Fine - Surgical- 9/64" X 9"- Threaded-Diamond - Point on one end- Non-Sterile	MicroAire/1636-609TNS		Disposable	6PK	5
E0-1383	Pin	Fine - Surgical- 9/64" X 9"-Diamond -Point on both ends	MicroAire/1636-709		Disposable	6PK	5
E0-1384	Pin	Fine - Surgical- 9/64" X 9"-Diamond -Point on both ends- Non-Sterile	MicroAire/1636-709NS		Disposable	6PK	5
E0-1385	Pin	Fine - Surgical- 9/64" X 9"- Threaded-Diamond -Point on both ends	MicroAire/1636-709T		Disposable	6PK	5
E0-1386	Pin	Fine - Surgical- 9/64" X 9"- Threaded-Diamond -Point on both ends- Non-Sterile	MicroAire/1636-709TNS		Disposable	6PK	5
E0-1387	Pin	Fine - Surgical- 5/32" X 9"- TrocarPoint on both ends	MicroAire/1640-109		Disposable	6PK	5
E0-1388	Pin	Fine - Surgical- 5/32" X 9"- TrocarPoint on both ends- Non-Sterile	MicroAire/1640-109NS		Disposable	6PK	5
E0-1389	Pin	Fine - Surgical- 5/32" X 9"- Threaded- TrocarPoint on both ends	MicroAire/1640-109T		Disposable	6PK	5
E0-1390	Pin	Fine - Surgical- 5/32" X 9"- Threaded- TrocarPoint on both ends- Non-Sterile	MicroAire/1640-109TNS		Disposable	6PK	5
E0-1391	Pin	Fine - Surgical- 5/32" X 9"- Trocar Points on one end	MicroAire/1640-509		Disposable	6PK	5
E0-1392	Pin	Fine - Surgical- 5/32" X 9"- Trocar Points on one end- Non-Sterile	MicroAire/1640-509NS		Disposable	6PK	5
E0-1393	Pin	Fine - Surgical- 5/32" X 9"- Threaded- Trocar Points on one end	MicroAire/1640-509T		Disposable	6PK	5
E0-1394	Pin	Fine - Surgical- 5/32" X 9"- Threaded- Trocar Points on one end- Non-Sterile	MicroAire/1640-509TNS		Disposable	6PK	5
E0-1395	Pin	Fine - Surgical- 5/32" X 9"-Diamond - Point on one end	MicroAire/1640-609		Disposable	6PK	5
E0-1396	Pin	Fine - Surgical- 5/32" X 9"-Diamond - Point on one end- Non-Sterile	MicroAire/1640-609NS		Disposable	6PK	5
E0-1397	Pin	Fine - Surgical- 5/32" X 9"- Threaded-Diamond - Point on one end	MicroAire/1640-609T		Disposable	6PK	5
E0-1398	Pin	Fine - Surgical- 5/32" X 9"- Threaded-Diamond - Point on one end- Non-Sterile	MicroAire/1640-609TNS		Disposable	6PK	5
E0-1399	Pin	Fine - Surgical- 5/32" X 9"-Diamond -Point on both ends	MicroAire/1640-709		Disposable	6PK	5
E0-1400	Pin	Fine - Surgical- 5/32" X 9"-Diamond -Point on both ends- Non-Sterile	MicroAire/1640-709NS		Disposable	6PK	5
E0-1401	Pin	Fine - Surgical- 5/32" X 9"- Threaded-Diamond -Point on both ends	MicroAire/1640-709T		Disposable	6PK	5
E0-1402	Pin	Fine - Surgical- 5/32" X 9"- Threaded-Diamond -Point on both ends- Non-Sterile	MicroAire/1640-709TNS		Disposable	6PK	5
E0-1403	Pin	Fine - Surgical- 3/16" X 9"- TrocarPoint on both ends	MicroAire/1648-109		Disposable	6PK	5
E0-1404	Pin	Fine - Surgical- 3/16" X 9"- TrocarPoint on both ends- Non-Sterile	MicroAire/1648-109NS		Disposable	6PK	5
E0-1405	Pin	Fine - Surgical- 3/16" X 9"- Threaded- TrocarPoint on both ends	MicroAire/1648-109T		Disposable	6PK	5

CLIN	Item Description	Salient Characteristics	Brand Name OR EQUAL	Installation Required	Disposable/Non-Disposable	Unit	Quantity
E0-1406	Pin	Fine - Surgical- 3/16" X 9"- Threaded- TrocarPoint on both ends- Non-Sterile	MicroAire/1648-109TNS		Disposable	6PK	5
E0-1407	Pin	Fine - Surgical	MicroAire/1648-508NS2		Disposable	2PK	5
E0-1408	Pin	Fine - Surgical- 3/16" X 9"- Trocar Points on one end	MicroAire/1648-509		Disposable	6PK	5
E0-1409	Pin	Fine - Surgical- 3/16" X 9"- Trocar Points on one end- Non-Sterile	MicroAire/1648-509NS		Disposable	6PK	5
E0-1410	Pin	Fine - Surgical- 3/16" X 9"- Threaded- Trocar Points on one end	MicroAire/1648-509T		Disposable	6PK	5
E0-1411	Pin	Fine - Surgical- 3/16" X 9"- Threaded- Trocar Points on one end- Non-Sterile	MicroAire/1648-509TNS		Disposable	6PK	5
E0-1412	Pin	Fine - Surgical- 3/16" X 9"-Diamond - Point on one end	MicroAire/1648-609		Disposable	6PK	5
E0-1413	Pin	Fine - Surgical- 3/16" X 9"-Diamond - Point on one end- Non-Sterile	MicroAire/1648-609NS		Disposable	6PK	5
E0-1414	Pin	Fine - Surgical- 3/16" X 9"- Threaded-Diamond - Point on one end	MicroAire/1648-609T		Disposable	6PK	5
E0-1415	Pin	Fine - Surgical- 3/16" X 9"- Threaded-Diamond - Point on one end- Non-Sterile	MicroAire/1648-609TNS		Disposable	6PK	5
E0-1416	Pin	Fine - Surgical- 3/16" X 9"-Diamond -Point on both ends	MicroAire/1648-709		Disposable	6PK	5
E0-1417	Pin	Fine - Surgical- 3/16" X 9"-Diamond -Point on both ends- Non-Sterile	MicroAire/1648-709NS		Disposable	6PK	5
E0-1418	Pin	Fine - Surgical- 3/16" X 9"- Threaded-Diamond -Point on both ends	MicroAire/1648-709T		Disposable	6PK	5
E0-1419	Pin	Fine - Surgical- 3/16" X 9"- Threaded-Diamond -Point on both ends- Non-Sterile	MicroAire/1648-709TNS		Disposable	6PK	5