

Wash D.C.

**McKESSON**  
Empowering Healthcare

VISN5

Quotation: H11VSN013-01  
Revision Date: September 17, 2010  
Project Name: HMI  
Facility Name: VA Washington  
Prepared By: Tony Douglas  
Contact Number: (804) 464-2851

---

### Horizon Medical Imaging Quotation Summary

Primary Quotation #H11VSN013-01

Installation and Professional Services

Professional Services Discount

Net Professional Services

**Subtotal Primary Quotation**

\$

**Net Total**

\$

The pricing set forth in this proposal represents McKesson's complete proposal for the products and or customer's facilities set forth herein (the "Pricing Proposal"), regardless of other proposals made by McKesson either simultaneously with this Pricing Proposal or otherwise regarding additional products or facilities that are not set forth herein.

McKesson has made no warranty or commitment with regard to any functionality not Generally Available as of the date of the Quotation or if the Quotation is included in a contract, then as of the Contract Supplement Effective Date

One time pricing concession based on McKesson's first VISN wide installation – establishing an HMI VA reference base for future sales.

This proposal requires some components from quotations H11VSN008 and H11VSN009 in order to provide a complete functional PACS system for the facility named herein.

Note that McKesson has a valid DIN PACS II Contract: Contract Number: SPM2D1-08-D-83402-13-08

Note that McKesson has responded to the DIN PACS III Solicitation with Solicitation Number: SPM2D1-10-R0011

- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION	UNIT PRICE <i>Includes Applicable Discounts</i>	TOTAL PRICE <i>Qty x Unit Price</i>
<b>Professional Services</b>					
1		SVC100 74006206	INBOUND ORDERS, ADT AND REPORT - INTEGRATION SERVICES  Professional services to interface Horizon Medical Imaging with an existing HL7 Compliant System for HMI to receive Orders, ADT Messages and textual reports		
2		SVC130 74010127	OUTBOUND REPORTS - INTEGRATION SERVICES  Professional services to interface Horizon Medical Imaging with an existing HL7 Compliant System for HMI to send textual reports		
3		SVC140 74008117	OUTBOUND STATUS UPDATES - INTEGRATION SERVICES  Professional services to interface Horizon Medical Imaging with an existing HL7 Compliant System for HMI to send updates when a study status changes.		
4		SER140 74010582	DICTATION INTERFACE SERVICES  Includes interface modifications by McKesson. Limited to the use of data which is already resident to the HMI database. Additional charges will apply for customization requests beyond this scope.		
5		CRX340 74006142	REPORT CONSULTING AND DESIGN SERVICES  Includes: Post-Installation Report Consulting Service (80 hours) to create custom modifications to the standard report forms. Report design projects exceeding 80 hours will be defined and quoted as custom software development projects, with pricing to be based on \$1600 per day (includes time and materials). Travel to customer site is not usually necessary, and is not included. If required, travel costs will be at the customer's expense.  All changes, customizations, modifications or improvements made or developed with regard to the Software pursuant to the Services provided as described herein shall remain the property of McKesson. Customer represents that any specifications or similar information provided to McKesson for purposes of performing the Services are not the proprietary information of any other party.		
6		TRV130 74-TBD	VTRIP APPLICATIONS TRAINING WITH DISPLAY PROTOCOL DESIGN AND BUILD - 3 DAY  Includes: Three consecutive days of on-site training, typically this would be sufficient for training of up to ten Radiologists.		
7		ASM100 75003140	IMPLEMENTATION AND APPLICATIONS COURSE  HMI Implementation and Applications is a pre-installation course intended to provide participants with the skills needed to implement HMI PACS. This includes workflow changes, infrastructure requirements, training and planning for PACS installations as well as intensive applications training.  The intended audience is PACS Administrators, IT personnel and Superusers who will be involved in the implementation and administration of the PACS system.  Available at McKesson's Richmond, British Columbia office. Travel expenses are not included.		

Quotation: H11VSN013 01  
Project: rHMI  
Site: VA Washington

Contact Name: Tony Douglas  
Contact Number: (804) 464-2851  
Revision Date: September 17, 2010

- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION	UNIT PRICE <i>Includes Applicable Discounts</i>	TOTAL PRICE <i>Qty x Unit Price</i>
8		ASM110 75003090	<p>ADVANCED APPLICATIONS COURSE</p> <p>HMI Advanced Applications is a post-installation advanced applications training course intended to reinforce and develop the skills necessary to support HMI PACS. Sessions will discuss and resolve site specific workflow and applications configuration concerns in a round-table format. The emphasis is on scenario-based learning and a large portion of the course will be devoted to hands on use of HMI PACS in a lab setting.</p> <p>The intended audience is PACS Administrators, IT personnel and Superusers.</p> <p>Prerequisites: Participants should have approximately 2 - 3 months of experience using HMI PACS as well as a strong understanding of their facility's clinical/radiology workflow.</p> <p>Available at McKesson's Richmond, British Columbia office. Travel expenses are not included.</p>		
9		ASM120 75003091	<p>ARCHITECTURE AND SUPPORT</p> <p>HMI Architecture and Support is a post-installation technical course intended to provide comprehensive introduction to HMI PACS architecture and the advanced skills required to streamline the PACS problem resolution process. The course will enable students to use system tools and logs to help diagnose problems.</p> <p>The intended audience is advanced technical users such as PACS Administrators, IT personnel and Superusers who will be responsible for the HMI PACS system internal processes and data flow and technical support.</p> <p>Prerequisites: This advanced course is specifically designed for the experienced PACS user and systems administrator. Students are expected to have the following prerequisites:</p> <ul style="list-style-type: none"> <li>- Completion of the HMI Implementation and Application training course</li> <li>- Intermediate level of Windows administration</li> <li>- Experience working with command line tools to administrate servers and workstations</li> <li>- LAN network experience</li> <li>- DICOM Basics</li> <li>- HL-7 Basics</li> </ul> <p>Available at McKesson's Richmond, British Columbia office. Travel expenses are not included.</p>		
10		INS100 74007667	<p>Installation and Training</p> <p>Includes:</p> <ul style="list-style-type: none"> <li>- Implementation and/or training services for implementation of this quotation</li> <li>- Any travel and onsite living expenses for McKesson personnel providing these services</li> <li>- Shipping and handling charges for shipment from McKesson to customer site</li> </ul>		
				<b>Grand Total</b>	

---

### McKesson Professional Consulting Services Quotation Summary

---

Primary Quotation #P11VSN002

Installation and Professional Services

---

<b>Subtotal Primary Quotation</b>	<b>\$</b>
-----------------------------------	-----------

---

---

<b>Net Total</b>	<b>\$</b>
------------------	-----------

---

The pricing set forth in this proposal represents McKesson's complete proposal for the products and or customer's facilities set forth herein (the "Pricing Proposal"), regardless of other proposals made by McKesson either simultaneously with this Pricing Proposal or otherwise regarding additional products or facilities that are not set forth herein.

McKesson has made no warranty or commitment with regard to any functionality not Generally Available as of the date of the Quotation or if the Quotation is included in a contract, then as of the Contract Supplement Effective Date.

- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION	UNIT PRICE <i>Includes Applicable Discounts</i>	TOTAL PRICE <i>Qty x Unit Price</i>
<b>Professional Services</b>					
1		PCS810 74023755	<p>MIGRATION SERVICES OF DICOM LEGACY PACS DATA            277,500 Studies to be migrated</p> <ul style="list-style-type: none"> <li>- Professional Services for migration of data from one(1) database to one(1) database</li> <li>- Includes basic data correction</li> <li>- Requires further investigation of existing legacy system's database and images to ensure successful migration.</li> <li>- Requires availability of DICOM Q/R SCU/SCP and DICOM Store SCU services on the Legacy system.</li> <li>- Migration of data above and beyond the specified scope will require additional charges at McKesson's standard time and materials rate.</li> </ul> <p>Customer is responsible for providing a complete list of data matching criteria, and for necessary manual matching/QC on studies identified as exceptions.</p> <p>Migration from the existing VistA PACS system.</p>		
<i>Subtotal Professional Services</i>					
<i>Grand Total</i>					

Quotation: H11VSN009-01  
Revision Date: September 17, 2010  
Project Name: HMI  
Facility Name: VA Washington  
Prepared By: Tony Douglas  
Contact Number: (804) 464-2851

## Horizon Medical Imaging Quotation Summary

Quotation #H11VSN009-01

Horizon Medical Imaging Software, Third Party Software and Equipment and Maintenance

### PACS Fee Per Study

93000 exams per year \$ \$

Option Year 2 Cost Per Study - PACS		
93000 exams per year	\$	\$
Option Year 3 Cost Per Study - PACS		
93000 exams per year	\$	\$
Option Year 4 Cost Per Study - PACS		
93000 exams per year	\$	\$
Option Year 5 Cost Per Study - PACS		
93000 exams per year	\$	\$

The pricing set forth in this proposal represents McKesson's complete proposal for the products and or customer's facilities set forth herein (the "Pricing Proposal"), regardless of other proposals made by McKesson either simultaneously with this Pricing Proposal or otherwise regarding additional products or facilities that are not set forth herein.

McKesson has made no warranty or commitment with regard to any functionality not Generally Available as of the date of the Quotation or if the Quotation is included in a contract, then as of the Contract Supplement Effective Date

One time pricing concession based on McKesson's first VISN wide installation – establishing an HMI VA reference base for future sales.

This proposal requires some components from quotations H11VSN008 and H11VSN009 in order to provide a complete functional PACS system for the facility named herein.

Note that McKesson has a valid DIN PACS II Contract: Contract Number: SPM2D1-08-D-83402-13-08

Note that McKesson has responded to the DIN PACS III Solicitation with Solicitation Number: SPM2D1-10-R0011



VISN5

Quotation: H11VSN009-01  
 Project: HMI  
 Site: VA Washington

Contact Name: Tony Douglas  
 Contact Number: (804) 464-2851  
 Revision Date: September 17, 2019

- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION		
<b>Horizon Medical Imaging Software</b>					
1		HMI100 72008731	<p>HORIZON MEDICAL IMAGING VOLUME LICENSE            Licensed for 93,000 Exams per year.</p> <p>Horizon Medical Imaging Software License, including the following functionality</p> <ul style="list-style-type: none"> <li>- Study Management: Image Acquisition, Online/Archive Storage and study tracking</li> <li>- Diagnostic Review: Horizon Rad Station - Advanced (HRS-A) and/or Horizon Rad Station - Distributed (HRS-D) for purposes of study interpretation. Includes integration to qualified dictation packages and 3rd party diagnostic tools (ie 3D)</li> <li>- Technologist Workflow: Horizon Rad Station - Advanced (HRS-A) and/or Horizon Rad Station - Distributed (HRS-D) for QA, Printing, Mark-up, Document Scanning, etc</li> <li>- Study Distribution: CD Publishing and Web-distribution tools including EMR integration</li> </ul> <p>Horizon Medical Imaging Software is licensed for the specified annual volume of Exams processed at the named facilities. Additional license fees are required to process higher exam volumes.</p>		
2		CRR100 72022392	<p>HORIZON MEDICAL IMAGING CRITICAL RESULTS REPORTING            Licensed for 93,000 Exams per year</p> <p>Critical Results Reporting provides the ability to:</p> <ul style="list-style-type: none"> <li>- Flag studies for critical result findings</li> <li>- Log communication acknowledgment</li> <li>- Generate aggregate statistical reports over an interval of time</li> </ul> <p>Critical Results Reporting Volume Software is licensed for the specified annual volume of exams processed at the named facilities. Additional license fees are required to process higher exam volumes</p> <p>Requires HMI 11.7 or higher</p>		



VISN5

Quotation: H11VSN009-01  
 Project HMI  
 Site VA Washington

Contact Name Tony Douglas  
 Contact Number (804) 464-2851  
 Revision Date September 17, 2010

- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION		
3		VTR005 72020388	<p>HORIZON MEDICAL IMAGING VTRIP SOFTWARE MODULE - 5 CONCURRENT USER LICENSES</p> <p>McKesson's Variable Thickness Regional Intensity Projection (VTRIP) software module allows for more efficient interpretation of thin-slice CT and MR studies by radiologists with clearer visual communication of key findings to referring physicians, by providing new visualization techniques embedded directly in the Horizon Rad Station user interface (both HRS-Advanced and HRS-Distributed).</p> <p>Features include:</p> <ul style="list-style-type: none"> <li>- Dynamic oblique Multi-Planar Reformats (MPR)</li> <li>- Dynamic sliding slab visualization, with interactive slab thickness adjustment</li> <li>- Multiple intensity projection options (MIP, AveIP, and minIP)</li> <li>- Display protocol support for the above</li> <li>- 3D Triangulation with instant zoom in/out toggle</li> <li>- Automatic documentation of key findings via intelligent derived views (HRS-Advanced only)</li> </ul> <p>Requires Horizon Medical Imaging 11.5.1 and workstation equipment that meets the minimum specification for HRS-D</p> <p>Please note that VTRIP will not be compatible with some future versions of Horizon Medical Imaging; please contact McKesson for details</p>		
<b>Third Party Software</b>					
4		SEG503 72019882	<p>SEGAMI OASIS - 3 CONCURRENT USER LICENSES</p> <p>Oasis Nuclear Medicine and PET/CT application suite</p> <p>Includes:</p> <ul style="list-style-type: none"> <li>Cardiology Applications</li> <li>Neurology Applications</li> <li>MultiModality Applications</li> <li>Planar Applications</li> <li>Pulmonology (Lung) Applications</li> <li>Generic Applications</li> <li>2 days of onsite training</li> <li>Integration with HMI</li> </ul> <p>Please note that Segami Oasis has specific workstation equipment requirements; configuration to be validated by McKesson.</p>		
5	1	SEG611 72019884	<p>CEDARS FOR SEGAMI OASIS - 1 CONCURRENT USER LICENSE</p> <p>Cedars QGS/QPS/ Companion - review and analysis of myocardial perfusion studies</p>		



- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION		
6		MMV301 72022170	MEDVIEW CLIENT (PACS PLUG-IN) CONCURRENT W/ CACHE SERVER SW 1ST LICENSE  Features include: -Display of dynamic, tomographic, MRP, MIP -Gated Cardiac review with CINE and rest/stress alignment -PET/CT review with fused displays and SUV/HU support with contrast presets -SPECT/PET/CT Current and Prior Volume Registration -MUGA review and analysis of planar gated blood pool studies -Save DICOM result back to HMI		
7		MMV302 72022171	MEDVIEW CLIENT (PACS PLUG-IN) CONCURRENT W/ CACHE SERVER SW 2ND AND SUB LICENSE  Features include: -Display of dynamic, tomographic, MRP, MIP -Gated Cardiac review with CINE and rest/stress alignment -PET/CT review with fused displays and SUV/HU support with contrast presets -SPECT/PET/CT Current and Prior Volume Registration -MUGA review and analysis of planar gated blood pool studies -Save DICOM result back to HMI		
8		MMV610 72022172	CORRIDOR 4DM MI (SPECT-PET) SINGLE SEAT LICENSE  Corridor 4DM MI by INVIA Medical Solutions, provides quantitative analysis of tomographic Nuclear SPECT and PETcardiac perfusion and blood pool studies. (The program does not acquire images or make diagnoses. Corridor 4DM-MI supports both PET and SPECT data. Data must be reconstructed & reoriented using the standards set by SNM, ACC and ASNC).		
9		ORC101 72008210	ORACLE SOFTWARE LICENSE  Oracle Standard Edition ONE license for 1x CPU		
<b>Equipment</b>					
10		WKS100HXP 80001506	HORIZON RAD STATION - HP  This chassis is not be used with 5MP Displays  Hardware: Workstation with a 100/1000BaseT NIC and Windows XP Software License.		
11		CDS610B 80004072	DISPLAY SYSTEM: BARCO CORONIS 6MP/ MDCC6130 DL  30" Color 6 MegaPixel Diagnostic Flat Panel with auto DICOM calibration, embedded front light sensor and protective glass.		

- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION		
12		ADI210BB 80004560	AUXILIARY DISPLAY SYSTEM BARCO MDRC-2120 20.1" Colour 1600 X 1200 Clinical Review Display with DICOM calibration		
13		WKS100HXP 80001506	HORIZON RAD STATION - HP This chassis is not be used with 5MP Displays Hardware: Workstation with a 100/1000BaseT NIC and Windows XP Software License		
14		CDS610B	DISPLAY SYSTEM: BARCO CORONIS 6MP/ MDCC6130 DL 30" Color 6 MegaPixel Diagnostic Flat Panel with auto DICOM calibration, embedded front light sensor and protective glass. Hardware to be provided by the customer, and must be an exact match to the standard McKesson Bill of Materials at time of implementation. This hardware is to be shipped (prepaid and insured) to McKesson's Canadian location for system configuration and testing by a date to be coordinated with McKesson's Production department		
15		ADI210BB	AUXILIARY DISPLAY SYSTEM BARCO MDRC-2120 20.1" Colour 1600 X 1200 Clinical Review Display with DICOM calibration Hardware to be provided by the customer, and must be an exact match to the standard McKesson Bill of Materials at time of implementation. This hardware is to be shipped (prepaid and insured) to McKesson's Canadian location for system configuration and testing by a date to be coordinated with McKesson's Production department		
16		WKS100HXP 80001506	HORIZON RAD STATION - HP This chassis is not be used with 5MP Displays Hardware: Workstation with a 100/1000BaseT NIC and Windows XP Software License.		

- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION		
17		CDS610B	<p>DISPLAY SYSTEM: BARCO CORONIS 6MP/ MDCC6130 DL</p> <p>30" Color 6 MegaPixel Diagnostic Flat Panel with auto DICOM calibration, embedded front light sensor and protective glass.</p> <p>Hardware to be provided by the customer, and must be an exact match to the standard McKesson Bill of Materials at time of implementation. This hardware is to be shipped (prepaid and insured) to McKesson's Canadian location for system configuration and testing by a date to be coordinated with McKesson's Production department.</p>		
18		ADI210BB	<p>AUXILIARY DISPLAY SYSTEM: BARCO MDRC-2120</p> <p>20.1" Colour 1600 X 1200 Clinical Review Display with DICOM calibration</p> <p>Hardware to be provided by the customer, and must be an exact match to the standard McKesson Bill of Materials at time of implementation. This hardware is to be shipped (prepaid and insured) to McKesson's Canadian location for system configuration and testing by a date to be coordinated with McKesson's Production department.</p>		
19		WKS100HXP 80001506	<p>HORIZON RAD STATION - HP</p> <p>This chassis is not be used with 5MP Displays</p> <p>Hardware: Workstation with a 100/1000BaseT NIC and Windows XP Software License.</p>		
20		GDS320B	<p>DISPLAY SYSTEM: DUAL BARCO CORONIS 3MP / MDCG 3120</p> <p>21" Greyscale 3 Mpixel Diagnostic Flat Panel with auto DICOM calibration, embedded front light sensor and protective glass.</p> <p>Hardware to be provided by the customer, and must be an exact match to the standard McKesson Bill of Materials at time of implementation. This hardware is to be shipped (prepaid and insured) to McKesson's Canadian location for system configuration and testing by a date to be coordinated with McKesson's Production department.</p>		

- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION		
21		ADI210BB	<p>AUXILIARY DISPLAY SYSTEM: BARCO MDRC-2120</p> <p>20.1" Colour 1600 X 1200 Clinical Review Display with DICOM calibration</p> <p>Hardware to be provided by the customer, and must be an exact match to the standard McKesson Bill of Materials at time of implementation. This hardware is to be shipped (prepaid and insured) to McKesson's Canadian location for system configuration and testing by a date to be coordinated with McKesson's Production department</p>		
22		WKS100HXP 80001506	<p>HORIZON RAD STATION - HP</p> <p>This chassis is not be used with 5MP Displays</p> <p>Hardware: Workstation with a 100/1000BaseT NIC and Windows XP Software License</p>		
23		GDS340B	<p>DISPLAY SYSTEM: QUAD BARCO CORONIS 3MP / MDCG 3120</p> <p>21" Greyscale 3 Mpixel Diagnostic Flat Panel with auto DICOM calibration, embedded front light sensor and protective glass</p> <p>Hardware to be provided by the customer, and must be an exact match to the standard McKesson Bill of Materials at time of implementation. This hardware is to be shipped (prepaid and insured) to McKesson's Canadian location for system configuration and testing by a date to be coordinated with McKesson's Production department</p>		
24		ADI210BB	<p>AUXILIARY DISPLAY SYSTEM: BARCO MDRC-2120</p> <p>20.1" Colour 1600 X 1200 Clinical Review Display with DICOM calibration</p> <p>Hardware to be provided by the customer, and must be an exact match to the standard McKesson Bill of Materials at time of implementation. This hardware is to be shipped (prepaid and insured) to McKesson's Canadian location for system configuration and testing by a date to be coordinated with McKesson's Production department</p>		
25		WKS100HXP 80001506	<p>TECHNOLOGIST RAD CASE VIEWING STATION - HP</p> <p>This chassis is not be used with 5MP Displays</p> <p>Hardware: Workstation with a 100/1000BaseT NIC and Windows XP Software License</p>		
26		ADI210BB 80004560	<p>DISPLAY SYSTEM: BARCO MDRC-2120</p> <p>20.1" Colour 1600 X 1200 Clinical Review Display with DICOM calibration</p>		

- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION		
27		WKS100HXP 80001506	PACS ADMINISTRATIVE WORKSTATION - HP  This chassis is not be used with 5MP Displays  Hardware: Workstation with a 100/1000BaseT NIC and Windows XP Software License		
28		ADI210BB 80004560	DISPLAY SYSTEM: BARCO MDRC-2 I20  20.1" Colour 1600 X 1200 Clinical Review Display with DICOM calibration		
29		WSA100 80000454	DOCUMENT SCANNING  Hardware: USB Document scanner with all required cabling and adapters. Includes an Automatic Sheet Feeder  Requires: Horizon Rad Station 3.3.2 or greater and Microsoft Windows XP		
30		CAL400 80001194	BARCO DISPLAY CALIBRATION KIT		
31		RVC101XP 80000387	BASIC VOICE CLIPS W/ HEADSET  Includes headset with integrated speakers and microphone  Requires Horizon Medical Imaging version 4.6 and Horizon Rad Station version 3.3 or greater		
32		NHA410H2K8x64 80004860	HIGH AVAILABILITY HORIZON MEDICAL IMAGING DB SERVER - UP TO 3,000,000 EXAMS PER YEAR - HP  Equipment: Dual servers in failover configuration to provide failover services in the event of single unit hardware failure. Windows Server license, Rack mountable servers, Gigabit ethernet adaptor, requires a copper Gigabit ethernet connection.  Requires 2x ORC101 Oracle Software License.		
33		HPA400 80001582	HORIZON MEDICAL IMAGING COMMUNICATIONS MANAGER - PAIR - UP TO 1,000,000 EXAMS  The Horizon Medical Imaging communications manager provides high-availability and single network IP access for DICOM, HTTP and HL-7 communications to an HMI Application Server farm		
34		HPA100 80003731	HORIZON MEDICAL IMAGING COMMUNICATIONS MANAGER - SINGLE - UP TO 1,000,000 EXAMS  The Horizon Medical Imaging communications manager provides a single network IP access for DICOM, HTTP and HL-7 communications to an HMI Application Server farm		
35		NSD710H 80004324	NETWORK ATTACHED STORAGE - 3,379 GB - SAS - HP  Includes: Network storage appliance with 3,379 GB usable space, 1000BaseT network connection.  Requires Horizon Medical Imaging version 4.4 or greater		

- BUDGETARY QUOTATION -

ITEM NO.	QTY	PART NO.	DESCRIPTION		
36		NSD620H 80004322	<p>NETWORK ATTACHED STORAGE - 7,578 GB - SATA - HP</p> <p>Includes: Network storage appliance with 7,578 GB usable space, 1000BaseT network connection.</p> <p>Requires Horizon Medical Imaging version 4.4 or greater.</p>		
37		UPS200H	<p>RACK MOUNTABLE UPS</p> <p>Hardware to be provided by the customer, and must be an exact match to the standard McKesson Bill of Materials at time of implementation. This hardware is to be shipped (prepaid and insured) to McKesson's Canadian location for system configuration and testing by a date to be coordinated with McKesson's Production department.</p>		

**System Configuration Settings**

(A) 12x WSA370

**WORKSTATION POWERSCRIBE INTEGRATION**

Includes: Integration of a Horizon Rad Station with an existing PowerScribe system for voice recognition. The PowerScribe client may be installed on the HRS-Advanced workstation or on a standalone PowerScribe station.

PowerScribe system and client licenses must be existing or purchased separately from PowerScribe, including PowerMics to be connected to Horizon Rad Station Advanced. The PowerScribe system must be interfaced to the site HIS/RIS via HL-7. All charges associated with the PowerScribe system in order to enable this interface are the responsibility of the customer.

McKesson recommends dual CPU workstations for VR processing. Desktop integration with VR software and 3D advanced visualization products may create performance issues, related to the competing requirements for CPU processing power. It is recommended that if both VR and Advanced Visualization are run on the same workstation that the work is done in a non-synchronous fashion. An alternative is to have a standalone 3D Advanced workstation with the standalone VR integration.

PLEASE NOTE: Requires confirmation by McKesson that the desired versions of PowerScribe and Horizon Rad Station have been successfully qualified together. The customer must provide McKesson with the specific version numbers to be integrated. Future software upgrades and/or charges may apply to validate integrations with other versions of PowerScribe or to make customizations to the qualified integration.

(B) 1x API100

**INTERFACE FOR 3RD PARTY EMR SYSTEM**

Provides McKesson standard API interface to enable image viewing of a 3rd party Electronic Medical Record [EMR] system. The Customer is responsible for any services or licensing required by the 3rd party EMR vendor. EMR100 must be purchased if customization of the standard API is required by the customer or 3rd party vendor.

(C) 1x DXS140

**HL7 OUTBOUND STATUS UPDATES**

Includes: HL7 interface to send updates to a RIS when the status of the study changes in Horizon Medical Imaging.

Requires Horizon Medical Imaging version 4.5 or greater.  
Requires DXS100 and SVC140.



VISN5

Quotation: H11VSN009-01PT  
Project HMI  
Site VA Washington

Contact Name Tony Douglas  
Contact Number (804) 464-7851  
Revision Date September 17, 2010

- BUDGETARY QUOTATION -  
3RD PARTY PASS THROUGH PACKAGES

Third Party Pass Through Packages are provided as a convenience to the customer. This purchase is as per the Third party vendor standard terms and conditions. As part of the terms and conditions of this project McKesson will act as the project manager and will coordinate installation. Customer will contract directly with Third Party Vendor for on-going service. Installation and service are provided by the Third Party Vendor. McKesson does not perform qualification of Third Party Pass Through components. Please refer to the Third Party Vendor for compatibility or operational issues.

ITEM NO.	QTY	PART NO.	DESCRIPTION		
Professional Services					
1		MMV310 TBD	MEDVIEW CACHE SERVER HW WIN 2003 SERVER Dell Server, Mini-Tower Chassis		