

DESCRIPTION OF REQUIREMENT

Salient Characteristics
For

The OER-PRO OLYMPUS REPROCESSOR (Brand Name or equal).

Systems must include the following features:

Applicable scopes: OER-PRO Olympus Reprocessor

Number of reprocessed endoscopes: Max. 2 (1 with certain models)

Cleaning method: Exterior surfaces – Ultrasonic cleaning, turbulent bath
Channel interiors – Fluid flushing
Valves – Ultrasonic cleaning, fluid flushing

Disinfection method: Exterior surfaces – Disinfection solution immersion
Channel interiors – Disinfection solution flushing and filling
Valves – Disinfection solution immersion

Cleaning time settings: 3 – 10 minutes (setting variable in 1 min. increments)

Disinfection time setting: Acecide-C – 7 minutes
ALDAHOL 1.8 – 10 minutes

Disinfectant solution temperature setting: 20 C (68 F)

Disinfectant solution heating method: Built-in heater in the reprocessing bath

Water discharge method: Forced draining using a pump (floor drain)

Disinfectant solution discharge method: Drains through disinfectant collection hose
Drains through drain hose (floor drain)

Reprocessing basin capacity: Approximately 14 L (3.7 gallons)

Disinfectant solution tank capacity: Approximately 17.5 L (4.6 gallons)

Disinfectant solution: Acecide-C
ALDAHOL 1.8

Detergent: EndoQuick

Visual leakage detection: Bubble detection during immersion

Alcohol flushing: Automatic flushing/draining using a pump and compressor

Dimensions: 450 (W) x 977 (H) x 765(D) mm (17.7 x 38.5 x 30.2 inch)

Weight: 120 kg dry condition (264.6 lbs)

Power supply: Voltage – 120 VAC
Frequency – 60 Hz
Input current – 5.5A
Voltage fluctuation - +/- 10%

Medical device classification: protection against electric impact – Class 1

OSHDP Seismic Anchorage: Preapproved (OPM-0137-13)

The item(s) are to be delivered on or before 30 days after receipt of the order at the expense of the offeror, FOB destination to: C.W. Bill Young VAMC, 10,000 Bay Pines Blvd, Bay Pines, FL. 33744.

PRICE/COST

ITEM INFORMATION

ITEM NUMBER	DESCRIPTION OF SUPPLIES/SERVICES	QUANTITY	UNIT	UNIT PRICE	AMOUNT
0001	OER-PRO: OER-PRO OLYMPUS REPROCESSOR MFG:OLYMPUS STOCK NUMBER: MF01-0033PL: MF01-	3.00	EA		
0002	0033PL EXTERNAL PRE-FILTRATION SYST MFG:OLYMPUS STOCK NUMBER: MF01-0014PL: MF01-	3.00	EA		
0003	0014PL 1 MICRON EXTERNAL FILTER MFG:OLYMPUS STOCK NUMBER: MAJ-865: MAJ-865 OER-	3.00	EA		
0004	PRO ENDOSCOPE HANGER	3.00	EA		

	MFG:OLYMPUS STOCK NUMBER: ACECIDE-C: ACECIDE-C				
0005	PERACETIC ACID 6 SETS/CASE MFG:OLYMPUS STOCK NUMBER: ACECIDE-C TEST:	6.00	EA		
0006	ACECIDE-C TEST: TEST STRIPS 100/btl MFG:OLYMPUS STOCK NUMBER: ENDOQUICK	6.00	EA		
0007	ENDOQUICK DETERGENT 3 BTLS/Case - 2L ea MFG:OLYMPUS STOCK NUMBER: MAJ-1511: MAJ-1511 OER-	6.00	EA		
0008	Pro CONNECTOR LF/ENT MFG:OLYMPUS STOCK NUMBER: MAJ-1513: MAJ-1513 OER-	3.00	EA		
0009	Pro CONNECTOR BF INSTRU/SX MFG:OLYMPUS STOCK NUMBER: MAJ-1904: MAJ-1904 OER-	3.00	EA		
0010	Pro CONNECTOR GYRUS SCOPES MFG:OLYMPUS STOCK NUMBER: MAJ-1515: MAJ-1515 OER-	3.00	EA		

0011	Pro CONNECTOR MAJ-855 &EBUS MFG:OLYMPUS STOCK NUMBER: MAJ-1503: MAJ-1503 OER-	3.00	EA	_____	_____
0012	Pro CONNECTOR DUAL CHANNEL MFG:OLYMPUS STOCK NUMBER: MAJ-840: MAJ-840 OER-	3.00	EA	_____	_____
0013	PRO EUS/EBUS SCOPE RACK MFG:OLYMPUS	1.00	EA	_____	_____
				GRAND TOTAL	_____