

Line #	Description	Qty
1	CX50 2D xMATRIX	1

Interface:

- 15.0 inch high resolution display with wide viewing angle
- Quick Keys and Active Mode
- Laptop style Alphanumeric QWERTY keyboard
- 8 TGCs and 2 LGCs
- Ergonomic carrying handle
- Includes AC adapter , power cord and system battery pack
- 2 USB flash drives on system
- 80 GB hard drive
- Internal DVD RW drive
- Windows Embedded Standard 7 Operating System

Architecture:

All-digital compact broadband beamformer, Microfine 2D focusing with Dynamic Focal Tuning that includes Advanced X-Res signal processing, 170 dB full time input dynamic range 18,432 digitally-processed channels, Continuously variable steering in 2D, color and Doppler modes 2D Opt signal processing with 4X multi-line parallel processing and frequency compounding.

Intelligent Controls:

The CX50 has been designed to make portable exams easy and efficient. With a single button, iSCAN technology automatically samples data for a new level of 2D and Doppler optimization iSCAN one-touch Intelligent Optimization, iSCAN one-touch Intelligent Color Optimization, iSCAN Doppler one-touch optimization.

Transducers:

Supports Compact family of transducers featuring PureWave imaging technology in the S5-1, CX7-2t, C5-1, D5CWC. Also supports the high resolution S12-4, S8-3, C8-5 and L12-3 transducers. All transducers provide breakthrough frequency bandwidths and array configurations. These transducers also have ergonomically designed lightweight flexible cables and compact connectors.

Modes:

- 2D
- M-mode
- Anatomical M-mode
- Color M-mode
- Pulsed Wave Doppler
- Color Power Angio (CPA)
- Continuous Wave Doppler

- Invert and Color Invert
- Color compare mode
- Dual mode
- Duplex for simultaneous 2D and Doppler
- 2D Optimization Signal Processing
- Live Compare
- Tissue Harmonic Imaging (THI)
- Reconstructed zoom with pan (read zoom)
- Write zoom
- Pulse Inversion Harmonic imaging
- Adaptive Doppler
- Adaptive Color Doppler
- Color Tissue Doppler imaging
- Pulsed Wave Tissue Doppler imaging
- Active Native Data - manipulation of image data
- DICOM Networking
 - Ethernet @100Mb/second Includes DICOM wired and wireless "G" and "N".
 - Provides DICOM 3.0 network print and store and storage commitment
 - Performed Procedure Step (pps)
 - Modality Worklist
- DICOM Structured Reporting

Cineloop review

- On-board workstation-class data management with thumbnail previews and storage of images, loops, and reports. Retrospective and prospective clip capture to internal drive or removable media
Integrated DVD/CD burning capability for storage of images or export in DICOM, JPEG and .avi for PC compatibility. Philips DICOM viewer option to imbed in media transfer for easy viewing of study on most PCs.
- Maintenance and Serviceability
- Remote Access for Expedient Clinical and Technical Support
- Flexible Service Agreements
- Clinical Application and Educational Support

Security-related features

Firewall policy blocks all unnecessary ports

Operating System hardening

- OS settings utilizing the DISA STIGS
- Disable unnecessary services
- Disable auto-run for removable media

Media export security

- Provides the ability to disable export of patient data to removable media

System Access Control

- No restrictions – users may perform exams and access all previously completed exams or MWL data
- Only patient data is locked – users may perform exams without requiring a login, but must successfully log in prior to accessing previously completed exams or MWL data
- Complete system is lockable – users and administrators must successfully log in prior to any system access

User management policy

- User management local
 - Local user management
 - Support for multiple unique user accounts
 - Support for multiple unique administrator accounts
- User management remote
 - Supports active directory authentication utilizing LDAP (system may not be joined to the domain)
 - Support for individual accounts or AD groups for users and administrators
 - May utilize LDAP or secure LDAP
 - Customer may configure the system to perform authenticated binding

Password policies

- Provides the ability to specify password policies for local accounts
 - Minimum password length (6-14 characters)
 - Maximum password length (6-63 characters)
 - Minimum password age (0-998 days)
 - Maximum password age (1-999 days)
 - Password complexity

Account lockout policies

- Lockout threshold (1-999 minutes)
- Lockout duration (1-999 minutes)
- Lockout counter reset (minutes)

Auto logoff

- Automatically logs off a user after the specified period of inactivity, user selectable

Hard drive encryption

- Encryption policies
 - 128 bit
 - 128 bit with diffuser

- 256 bit
- 256 bit with diffuser

Audit log export

- Audit logs may be exported
- Available protocols are UDP or TLS

Clinical Education

1 Day PAS Onsite - Ultrasound system or upgrade onsite training provided by a PAS (Product Applications Specialist) for specific system applications or upgrades; not per modality. *Education is provided Monday - Friday during normal business hours.* Note: Philips Healthcare personnel are not responsible for actual patient contact or operation of equipment during education sessions except to demonstrate proper equipment operation. The training sessions should be attended by the appropriate healthcare professional as identified by the department director. *Repeat training for staff non-attendance will not be accepted.* Site must be patient-ready to meet training expectations. All onsite training day expires within 90 days from system or upgrade installation date. Exceptions are for 3D Stress onsite training (which expires 9 months from system or upgrade installation date) and Fusion & Needle Navigation onsite training (which expires 180 days from system or upgrade installation date).

THE NUMBER OF ONSITE TUITIONS YOU RECIEVE MAY VARY BASED ON PURCHASED OPTIONS. PLEASE CONSULT YOUR SALES REPRESENTIVE FOR FURTHER DETAILS

All Tuitions must be registered prior to the expiration date. The course chosen must be taken within 90 days of expiration.

2	Adult Echo, Clinical Option	1
	Tissue Specific Imaging software for adult cardiac ultrasound applications. Display optimization software with Tissue Specific presets for adult cardiac imaging and Doppler applications. Analysis software package includes cardiac imaging protocol measurements and configurable reports and finding codes. Active native data for post-process optimization and advanced XRES adaptive image processing for improved tissue conspicuity. iSCAN intelligent one-button optimization for adaptive gain compensation in 2D, Doppler, Tissue Doppler Imaging and LVO contrast functions. Includes Live compare mode, cardiac High-Q Automatic Doppler Analysis and respiration waveform from chest impedance. Allows operation of S5-1, CX X7-2t and D2cwc transducers.	
3	Cart	1
	Highly mobile cart that includes: 4 swivel wheels with 4 locking casters, rear handle, micropositioning grips, quick-connect tray, utility drawer, storage shelf, footrest, internal isolation transformer, B&W printer brackets, integrated transducer connector holder, gel holders and cable management. Includes USB hub for additional connectivity.	
4	USA Power Cord	1
5	Government Security	1
	Required by all DoD customers. This option disables VNC capabilities (which if enabled would provide remote desktop support) for increased security of data.	
6	S5-1 Transducer	1

Features compact connector designed for reliability and improved ergonomics. Manufactured in accordance with the European Union's Restriction of Hazardous Substances (RoHS) directive.

PureWave crystal Sector array transducer with 5 to 1 MHz extended operating frequency range for adult cardiology, abdominal, vascular, TCD and Acute Care.

7	D2CWC Transducer	1
	Non-imaging 2 MHz PW/CW Doppler transducer for cardiac applications	
8	English Manual	1
	Operation Manual	
9	1st SVC Manual for Gov	1