

630-B60017
VAMC BROOKLYN, NY

NM BONE DEN

Qty

Product Description

1

HORIZON-C

Advanced Fan-Beam DXA Technology

One Pass™ Acquisition Technique; Multi-Detector Array Scanning Method
New 128 high-resolution multi element ceramic digital detector array
New High Frequency X-ray Generator
Indexing Scan Table with Positioning Accessories
Motorized Table and C-arm
Dynamic Internal Reference System for Continuous Calibration
Computer Console
QDR Anthropomorphic Spine Phantom

Standard Computer Package

Computer Workstation with Dual Core 3 GHz
Windows® 7 Professional
500 GB hard drive
4 GB RAM
19" Widescreen LCD Monitor
HP Professional Series Color DeskJet® printer
DVD RAM drive

Standard Configuration

Hologic APEX 5™ Operating System
Automatic PASS/FAIL Quality Control
Express BMD™ 10 Second Acquisition
Single Energy Scan Display Capability
Window/Level Control for Image Optimization

APEX Productivity Tools

Express Exam™ Workflow Management
OneTime™ Auto Analysis with Histogram
ProTech with DXApro™
Auto Hip Positioning
Reposition/Rescan Feature
Automatic Scan Comparison for Serial Exams
Least Significant Change Configuration

		<p>Horizon Advance Reporting Solutions</p> <p>QDR OnePage™ Report with Rate of Change Assessment</p> <p>FRAX® 10 Year Fracture Assessment</p> <p>New Dual Hip™ Report</p> <p>New integrated Report Writer DX™</p> <p>Horizon Scan and Analysis Protocols</p> <p>AP Lumbar Spine with Automatic Low Density Analysis and Scoliosis Analysis</p> <p>Proximal Femur, Automatic Low Density Analysis and Hip Structure Analysis (HSA)</p> <p>Dual Hip™</p> <p>Forearm</p> <p>IVA HD with Image Pro High Resolution Imaging Capability</p> <p>Quantitative Morphometry</p> <p>Integrated Physician's Viewer with MXApro™</p> <p>New Atypical Femur Fracture Assessment (AFF) High Resolution Imaging Capability</p> <p>Pediatric Analysis for Spine, Femur, Infant Spine and Forearm</p> <p>DXA Practice Manager</p> <p>Patient Database Manager</p> <p>On-line Normative Reference Values and User-Definable Values</p> <p>NHANES, Ethnic and Pediatric/BMDCS Reference Data</p> <p>Context Sensitive Help</p> <p>General Region of Interest</p> <p>DAP Reporting</p> <p>QDR Reference Manual CD</p> <p>Installation and a Twelve (12) Month Comprehensive Warranty</p> <p>Hologic Connect™ Remote service tool</p>	
1		18.5" FLAT PANEL MONITOR	
1		ASSY, P4 COMPUTER, PCIE, W7	
1		FRAX OPTION FOR APEX	
1		ATYPICAL FEMUR FRACTURE OPTION	
1		<p>DXA TRAINING, INITIAL, 1.5 DAYS (12 HOURS), 1SITE, MAX 4 STAFF</p> <ul style="list-style-type: none"> • One and a half (1.5) consecutive days (total of 12 hours) of Applications Training for up to (4) Operators per site. • Applications must be completed within 12 months of equipment shipment. • Please note: Cancellation must be made 6 days prior to the confirmed and scheduled applications training start date. • Initial Applications 	
1		TRADE IN - HOLOGIC - BONE	
1		<p>INKJET-COLOR</p> <p>HP Office Inkjet-Color Printer (CMP-00993) included in all APEX 4.XX Windows 7 Upgrades</p>	
1		<p>KIT, PCIE CABLE, INTERFACE 30'</p> <p>30 FT CABLE INTERFACE Kit</p>	
1		POWER KIT FOR 120VAC	
1		<p>APEX IRIS CONNECTIVITY SUITE W/ENTERPRISE DATA MANAGEMENT</p> <p>IRIS Enterprise-APEX for systems on APEX 4.xx Includes</p> <p>DICOM Storage</p> <p>DICOM Modality Worklist</p> <p>Remote Physician's Viewer</p> <p>Physician's Report Writer</p> <p>HL7</p> <p>Enterprise Data Management</p> <p>Installation and Twelve (12) Month Comprehensive Warranty for IRIS Enterprise</p>	