

The Lake Nona Orlando VA Medical Center is seeking the below items to network with existing equipment purchased on July 14, 2016 for the facility. Spectralis Tracking OCT System-OCT w/OCT2 advanced-GMPE Software: ophthalmic diagnostic imaging platform is used to examine the structures of the retina, macula and optic nerve. It combines two technologies confocal scanning laser ophthalmoscopy (cSLO) and spectral domain optical coherence tomography (SD-OCT).

0226-OCT2-ADV	SPECTRALIS TRACKING OCT SYSTEM-OCT W/OCT2- ADVANCED	1 EA
0226-680-2	SPECTRALIS CAMERA-OCT FIXED MOUNT W/OCT2-ADVANCED	1 EA
0231-158	COMPUTER-SPEC OCT 2, WIN 10, HP Z440 6.6	1 EA
0235-122	MONITOR, NEC NULISYNC, 24" EA244 WMI-BK	1 EA
0236-142	PRINTER HP OFFICEJET PRO 8210	1 EA
0237-257	TABLE SPECTRALIS OCT FORMICA (REV B)	1 EA
0238-281	HEYEX NETWORK LICENSE (INCLUDED W/SYSTEM)	1 EA
0238-479	ISOLATION TRANSFORMER 3RD EDITION REV 2	1 EA
0238-483	HARD DRIVE EXTERNAL, 4TB ENTERPRISE MINISTACK	1 EA
0998-OCT	1 YR WARRANTY INCLUDES INSTALLATION & TRAINING	1 EA
0226-ASM	SPECTRALIS ANTERIOR SEGMENT MODULE	1 EA
0226-093	SPECTRALIS ANTERIOR SEGMENT LENDS	1 EA
0226-960	SPECTRALIS ANTERIOR SEGMENT SOFTWARE VER 6.3	1 EA
0226-963	SPECTRALIS SOFTWARE VERSION 6.6	1 EA
0226-962B	GMPE SOFTWARE	1 EA
0226-955	SPECTRALIS SOFTWARE-GMPE ACQUISITION	1 EA
0226-956	SPECTRALIS SOFTWARE-GMPE VIEWER	1 EA
14208	GMPE KIT, US	1 EA
0226-962T	GMPE SOFTWARE-INSTALLATION & ON-SITE TRAINING	1 EA

NOTE: Questions to this request must be received by March 20, 2018 by email to EarlCi.Johnson@va.gov. The following must be listed in the subject title "36C24818Q9083 Spectralis Tracking OCT System-OCT w/OCT2" for all responses.

NOTE: The end date for responses is 12:00 p.m. on March 26, 2018