

XR RAD, VAMC MUSKOGEE, OK
PO# 623-B68005

Item Qty		Description
No.		
1	1	<p>XR656P 2D BASE W AUTORAD</p> <p>Discovery XR656 Plus Base System with TWO FlashPad Detectors</p> <p>The Discovery XR656 Plus is an advanced digital radiographic imaging system that provides excellent image quality, a variety of image manipulation and post-processing tools as well as the option to utilize GE Healthcare's exclusive advanced clinical applications. Powered by FlashPad, GE's wireless digital detector, workflow is streamlined with TWO portable detectors that can be used in the table, wall-stand or freely to best accommodate most 2D exam requirements as well as advanced radiographic requirements.</p> <p>At its core, the Discovery XR656 Plus delivers improved productivity through an efficient and intuitive workflow. The base system comes</p>

equipped with a systems cabinet,
acquisition
review workstation, image processing
tools,
short-term storage and quick in-room
viewing of
images. Also included is a host of
networking
and connectivity options, dose reporting
and
system maintenance. These features are
designed
to make this system easy to use and
reliable
while providing high quality radiographic
images
in a digital environment.
The Discovery XR656 Plus's core feature
set can
help streamline even the toughest exam.
Consistent image quality helps reduce
retakes,
and unique image processing helps save
time.
With DICOM connectivity, the Discovery
XR656
Plus connects easily to HIS, RIS or PACS
for
efficient data transfer.
The Discovery XR656 Plus base system
includes:
Two wireless Digital Detectors with
batteries
and optional tether.
o Single panel (non-tiled) amorphous
silicon
detector with a Cesium Iodide scintillator
and two handgrips that facilitate easy
positioning and a secure grip.
o Wireless connectivity through
Ultrawide

Band Technology for fast, efficient and secure communications between the detector and the base system.

- o Image area 40.4 cm x 40.4 cm (15.9 in. x 15.9 in.)
- o Active matrix 2022 x 2022 pixels.
- o Pixel pitch 200 microns.
- o Typical upper dynamic range 7.8 mR @ RQA5.
- o Typical DQE 68 @ 0lp @ RQA5.
- o Weighs 4.32 kg (9.52 lb.) without battery, 0.18 kg (0.40 lb.) battery weight.
- o Dimensions: L 580 mm, H 452 mm, T 24 mm.
- o Can support up to 160 kg (352 lb.) of distributed load.
- o Battery or tethered operated. Includes two rechargeable and exchangeable batteries and 7 m cable (4 m or 10 m optional) for optimal connectivity and power.
- o One 6:1 clip-on grid for FlashPad detectors.
- o Battery allows for 150 images to be taken in 3 hours.
- o QAP (Quality Assurance Procedure)
- o FlashPad detector can support Advanced Applications that are options available on the Discovery XR656 Acquisition Workstation

The Acquisition Workstation is the primary interface to the network and provides

image post-processing capabilities. The System Controller Module provides single point control, directing and coordinating overall system operation, while monitoring all system modules automatically through software. The Acquisition Workstation includes:

- o Two 19-inch Flat Panel Monitors that help minimize desktop space required.
- o Keyboard and Mouse
- o CPU Tower with 6GB RAM, 320 GB of hard drive storage and capacity for over 22,000 images.
- o 120/140 VAC, 50/60 Hz.
- o Easy Image Manipulation and
- o Image Display Tools which include:
 - o Window width and level
 - o Gray scale/contrast invert
 - o Interpolated zoom and roam
 - o Image flips (horizontal, vertical) with automatic indicator
 - o Image Rotate - 90 degree increments
 - o Free rotation - 360 degrees
 - o Image orientation management
 - o Electronic Left/Right Markers
 - o Free text annotation
 - o Manual shuttering
 - o Image Annotations and Measurement Tools
- o Multi-Resolution Post Processing
- o Customizable Image Processing to Match
- o User Preferences
- o CD-RW and DVD Drive for Image Archive

- o Image Viewer on Archive CD's and DVD's
- Acquisition Control System
 - o Single Point System Control and Monitoring.
 - o Auto-Protocol Programming: comes with default set of exam and view protocols and the ability to build an infinite number of exams or views through the editor function.
- User
 - also has the ability to attach acquisition protocols.
- o System manager allows equipment error logging and provides resident power-up diagnostics.
- o kVp selections range from 40-150 in 1 kVp increments
- o Focal Spot Sizes of 0.6 and 1.2 mm - mA
- Patient Management
 - o DICOM Modality Worklist for HIS/RIS - SCU (with programmable auto refresh)
 - o Patient edit/auto-folding (Copy exam)
 - o Patient Dose Reporting Calculated Dose Monitoring (mGy)
 - o "Patient Directory" provides fast access to the image and exam database for case reviews and file management.
 - o Detector Exposure Indicator: tool for detector dose feedback to ensure exposure was within normal limits.
 - o Emergency Patient Exam Capability

Emergency patient feature - allows user to open exam and acquire images without a worklist entry.

- o The Copy Patient function allows merging of the patient information with the exam images.
- o Fast Image Display - Average time for a partially processed image is approximately three seconds and less than eight seconds for a fully processed image. Times may vary based on how the detector is connected to the system (i.e. docked, tethered or wireless)
- o Set of default adult and pediatric protocols allows quick selection of the appropriate techniques for common procedures/exams with the ability to define unlimited number of custom protocols.
- o Set of 4 Factory (GE pre-set) image processing selections (looks) optimized for each anatomical view with the ability to define multiple Custom looks for each anatomical view/ patient size combination.
- o Automatic image storage and print with DICOM 3.0 and IHE Compliant networking, further increasing exam throughput and

decreasing examination time for patients.

Image Quality and Dose

- o Multi-resolution image processing capability

- o Tissue Equalization used to correct over-penetrated and under-penetrated areas

within the image.

- o Auto and manual image shuttering cropping

tool.

- o Automated brightness/contrast setting (Smart

Windowing)

- o Orthopedic Magnification/Print.

- o Detector Exposure Index (DEI) - dose tracking and QC metric.

- o Dose Area Product (DAP) - entrance dose

metric.

- o Grid Line Reduction (selectable)

- o Intelligent Collimator Edge Detection

automated, image based

cropping/shuttering

tool that relies solely on image

information

to locate the collimator edges present in the image.

- o Orthopedic Magnification/Print

DICOM and Standard Networking

Capabilities

Images may be transmitted manually or automatically through the DICOM interface to

printers, archival devices, servers or review

workstations.

Please refer to the DICOM Conformance Statement

for complete definition of supported

DICOM
connectivity services.
DICOM and Standard Networking
Capabilities
include:
o Ethernet Network Link - DICOM 3.0.
o DICOM Storage (with auto-send to
different
locations)
o DICOM Storage Commitment (with
programmable
auto delete function)
o DICOM Modality Worklist for HIS/RIS
(with
auto refresh)
o DX/CR Worklist Filtering.
o DICOM Media Interchange on DVD-R.
o DICOM Modality Perform Procedure
Step (MPPS)
feedback to the HIS/RIS (SPS PPS)
o DICOM Grayscale Print (with print
layout at
the console)
o DICOM Query/Retrieve (retrieves
images back
from PACS)
o System Access & Authorization Control
to
support HIPAA Compliance.
o Full Range of Printing Options.
o Numerous Layout and Format Options.
Systems Cabinet
o Built-in System Distribution Power
Module
and Circuit Breaker for Single Point
Power
Feed to Room Subsystems
o Modular Designed X-ray Systems
Based on a
Digital Communications Network for
Improved

Reliability and Image Quality.
Quality Control/System Reliability
Features
Preventing customer experienced
system failures
and reducing unplanned system
downtime are
critical.
The following features help to achieve
these
goals:

- o Using the integrated system Quality Assurance Procedure (QAP), image quality checks can be easily performed by the customer.
- o The QAP includes a phantom, optimized for Digital Image Quality testing and is included with the system.
- o System changes are highlighted and can be corrected before they become a problem.

Discovery XR656 Overhead Tube Suspension
The Overhead Tube Suspension (OTS) system with motorized movement delivers excellent levels of operational support for efficient operation and precise positioning.

The Overhead Ceiling Tube Suspension Package includes:

- o Auto-Positioning Package (included in base)
- o Auto-Positioning enables the users to select a predefined system position from the system console and automatically

move
 the equipment by simply holding the
 "Auto
 Positioning" buttons. This feature
 reduces
 user fatigue and increases the
 productivity of the operator.

- o Auto-Positioning is controlled at the acquisition workstation or with the IR remote control, allowing the user to remain in the room while moving the system.
- o Pre-set positions at the table, wallstand and park position at various SIDs and vertical and horizontal orientations.
- o Auto-Positioning will incorporate angulation of the tube, longitudinal, lateral, rotational and vertical positioning of OTS, table detector longitudinal positioning, wallstand detector vertical positioning.
- o Auto-Positioning comes with 7 default positions and up to 10 additional user defined positions can be added to the system.
- o Tube and Detector Tracking
 Synchronized tube and detector tracking enables convenient workflow by keeping the tube correctly aligned with the wall-stand or table detector.
- o Motorized 5-axis movements
 Supported
 Positions: park, table 100 cm SID (head, center, foot), wallstand 100 cm and 180 cm CID, and horizontal wallstand 100 cm SID.
- o IR Remote - The infrared remote is an

in-room control allowing the technologist greater flexibility & ability to pre-position the system automatically in preparation for the next clinical exam/view.

- o Auto or manual positioning with single lock release and auto detents for assisted manual positioning.
- o Touch-screen user interface with LCD screen display helps confirm patient data, review techniques, receptor selection and modify/confirm wireless detector association. The easy to Read, Auto Rotating user screen also includes a message readout line and easy to see light indicators.
- o Patient Side Touch Screen User provides the following functions to the user:
 - o Lock, Detent Control
 - o Field of View Image Size Selection
 - o Collimator Field Light Selection
 - o Technique Adjust (kVp, mAs)
 - o Receptor Selection (table, wallstand, wireless or cassette)
 - o Exam Inhibit Display
 - o Collimator Manual Override
 - o Position Display (Source-to-Image Distance, X-ray Tube Angle, Column Rotation Angle)
 - o Display of Patient Name for In-Room Verification (this feature can be disabled)

Repeat and Reject Analysis for the Discovery XR656

The Repeat and Reject Analysis (RRA) software package is a quality control (QC) application available for the Discovery XR656 that allows for repeat or reject images to be captured and categorized by technologist. It is designed to help track and analyze the X-ray repeat rate

- o The easy-to-use operator interface helps technologists classify each image they reject and to select a specific reason for the repeat /image rejection.
- o RRA can be a helpful teaching tool because it includes links to actual JPEGs of the rejected images to help the user analyze why the image was rejected.
- o The RRA application tracks the rejected image data by operator, exam type, date and reason code.
- o Reports can be exported in DVD, CD or USB format for ease of use.

Auto Field of View for the Discovery XR656

Auto Field of View enables the user to pre define the collimation size on an individual view basis. Helps improve workflow by reducing the time it takes to make the final collimation adjustment for a particular

exam.
 Auto Protocol Assist for the Discovery
 XR656
 Auto-Protocol Assist - the Discovery
 XR656
 system will automatically transition
 directly
 to the Acquire screen when the protocol
 code
 downloaded from the HIS/RIS
 (automatically
 performed with worklist refresh)
 matches the
 exam code contained in the protocol
 database.
 This tool eliminates the user steps
 required
 to select patient exam types and initiate
 an exam.

2	1	2, 3, or 4 Meter Longitudinal Drive Belt Kit 2, 3 or 4 Meter Longitudinal Drive Belt Kit
3	1	Inboard Rail Select 2, 3 or 4 Meter Longitudinal Rail Select (Dependent on Room Size)
4	1	2, 3, or 4 Meter Bridge Select 2, 3 or 4 Meter Bridge Select (Dependent on Room Size)
5	1	2, 3, or 4 Meter Bridge Cable Select 2, 3 or 4 Meter Bridge Cable Select
6	1	2, 3, or 4 Meter Bridge Cable Drape

Select
2, 3 or 4 Meter Bridge Cable Drape Select

7	1	Required Language Labels Required System Language Labels - English
8	1	80kW High Frequency Generator 80kW High Frequency Generator, 50 or 60Hz The high frequency 80kW power unit is designed for radiographic applications and utilizes microprocessor controlled power and parameter adjustment Specifications: <ul style="list-style-type: none">• 1000mA at 80kVp• 800mA at 100kVp Input Power: 360-480VAC, 3 Phase, 60Hz kVp Range: 40-150kVp, 1kVp increments <ul style="list-style-type: none">• kVp Accuracy: 3% +/-2kVp mA Range: Small Focal Spot: 10-320mA <ul style="list-style-type: none">• Large Focal Spot: 160-1000mA
9	1	XR656 PLUS OD TBL W 100 Elevating table with floating tabletop and 3 AEC ion-chamber. Includes cross lateral detector holder. 100 cm Focus Table Grid - 100 cm (40 inch) focus - 12:1 ratio, 70 lines/cm
10	1	Table Top FlashPad Lateral Detector Holder Table Top Lateral Detector Holder Wireless DR detector holder, designed

specifically for GE, secures the detector in a vertical position on the tabletop for cross-table imaging.

11	1	Comprsn Band XR656 Plus Compression Band
12	1	TABLE REAR FOOT PEDAL Table Rear Foot Pedal
13	1	<p>656 Plus Ext WS w 3 Grids Discovery XR656 Plus Tilting Wallstand with Extended Length Arm Discovery XR656 Plus Digital Tilting Wallstand is designed for radiography applications with the patient standing, sitting or lying on a gurney. The Discovery XR656 Plus wallstand is designed for use with GE's exclusive FlashPad wireless digital detector, overhead tube suspension, ion chamber and removable non-reciprocating grid. The FlashPad wireless digital detector can be operated docked, tethered, or in a wireless digital cassette mode</p> <ul style="list-style-type: none">• The wallstand is motorized Electromagnetic braking secures vertical motion.• Motorized receptor tilting controlled with either IR remote control or hand switches located on the arm.• Vertical motorization of the wallstand with foot switch or IR remote control.• Auto-tracking from the overhead tube suspension.

- Graphic outlines of image sizes and ion chamber scan areas on the front panel enhance accuracy and safety Preparation is fast and simple for better patient throughput.
 - The wallstand tilts from -20 degrees to 90 degrees.
 - The wallstand grids are removable from the side.
 - It is configurable for either left or right side insertion.
 - The wallstand comes with one removable grid
 - 130 cm Focus Grid
 - 180 cm (72 in.) focus grid with a SID range of 145 cm - 245 cm (70 lines/cm, 13:1 ratio)
 - Automatic Exposure Control (AEC) utilizes three ion chamber sensors, which are mounted between the patient and digital detector.
 - Automatic Exposure Control (AEC) on this extended arm utilizes three ion chamber sensors which are mounted between the patient and digital detector on the wallstand. The addition of the 4th ion chamber provides for supine chest and abdominal imaging when a stretcher is used (anytime the wallstand housing is horizontal). The standard right, left and center ion chambers are available when the wallstand is in the vertical position.
 - Includes patient handgrips and a lateral support bar.
- 100 cm Focus Wallstand Grid - 100 cm (40 inch) focus - 13:1 ratio, 70 lines/cm

- 14 1
- Dual Energy - Chest and Abdomen
Dual Energy Subtraction Chest and Abdomen
Dual Energy Subtraction is an Advanced Clinical Application.
- Provides three diagnostic images:
Standard
 - Soft tissue only
 - Bone only
 - Technique is available for PA/AP chest and abdomen imaging.

- 15 1
- VolumeRAD Table and Wallstand - 4 Sweeps
VolumeRAD at Wallstand and Table with Extended Arm (4 sweeps)
VolumeRAD is an advanced application for radiographic X-ray imaging. It allows retrospective reconstruction of a number of planes (slices) from a series of low dose exposures acquired within a limited angular range. These slices show anatomical structures at different depths
- This application supports the following four sweeps at wallstand:
Table Horizontal sweep of the tube over the table for prone or supine imaging
 - Vertical tube sweep at the wallstand for upright imaging
 - Horizontal sweep of the tube over the stretcher for prone or supine imaging (requires the selection of the extended arm wallstand)
 - Lateral tube sweep at the wallstand for cross table imaging over the stretcher for prone or supine

imaging

16 1

Auto Image Paste for the Discovery
XR656 Wallstand and Table
Wallstand Auto Image Paste for Spine
and Long Bone Imaging

- Fully automated acquisition and processing of a series of images with user defined start and stop locations on the anatomical regions of interest.
- Average acquisition time for a 3-image exam (90 cm coverage) is less than 10 seconds. Image pasting and processing time for a 3-image exam is less than 15 seconds from the last exposure.
- Allows 2 to 5 images to be pasted together with a maximum range of 150 cm.
- Includes imaging of the spine for scoliosis evaluation and imaging of the legs for orthopedic evaluations.
- Supports anatomies/view combinations of Spine Antero-posterior, Spine Postero anterior, Spine Lateral, Leg Antero-posterior, Leg Postero-anterior.
- Includes a patient stand with screen to keep the patient comfortable during acquisition.
- Table Auto Image Paste for Spine and Long Bone Imaging Fully automated acquisition and processing of a series of images with user defined start and stop locations on the anatomical regions of interest.

- Average acquisition time for a 3-image exam (90 cm coverage) is less than 10 seconds. Image pasting and processing time for a 3-image exam is less than 15 seconds from last exposure.
- Allows 2 to 3 images to be pasted together with a maximum range of 100 cm.
- Includes imaging of the spine for scoliosis evaluation and imaging of the legs for orthopedic evaluations.
- Supports anatomies/view combinations of Spine Antero-posterior, Spine Postero anterior, Spine Lateral, Leg Antero-posterior.

17 1

Infrared Remote
Infrared Remote Control
The infrared remote optionally allows remote control of several system features (e.g., collimator light, FOV selection, wallstand and OTS motorization).

18 1

Grid Holder Assembly

- Grid Holder Assembly Extra Wall Mounting Grid Holder for up to 4 Grids

19 1

PATIENT ORI-TAG PERMANENT
PO TAG REMOVE PERMANENT

20	1	UPS Uninterrupted Power Supply
21	1	<p>25 KAIC X-Ray Main Disconnect Panel 80 Amp, 480 V / 208 V</p> <p>25 KAIC X-Ray Main Disconnect Panel 80 Amp, 480 V / 208 V</p> <p>FEATURES/BENEFITS</p> <ul style="list-style-type: none"> • Serves as the main power disconnect between the X-Ray system and the facility 480V or 208V power source • Provides emergency shut down, undervoltage protection and overcurrent protection for the X-Ray power distribution cabinet • Standardized design provides a platform for future upgrades of the system • Offers a number of advantages by combining a variety of individual components into a single pre-engineered and factory tested panel • UL and cUL listed for compliance with NEC Article 100 and Article 110-3 • Remote emergency off pushbutton located by X-Ray control provides immediate shut down of the entire system to comply with NEC required disconnecting means • Surface or semi-flush mounting <p>SPECIFICATIONS</p> <ul style="list-style-type: none"> • Dimensions (H x W x D): 48" x 20" x 6.68" • Weight: 80 lbs. • Mounting: via keyhole slots; Width is 16" on centers, Height is 45.5" on

centers

COMPATIBILITY

- GE Three Phase X-Ray generators
- NOTES: Customer is responsible for rigging and arranging for installation with a certified electrician
- ITEM IS NON-RETURNABLE AND NON-REFUNDABLE

22 1

VolumeRAD TiP Training: 4 Days Onsite (3 Days + 1 Day)
VolumeRAD TiP Training: 4 Days Onsite (3 Days + 1 Day)
One 3 day and one 1 day TiP onsite applications training visit covering the VolumeRAD software functionality. Onsite training provided from 8AM to 5PM, Monday through Friday. Includes T&L expenses.
This training program must be scheduled and completed within 12 months after the date of product delivery.

23 1

Auto Image Paste TiP Training: 2 Consecutive Days Onsite
Auto Image Paste TiP Training: 2 Consecutive Days Onsite
Two consecutive days of TiP Applications Onsite X-ray training covering the Auto Image Paste function. Onsite training provided from 8AM to 5PM, Monday through Friday. Includes T&L expenses.
This training program must be scheduled and completed within 12 months after the date of product delivery.

24	1	X-ray Training, 4 Consecutive Days Onsite X-ray TiP Training, 4 Consecutive Days Onsite Four Days X-ray Onsite Training provided from 8AM to 5PM, Monday through Friday. Includes T&L expenses. Days provided consecutively. This training program must be scheduled and completed within 12 months after the date of product delivery.
----	---	---

25	1	WALLSTAND CABLE SELECT Connection cables for the wallstand
----	---	---

26	1	Table cable select Table cable select
----	---	--

27	1	Rigging for deinstallation of Philips system
----	---	---

Philips Easy Diagnost Trade-In

Options

(These items are not included in the total quotation amount)

Item No.	Qty	Catalog No.	Description
28	1		X-ray TiP Training, 2 Consecutive Days Onsite X-ray TiP Training, 2 Consecutive Days Onsite Two Day X-ray Onsite Training provided from 8AM to 5PM, Monday through Friday. Includes T&L expenses. Days provided consecutively. This training program must be scheduled and completed within 12 months after the date of product delivery.
29	1		REV DIGITAL RAD SYSTEM IR Revolution Digital Rad Systems (Class/Lab) The Revolution Digital Rad Systems service training consists of an online course & 1 week of in-residence class and labs. The systems taught in the course include: Revolution XQI, Revolution XRd (Gipeto), Revolution XRd-2x (Lightning), Definium 8000 (Thunder), Optima XR640, and Discovery XR650. This course