

Identification		Tracking and Location								Dimensions / Mntg.				Electrical								Piped Utilities								Env	
Eq.Tag	Equipment Description	Quantity	Existing	New	Future	Existing Room	Existing Room Number	New Room	New Room Number	Width (in)	Depth (in)	Height (in)	Mounting Location	Number of Outlets	Hard Wire / Plug	Volt	Amp	Ded. Circuit	Watts	Emergency	UPS	Cold Water	Hot Water	Drain	Steam	Comp. Air	Vacuum	Gases	Exhaust	Remarks	
<b>NOTES:</b> 1    Equipment tags not shown were not used. 2    Loose equipment is shown in survey for information only; future locations not tracked. 3    Equipment tag based on survey location																															
Chemistry-Hematology																															
CH1	Coagulation Centrifuge		X			Hematology	3B136	Chemistry	3B134	20	24	16	B - Bench	1	P	120	7	TBD		Yes										CH1 or CH2 to be on emergency power	
CH2	Coagulation Centrifuge		X			Hematology	3B136	Chemistry	3B134	20	24	16	B - Bench	1		120	7	TBD		Yes											
CH3	Cell Video Camera System		X			Hematology	3B136	Chemistry	3B134	21	24	24	B - Bench	1	P	120															
CH4	Cellular Analysis System		X			Hematology	3B136	Chemistry	3B134	72	36	60	F - Floor		P	120	15	X	720	Yes				Yes						Width listed includes a computer system	
CH5	Cellular Analysis System		X			Hematology	3B136	Chemistry	3B134	72	36	60	F - Floor		P	120	15	X	720	Yes				Yes						Width listed includes a computer system	
CH6	Not used																														
CH7	Not used																														
CH8	Not used																														
CH9	Cytocentrifuge		X			Hematology	3B136	Chemistry	3B134	16	24	24	B - Bench	1	P	120															
CH10	Refrigerator		X			Hematology	3B136	Chemistry	3B134	80	36	74	F - Floor	1	P	120	12	X		Yes											
CH11	Freezer		X			Hematology	3B136	Serology	3B135	30	36	72	F - Floor	1	P	120				Yes											
CH12	Freezer		X			Hematology	3B136	Chemistry	3B134	36	40	72	F - Floor	1	P	120				Yes											
CH13	Evolution Coagulation		X			Hematology	3B136	Chemistry	3B134	36	26	56	F - Floor	1	P	120				Yes										With a workstation 34w 24d 30h	
CH14	Compact Coagulation		X			Hematology	3B136	Chemistry	3B134	45	29	48	F - Floor	1	P	120				Yes										With a workstation 34w 24d 30h	
CH15	Chemistry Analyzer		X			Chemistry	3B134			70	40	60	F - Floor																	Equipment to be deleted	
CH16	Chemistry Analyzer		X			Chemistry	3B134			70	40	60	F - Floor																	Equipment to be deleted	
CH17	Pure Water Systems with tank		X			Chemistry	3B134	Serology	3B135	60	30	24 -	BenchW - Wall		P							Yes		Yes						Two water purifiers and a tank.	
CH18	Not used																														
CH19	Centrifuge		X			Chemistry	3B134			20	24	16	B - Bench	1	P	120	7	TBD												replaced by the automated centrifuge	
CH20	Auromated Fecal Occult Analyzer		X			Chemistry	3B134	Chemistry	3B134	12	22	15		1	P	100-240	2.1-1.3							yes						cup sink; Loose; 50/60 Hz	
CH21	Chemistry Analyzer		X			Chemistry	3B134	Chemistry	3B134	92	32	64		1	P	See Remarks		X												200-240 volts, 20 amp, 60 Hz	
CH22	Chemistry Analyzer		X			Chemistry	3B134	Chemistry	3B134	92	32	64		1	P	See Remarks		X												200-240 volts, 20 amp, 60 Hz	
CH23	Centrifuge		X			Chemistry	3B134	Serology	3B135	20	24	16	B - Bench		P	120	7	TBD												replaced by the automated centrifuge	
CH24	Centrifuge		X			Chemistry	3B134	Serology	3B135	20	24	16	B - Bench		P	120	7	TBD												replaced by the automated centrifuge	
CH25	Centrifuge		X			Chemistry	3B134	Serology	3B135	19	22	14		1	P	120	4.5													replaced by the automated centrifuge	
CH26	Refrigerator		X			Chemistry	3B134	Dry Goods	3B134B	80	36	74	F - Floor	1	P	120	12	X		Yes											
CH27	Refrigerator / Freezer		X			Chemistry	3B134	Dry Goods	3B134B	57	36	76	F - Floor	1	P	115	15.9			Yes											
CH28	Analyzer HDAIC		X			Chemistry	3B134	Chemistry	3B134	22	22	19	B - Bench	1	P	100-240								yes						cup sink; 180 VA 50/60 Hz; Analyzer Hemoglobin G8, Tosoh Bioscience, Dennis Smith 1-800-695-6550 ext. 4342	
CH29	Freezer		X			Chemistry	3B134	Serology	3B135	32	30	66	F - Floor							Yes											
CH30	-70 Freezer		X			Chemistry	3B134	Serology	3B135	36	36	74	F - Floor	1	P	115	16			Yes										60 Hz, 1 phase	
CH31	-20 Freezer		X			Chemistry	3B134	Serology	3B135	52	39	76	F - Floor	1	P	120	12	X		Yes											
CH32	Undercounter Refrigerator		X			Chemistry	3B134	Serology	3B135	28	27	33	C - Undercounter							Yes										info not available. See photo	
CH33	Electrophoresis Instrument, Densitometer @ computer 30" x 20"		X			Chemistry	3B134	Serology	3B135	56	24	14	B - Bench		P															Dimensions include adjacent water baths; Spife 3000 Helena Industries, customer #94740537, serial #580226802, 1-800-231-5663, Kathleen Morris 540-336-3778	
CH34	Water Purification Assembly		X			Chemistry	3B134	Chemistry	3B134	40	24	80	W - Wall	4	P	120						Yes	Yes							(2) Elga DV25, 1 additional Elga	
CH35	Cellular Analysis System		X			Hematology	3B136	Chemistry	3B134	36	31	72																		780 LB	
CH36	Chemistry Analyzer	3	X			Chemistry	3B134	AUTO LINE	3B134							See Remarks		X												200-240 volts, 20 amp, 50-60 Hz	



[illegible]

Frozen Sectioning							3B131																						
FS2	Flammables Storage Cabinet		X			Frozen Section	3B131	Frozen Section	3B131	36	26	36	F - Floor																
FS4	Cryostat		X			Frozen Section	3B131	Frozen Section	3B131	30	30	44	F - Floor	1	P	120		1600 VA											45 A eff. Max. start-up current for 5 sec
FS5	Cryostat		X			Frozen Section	3B131	Frozen Section	3B131	31	30	42	F - Floor	1	P	120		1800 VA											
FS6	Carbon Filter Ductless Hood		X			Frozen Section	3B131			27	14	18	B - Bench	1	P	120													to be discarded after renovation
FS7	BSC (4')			X				Frozen Section	3B131	53"												X		X			X	X	sink
FS8	Fume Hood			X				Frozen Section	3B131													X		X				X	sink for formalyn xylene
Histology							3B132																						
ST1	Autoclave		X			Sterilization	3B134A																						
ST2	Flammables Storage Cabinet		X			Sterilization	3B134A	Histology	3B132	36	24	36	F - Floor																note: short (3' high) cabinet
ST3	Solvent Recycler		X			Sterilization	3B134A	Histology	3B132	33	24	68	F - Floor																Adjacent storage for liquid containers
FS1	Stainer		X			Frozen Section	3B131	Histology	3B132	84	30	18	B - Bench		P														Width includes +/- 36" wide stainer +
FS3	Undercounter Refrigerator		X			Frozen Section	3B131	Histology	3B132	24	27	35	C - Undercount	1	P	115	4.3												
H1	Refrigerator		X			Histology	3B132	Histology	3B132	28	33	78	F - Floor	1	P	120													
H2	Microtome		X			Histology	3B132	Histology	3B132	18	22	14	B - Bench	1	P	120													
H3	Microtome		X			Histology	3B132	Histology	3B132	18	22	14	B - Bench	1	P	120													
H4	Microtome		X			Histology	3B132	Histology	3B132	18	22	14	B - Bench	1	P	120													
H5	Undercounter Incubator (oven)		X			Histology	3B132	Histology	3B132	21	24	30	C - Undercount	1	P	120													
H6	Biological Safety Cabinet, Class II, B2		X			Histology	3B132	Histology	3B132	78	32	80	F - Floor			120	14	X				X		X			X	X	has a gas fixture, not connected
H7	Undercounter Freezer		X			Histology	3B132	Histology	3B132	29	24	35	F - Floor	1	P	120													
H8	Tissue Tek		X			Histology	3B132	Histology	3B132	48	24	63	B - Bench	1	P	120	11	X											
H9	Tissue Tek		X			Histology	3B132	Histology	3B132	24	24	55	F - Floor	1	P	120	11	X											Stacked
H10	Flammables Storage Cabinet		X			Histology	3B132	Histology	3B132	44	20	67	F - Floor																
H11	Tissue tek TEC Embedder		X			Histology	3B132	Histology	3B132	45	32	18	B - Bench	1	P	120													
H12	Coverslipper		X			Histology	3B132	Histology	3B132	26	24	26	B - Bench	1	P	120													
H13	Autostainer		X			Histology	3B132	Histology	3B132	43	24	20	B - Bench	1	P	120							Yes		Yes				
H14	Electric Oven		X			Histology	3B132	Histology	3B132	21	20	32	C - Undercount	1	P	120													
H15	Auto Slide Stainer & Coverslipper			X				Histology	3B132																X				To floor drain, replaces H12 & H13
H16	Embedding Center				X			Histology	3B132																				
Cytology							3B132A																						
C1	Stainer / Slide maker		X			Cytology	3B132A	Histology	3B132	44	24	24	B - Bench	1	P	120													connected to sink; alchocol, fresh
C2	Refrigerator		X			Cytology	3B132A	Histology	3B132	28	33	78	F - Floor	1	P	120													
C3	Centrifuge		X			Cytology	3B133	Histology	3B132	15	22	8	B - Bench	1	P	120													
C4	Centrifuge		X			Cytology	3B133	Histology	3B132	24	24	18	B - Bench	1	P	120													
C5	Thin Prep 2000		X			Cytology	3B133	Histology	3B132	16	16	20	- BenchF - Floo	2	P	120													with vacuum pump below
C6	Tripeth Imaging		X			Cytology	3B133	Histology	3B132	15	24	16	B - Bench	1	P	120													
C7	Tripath Respstain		X			Cytology	3B133	Histology	3B132	72	30	51	B - Bench	4	P	120													With vacuum pump, with PC, on
C8	Biological Safety Cabinet, Class II, B2		X			Cytology	3B133	Histology	3B132	77	30	92	F - Floor			120	14	X										X	
C9	Auto Slide Stainer & Coverslipper				X			Histology	3B132																X				To floor drain, replaces C1

