

SALIENT CHARACTERISTICS FOR UROLOGY C-ARM TABLE

Specifications/Salient Characteristics:

Form:

- The Urology table must be compatible with the Biodex 800 Urology C-Arm. This table is to be used directly with the C-ARM.
- Urology table acquisition must include a rail mounted Phenolic Armboard and Pneumatic, Boot-Style Stirrups.

Fit:

- Compatible with the Biodex 800 Urology C-Arm
-
- Compatible with existing supporting supplies and disposables. Extremely low and high working positions.
- Must be able to support a 500 lb. patient
- Joint adjustment of headrest
- Flexible patient preparation by means of removable components
- Permits conventional X-raying during OR
- Ergonomic Y-shaped base for optimal access to patient in gynecology/urology
- Optimal docking capabilities for robot-assisted surgery

Function:

The Urology table is required to operate with the C-ARM mentioned above.

Installation: Vendor shall install the C-Arm Urology table and all associated components.