

## **SALIENT CHARACTERISTICS FOR UROLOGY C-ARM TABLE**

### **Specifications/Salient Characteristics:**

#### **Form:**

- The Urology table must be compatible with the OEC 9900 C-ARM. This table is to be used directly with the C-ARM.
- Urology table acquisition must include a rail mounted Phenolic Armboard and Pneumatic, Boot-Style Stirrups.

#### **Fit:**

- Compatible with the OEC 9900 C-ARM.
- 
- Compatible with existing supporting supplies and disposables. Extremely low and high working positions.
- Must be able to support a 500 lb. patient
- Joint adjustment of headrest
- Flexible patient preparation by means of removable components
- Permits conventional X-raying during OR
- Ergonomic Y-shaped base for optimal access to patient in gynecology/urology
- Optimal docking capabilities for robot-assisted surgery

#### **Function:**

The Urology table is required to operate with the C-ARM mentioned above.

**Installation:** Vendor shall install the C-Arm Urology table and all associated components.