

SHIP TO:
 CHIEF,A&MM B82067
 V.A. Medical Center
 V. A. WESTSIDE (537)
 2030 W. TAYLOR ST.
 CHICAGO, IL 60612

P.O.# 537-B82067

Item	Part Number - Description	Qty
1	L15600 - X-Porte Ultrasound Kiosk System with Stand	1
<p>The SonoSite X-Porte is a high resolution, all digital ultrasound kiosk designed to meet the needs of the demanding clinical environment. The X-Porte Kiosk includes: a capacitive multi-touch panel user interface, adjustable high resolution 19" LED monitor for image display and review, a system stand featuring a triple transducer connection, battery pack and a storage basket.</p> <p>In addition to high resolution 2D imaging the system includes the following imaging modes: M-mode, Pulsed Wave Doppler (update and simultaneous), Continuous Wave Doppler, Velocity Color Doppler, Color Power Doppler and Tissue Doppler Imaging.</p> <p>The X-Porte kiosk features proprietary software technologies: XDI – Extreme Definition Imaging, Steep Needle Profiling, SonoHD2 Imaging Technology, SonoMB® Multibeam Technology, SonoADAPTTM Tissue Optimization Technology, Pulse Inversion Tissue Harmonic Imaging, Color HDTM Technology, AutoGain (see key for Technology descriptions) ,</p> <p>For maximum imaging performance the X-Porte visualization tool features proprietary software technologies: XDI – Extreme Definition Imaging, Steep Needle Profiling, SonoHD2 Imaging Technology, SonoMB Multibeam Technology, SonoADAPT Tissue Optimization Technology, Pulse Inversion Tissue Harmonic Imaging, Color HD Technology, AutoGain</p>		

Item	Part Number - Description	Qty
------	---------------------------	-----

- XDI Extreme Definition Imaging -shapes the ultrasound beam to pinpoint precision. As a result, tissue resolution is improved and side lobe image artifacts are minimized to produce images with less haze and increased clarity.
- Steep Needle Profiling Technology - enhances needle visualization while maintaining image quality of the target and surrounding anatomy.
- SonoHD2 Imaging Technology - reduces speckle noise and image artifacts while enhancing tissue margins, thereby dramatically improving contrast and structural resolution.
- SonoMB MultiBeam Technology - improves contrast resolution; tissue differentiation of small structures and enhances border delineation.
- SonoADAPT Tissue Optimization Technology – provides continuous automatic adaptive manipulation of multiple imaging parameters to reduce manual time consuming system adjustments.
- Pulse Inversion Tissue Harmonic Imaging, utilizes higher frequencies from the ultrasound signal to create high resolution images with dramatically cleaner and sharper contrast between structures.
- Color HD Technology- increasing color performance, sensitivity and frame rates for more diagnostic information.
- AutoGain – provides automatic 2D gain image optimization with the push of a button.

The Triple Transducers Quick Connect allows 3 transducers to be simultaneously connected to the system, for efficient selection of transducers via the touchscreen control panel.

Image Storage and Connectivity is accomplished through: - The kiosk has 64GB of on-board Flash memory to retain enough data for 320 exams and features thumbnail review of saved images and video clips. - 6 USB ports for sharing of images (jpeg and bmp) and video clips (AVI) to a PC or Mac. - Embedded Digital Video Recorder for customized digital video capture. - DICOM package with the capability to store image and clips and interface within your PACS network with the complete offering of DICOM Print, Store, Modality Worklist, Perform Procedure Step (PPS), and Storage Commitment. - Wireless capabilities (B and G networking) for transmission of images and information wirelessly. (Subject to Country Certification) - HDMI and VGA out capability

The Kiosk features SonoSite's Visual Guides, a library of animated learning tutorials are featured on-board the kiosk for just in time learning and reference. The learning tutorials are displayed on the control panel so users can view the tutorials and scan simultaneously. Visual Guide packages included are: Ultrasound Basics, Acute Care Procedures, MSK and Anesthesia.

X-Porte's reliability and durability are backed by a five-year standard warranty, which is included in the price of the system. Software upgrades to the main system are accomplished through a USB ports on the system.

(1) One Transducer is required with each system purchase. Transducer selection may be changed, and transducers may be added or removed as required to support specific clinical applications.

2	P16038 - HFL38XP / 13-6 MHZ Transducer, ROHS	1
	HFL38xp/13-6 MHz multi-frequency, broadband, 38 mm linear array transducer. This transducer is biopsy compatible.	
3	P14566 - L25XP / 13-6 MHZ Transducer	1
	L25xp/13-6 MHz multi-frequency, broadband, 25 mm linear array transducer. This transducer is biopsy compatible.	

Item	Part Number - Description	Qty
4	P19620 - TEExp/ 8-3 MHZ Transducer TEExp / 8-3 MHz transducer with 2D, M-mode, color flow Doppler, Pulsed Wave and Continuous Wave Doppler. Supports cardiac exams. Maximum depth is 18 cm	1
5	P22557 - rP19XP / 5-1 MHZ Transducer rP19x / 5-1 MHz multi-frequency, broadband, 19 mm phased array transducer featuring DirectClear technology (patent-pending transducer design that increases resolution and penetration without compromising SonoSite's renowned durability).	1
6	P20870 - TEEp/xi Transducer Bite Hole Detector Kit The purpose of the bite hole detecto is to test the endoscopic shaft and the transducer tip for bite holes after procedures. Kit includes: - Bite Hole Indicator meter with leads - Endoscope Puncture Test Fork - A conductor plat (one is included with bite hold indicator kit)	1
7	P19617 - C35XP 8-3MHZ Transducer C35xp/8-3 MHz multi-frequency, broadband, 35 mm curved array transducer. This transducer is biopsy compatible.	1
8	P14567 - HFL50XP / 15-6 MHZ Transducer HFL50xp/15-6 MHz multi-frequency, broadband, 50 mm linear array transducer. This transducer is biopsy compatible.	1
9	P16038 - HFL38XP / 13-6 MHZ Transducer, ROHS HFL38xp/13-6 MHz multi-frequency, broadband, 38 mm linear array transducer. This transducer is biopsy compatible.	1
10	P18657 - HSL25XP / 13-6 MHZ Transducer, ROHS HSL25xp/13-6 MHz multi-frequency, broadband, 25 mm linear array transducer for nerve, musculoskeletal, arterial, venous, lung, superficial and ophthalmic exams. This transducer has a maximum depth of penetration of 6 cm. Transducers are designed for 1 meter drop test durability.	1
11	P14561 - C60XP / 5-2 MHZ Transducer C60xp/5-2 MHz multi-frequency, broadband, 60 mm curved array transducer for abdominal, obstetrics, gynecology, nerve and musculoskeletal exams. This transducer is biopsy compatible with a maximum depth of penetration of 30 cm. Transducers are designed for 1 meter drop test durability.	1
12	L12830 - Powerpark System	1

Item	Part Number - Description	Qty
	PowerPark docking station for the Universal stand. This stand module and docking station allows the ultrasound systems to charge without plugging in the power cord.	
13	P22511 - Sony UP-D711MD B&W Printer	1
	Sony UP-D711MD compact, medical grade digital black and white printer with USB Interface. Powered directly by the connected ultrasound system battery, eliminating the need for additional AC outlets.	
14	P14689 - USB Footswitch	1
	Hands-free operation with a customizable two pedal footswitch with USB connection. Allows user customization for high use functions such as Freeze, Save, Save Clip, Record, or Print.	
15	P16620 - Service Manual, X-Porte	1
	X-Porte Ultrasound System Service Manual	
16	P14645 - User Guide, X-Porte, English (Paper)	1
	X-Porte User Guide (CD Only), English	
17	P16620 - Service Manual, X-Porte	1
	X-Porte Ultrasound System Service Manual	
18	P20262 - Kit, ECG Slave Cable and Adapter	1
	The ECG slave cable is an adapter for the existing SonoSite ECG Module, allowing an ECG signal to be slaved directly from a patient monitor into the SonoSite Ultrasound system. With the ECG slave cable, only one set of ECG leads is required.	
	This item requires the SonoSite ECG Option (Sold Separately)	
19	P15169 - ECG Option	1
	ECG option includes 6 ft. (1.8 m) cable connects to an ECG port on the back of X-Porte system. Allows for display of ECG trace on clinical display. Includes a 2 ft. 3-leadwire set, for a total cable length of 8 ft.	
20	P14145 - SonoSite Ultrasound Phantom	1
	Custom ultrasound phantom manufactured by Blue Phantom that contains one simulated vessel and one simulated nerve. An excellent training tool to help clinicians develop, practice and maintain the skills necessary to use ultrasound for guiding procedures.	
21	L12920 - 18 Gauge Needle Guide Starter Kit For The L25 Linear Array Transducer	1
	18 gauge Needle Guide Starter Kit for the L25/13-6 linear array transducer. Each starter kit contains 5 sterile procedure kits including: needle guide, needle guide bracket, telescopically-folded CIV-Flex transducer cover, gel packet and colored elastic bands. Sheath Size: 10.2 x 147cm.	

Item	Part Number - Description	Qty
22	MT_XP_TRADE_IN - M-Turbo to X-Porte Trade In Credit Serial No:SN 03BCNV	1
23	5-YR-WARRANTY - 60 MONTHS - STANDARD COVERAGE WARRANTY Please refer to Warranty Schedule. 60 MONTHS - STANDARD COVERAGE WARRANTY	1
24	5-YR_SW_MAINT - X-PORTE SOFTWARE MAINTENANCE 5-YR Warranty and FUJIFILM SonoSite Software and Firmware Maintenance - 60 Months - Standard Warranty and Maintenance for FUJIFILM SonoSite Software and Firmware (excludes optional third party software).	1
25	5-YR_SW_MAINT - X-PORTE SOFTWARE MAINTENANCE 5-YR Warranty and FUJIFILM SonoSite Software and Firmware Maintenance - 60 Months - Standard Warranty and Maintenance for FUJIFILM SonoSite Software and Firmware (excludes optional third party software).	1