

<div><div><div><div><div></div><div>8 7/8"</div></div><div>PLAN</div><div><div><div><div></div><div>1'-7 7/8"</div></div><div>FRONT</div><div><div><div><div></div><div>1'-7 13/16"</div></div><div>SIDE</div></div></div></div><div><div><div><div>VIT (CPU)</div><div>HEAT OUTPUT (BTU/HR)</div><div>1,000.00</div><div>WEIGHT (LBS)</div><div>42.00</div></div></div></div></div><div><div><div><div></div><div>1'-8 1/16"</div></div><div>PLAN</div><div><div><div><div></div><div>1'-9 11/16"</div></div><div>FRONT</div><div><div><div><div></div><div>9 5/16"</div></div><div>SIDE</div></div></div></div></div></div></div></div></div></div>						<div><div><div>TOSHIBA</div><div>Leading Innovation >>></div><div>TOSHIBA AMERICA MEDICAL SYSTEMS INC.</div><div>www.medical.toshiba.com</div></div></div>	
<div><div><div>1</div><div>VITREA WORKSTATION</div></div><div>SCALE: 1" = 1' - 0"</div><div>REVISION: 04-01-08</div></div>							
						<div><div><div>DEPT. OF VETERAN AFFAIRS</div><div>NY HARBOR HEALTH CARE</div><div>AQUILION - 64</div><div>800 POLY PL.</div><div>BROOKLYN, NY 11206</div></div><div><div>THESE TOSHIBA PLANS ARE FOR INFORMATIONAL PURPOSES ONLY AND SHALL NOT BE USED FOR ANY PURPOSE OTHER THAN THAT AGREED UPON BETWEEN TOSHIBA AND THE CUSTOMER. THESE SITE PLANS ARE NOT TO BE USED FOR CONSTRUCTION PURPOSES.</div><div>DATE: 10/22/10</div><div>SCALE: AS NOTED</div><div>DRAWN: J.A.D.</div><div>QUOTE: 98050</div><div>PROJECT: 100010904CTP</div></div></div>	
						A4	