

Appendix I

The following laboratory pictures, list of equipment and list of chemicals to be moved are provided to support the solicitation of bids for relocating the Laboratory of Proteomics and Protein Science group from the Research program (688/151) at the Washington DC VAMC, 50 Irving St. NW, Washington DC to the Research program (512/151) at the Baltimore VAMC, 10 N. Green St., Baltimore MD.

Table of contents:

Pictures of currently occupied space at the Washington DC VAMC with equipment	2-7
List of equipment to be moved from Washington DC to Baltimore	8-11
List of chemicals to be moved from Washington DC to Baltimore	12-20
(DOT shipping information included for regulated chemicals)	

Laboratory pictures

Current space occupied by the Laboratory of Proteomics and Protein Science

Washington DC VA Medical Center

50 Irving St, NW

Washington DC 20422

Notes: Offices 1 and two and labs 1, 2, 3, 4 are all on the main (upper level) of the Research Building
Office 2 is on the lower level of the Research Building at the DC VA Medical Center

Access to the outside is either (1) down a freight elevator to the lower level and out a double door to the delivery area outside the research building. The pavement is a little uneven and there is no loading dock, so trucks would need a lift gate.

Or 2) across the pedestrian bridge that connects the upper level of the Research building to the hospital, and through the main hospital to the hospital loading dock

Laboratory I
400 square feet

Pictures taken from left
Of the door,
sweeping right



Laboratory 2
400 square feet

Pictures taken from left
Of the door,
sweeping right



Laboratory 3
400 Square feet

Pictures taken from left
Of the door,
sweeping right



Laboratory 4
400 Square feet

Pictures taken from left
Of the door,
sweeping right





Office 1
Approx. 100 sq ft

Office 2
100 sq ft



Office 3
100 sq ft



List of Equipment to be moved

Description	DC VA location	Baltimore location	Power Requirements	Footprint	manufacturer
Cryostat	1F-125	3C-105	120V/15A	3Ft X 3 Ft	LEICA MICROSYSTEMS INC
VACUOBOX	1F-125	3C-105	120V/15A	1Ft x 1Ft	BRINKMANN INSTRUMENTS INC
CENTRIFUGE, TABLETOP, LOW-SPEED, REFRIGERATED	1F-127	3C-105	120V/15A	1Ft x 1Ft	EPPENDORF SCIENTIFIC INC
TISSUE PUNCHER	1F-129	3C-105	N/A (Shelf)	1Ft x 1Ft	BEECHER INSTRUMENTS INC
MICROSCOPE, LIGHT, LABORATORY	1F-129	3C-105	120V/15A	3Ft W x 2 Ft D	OLYMPUS AMERICA INC
BALANCE, ELECTRONIC	1F-129	3C-105	120V/15A	1ft X 2.5 ft	SARTORIUS CORP
PHOTOCHEMICAL REACTOR	1F-129	3C-105	120V/15A	1Ft x 1Ft	SOUTHERN NEW ENGLAND ULTRA VIOLET C
Water bath	1F-129	3C-105	120V/15A	1Ft x 1Ft	ORGANOMATION ASSOC INC
CENTRIFUGE, TABLETOP	1F-129	3C-105	120V/15A	1ft X 2.5 ftD	EPPENDORF SCIENTIFIC INC
MICROSCOPE, LIGHT, LABORATORY	1F-129	3C-105	120V/15A	2Ft W x 2 Ft D	OLYMPUS AMERICA INC
BATH, TISSUE FLOTATION	1F-129	3C-105	120V/15A	1Ft x 1Ft	NESLAB INSTRUMENTS INC
Water bath	1F-129	3C-105	120V/15A	1Ft x 1Ft	FISHER SCIENTIFIC CO
VORTEX GENIE	1F-129	3C-105	120V/15A	1Ft X 1Ft	SCIENTIFIC SUPPLIES AND TECHNOLOGY
HOT PLATE	1F-129	3C-105	120V/15A	1Ft x 1Ft	BARNSTEAD/THERMOLYNE CORP, DIV SYBR
vortexer	1F-129	3C-105	120V/15A	0.5 Ft x 0.5 Ft	VWR
MICROTOME, ROTARY	1F-129	3C-105	N/A	2.5 ft x 2.5 ft	LEICA MICROSYSTEMS INC
waterbath, circulating	1F-129	3C-105	120V/15A	1Ft X 2.5 Ft	POLYSCIENCES INC
MALDI Mass Spectrometer	1F-129	3C-105	208V/20A/UPS (NEMA L8-20R)	5ft W x 4 ft deep (floor)	Bruker Daltonics
Image Prep System (slide coater)	1F-129	3C-105	120V/15A	2ft W x 2ft d (TC hood)	Bruker Daltonics
Computer workstation	1F-129	3C-105	120V/15A/UPS	3Ft X 2.5 ft	HEWLETT-PACKARD CO
Monitor	1F-129	3C-105	120V/15A/UPS		Hewlett-Packard
INCUBATOR, AEROBIC	1F-136	3C-105	120V/15A (needs gas tank)	26" W X 2Fd (Stacks)	THERMO SCIENTIFIC LABORATORY SAMPLE
INCUBATOR, AEROBIC	1F-136	3C-105	120V/15A (needs gas tank)	26" W X 2Fd (Stacks)	THERMO SCIENTIFIC LABORATORY SAMPLE
THERMOSTAT, ALS	1F-137	3C-105	120V/15A/UPS	2Ft x 2 Ft (stacks)	AGILENT TECHNOLOGIES
PUMP, LABORATORY	1F-137	3C-105	120V/15A/UPS	2Ft x 2 Ft (stacks)	AGILENT TECHNOLOGIES
VACUUM DEGASSER	1F-137	3C-105	120V/15A/UPS	2Ft x 2 Ft (stacks)	AGILENT TECHNOLOGIES
Sample Injector, model G1377A	1F-137	3C-105	120V/15A/UPS	2Ft x 2 Ft (stacks)	AGILENT TECHNOLOGIES
SPECTROMETER, MASS, LABORATORY	1F-137	3C-105	208V/20A/UPS (NEMA L8-20R)	2Ft x 2 Ft (stacks)	AGILENT TECHNOLOGIES
GENERATOR GAS	1F-137	3C-105	120V/15A/UPS	2Ft X 4ft (has quiet cover)	PARKER HANNIFIN CORP
PUMP, LABORATORY	1F-137	3C-105	Part of Mass spectrometer	2ft X 4ft (has quiet cover)	AGILENT TECHNOLOGIES
COMPUTER, MEDICAL SYSTEM	1F-137	3C-105	120V/15A/UPS	2Ft x 2 Ft	HEWLETT-PACKARD CO (COMPUTERS)
COMPUTER, MEDICAL SYSTEM	1F-137	3C-105	120V/15A/UPS	2Ft x 2 Ft (stacks)	HEWLETT-PACKARD CO (COMPUTERS)
MONITOR, MEDICAL SYSTEM	1F-137	3C-105	120V/15A/UPS		HEWLETT-PACKARD CO (COMPUTERS)
MONITOR, MEDICAL SYSTEM	1F-137	3C-105	120V/15A/UPS		HEWLETT-PACKARD CO (COMPUTERS)
DIODE-ARRAY DETECTOR	1F-137	3C-105	120V/15A/UPS	2Ft x 2 Ft (stacks)	AGILENT TECHNOLOGIES
WARMING UNIT, MICROSCOPE SLIDE	1F-137	3C-105	120V/15A	1 ft X 2 Ft	ONCOR INC
FREEZER, LABORATORY	1F-137	3C-105	120V/15A	42" X 3 Ft deep	REVCO SCIENTIFIC INC
Drying oven	storage	3C-105	120V/15A	2 Ft X 2 Ft	BOEKEL SCIENTIFIC INC
PUMP, CIRCULATING-FLUID	1F125	3C-117	120V/15A	1Ft X 1Ft	GARDNER DENVER
Fluorescence Detector	1F-125	3C-117	120V/15A	2Ft x 2 Ft (stacks)	AGILENT TECHNOLOGIES
REFRIGERATOR, LABORATORY	1F-125	3C-117	120V/15A	3 Ft X 3 Ft	REVCO SCIENTIFIC INC
CELL DISRUPTION BOMB	1F-125	3C-117	N/A	1 ft X 1 Ft	PARR INSTRUMENT COMPANY
WATER BATH CIRCULATOR	1F-125	3C-117	120V/15A	1 ft X 2 Ft	NESLAB INSTRUMENTS INC
VESICLE EXTRUDER	1F-125	3C-117	120V/15A	1 ft X 1 Ft	SCIEMA TECHNICAL SERVICES LTD
BATH, TISSUE FLOTATION	1F-125	3C-117	120V/15A	1 ft X 1 Ft	PRECISION MEDICAL INC
BATH, TISSUE FLOTATION	1F-125	3C-117	120V/15A	1 ft X 1 Ft	PRECISION MEDICAL INC
CENTRIFUGE, TABLETOP	1F-125	3C-117	120V/15A	1 ft X 1 Ft	EPPENDORF SCIENTIFIC INC
SONIC DISMEMBRATOR	1F-125	3C-117	120V/15A	1 ft x 2 ft	FISCHER INDUSTRIES INC
COMPUTER, MEDICAL SYSTEM	1F-125	3C-117	120V/15A	1 ft X 2 Ft	COMPAQ COMPUTER CORP
Light Scattering Spectrometer	1F-125	3C-117	120V/15A	2 Ft X 2 Ft	NICOLET BIOMEDICAL INC
CENTRIFUGE, TABLETOP	1F-125	3C-117	120V/15A	0.5 Ft x 0.5 Ft	FISHER SCIENTIFIC CO
PUMP, CIRCULATING-FLUID	1F-125	3C-117	120V/15A	1 ft X 1 Ft	SAVANT INSTRUMENTS INC
REFRIGERATOR, LABORATORY	1F-125	3C-117	120V/15A	3 Ft X 3 Ft	PRECISION MEDICAL INC
Circulating water bath	1F-125	3C-117	120V/15A	1 ft x 2 ft	NESLAB INSTRUMENTS INC

FLUOROMETER	1F-125	3C-117	120V/15A	1 Ft X 1 ft	BIO-RAD LABORATORIES, LIFE SCIENCE GR
CALORIMETER	1F-125	3C-117	120V/15A	2 Ft X 2 Ft	CALORIMETRY SCIENCES CORP
CENTRIFUGE, TABLETOP, LOW-SPEED, REFRIGERATED	1F127	3C-117	120V/15 A	3Ft X 2.5 ft	EPPENDORF SCIENTIFIC INC
ROTATOR	1F127	3C-117	120V/15A	1ft X 0.5 ft	JANKE&KUNKEL
PUMP, CIRCULATING-FLUID	1F127	3C-117	120V/15A	1ft X 1 Ft	GARDNER DENVER
HOT PLATE	1F127	3C-117	120V/15A	1ft x 1ft	BARNSTEAD/THERMOLYNE CORP
MICROPLATE READER	1F-127	3C-117	120V/15A	2.5 ft x 2.5 ft	MOLECULAR DEVICES CORP
SPEED VAC CONCENTRATOR	1F-127	3C-117	120V/15A	2.5 ft x 2.5 ft	THERMO ELECTRON CORP
CENTRIFUGE, TABLETOP, LOW-SPEED, REFRIGERATED	1F-127	3C-117	120V/15A	2 ft x 2.5 ft	THERMO ELECTRON CORP
PUMP, LABORATORY	1F-127	3C-117	120V/15A	1 ft x 1ft	THERMO ELECTRON CORP
MICROPLATE WASHER	1F-127	3C-117	120V/15A	2 ft x 1 ft	BIO-TEK INSTRUMENTS INC
DNA SEQUENCING SYSTEM	1F-127	3C-117	120V/15A	2 ft x 1.5 ft	BIO-RAD LABORATORIES, LIFE SCIENCE GR
COMPUTER, MEDICAL SYSTEM	1F-127	3C-117	120V/15A	1 ft x 2ft	HEWLETT-PACKARD CO (COMPUTERS)
MONITOR, MEDICAL SYSTEM	1F-127	3C-117	120V/15A		VIEWSONICS INC
BATH, TISSUE FLOTATION	1F-127	3C-117	120V/15A	1ft X 1 Ft	BRINKMANN INSTRUMENTS INC
BALANCE, ELECTRONIC	1F-127	3C-117	120V/15A	1 ft x 2 ft	METTLER ELECTRONICS CORP
CENTRIFUGE, TABLETOP	1F-127	3C-117	120V/15A	1ft x 1ft	EPPENDORF SCIENTIFIC INC
CENTRIFUGE, TABLETOP	1F-127	3C-117	120V/15A	0.5 Ft x 0.5 Ft	FISHER SCIENTIFIC CO
VORTEX GENIE	1F-127	3C-117	120V/15A	0.5 ft x 0.5 ft	FISHER INDUSTRIES INC
DNA SEQUENCING SYSTEM	1F-127	3C-117	120V/15A	1 ft x 2 ft	PERKIN ELMER ANALYTICAL INSTRUMENTS
BATH, TISSUE FLOTATION	1F-127	3C-117	120V/15A	1ft X 1 Ft	BRINKMANN INSTRUMENTS INC
VORTEX GENIE	1F-127	3C-117	120V/15A	0.5 Ft x 0.5 Ft	GENIE INDUSTRIES
POWER SUPPLY, REGULATED, ELECTROPHORESIS	1F-127	3C-117	120V/15A	0.5 Ft x 0.5 Ft	LIFE TECHNOLOGIES INC
BATH, TISSUE FLOTATION	1F-127	3C-117	120V/15A	1 ft x 1 ft	PRECISION MEDICAL INC
POWER SUPPLY, REGULATED, ELECTROPHORESIS	1F-127	3C-117	120V/15A	0.5 Ft x 0.5 Ft	LIFE TECHNOLOGIES INC
POWER SUPPLY, REGULATED, ELECTROPHORESIS	1F-127	3C-117	120V/15A	0.5 Ft x 0.5 Ft	LIFE TECHNOLOGIES INC
Pump, circulating	1F-127	3C-117	120V/15A	1 ft x 1 ft	SAVANT INSTRUMENTS INC
vortexer	1F-127	3C-117	120V/15A	0.5 Ft x 0.5 Ft	VWR
vortexer	1F-127	3C-117	120V/15A	0.5 Ft x 0.5 Ft	VWR
OVEN, DRYING	1F-127	3C-117	120V/15A	2 ft X 2 ft	Blue M
CENTRIFUGE, TABLETOP, LOW-SPEED, REFRIGERATED	1F-129	3C-117	120V/15A	1 Ft X 1 Ft	EPPENDORF SCIENTIFIC INC
CENTRIFUGE, TABLETOP	1F-129	3C-117	120V/15A	0.5 Ft x 0.5 Ft	FISHER SCIENTIFIC CO
MIXER, CLINICAL LABORATORY	1F137	3C-117	120V/15A	1 ft x 2 ft	PHARMACIA CORP
HPLC Pump, model 306	1F137	3C-117	120V/15A	1 ft X 1 Ft	GILSON INC
PUMP, CIRCULATING-FLUID	1F137	3C-117	120V/15A	1 ft x 1 ft	GARDNER DENVER
LAMINAR AIR FLOW UNIT, MOBILE	1F-137	3C-117	120 V/15A	4 ft x 2.5 ft	AIRCLEAN SYSTEMS
WARMING UNIT, MICROSCOPE SLIDE	1F-137	3C-117	120 V/15A	2 Ft X 0.5 ft	FISHER SCIENTIFIC CO
PUMP, CIRCULATING-FLUID	1F-137	3C-117	120 V/15A	1ft x 1 ft	PHARMACIA CORP
COMPUTER WORKSTATION, CPU/MONITOR/KEYBOARD	GF147/or 1F137	3C-117	120 V/15A	1 ft x 2 ft	DELL INC
COMPUTER WORKSTATION, CPU/MONITOR/KEYBOARD	GF147/or 1F137	3C-117	120 V/15A	1 ft x 2 ft	DELL INC
HEATER TIMER	storage	3C-117	120 V/15A	1 Ft x1 ft	J-KEM SCIENTIFIC INC
POWER SUPPLY, REGULATED	storage	3C-117	120 V/15A	1 ft x 1 ft	TEKTRONIX INC
VISCOTIMER	storage	3C-117	120 V/15A	1 ft x 2 ft	BRINKMANN INSTRUMENTS INC
CONTROLLER, TEMPERATURE	storage	3C-117	120 V/15A	1 ft x 1ft	ALPHA OMEGA
Diaphragm Pump, Model 31	storage	3C-117	120 V/15A	1 ft x 1 ft	ALCATEL DRYTEL
CENTRIFUGE, TABLETOP	storage	3C-117	120 V/15A	1 ft x 1 ft	WESCOR INC
OSCILLOSCOPE	storage	3C-117	120 V/15A	1.5 ft x 1 ft	TEKTRONIX INC
Drying oven, vacuum	1F-125	3C-117	120 V/15A	1.5 ft x 2 ft	National lab
PCR thermo cycler	1F-127	3C-117	120V/15A	1 ft x 1.5 ft	PERKIN ELMER ANALYTICAL INSTRUMENTS
PCR preparation robot	storage	3C-117?	120 V/15A	2 ft x 2 ft	CORBETT ROBOTICS INC
COMPUTER, MEDICAL SYSTEM	storage	3C-117?	120 V/15A	1 ft x 2 ft	DELL INC
MONITOR, MEDICAL SYSTEM	storage	3C-117?	120 V/15A		DELL INC
FREEZER, LABORATORY	storage	3C-117?	120 V/15A		FISHER SCIENTIFIC CO
minus 80 FREEZER, LABORATORY	1F-120	basement freezer rm	120 V/15A	3 ft x 3 ft	REVCO SCIENTIFIC INC
COMPUTER, MEDICAL SYSTEM	1F-125	Bressler	120 V/15A	42 " X 3 ft	HEWLETT-PACKARD CO (COMPUTERS)
				1 ft x 2 ft	

MONITOR, MEDICAL SYSTEM	1F-125	Bressler	120 V/15A		HEWLETT-PACKARD CO (COMPUTERS)
CALORIMETER	1F-125	Bressler	120 V/15A	1 ft x 2 ft	TA INSTRUMENTS
ULTRACENTRIFUGE	1F-125	Bressler	120V/15A	3 ft x 2 ft	BECKMAN COULTER INC
BALANCE, ELECTRONIC	1F-125	Bressler	120 V/15A	1 ft x 1 ft	METTLER ELECTRONICS CORP
Ink Jet printer	1F-125	Bressler	120 V/15A	2 ft x 1 ft	H-P PRODUCTS INC
ANALYZER, PROTEIN	1F-137	Bressler	120 V/15A	3 ft x 2.5 ft	BIO-RAD LABORATORIES, LIFE SCIENCE GR
COMPUTER, MEDICAL SYSTEM	1F-137	Bressler	120 V/15A	1 ft x 2 ft	HEWLETT-PACKARD CO (COMPUTERS)
MONITOR, MEDICAL SYSTEM	1F-137	Bressler	120 V/15A		NEC CORP
CD Spectropolarimeter, J-715	1F-137	Bressler	120 V/15A	3 ft x 2.5 ft	JASCO INC
Circulating water bath	1F-137	Bressler	120 V/15A	1 ft x 2 ft	NESLAB INSTRUMENTS INC
ULTRACENTRIFUGE	1F-137	Bressler	208V/20 A (NEMA 6-30R)	3 ft x 3 ft (floor)	BECKMAN COULTER INC
VOLUMETRIC IMAGER SYSTEM	Storage/stan	Bressler	120 V/15A	2.5 ft x 2.5 ft	ACTUALITY SYSTEMS INC
COMPUTER, MEDICAL SYSTEM	Storage/stan	Bressler	120V/15A	1ft x 2 ft	DELL INC
MONITOR, MEDICAL SYSTEM	Storage/stan	Bressler			DELL INC
PRINTER, COMPUTER/NETWORK, LASER	offices	offices	120V/15A	1 ft x 1.5 ft	HEWLETT-PACKARD CO (COMPUTERS)
PRINTER, COMPUTER/NETWORK, LASER	offices	offices	120V/15A	1 ft x 1.5 ft	HEWLETT-PACKARD CO (COMPUTERS)
PRINTER, COMPUTER/NETWORK, LASER	offices	offices	120V/15A	1 ft x 1.5 ft	HEWLETT-PACKARD CO (COMPUTERS)
PRINTER, COMPUTER/NETWORK, LASER	offices	offices	120V/15A	1 ft x 1.5 ft	HEWLETT-PACKARD CO (COMPUTERS)
PRINTER, COMPUTER/NETWORK, LASER	offices	offices	120V/15A	1 ft x 1.5 ft	HEWLETT-PACKARD CO (COMPUTERS)
MICROWAVE OVEN, Antigen retrieval	storage	common	120 V/15A	2 ft x 2 ft	CEM CORPORATION
Edge Stereo Microscope)	storage	common	120 V/15A	3 ft x 3 ft	NIKON INC, INSTRUMENT GROUP
CARS system (3 boxes?)	Storage	common	120 V/15A	5 X 10 ft	
Peptide Synthesizer	storage	common or give-away?	120 V/15A	2.5 ft x 2.5 ft	APPLIED BIOSYSTEMS DIV PERKIN ELMER IN
COMPUTER, MEDICAL SYSTEM	storage	common or give-away?	120 V/15A	1 ft x 2 ft	DELL INC
MONITOR, MEDICAL SYSTEM	storage	common or give-away?	120 V/15A		DELL INC
COMPRESSOR, MEDICAL-AIR	1F107	Common?	120 V/15A	1 ft x 1 ft	ULVAC KIKO INC
COMPUTER, MEDICAL SYSTEM	1F107	Common?	120 V/15A	1 ft x 2 ft	DELL INC
MICROSCOPE, ELECTRON, SCANNING	1F-107	Common?	120 V/15A	2 ft x 2 ft	JEOL USA INC
MONITOR, MEDICAL SYSTEM	1F-107	Common?	120 V/15A		DELL INC
MICROSCOPE, LIGHT, LABORATORY	1F-123B	Common?	120 V/15A	1 ft x 2.5 ft	LEICA MICROSYSTEMS INC
CAMERA, CCD HIGH RESOLUTION	1F-123B	Common?	120 V/15A	goes with microscope	QUANTITATIVE IMAGING CORPORATION
Atomic Force Microscope	1F125	common?	120 V/15A	2.5 ft x 2.5 ft	DIGITAL EQUIPMENT CORP
COMPUTER, MEDICAL SYSTEM	1F-125	Common?	120 V/15A	1 ft X 2 ft	HEWLETT-PACKARD CO (COMPUTERS)
MONITOR, MEDICAL SYSTEM	1F-125	common?	120 V/15A		HEWLETT-PACKARD CO (COMPUTERS)
PRINTER, COMPUTER/NETWORK, LASER	1F-125	Common?	120 V/15A	1.5 ft x 1 ft	HEWLETT-PACKARD CO (COMPUTERS)
FRACTION, COLLECTOR	1F-125	Common?	120 V/15A	2 ft x 2 ft (stackable)	AGILENT TECHNOLOGIES
THERMOSTAT, ALS	1F-125	Common?	120 V/15A	2 ft x 2 ft (stackable)	AGILENT TECHNOLOGIES
THERMOSTAT, ALS	1F-125	Common?	120 V/15A	2 ft x 2 ft (stackable)	AGILENT TECHNOLOGIES
AUTOSAMPLER	1F-125	Common?	120 V/15A	2 ft x 2 ft (stackable)	AGILENT TECHNOLOGIES
FRACTION, COLLECTOR	1F-125	Common?	120 V/15A	2 ft x 2 ft (stackable)	AGILENT TECHNOLOGIES
THERMOSTATED COLUMN COMPARTMENT	1F-125	Common?	120 V/15A	2 ft x 2 ft (stackable)	AGILENT TECHNOLOGIES
PUMP, LABORATORY	1F-125	Common?	120 V/15A	2 ft x 2 ft (stackable)	AGILENT TECHNOLOGIES
DIODE-ARRAY DETECTOR	1F-125	Common?	120 V/15A	2 ft x 2 ft (stackable)	AGILENT TECHNOLOGIES
VACUUM DEGASSER	1F-125	Common?	120 V/15A	2 ft x 2 ft (stackable)	AGILENT TECHNOLOGIES
COMPUTER, MEDICAL SYSTEM	1F-125	Common?	120 V/15A	2 ft x 2 ft (stackable)	DELL INC
MONITOR, MEDICAL SYSTEM	1F-125	Common?	120 V/15A		DELL INC
CENTRIFUGE, TABLETOP, LOW-SPEED, REFRIGERATED	storage	common?	120 V/15A	3 ft x 2.5 ft	EPPENDORF SCIENTIFIC INC
SPECTROPHOTOMETER, ULTRAVIOLET/VISIBLE (nanodrop	1F-137	common?	120 V/15A	1 ft x 1 ft	THERMO SCIENTIFIC LABORATORY SAMPLE
BioRAD Gel documentation system	1F-137	common?	120 V/15A	2.5 ft x 2.5 ft	BIO-RAD LABORATORIES, LIFE SCIENCE GR
COMPUTER, MEDICAL SYSTEM	1F-137	common?	120 V/15A	1 ft x 2 ft	HEWLETT-PACKARD CO (COMPUTERS)
MONITOR, MEDICAL SYSTEM	1F-137	common?	120 V/15A		HEWLETT-PACKARD CO (COMPUTERS)
PRINTER, THERMAL	1F-137	common?	120 V/15A	1ft x 1 ft	MITSUBISHI ELECTRONICS AMERICA, MEDIC
REFRIGERATOR, LABORATORY	1F-137	common?	120 V/15A	42" x 3 ft	Fisher
OVEN, DRYING	1F-137	common?	120V/15A	2ft x 2 ft	LINAK
LIGHT CONTROL SYSTEM	common	common?	120V/15A		NEWPORT ELECTRONICS INC

MIXER, CLINICAL LABORATORY	storage	common?	120V/15A		THERMOLYNE CORP
CABINET, BIOLOGICAL SAFETY, CLASS II	1F-127		FISHER HAMILTON LLC	54L933	10071
tissue embedding station	1F-129	Histology lab?	120 V/15A	3FtX 2.5 ft	LEICA MICROSYSTEMS INC
Cryo console, embedding station	1F-129	Histology lab?	120V/15A	2.5 ft x 2.5 ft	LEICA MICROSYSTEMS INC
tissue processor, automatic	1F-137	Histology lab?	120V/15A	4Ft x 4ft	LEICA MICROSYSTEMS INC

List of Chemicals to be moved

Product Name	Manufacturer	CAS Number	Quantity	MSDS on file?	Room	Storage Method ^a	Physical State ^b	Hazardous Material Class ^c				DOT HAZARD Class (section 14 of MSDS)
								H	F	R	SH	
11-Aminoundecanoic acid	Sigma	28164-41-8	100G	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous Goods
1-Octanol	Sigma	111-87-5	100 ml	Yes	1F125	bottle/shelf	L	2	2	0		UN 1993 Class: CBL PG: III
2,6-DI-TERT-BUTYL-P-CRESOL	Aldrich	128-37-0	50 G	Yes	1F125	bottle/shelf	S	1	1	0		
2,6-PYRIDINEDICARBOXYLIC ACID	Aldrich	100-26-5	100g	Yes	1F125	bottle/shelf	S	1	0	0		
5-SULFOSALICYLIC ACID	Sigma	5965-83-3	100 Gm	Yes	1F125	bottle/shelf	S	2	0	0		Not Dangerous Goods
8-Anilo-1-naphthanesulfonic Acid	Sigma	10417-5G-F	5 G	Yes	1F125	bottle/shelf	S	2	0	0		Not Dangerous goods
ADA	Sigma	26239-55-4	25 G	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous Goods
AG 11A8, BioRad	BioRad	68955-09-9	25 G	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
AG 501 X8(D), BioRad	BioRad	68955-09-9	25 G	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
AGAROSE, HMP	Gibco	9012-36-6	500 G	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
AGAROSE, LMP	Gibco	9012-36-6	100 G	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
alpha-cyano-4 hydroxycinnamic acid matrix	bruker	28166-41-8	1 g	Yes	1F125	refrigerator	S	2	0	0		Not Dangerous goods
ALUMINA	Sigma	1344-28-1	1 Kg	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
AMBERLITE, MB-1	Sigma	100915-96-6	25 G	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
AMMONIUM ACETATE	Sigma	631-81-8	100 G	Yes	1F125	bottle/shelf	S	1	1	0		Not Dangerous goods
AMMONIUM SULFATE	Sigma	7783-20-2	1Kg	Yes	1F125	bottle/shelf	S	1	1	0		Not Dangerous goods
ASPARAGINE	Sigma	70-47-3	500 G	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
Aurintricarboxylic acid ammonium salts	Sigma	569-58-4	50 g	Yes	1F125	bottle/shelf	S	2	0	0		Not Dangerous goods
BCA PROTEIN ASSAY KIT	Pierce	7758-99-8	1 kit	Yes	1F125	bottle/shelf	L	2	0	0		Not Dangerous goods
BenchMark Protein Ladder	Life Technologies	151-21-3	0.5 mL	Yes	1F125	freezer	L	1	0	0		
BIO-BEADS, SM-2	BioRad	9052-95-3	100 G	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
BioGel P10	BioRad	25034-58-6	100G	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
BioGel P100	BioRad	25034-58-6	100G	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
BioGel P6DG	BioRad	25034-58-6	100 G	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
B-Lactoglobulin B	Sigma	L8005	25 mg	Yes	1F125	refrigerator	S	0	0	0		Not Dangerous goods
BOILING CHIPS	Aldon Corp	471-34-1	100 G	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
BORAX, anhydrous	Sigma	1330-43-4	500 G	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
BORIC ACID	Sigma	10043-35-3	1 Kg	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
Bovine serum albumin	Sigma	9048-48-8	200 g	Yes	1F125	refrigerator	S	0	0	0		Not Dangerous goods
CALCIUM CHLORIDE	Sigma	10043-52-4	500 G	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
CALCIUM IONOPHORE A23187	Sigma	52665-69-7	10 mg	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
(-) Camphor	Sigma	464-48-2	1 G	Yes	1F125	bottle/shelf	S	2	2	2		UN2717/ class 4.1/ PG III
Camphorsulfonic acid	Sigma	82509-30-6	25 g	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
CAPRIC ACID	Sigma	334-48-5	100 G	Yes	1F125	bottle/shelf	S	2	1	0		Not Dangerous goods
CAPS	Sigma	C2632	25 G	Yes	1F125	bottle/shelf	S	2	0	0		Not Dangerous goods
CAPSO	Sigma	734636-39-5	25 G	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
Carboxypeptidase A	Sigma	108-88-3, 11075-17-5	500 U	Yes	1F125	refrigerator	S	2	3	0		UN 1294/ Class 3/ PG II
casein blocking buffer	ThermoFisher	10049-21-5	1 L	Yes	1F125	refrigerator	L	2	0	0		Not Dangerous goods
CHARCOL, activated 100-400 mesh	Sigma	64365-11-3	1 Kg	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
CHELEX 100 RESIN	BioRad	68954-42-7	500 G	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
CHOLIC ACID, Sodium salt	Sigma	73163-53-8	100 G	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
CITRIC ACID	Fisher	77-92-9	500 G	Yes	1F125	bottle/shelf	S	2	0	0		Not Dangerous goods
CITRIC ACID, anhydrous	Fisher	77-92-9	500 G	Yes	1F125	bottle/shelf	S	2	0	0		Not Dangerous goods
CM Sephadex C50	Sigma	9047-08-9	10 g	Yes	1F125	bottle/shelf	S	0	0	0		not dangerous goods
Congo Red	Sigma	573-58-0	25 g	Yes	1F125	bottle/shelf	S	2	0	0		Not Dangerous goods
COPPER (I) CHLORIDE	Sigma	7758-89-6	1 G	Yes	1F125	bottle/shelf	S	2	0	0		UN 2802/Class 8/ PG III

Cresyl violet acetate	Sigma	10510-54-0	10 Gm	Yes	1F125	bottle/shelf	S	0	0	0	Not Dangerous goods
CYSTAMINE	Sigma	56-17-7	25g	Yes	1F125	bottle/shelf	S	1	0	0	Not Dangerous goods
CYSTINE	Sigma	56-89-3	25 Gm	Yes	1F125	bottle/shelf	S	0	0	0	Not Dangerous goods
DDPE didodecylphosphocholine	Avanti	59752-57-7	25 mg	Yes	1F125	freezer	S	0	0	0	Not Dangerous goods
DEAE-Dextran	Sigma	9015-73-0	50 g	Yes	1F125	bottle/shelf	S	0	0	0	Not Dangerous goods
DEAE-sepharose	Sigma	7732-15-5, 64-17-5	300 mL	Yes	1F125	refrigerator	L	2	2	0	Not Dangerous goods
Decanoic-d19 acid	Sigma	88170-22-3	1 G	Yes	1f125	bottle/shelf	S	2	1	0	Not Dangerous goods
Deoxycholic Acid, sodium salt	Sigma	302-95-4	25 g	Yes	1F125	bottle/shelf	S	2	0	0	Not Dangerous goods
DEUTERIUM OXIDE	Sigma	7789-20-0	10 X 1 mL vial	Yes	1F125	bottle/shelf	L	0	0	0	Not Dangerous goods
Dextran Sulfate Sodium Salt	Sigma	9011-18-1	100 g	Yes	1F125	bottle/shelf	S	0	0	0	Not Dangerous goods
Dithiothreitol	Sigma	3483-12-3	2 g	Yes	1F125	refrigerator	S	2	1	0	UN 3335/Class 9/PG III
DMPC	Avanti	18194-24-6	200 mg	Yes	1F125	freezer	S	0	0	0	Not Dangerous goods
DMPE PEG2000 M	Avanti	474922-82-2	25 mg	Yes	1F125	freezer	S	0	0	0	Not dangerous goods
DMPG	Avanti	200880-40-6	25 mg	Yes	1F125	freezer	S	0	0	0	Not Dangerous goods
DNA ladder FX174/Hae III	Gibco BRL	77-86-1, 60-00-4	2 mL	Yes	1F125	freezer, box	S	0	0	0	Not Dangerous goods
DNA low mass ladder	Gibco BRL	na	2 mL	Yes	1F125	freezer, box	S	0	0	0	Not Dangerous goods
DNA, herring sperm	Sigma	100403-24-5	10 g	Yes	1F125	refrigerator	S	0	0	0	Not Dangerous goods
DOPC, CHCl3 solution	Avanti	67-66-3, 64-17-5, 4235-95-4	1 g	Yes	1F125	refrigerator	L	2	0	0	UN 1888/Class 6.1/PG III
DPPC	Avanti	63-89-8	25 mg	Yes	1F125	freezer	S	0	0	0	Not Dangerous goods
DPPE	Avanti	923-61-5	25 mg	Yes	1F125	freezer	S	0	0	0	Not Dangerous goods
DRIERITE	Fisher	7778-18-9	5 Lb	Yes	1F125	bottle/shelf	S	1	0	0	Not Dangerous goods
DRIERITE, in desiccant cylinder	Fisher	7778-18-9	1 Lb	Yes	1F125	bottle/shelf	S	1	0	0	Not Dangerous goods
DSPC	Avanti	816-94-4	25 mg	Yes	1F125	freezer	S	0	0	0	Not Dangerous goods
DSPE	Avanti	1069-79-0	25 mg	Yes	1F125	freezer	S	0	0	0	Not Dangerous goods
DSPE-PEG(2000) Amine	Avanti		10 mg	Yes	1F125	freezer	S				
DSPE-PEG(2000) Biotin	Avanti	385437-57-0	11 mg	Yes	1F125	freezer	S	0	0	0	Not Dangerous goods
DSPE-PEG(2000) Carboxylic Acid	Avanti	474922-20-8	10 mg	Yes	1F125	freezer	S	0	0	0	Not Dangerous goods
Dulbecco's modified Eagle's media, high glucose	Sigma	na	100 g	Yes	1F125	refrigerator	S/box	0	0	0	Not Dangerous goods
Dulbecco's modified Eagle's media, low glucose	Sigma	na	100 g	Yes	1F125	refrigerator	S/box	0	0	0	Not Dangerous goods
Dynabeads M270 COOH	Invitrogen	7732-18-5	2 mL	Yes	1F125	refrigerator	L	0	0	0	Not Dangerous goods
EDTA	Sigma	60-00-4	100 g	Yes	1F125	bottle/shelf	S	2	0	0	UN 3077/Class 9/ PG III
EGTA	Sigma	67-42-5	25 g	Yes	1F125	bottle/shelf	S	0	0	0	Not Dangerous goods
ETHANOLAMINE	Sigma	141-43-5	25 mL	Yes	1F125	bottle/shelf	L	3	2	0	UN2491/ Class 8/ PG III
Ferrocenium hexafluorophosphate	Sigma	11077-24-0	5 g	Yes	1F125	bottle/shelf	S	0	0	0	UN 3261/Class 8/ PG II
FICOLL PM 400	Sigma	26873-85-8	25 Gm	Yes	1F125	bottle/shelf	S	0	0	0	Not Dangerous goods
FICOLL PM 70	Sigma	72146-89-5	25 Gm x 2	Yes	1F125	bottle/shelf	S	0	0	0	Not Dangerous goods
Florisil	Sigma	1343-88-0	250 G	Yes	1F125	bottle/shelf	S	0	0	0	Not Dangerous goods
fluorescent biotin quant. Kit	ThermoFisher	10049-21-5, 7647-14-5, 7558-80-7	box	Yes	1F125	refrigerator	box	2	0	0	Not Dangerous goods
GALACTOSE, D	Sigma	59-23-4	10 G	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
Gelatin	Polysciences	009000708	100 g	Yes	1F125	bottle/shelf	S	0	1	0	? (not dangerous - sigma)
GLASS BEADS, 100 mesh	Analabs	na	1 lb	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
GLASS BEADS, 150 mesh	Sigma	na	500 G	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
Glass Wool	Sigma	65997-17-3	1 lb	Yes	1F125	bottle/shelf	S	2	0	0	Not dangerous goods
GLUCOSE, d	Mallinckrodt	50-99-7	500 G	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
Glutamine, L	Sigma	56-85-9	100 g	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
GLYCINE HYDROCHLORIDE	Sigma	56-40-6	100 g	Yes	1F125	bottle/shelf	S	1	0	0	Not dangerous goods

Go Taq hot start polymerase colorless	Promega	57-50-1, 1185-53-1	25 mL	Yes	1F125	freezer	L	0	0	0		Not dangerous goods
GUANIDINE HYDROCHLORIDE	Sigma	50-01-1	1 Kg	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods
Hematoxylin	SDI Chem	517-28-2	25 g	Yes	1F125	bottle/shelf	S	2	1	0		Not dangerous goods
HEPES BUFFER	Sigma	7365-45-9	500G	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
HEPTADECANOIC ACID	Sigma	506-12-7	5 G	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods
		288-32-4, 60-00-4, 77-86-1, 288-32-4, 10101-97-0										
His-Bind purification kit	Novagen	10101-97-0	box/kit	Yes	1F125	refrigerator	box	3	0	0		UN 3316/Class 9/PG III
HPLC SORBENT	Sigma	71889-02-6	50 G	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
IEF buffer, anode	Novex	7664-38-2	100 mL	Yes	1F125	refrigerator	L	0	0	0		not dangerous goods
IEF buffer, cathode	Novex	56-81-5, 26628-22-8	100 mL	Yes	1F125	refrigerator	L	2	0	1		not dangerous goods
IEF gel	invitrogen	79-06-1	box	Yes	1F125	refrigerator	S	3	0	0		not dangerous goods
IEF Marker 3-10	Serva	na	0.5 mL	Yes	1F125	refrigerator	L	0	0	0		not dangerous goods
IEF sample buffer, pH 3-10	Novex	56-81-5, 26628-22-8	25 mL	Yes	1F125	refrigerator	L	0	0	0		Not dangerous goods
Lactose monohydrate	Sigma	5989-81-1	1 Kg	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
LAURYL ALCOHOL	Sigma	112-53-8	100 G	Yes	1F125	bottle/shelf	S	2	1	0		Not dangerous goods
Lauryl Sulfate (see Sodium Dodecyl Sulfate)	Sigma	151-21-3	1 Kg	Yes	1F125	bottle/shelf	S	2	3	3		UN 2926/Class 4.1/PG II
L-Broth Agar	Bio 101	7647-14-5	m,	Yes	1f125	bottle/shelf	S	1	0	0		Not dangerous goods
LIPOPHILIC SePHADEX LH-20	Sigma	9041-37-6	100G	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
Lithium Carbonate	Sigma	554-13-2	100 G	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods
Lithium Chloride	Sigma	7447-41-8	25 mg	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods
Luminol	Biochimica	208-309-4	1 Gm	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods
LYSINE	Sigma	657-27-2	100 Gm	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
MagnaBind amine beads	ThermoFisher	1317-61-9	50 mL	Yes	1f125	refrigerator	L	1	0	0		not dangerous goods
MagnaBind COOH beads	ThermoFisher	1317-61-9	10 mL	Yes	1f125	refrigerator	L	0	0	0		not dangerous goods
MagnaBind go anti mo IgG	ThermoFisher	1317-61-9	50 mL	Yes	1F125	refrigerator	L	1	0	0		Not dangerous goods
MagnaBind go anti RA IgG	ThermoFisher	1317-61-9	50 mL	Yes	1F125	refrigerator	L	1	0	0		Not dangerous goods
MAGNESIUM CHLORIDE HEXAHYDRATE	Sigma	7791-18-8	500g	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
MAGNESIUM SULFATE	Sigma	7487-88-9	500g	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
		1309-38-2;										
Magnosphere MK230 COOH	JSR/Cosmo	7732-18-5	10 mL	Yes	1F125	refrigerator	L	1	0	0		Not Dangerous goods
MALONAMIDE	Aldrich	108-13-4	5 Gm	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
MES	Sigma	4432-31-9	25 Gm	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
Methylene blue	Sigma	7220-79-3	25 Gm	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods
Mineral oil	Sigma	8042-47-5	1 L	Yes	1F125	bottle/shelf	L	0	1	0		Not dangerous goods
		12401-88-4; 1344-28-1; 14808-60-7										
MOLECULAR SIEVES	Sigma	14808-60-7	250 Gm	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
MOLYBDIC ACID	Sigma	7782-91-4	500g	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
MOPS	Sigma	7119-22-7	100 G	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
Mowiol® 40-88	Sigma	9002-89-5	100 G	Yes	1F125	bottle/shelf	S	0	1	0		Not Dangerous goods
neutravidin	ThermoFisher	na	10 mg	Yes	1F125	refrigerator	S	0	0	0		Not Dangerous goods
Neutravidin agarose resin	ThermoFisher	na	50 mg	Yes	1F125	refrigerator	S	0	0	0		Not Dangerous goods
Ninhydrin	Sigma	485-47-2	10 G	Yes	1F125	bottle/shelf	S	2	0	0		not dangerous goods
N-Lauroylsarcosine sodium salt	Sigma	137-16-6	100 g	Yes	1F125	bottle/shelf	S	4	1	0		UN 2811/ Class/ 6.1 PG: II
		26628-22-8, 79-06-1										
NuPAGE Bis-Tris protein gels	Novex	26628-22-8, 79-06-1	2 box	Yes	1F125	refrigerator	S	0	0	0		not dangerous goods
Octyl-B-D-n-glucopyranoside	Sigma	29836-26-8	1.1 g	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
ORCINITOL MONOHYDRATE	Aldrich	6153-39-5	25 gM	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods
PALMITIC ACID	Sigma	57-10-3	5 Gm	Yes	1F125	bottle/shelf	S	2	1	0		Not dangerous goods

Particle size standards	ThermoFisher	9003-53-6, 9003-70-7, 7732-18-5	10 x 25 mL	Yes	1F125	refrigerator	L	1	0	0		Not dangerous goods
Peptide calibration standard II	Bruker	222570	1 mg	Yes	1F125	Freezer	S	0	0	0		not regulated as dangerous goods
PHENOL RED	Sigma	34487-61-1	5g	Yes	1F125	bottle/shelf	S	2	0	1		Not dangerous goods
PHENYLALANINE	Sigma	63-91-2	25 Gm	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
phosphocholine (PC)	Sigma	8002-43-5	100 mg	Yes	1F125	freezer	S	0	0	0		Not Dangerous goods
PIPES BUFFER	Sigma	76836-02-7	25 Gm	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
Poly(vinyl alcohol)	Sigma	9002-89-5	250 g	Yes	1F125	bottle/shelf	S	0	1	0		Not dangerous goods
POLYETHYLENE GLYCOL 8000	Sigma	25322-68-3	1 Kg	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
Polyvinylpyrrolidone	Sigma	9003-39-8	100 g	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
POTASSIUM ACETATE	Sigma	127-08-2	500 Gm	Yes	1F125	bottle/shelf	S	1	0	0		Not dangerous goods
Potassium Bromide	Fisher	7758-02-3	25 g	Yes	1F125	bottle/shelf	S	1	0	0		Not Dangerous goods
POTASSIUM CARBONATE	Sigma	584-08-7	500g	Yes	1F125	bottle/shelf	S	2	0	0		Not Dangerous goods
POTASSIUM CHLORIDE	Sigma	7447-40-7	500 Gm	Yes	1F125	bottle/shelf	S	1	0	0		Not dangerous goods
Potassium dichromate	Sigma	778-50-9	50 g	Yes	1F125	bottle/shelf	S	4	0	3	ox	UN 3086/Class 6.1/PG I
Potassium Ferricyanide (III)	Sigma	13746-68-2	5 g	Yes	1F125	bottle/shelf	S	0	0	0		UN3077/ Class: 9/ PG III
Potassium Iodide	Fisher	7681-11-0	500 g	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods
Potassium Metabisulfite	Sigma	16731-55-8	500 g	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods
POTASSIUM PHOSPHATE, monobasic	Sigma	7778-77-0	1 Kg	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
pre-coated plates	Thermo	na	box	Yes	1F125	refrigerator	box	0	0	0		Not dangerous goods
Protein A/G agarose	Santa Cruz	64-17-5	2 mL	tested from	1F125	refrigerator	L	1	0	0		Not dangerous goods
PROTEIN ASSAY RGT S	BioRad	na	5 mL	Yes	1F125	bottle/shelf	L	1	0	0		
PROTEIN ASSAY dye reagent concentrate	BioRad	7664-38-2, 7732-18-5, 67-56-1	500 mL	Yes	1F125	refrigerator	L	3	3	0		UN2924/Class3+8/PG III
PROTEIN ASSAY RGT A	BioRad	na	250mL	Yes	1F125	bottle/shelf	L	1	0	0		
PROTEIN ASSAY RGT B	BioRad	na	150 mL	Yes	1F125	bottle/shelf	L	1	0	0		
Q sepharose fast flow	Sigma	64-17-5	50 mL	Yes	1F125	refrigerator	L	2	2	0		UN1993/Class CBL/PG III
Quercetin	sigma	117-39-5	10 g	Yes	1F125	bottle/shelf	S	2	0	0		UN2281/Class 6.1/PG III
RIBOFLAVIN, Vitamin B12	Sigma	83-88-5	50 G	Yes	1F125	bottle/shelf	S	1	0	0		
RNALATER	Ambion	N/A	500 ml	Yes	1F125	bottle/shelf	L	0	0	0		Not Dangerous goods
SAPONIN	C & B	8047-15-2	0.25 Lb	Yes	1F125	bottle/shelf	S	2	0	0		Not Dangerous goods
SDSOut precipitation reagent	Pierce	16788-57-1, 7778-77-0	kit	Yes	1F125	refrigerator	S	2	0	0		Not Dangerous goods
SEBACIC ACID	Sigma	111-20-6	100g	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
Select Agar	Gibco	Na	500 g	Yes	1F125	bottle/shelf	S	0	0	0		Not Dangerous goods
SEPHADEX LH-20	Sigma	9041-37-6	50 Gm	Yes	1F125	bottle/shelf	S	0	0	0		not dangerous goods
SEPHADEX G-100	Sigma	9050-94-6	50 Gm	Yes	1F125	bottle/shelf	S	0	0	0		not dangerous goods
SEPHADEX G-25	Sigma	9041-35-4	4 X 50 Gm	Yes	1F125	bottle/shelf	S	0	0	0		not dangerous goods
SEPHADEX G-50	Sigma	9048-71-9	3 X 50 Gm	Yes	1F125	bottle/shelf	S	0	0	0		not dangerous goods
SEPHADEX G-75	Sigma	37224-29-6	3 X 50 Gm	Yes	1F125	bottle/shelf	S	0	0	0		not dangerous goods
sepharose CL-2B	Sigma	64-17-5	100 mL	Yes	1F125	refrigerator	L	2	2	0		UN1993/Class CBL/PG III
sepharose CL-4B	Sigma	64-17-5	600 mL	Yes	1F125	refrigerator	L	2	2	0		UN1993/Class CBL/PG III
SILICA GEL	BioRad	7631-88-9	1 Kg	Yes	1F125	bottle/shelf	S	1	0	0		
SILICA GEL, 60	EM Sciences	7631-88-9	1 Kg	Yes	1F125	bottle/shelf	S	1	0	0		
Slide-a-lyzer	Thermo	na	box	Yes	1F125	refrigerator	S	0	0	0		not dangerous goods
SODIUM ACETATE	Sigma	6131-90-4	1kg	Yes	1F125	bottle/shelf	S	1	0	0		Not dangerous goods
SODIUM BICARBONATE	Mallinkrodt	144-55-8	1kg, 500 g	Yes	1F125	bottle/shelf	S	1	0	0		Not dangerous goods
SODIUM BISULFITE	Aldrich	7631-90-5	100g	Yes	1F125	bottle/shelf	S	2	0	1		not hazardous for transport
SODIUM BORATE	Sigma	1303-96-4	500g	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods
SODIUM CARBONATE	Sigma	497-19-8	1 Kg	Yes	1F125	bottle/shelf	S	1	0	1		Not dangerous goods
SODIUM CHLORIDE	Sigma	7647-14-5	1500 G	Yes	1F125	bottle/shelf	S	1	0	0		Not dangerous goods
SODIUM CITRATE	Sigma	6132-04-3	500g/1kg	Yes	1F125	bottle/shelf	S	0	0	0		Not dangerous goods
Sodium Cynaoborohydride	Sigma	25895-60-7	25 g	Yes	1F125	bottle/shelf	S	4	3	3		UN 3179/ Class 4.1/PG II
Sodium Deoxycholate	Sigma	302-95-4	10 g	Yes	1F125	bottle/shelf	S	2	0	0		Not dangerous goods

Sodium Dodecyl sulfate	Sigma	151-21-3	1 Kg	Yes	1F125	bottle/shelf	S	2	3	3	UN 2926/Class 4.1/PG II
Sodium Dodecyl Sulfate	Sigma	151-21-3	500 g	Yes	1F125	bottle/shelf	S	2	3	3	UN 2926/Class 4.1/PG II
SODIUM PHOSPHATE, dibasic	Sigma	7558-79-4	500 Gm	Yes	1F125	bottle/shelf	S	1	0	1	Not dangerous goods
SODIUM PHOSPHATE, monobasic	Sigma	7758-80-7	500 Gm	Yes	1F125	bottle/shelf	S	1	0	1	Not dangerous goods
Sodium Sulfate decahydrate	Sigma	7727-73-3	100 g	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
SODIUM SULFITE	Sigma	7757-83-7	50Gm/500g	Yes	1F125	bottle/shelf	S	1	0	0	Not dangerous goods
SODIUM THIOSULFATE ANHYDROUS	Mallinkrodt	7772-98-7	1 Lb	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
Solvent Blue 38 Dye	Sigma	1328-51-4	25 g	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
Spectral calibration kit I and II	Applied Biosystems	56-81-5, 26628-22-8	2 boxes	Yes	1F125	freezer	S	0	0	0	Not dangerous goods
SSPE Buffer	Sigma	7558-79-4	100 Gm	Yes	1F125	bottle/shelf	S	1	0	0	Not dangerous goods
STEAROYLAMINE	Sigma	124-30-1	25 Gm	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
SUCROSE	Sigma, Fisher	57-50-1	250 G X 2/500G	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
SUDAN BLACK B DYE	Fisher	4197-25-5	10 g	Yes	1F125	bottle/shelf	S	0	2	0	UN1325/ Clas 4.1/ PG III
Tagman universal master mix II	Applied Biosystems	56-81-5, 26628-22-8, 77-86-1	25 mL	Yes	1F125	refrigerator	L	1	0	0	Not dangerous goods
TAUROCHOLIC ACID	Sigma	345909-26-4	1 Gm	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
TCEP HCl	Thermo	51805-45-9	1 Gm	Yes	1F125	bottle/shelf	S	2	0	0	UN3261/ Class 8/ PG II
THAM BUFFER	Fisher	77-86-1	1 Kg	Yes	1F125	bottle/shelf	S	2	0	0	Not dangerous goods
Thioflavin-T	Sigma	2390-54-7	25 g	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
Trichloroacetic acid	Sigma	T6399	100 G	Yes	1F125	acid cabinet	S	3	1	0	UN1839/Class 8/PG II
TRIS BASE	Sigma	77-86-1	1 Kg	Yes	1F125	bottle/shelf	S	2	0	0	Not dangerous goods
TRYPTOPHAN	Sigma	73-22-3	25 g	p	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
TYROSINE	Sigma	60-18-4	25 g	Yes	1F125	bottle/shelf	S	0	0	0	Not dangerous goods
UREA	Sigma	57-13-6	500g	Yes	1F125	bottle/shelf	S	2	0	0	Not dangerous goods
ZINC SULFATE	Fisher	7446-20-0	500g	Yes	1F125	bottle/shelf	S	2	0	0	UN3077/ class 9/ PG III
ZINC SULFIDE	Fisher	1314-98-3	25 g	Yes	1F125	bottle/shelf	S	1	0	0	Not dangerous goods
AG 50 WX8	BioRad, supelco	69011-20-7	2 X100 G	Yes	1F-125	bottle/shelf	S	0	0	0	Not Dangerous goods
Alcian Blue	United States Biochemical Corporation	12040-44-7	25 g	Yes	1F-125	bottle/shelf	S	0	0	0	Not Dangerous goods
Alizeran-3-methylaminodiacetic acid	Sigma	3952-78-1	25 g	Yes	1F-125	bottle/shelf	S	0	0	0	Not Dangerous goods
Bromophenol blue	Sigma	62625-28-9	10 g	Yes	1F-125	bottle/shelf	S	1	1	0	Not Dangerous goods
Calcein	Biochemica	1461-15-0	500 mg	Yes	1F-125	bottle/shelf	S	0	0	0	Not Dangerous goods
Carbowax PEG 8000	Fisher	25322-68-3	500 g	Yes	1F-125	bottle/shelf	S	1	1	0	Not Dangerous goods
DEAE-Sephadex A-25 chloride	Sigma	39455-71-7	50 g	Yes	1F-125	bottle/shelf	S	0	0	0	Not Dangerous goods
Dimedone (5,5-Dimethyl-1,3-cyclohexanedione)	Acros	128-81-8	25 g	Yes	1F-125	bottle/shelf	S	0	0	0	Not Dangerous goods
Erioglaucine, Disodium salt	Sigma	3844-45-9	25 g	Yes	1F-125	bottle/shelf	S	0	0	0	Not Dangerous goods
L-Arginine	Sigma	1119-34-2	100 g	Yes	1F-125	bottle/shelf	S	0	0	0	Not Dangerous goods
Lysozyme	Sigma	12650-88-3	1 g	Yes	1F125, 1F137	freezer	S	0	0	0	not dangerous goods
pepsin	Sigma	9001-75-6	25 g	Yes	1F125, 1F137	freezer	S	3	0	0	Not dangerous goods
AMMONIUM BICARBONATE	Sigma	1066-33-7	500 G	Yes	1F125/129	bottle/shelf	S	1	0	0	UN 3077/class 9/ PG III
Agarose, ultrapure	Invitrogen	9012-36-6	200 g	Yes	1F127	bottle/shelf	S	0	0	0	Not Dangerous Goods
Buffer solution pH 10	Fisher	7732-18-5, 1332-57-6, 144-55-8	500 mL	Yes	1F127	bottle/shelf	L	1	0	0	Not Dangerous goods
Buffer solution pH 4	Fisher	7732-18-5, 877-24-7, 50-00-0, 67-56-1	500 mL	Yes	1F127	bottle/shelf	L	1	0	0	Not Dangerous goods

Buffer solution pH 7	Fisher	7732-18-5, 7558-79-4, 7778-77-0	500 mL	Yes	1F127	bottle/shelf	L	1	0	0	Not Dangerous goods
Coomassie blue	Bio-Rad	6104-58-1	10 g	Yes	1F127	bottle/shelf	S	0	0	0	Not Dangerous goods
crystal violet	Sigma	548-62-9	25 g	Yes	1F127	bottle/shelf	S	2	0	0	UN 3077/Class 9/ PG III
DNeasy blood & tissue kit	Qiagen	50-01-1, 151- 21-3, 39450-01- 6	500 g	Yes	1F127	box/shelf	S/L	2	0	0	Not Dangerous goods
Dulbecco's PBS (10X)	Cambrex	7558-79-4, 7558-80-7, 7647-14-5	3 L	Yes	1F127	bottle/shelf	L	0	0	0	Not Dangerous goods
Electrode storage solution	Ross/Fisher	877-24-7, 7447- 40-7, 7732-18- 5	1 L	Yes	1F127	bottle/shelf	L	2	0	0	Not Dangerous goods
Ficoll	Fisher	72146-89-5	100 g	Yes	1F127	bottle/shelf	S	0	0	0	Not Dangerous goods
Glycerol	Fisher	56-81-5	500 mL	Yes	1F127	bottle/shelf	L	0	1	0	Not dangerous goods
Heat transfer fluid (Cryocool)	Dow/Fisher	000107-51-7, 000141-62-8, 000141-63-9, 083148-62-9	500 mL	Yes	1F127	bottle/shelf	L	1	2	0	Not dangerous goods
methyl beta cyclodextrin	Sigma	128446-36-8	1 g	Yes	1F127	bottle/shelf	S	0	1	0	Not dangerous goods
Nonidet P-40	Sigma	9036-19-5	50 ml	Yes	1F127	bottle/shelf	L	0	1	0	not dangerous goods
PHOPHATE BUFFERed saline PACKS + TWEEN	Sigma	7447-40-7	500 g	Yes	1F127	bottle/shelf	S	1	0	0	Not dangerous goods
PHOSPHATE BUFFERed saline PACKS	Sigma	7447-40-7	500 g	Yes	1F127	bottle/shelf	S	1	0	0	Not dangerous goods
Polyethyleneglycol Bisphenol A Epichlorohydrin Copolymer	Sigma	42617-82-3	500 g	Yes	1F127	bottle/shelf	S	0	0	0	Not dangerous goods
RNASE Away	MBP	7732-18-5	250 ml	Yes	1F127	bottle/shelf	L	0	0	0	Not Dangerous goods
Saline Sodium Citrate (20X)	Fisher	6132-04-3, 7647-14-5, 7732-18-5	1 Lb	Yes	1F127	bottle/shelf	L	1	0	0	Not Dangerous goods
Sodium acetate, 3 M solution	Q Biol	127-09-3	100 mL	Yes	1F127	bottle/shelf	L	1	1	0	non hazardous for transport
Sodium azide	Sigma	26628-22-8	25 g	Yes	1F127	bottle/shelf	S	4	0	0	UN 1687/Class 6.1/PG: II
Sodium dodecyl sulfate, 10% solution	Q Biol	151-21-3	500 mL	Yes	1F127	bottle/shelf	L	1	1	0	UN 2926/Class 4.1/PG II
Tris buffer	Sigma	77-86-1	1.8 L	Yes	1F127	bottle/shelf	L	2	0	0	Not dangerous goods
Tris-EDTA, 10X solution	Research Genetics	1185-53-1, 60- 00-4	1 L	Yes	1F127	bottle/shelf	L	2	0	0	Not dangerous goods
Triton X-100	Sigma	9002-93-1	100 mL	Yes	1F127	bottle/shelf	L	2	1	0	UN 3082/ Class 9/ PG: III
Trypan Blue	Mallingkrodt	72-57-1	25 g	Yes	1F127	bottle/shelf	S	1	2	1	Not dangerous goods
Tween-20	Sigma	9005-64-5	200 ml	Yes	1F127	bottle/shelf	L	1	0	0	
vacuum grease	Beckman	7631-86-9, 63148629, 70131678	100 g	Yes	1F127	bottle/shelf	S	1	1	0	Not dangerous goods
Water	Quality Biologicals and others	7732-18-5	7.5 L	Yes	1F127	bottle/shelf	L	0	0	0	Not dangerous goods
Xylene cyanol FF	Sigma	2650-17-1	10 g	Yes	1F127	bottle/shelf	S	2	0	0	Not dangerous goods
EDTA 0.5 M solution	Research Genetics	7732-18-5, 139- 33-3	1.1 L	Yes	1F127, 1F129	bottle/shelf	L	0	0	0	Not Dangerous goods
Eosin Y solution, alcoholic, with phloxine	Sigma	64-17-5, 67-56- 1	1 L	Yes	1F129	bottle/shelf	L	2	3	0	UN1170/class3/PGII
Fmoc-Ala-OH amino acid	Advanced Chem Tech	35661-39-3	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Arg(Pmc)-OH amino acid	Advanced Chem Tech	119831-72-0	5 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport

Fmoc-Asn(Trt)-OH amino acid	Advanced Chem Tech	132388-59-1	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Asp(Obut)-OH amino acid	Advanced Chem Tech	71989-14-5	5 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Gln(Trt)-OH amino acid	Advanced Chem Tech	132327-80-1	5 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Glu(Obut)-OH amino acid	Advanced Chem Tech	71989-18-9	5 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Gly-OH amino acid	Advanced Chem Tech	29022-11-5	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-His(Trt)-OH amino acid	Advanced Chem Tech	109425-51-6	5 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Ile-OH-amino acid	Advanced Chem Tech	71989-23-8	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Leu-OH amino acid	Advanced Chem Tech	35861-60-0	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Lys(Boc)-OH amino acid	Advanced Chem Tech	71989-26-9	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-lys-OH amino acid	Advanced Chem Tech	105047-45-8	5 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Met-OH amino acid	Advanced Chem Tech	71989-28-1	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Phe-OH amino acid	Advanced Chem Tech	35861-40-8	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Ser(But)-OH amino acid	Advanced Chem Tech	71989-33-8	5 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Thr(But)-OH amino acid	Advanced Chem Tech	71989-35-0	5 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Trp(Boc)-OH	Advanced Chem Tech	143824-78-8	5 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Tyr(But)-OH amino acid	Advanced Chem Tech	71989-38-3	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Fmoc-Val-OH amino acid	Advanced Chem Tech	68858-20-8	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Mayer's hematoxylin solution	Sigma	302-17-0, 7784-26-1	500 mL	Yes	1F129	bottle/shelf	L	2	0	0	Not dangerous goods
N-Acetyl-L-Cysteine	Sigma	616-91-1	25 g	yes	1F129	dessicator	S	0	0	0	Not regulated for Transport
NuPAGE LDS 4X sample buffer	Invitrogen	56-81-5, 2044-56-8	250 mL	Yes	1F129	bottle/shelf	L	0	0	0	Not dangerous goods
NuPAGE MES SDS running buffer	Invitrogen	151-21-3	500 mL	Yes	1F129	bottle/shelf	L	1	0	0	Not dangerous goods
OCT compound	Tissue-Tek	9002-89-5, 25322-68-3	12X118 mL	Yes	1F129	box	L	0	0	0	
Paraffin (Polyfin wax)	Polysciences	009017270, 064742514	1.5 kg	Yes	1F129	bottle/shelf	S	0	1	0	Not dangerous goods
ParaMaX Wang Resin	Advanced Chem Tech	201058-08-4	25 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
Protein Extraction Reagent Type 4	Sigma	57-13-6, 62-56-6, 77-88-1	4 BTL	Yes	1F129	bottle/shelf	L	2	0	0	
Rink Resin HS	Advanced Chem Tech	N/A	5 g	yes	1F129	dessicator	S	1	0	0	Not regulated for Transport
SilverQuest Destainer A	Life Technologies	none given	2X60 mL	Yes	1F129	bottle/shelf	L	0	0	0	Not dangerous goods
SilverQuest Destainer B	Life Technologies	none given	60 mL	Yes	1F129	bottle/shelf	L	0	0	0	Not dangerous goods
SilverQuest Developer	Life Technologies	none given	250 mL	Yes	1F129	bottle/shelf	L	0	0	0	Not dangerous goods

SilverQuest Developer Enhancer	Life Technologies	50-00-0	2 mL	Yes	1F129	bottle/shelf	L	3	2	0		UN 2209/ Class 8/ PG: III
SilverQuest Sensitizer	Life Technologies	4432-31-9, 68-12-2	25 mL	Yes	1F129	bottle/shelf	L	1	0	0		Not dangerous goods
SilverQuest Stainer	Life Technologies	7761-88-8	60 mL	Yes	1F129	bottle/shelf	L	3	0	2		UN 1493/ Class none/ PG: II
SilverQuest Stopper	Life Technologies	none given	250 mL	Yes	1F129	bottle/shelf	L	0	0	0		Not dangerous goods
SimplyBlue SafeStain	Life Technologies	none given	7 L	Yes	1F129	bottle/shelf	L	0	2	0		Not dangerous goods
ammonium phosphate	sigma	218003	100 g	yes	1F-129	bottle/shelf	S	0	0	0		
Methyl cellulose	Sigma	9004-67-5	100 G	Yes	1F-129	bottle/shelf	S	0	0	0		Not dangerous goods
sodium dodecyl sulfate solution	sigma	151-21-3	100 mL	yes	1F-129	bottle/shelf	L	2	3	3		UN 2926/Class 4.1/PG II
1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-(carboxyfluorescein) (ammonium salt)	Avanti Polar Lipids	799268-49-8	10 mg	Yes	1F137	Freezer	S	0	0	0		Not dangerous goods
1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-(dodecylamine)	Avant Polar Lipids	474944-10-0	10 mg	Yes	1F137	Freezer	S	0	0	0		Not dangerous goods
1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-(lissamine rhodamine B sulfonyl) (ammonium salt)	Avanti Polar Lipids	384833-00-5	2 x 10 mg	Yes	1F137	freezer	S	0	0	0		Not dangerous goods
1,2-dipalmitoyl-sn-glycero-3-phosphoethanolamine-N-[4-(p-maleimidophenyl)butyramide] (sodium salt)	Avanti Polar Lipids	384835-50-1	10 mg	Yes	1F137	freezer	S	0	0	0		Not dangerous goods
1,2-distearoyl-sn-glycero-3-phospho-L-serine (sodium salt)	Avant Polar Lipids	321595-13-5	10 mg	Yes	1F137	Freezer	S	0	0	0		Not dangerous goods
Actin	Sigma	51005-14-2	5 mg	Yes	1F137	freezer	S	0	0	0		Not Dangerous Goods
Alcohol dehydrogenase	Sigma	9031-72-5	20 U	Yes	1F137	freezer	S	0	0	0		Not Dangerous Goods
Amyloglucosidase	Sigma	9032-08-0	25 mg	Yes	1F137	freezer	S	0	0	0		Not Dangerous goods
Carbonic Anhydrase	Sigma	9001-03-0	100 mg	Yes	1F137	freezer	S	0	0	0		Not Dangerous goods
Cytochrome C	Sigma	9007-43-6	50 mg	Yes	1F137	freezer	S	0	0	0		Not Dangerous goods
DSPE-PEG(2000) Maleimide	Avanti	474922-22-0	10 mg	Yes	1F137	freezer	S	0	0	0		Not Dangerous goods
Fetal bovine serum	Atlanta	N/A	6 L	Yes	1F137	freezer	L	0	0	0		Not Dangerous goods
Hemoglobin	Sigma	9008-02-0	1 g	Yes	1F137	freezer	S	0	0	0		Not dangerous goods
Lectin from Lens culinaris	Sigma	N/A	5 mg	Yes	1F137	freezer	S	0	0	0		not dangerous goods
L-Lactic Dehydrogenase	Sigma	9001-60-9	25 kU	Yes	1F137	freezer	S	0	0	0		not dangerous goods
Melittin	Sigma	37231-28-0	10 mg	Yes	1F137	freezer	S	1	1	0		Not dangerous goods
Myoglobin	Sigma	100684-32-0	1.25 g	Yes	1F137	freezer	S	0	0	0		Not dangerous goods
N-Ethylmaleimide	Sigma	128-53-0	1 g	Yes	1F137	freezer	S	3	2	1		UN2928/Class 6.1(8)/PG II
Papain, from papsys latex	Sigma	9001-73-4	50 mg	Yes	1F137	freezer	S	2	0	0		Not dangerous goods
penicillin/streptomycin solution	Gibco BRL	3810-74-0, 69-57-8	1 L	Yes	1F137	freezer	box	1	0	0		Not dangerous goods
polylactide-co-glycolide	Sigma	26780-50-7	10 g	Yes	1F137	freezer	S	0	0	0		Not Dangerous goods
Poly-L-aspartic acid	Sigma	31871-95-1	100 mg	Yes	1F137	freezer	S	0	0	0		Not dangerous goods
Poly-L-lysine hydrobromide	sigma	25988-63-0	25 mg	Yes	1F137	freezer	S	0	0	0		Not dangerous goods
Protease from Streptomyces griseus	Sigma	9036-08-0	1 g	Yes	1F137	freezer	S	2	0	0		Not dangerous goods
Protease Inhibitor Cocktail	Calbiochem	6381-92-6	2 vials	Yes	1F137	Freezer	L	0	0	0		Not dangerous goods
Protease, from Bacillus licheniformis	Sigma	9014-01-1	50 mg	Yes	1F137	Freezer	S	2	0	0		Not dangerous goods
Proteinase K	Acros	39450-01-6	100 mg	Yes	1F137	freezer	S	1	1	0		not dangerous goods
purified recombinant HIV-1 p24 core Ag	Feldan	N/A	4 mg	Yes	1F137	Freezer	S	0	0	0		Not classified as dangerous goods
Ribonuclease A	Sigma	9001-99-4	100 mg	Yes	1F137	freezer	S	0	0	0		Not Dangerous goods

Transferrin	Sigma	11096-37-0	100 mg	Yes	1F137	freezer	S	0	0	0		Not dangerous goods
	promega, protea, calbiochem, sigma											
trypsin		9002-07-7	1 g	Yes	1F137	freezer	S	2	0	0		Not dangerous goods
trypsin-EDTA solution	Atlanta Biolog	na	1 L	Yes	1F137	freezer	L	0	0	0		Not dangerous goods
α -Chymotrypsin	Sigma	9004-07-3	350 mg	Yes	1F137	freezer	S	3	0	0		Not Dangerous goods
α -Lactalbumin	Sigma	N/A	1 g	Yes	1F137	freezer	S	0	0	0		not dangerous goods
β -Lactoglobulin B	Sigma	9066-45-9	25 mg	Yes	1F137	freezer	S	0	0	0		not dangerous goods