

Item No.	Qty	Description
1	1	<p>FIPS 140-2 Wireless Security The Optima XR220 includes FIPS 140-2 validated cryptography module for secure wireless connection to the hospital network.</p> <ul style="list-style-type: none"><li>• DHCP &amp; static IP address are supported</li><li>• Authentication Methods: EAP TLS &amp; WPA2 PSK are supported</li></ul>
2	1	<p>Wireless Connectivity for Optima XR220amx and Optima XR200amx</p> <p>802.11 a/b/g n-compatible wireless connectivity to hospital network</p> <p>Wi-Fi Certified</p> <p>Compatible with:</p> <ul style="list-style-type: none"><li>• 802.11i, Wi-Fi Protected Access 2 (WPA2), WPA 802.1X</li><li>• AES - TKIP</li><li>• 64-, 128-WEP</li><li>• VPN: IPSec - IKE</li><li>• Management Frame Protection (MFP) EAP Types:<ul style="list-style-type: none"><li>- LEAP</li><li>- LEAP + 128-WEP</li><li>- LEAP + WPA</li><li>- EAP - TLS</li><li>- EAP-TTLS/MSCHAPv2</li><li>- EAP-FAST</li><li>- PEAP-GTC</li><li>- PEAP/MSCHAPV2</li></ul></li></ul>
3	1	<p>Repeat/Reject Analysis for Optima XR220amx/Upgraded Optima</p>

Item No.	Qty	Description
4	2	<p>XR200amx</p> <p>Optima XR200amx Training</p> <p>One 1 day TiP onsite training visit for Optima XR200amx.</p> <p>Includes T&amp;L expenses. Days provided in two customer visits.</p> <p>This training program must be scheduled and completed within 12 months after the date of product delivery.</p>