

# 6154

## Flush-Mounted In-Floor Scale

for Stretchers, Wheelchairs,  
Chairs, and Stand-On Weighing

**In-Floor  
Scales:**

**6154**

[6102](#)

[5152](#)

- Scale is recessed into the floor.
- 1000 lbs/445 kg capacity.
- Heavy duty handrail for patient safety.



### SPECIFICATIONS

#### 6154 Flush-Mounted In-Floor Scale (computer enhanced)

|                     |   |
|---------------------|---|
| Accuracy:           | 1/10 lb-100gm.  |
| Weighing Range:     | 1000 lb-445 kg.   |
| Readout:            | Digital LED display in pounds and kilograms, or kilograms only. Surface or flush mounted. |
| Power Source:       | 110 VAC (220v available) with rechargeable battery back up.                               |
| Platform (deck)     | Anodized aluminum, 35" X 54".   |
| Weighing Mechanism: | Electronic load cells.  |
| Rough-In Frame Kit  | Allows for a uniform, finished-looking installation.                                      |
| Automatic Zero:     | Returns to zero, ready for the next patient.  |
| Reweigh:            | Recomputes weight while patient is on scale.  |
| Weight Recall:      | Recalls last weight.  |

[Need more  
info?](#)

Click [here](#) to request a brochure, demo, prices or to have a Scale-Tronix representative contact you.

Or call:  
1 800 873  
2001

|                      |  |
|----------------------|--|
| Successive Weighing: | Weigh a continuous stream of patients without re-zeroing.      |
| Audible Signal:      | Indicates weight is in memory.                                 |
| Printer :            | (accessory) Built-in, paper-tape, uses standard thermal paper. |
| Computer Capability: | RS-232 output (optional).                                      |

---

The 6154 platform is level with the floor for easy-access, multi-purpose weighing. The patient simply walks, or is wheeled onto the weighing platform and is weighed.

Weight is displayed to an accuracy of 100 gm or 1/10 lb on a recessed or surface-mounted digital readout. It's line-cord powered for continual use, with a rechargeable battery backup providing operation during power outages. The scale requires no zeroing between patients.

The optional heavy-duty handrail provides patient safety and convenience, and is positioned to remain clear of wheelchairs and beds.

The 6154 is a computer enhanced scale. Computer enhancement provides:

- Digital display in either pounds and kilos.
- Automatic zero.
- Reweigh (recomputes weight while patient is on scale).
- Weight recall (recalls last weight).
- Audible signal indicating when weight is in memory.
- Line-cord power with battery backup.
- Flush or surface-mounted readout.
- Weight printer or RS-232 capability.

We're



and  
proud  
of it!



---

[Scale-Tronix Home Page](#) / [Products](#) / [Frequently Asked Questions](#)  
[Company History and Information](#) / [Contact Scale-Tronix](#) / [Info Request Form](#)

Questions or Problems concerning this site can be sent to our [webmaster](#).