

ATTACHMENT #3 Annual Compliance Testing for Ionizing x-ray Equipment Option Year 2

Price-Computed
Tomography (CT)
Protocol Optimization

| MANUFACTERER | MODEL/NAME | LOCATIO N | VA NUM | SITE NUM | Serial Number | Price | |
|--------------------------------|-------------------------|--------------|--------|---------------|---------------|-------|----|
| Battle Creek VA Medical Center | | | | | | | |
| Hologic | Discovery W ver 4.0.2 | 2231-2 | 132206 | | 87371 | \$ | |
| Siemens | C-cam | 2270-2 | 132488 | 400-429726 | C.Cam - 1745 | \$ | |
| | | 2229-2 | | | | | |
| | Multix 475517 | Room 3 | 121388 | 400-276928 | 2757 | \$ | |
| Toshiba | AQUILION PRIME TSX-303a | 2269-2 | 143607 | 30002699 | FOB16x2133 | \$ | \$ |
| | Radrex DRAD-3000E | 2233-2 | 143009 | 30008139 | U8B1662030 | \$ | |
| Philips | Easy Diagnost Elva | 2237-2 | | | | \$ | |
| General Electric | AMX 4 plus | 0017-2 | 105769 | 616966AMX41 | 521167wk7 | | |
| | AMX 4 plus | 150-82 | 118097 | 269966amx4 | 1011791wk9 | \$ | |
| Toshiba MRI | MRI-1504 | 100-82MRI | 132701 | 362189 | U4A1342030 | \$ | |
| LUDLUM MEASUREMENTS | Survey Meter 14c | 2271B-2 | 119790 | | 247219 | \$ | |
| | Survey Meter 14c | 2271B-2 | 126973 | | 267278 | \$ | |
| | Survey Meter 14c | 2271B-2 | 126974 | | 267284 | \$ | |
| Bicron | Radiac Survey 2000 | 2271b-2 | 7933 | | A208M | \$ | |
| GE Healthcare | Infinia xw4600 | 2271-2 | n/a | 269223VABCINF | 18063 | \$ | |

| | | | | | | |
|------------|--------------------------|---------|--------|-------------|------------|----|
| | Captus 3000 | 2271-2 | 116874 | | 2545sb1 | \$ |
| Capentec | Dose Calibrator CRC-25w | 2271-2 | 126525 | | 152070 | \$ |
| Pinestar | Dose Management computer | 2271-2 | 132594 | VHABACNMIS1 | G26TFX1 | \$ |
| ARIBEX INC | NOMAD PRO | 2249A-2 | 119698 | | AX01230812 | \$ |
| Planmeca | IXW65 Intraoral xray | 2247-2 | 117001 | | itha55274 | \$ |
| | | 2242-2 | 117002 | | itha54709 | \$ |
| | | 2243-2 | 117003 | | itha55302 | \$ |
| | | 2246-2 | 117004 | | itha55295 | \$ |
| | | 2241-2 | 117005 | | itha55305 | \$ |
| | | 2245-2 | 117006 | | itha55223 | \$ |
| | | 2256-2 | 127866 | | ITHc82250 | \$ |
| Planmeca | Promax 2d Panoramic | 2240-2 | 118111 | | RPX236533 | \$ |

Wyoming CBOC

| | | | | | | |
|------------|-------------------------------|--------|--------|--------------|------------|----|
| | | 1J129- | | | | |
| GE Medical | C-arm OEC 9900elite | HCC | 142225 | E2-5082 | rm3394 | \$ |
| | | 2H114- | | | | |
| Philips | Easy Diagnost Eleva Release 5 | HCC | | 62634084 | 363 | \$ |
| | | 2H111- | | | | |
| Toshiba | AQUILION PRIME TSX-303A/5L | HCC | 138505 | 30009227 sid | 5LB14X2091 | \$ |
| | | 2H109- | | | | |
| Hologic | DEXA | HCC | | | 200131 | \$ |
| | | 1F128- | | | | |
| Planmeca | Pro X | HCC | 140641 | | IPX012871 | \$ |
| | | 1F127- | | | | |
| | Pro X | HCC | 140642 | | IPX012876 | \$ |
| | | 1F125- | | | | |
| | Pro X | HCC | 140644 | | IPX012872 | \$ |
| | | 1F123- | | | | |
| | Pro X | HCC | 140645 | | IPX012868 | \$ |

| | | | | | |
|------------------|--------|--------|--------------|----|-----------------------------|
| Pro X | 1F121- | | | | |
| | HCC | 140646 | IPX012866 | \$ | |
| | 1F118- | | | | |
| | HCC | 140647 | IPX012889 | \$ | |
| | 1F117- | | | | |
| | HCC | 140648 | IPX012881 | \$ | |
| | 1F114- | | | | |
| | HCC | 140649 | IPX012867 | \$ | |
| Pro X | 1F113- | | | | |
| | HCC | 140650 | IPX012870 | \$ | |
| Pro X | 1F112- | | | | |
| | HCC | 140651 | IPX012866 | \$ | |
| Pro X | 1F111- | | | | |
| | HCC | 140652 | IPX012869 | \$ | |
| Pro X | 1F107- | | | | |
| | HCC | 140653 | IPX012865 | \$ | |
| Promax 2d and 3d | 1F126 | 140483 | TPX358632 | \$ | |
| | | | Total | \$ | \$ |
| | | | | | Aggregate Total Base |
| | | | | | Period ACR Testing |
| | | | | | and CT Optimization |
| | | | | | \$ |

