

Attachment  
Sole Source Synopsis - VA246-17-N-1941

**EA=each, VI=vial, UD=unit dose, KT=kit, CS=case**

<b>Item No.</b>	<b>Description</b>	<b>Est. Qty</b>	<b>Unit</b>
<b>1</b>	99mTc Pertechnetate (Point Flood Source), 500 uCi - 3 mCi HEU	13	EA
<b>2</b>	99mTc Pertechnetate (Linearity Source), 50-200 mCi HEU	5	VI
<b>3</b>	99mTc Pertechnetate, 0 to 25 mCi HEU	25	UD
<b>4</b>	99mTc Sulfur Colloid, 0 to 6 mCi HEU	5	UD
<b>5</b>	99mTc Sulfur Colloid, (filtered) HEU	12	UD
<b>6</b>	99mTc MAA, 0 to 5 mCi HEU	30	UD
<b>7</b>	99mTc MDP, 0 to 25 mCi HEU	01	UD
<b>8</b>	99mTc HDP, 0 to 25 mCi	125	UD
<b>9</b>	99mTc DTPA, 0 to 20 mCi HEU	01	UD
<b>10</b>	99mTc DTPA, 21 to 75 mCi HEU	30	UD

<b>Item No.</b>	<b>Description</b>	<b>Est. Qty</b>	<b>Unit</b>
<b>11</b>	99mTc MAG 3, 0 to 10 mCi HEU	10	UD
<b>12</b>	99mTc Choletec, 0 to 10 mCi HEU	25	UD
<b>13</b>	99mTc Hepatolite, 0 to 10 mCi HEU	01	UD
<b>14</b>	99mTc Neurolite, 0 to 25 mCi –Exametazime HEU	01	UD
<b>15</b>	Indium-111 Capromab (Prostascint)	01	UD
<b>16</b>	Indium-111Pentetrotide (Neuroendocrine)	01	UD
<b>17</b>	67 Gallium Citrate, per mCi	03	EA
<b>18</b>	201Thallium Chloride, per mCi	5	EA
<b>19</b>	123 Iodine Diagnostic Capsule, 100 uCi	01	EA
<b>20</b>	123 Iodine Diagnostic Capsule, 200uCi	15	EA

<b>Item No.</b>	<b>Description</b>	<b>Est. Qty</b>	<b>Unit</b>
<b>21</b>	123 Iodine MIBG 10 mCi	1	EA
<b>22</b>	123 Iodine-WB-Diagnostic Capsule 1-3 mCi	1	EA
<b>23</b>	Hetastarch WBC Set-up Syringe	1	KT
<b>24</b>	WBC Set-up Syringe, 1000 units Heparin	1	EA
<b>25</b>	Ceretec (Tc99m Labeled) WBC HEU	1	UD
<b>26</b>	(In111 Labeled) WBC	1	UD
<b>27</b>	Pyrophosphate, per vial	30	VI
<b>28</b>	Sincalide, (generic Kinevac), 0-2 mcg unit dose	15	EA
<b>29</b>	Sincalide (generic Kinevac), each additional mcg	5	EA
<b>30</b>	RBC Tagging Kit, Ultratag per kit	1	KT

<b>Item No.</b>	<b>Description</b>	<b>Est. Qty</b>	<b>Unit</b>
<b>31</b>	AeroVent Jr-MediNuclear-cat.# AV-100H 20/case	02	CS
<b>32</b>	AeroVent Jr-MediNuclear-cat.# AV-100HM 20/case	02	CS
<b>33</b>	99mTc Cardiolite, 10 mCi HEU	650	UD
<b>34</b>	99m Tc Cardiolite, 30 mCi HEU	650	UD