

**DEPARTMENT OF VETERANS AFFAIRS**

**Justification for Single Source Awards IAW [FAR 13.106-1](#)  
For  
Over Micro-Purchase Threshold but Not Exceeding the SAT (\$150K)**

**Acquisition Plan Action ID: VA770-18-AP-0303**

**1. Contracting Activity:**

Department of Veterans Affairs  
Consolidated Mail Outpatient Pharmacy (CMOP) Acquisitions Division.  
Department of Veterans Affairs  
Network 15 Contracting Office (NCO 15)  
CMOP Division  
3450 S. 4th St Trafficway  
Leavenworth, KS 66048  
**2237# 760-18-2-050-0476**

**2. Brief Description of Supplies/ Services required and the intended use/Estimated Amount:**

One (1) Med-Surg Line Items 1.

Estimated cost: [REDACTED]

Item #	Description/Part Number*	Qty.	Price	Total extended price	Brand name
1	SENSOR, GLUCOSE, CONTINUOUS, DEXCOM G4 PLATINUM (XI43)	312	[REDACTED]	[REDACTED]	DEXCOM

Delivery Date: **2/5/2018**

**3. Unique characteristics that limit availability to only one source, with the reason no other supplies or services can be used:**

Within the VHA brand name items may be prescribed by physicians when in their clinical judgment, the brand name item will provide a better therapeutic outcome for patients. The brand name of the product is in the CMOP print name. Pharmacists legally must fill the prescription as it was written for Line Item 1.

**Line Item 1**

SENSOR, GLUCOSE, CONTINUOUS, DEXCOM G4 PLATINUM (XI43)

Dexcom G4 Platinum Continuous Glucose Monitoring Sensor, 40 to 400mg/dL Glucose Range, Stays Up to 7 Days, Water Resistant, Disposable, Sterile, 4/Pack, Retail

VHAPM Part 813.106 Simplified Acquisition Procedures: Soliciting from a Single Source  
Attachment 1: Single Source Justification for SAP under \$150K

These Dexcom G4 and G5 platinum CGM sensors are made up of an applicator, a sensor pod, and a sensor wire. You remove the applicator after insertion. The sensor pod stays on your skin for the entire sensor session, up to 7 days. The pod is made of plastic and an adhesive patch. The sensor wire is thin and flexible, and inserts just under your skin. It is attached to the sensor pod, and is made of silver and platinum metal with polymer membranes. You discard the sensor at the end of the session.

- With a sensor needle that's more than 1/3 smaller than the next available brand, virtually painless, simple to use, and easy to wear.
- Designed to deliver exceptional sensing and outstanding performance to help drive decisions that empower you.
- Improved accuracy in the low and high glucose ranges means better detection for diabetes decisions.
- The only sensor on the market approved for up to 7 days use for added value and greater convenience for you and your schedule.
- This sensor works with both Dexcom G4 and G5 devices.

**4. Description of market research conducted and results or statement why it was not conducted:**

This is available through multiple sources on the open market but there is only one manufacturer of this insulin pump component per physician prescriptions only Brand Name is acceptable.

NAC, FSS, and GSA ADVANTAGE were searched but no contract items were found but only a WOSB has been found on GSA ADVANTAGE.

**5. Contracting Officer's Certification:** *Purchase is approved in accordance with FAR13.106-1(b). I certify that the foregoing justification is accurate and complete to the best of my knowledge and belief.*

James D. Justus 1030208

Digitally signed by James D. Justus

1030208

Date: 2018.01.29 07:58:26 -06'00'

\_\_\_\_\_  
Name

Title

\_\_\_\_\_  
Date