

Base Year					
Line Item	Description	Qty	Unit	Price	Cost
0001	<p>Must have: one (1) squeeze grip handle that operates all movements into all positions of the stirrups: provides easy intraoperative adjustment without compromising the sterile field. Simply release the handle to secure the leg holder/stirrup in all directions</p> <p>Must accommodate and be certified to manage/hold a minimum of 500 lbs.</p> <p>Must have: extended "FIN" boot design – support of outer leg well beyond the mid-calf level: extended lateral fin is designed to cover and protect the head of the fibula and the peroneal nerve</p> <p>Must have: ball and socket joint design which mimics the hip joint and movements</p> <p>Must have: a lithotomy range (from sterile field) of +84 degrees to -33 degrees</p> <p>Must have: an abduction range +25 to -9 degrees</p> <p>Must have: Lift-Assist technology: allows easy movement of the leg when placing it in the desired position</p> <p>Must have: a floating boot – self -adjusting boot minimizes pressure on the calf when moving the stirrup</p> <p>Must have: lithotomy and length indicators – visual indicators that help ensure precise positioning</p> <p>Must have: degree markings on stirrup – somewhere located on the hydraulic piston of the mechanism of the Lift-assist.</p>	4	EA		
0002	Reusable Pad Set	4	EA		
0003	Easy Lock Blade Clamp US Rail	8	EA		