

## Attachment C

### Brand Name or Equal:

MFR: Audioscan Verifit

### ITEM INFORMATION

ITEM NO.	DESCRIPTION OF SUPPLIES/SERVICES	QTY	UNIT	UNIT PRICE	AMOUNT
0001	Audioscan Verifit 2 Clinical REM/HIT (or EQUAL) Part No: 3547	3.00	EA		
0002	Audiosan NOAH Module (or EQUAL) Part No: VA-201	3.00	EA		

### Salient Characteristics:

#### 1. Essential/significant physical, functional, or performance characteristics.

##### System must have/be:

- Able to perform verify binaurally — both on the ear, and in the test box.
- Positioning hearing devices as-worn for an accurate real-time binaural
- Puttyless coupling options
- Has the ability to watch as paired instruments coordinate volume and programming
- Teletest handset provides visually confirmation streaming between instruments.
- Ability to verify whether wideband benefit is being delivered.
- System has 16kHz acoustic bandwidth providing the ability to see accurate 1/3rd octave results up to 12.5kHz.
- Speechmap include graphical SII display for easy counseling.
- Unified controls for most binaural testing.